



*Invacare*® | **Matrix® Flo-tech™ Solution**

Slimline cushion with unique dual fluid sac system

**matrix**®  
SEATING SERIES

**INVACARE**  
Yes, you can.™

The *Invacare **Matrx Flo-tech Solution*** cushion uses a combination of a slim, shaped foam cushion base and a dual layer of fluid sacs over the ischial and sacral area. The dual fluid sacs and the unique moulded foam shape of the solution increase the support surface area. This delivers a superior level of protection for patients identified as being at Very High Risk of developing pressure ulcers, or for those with already established tissue damage.

## Features and Options



### Ingress resistant cover

Cover is two-way stretch, water-resistant and vapour-permeable. Cover features an anti-slip base with carry handle for easy transportation.



### Available in 80 sizes

The *Invacare **Matrx Flo-tech Solution*** is available in a range of 80 sizes, from 14-24" widths and 13-22" depths. (See prescription chart for details).



### *Invacare **Matrx Flo-tech Deep Solution** option*

Option of a 10 cm (4") "Deep" foam base is also available in 12 sizes. (See prescription chart for details).



### Ultra-fresh™ added to all cushions for antimicrobial and odour protection

Antimicrobial treatment has been added to all **Flo-tech** foam to help make it resistant to the growth of odour causing bacteria. Helps provide longer lasting freshness.

### Slimline contoured foam base with large dual fluid sac covering the entire top surface of the cushion

The contact area between the individual and the foam cushion is protected, as the fluid displaces to redistribute pressure away from bony prominences. The additional layer of fluid is incorporated under the ischial and sacral area to further enhance the pressure reducing qualities of this cushion, making this cushion suitable for those at very high risk.

### Slimline contoured cushion

Some individuals prefer a slimline cushion, finding it more convenient with their daily living activities. The low profile is particularly beneficial for the facilitation of side transfers.

## Technical data


For more comprehensive pre-sales information about this product, including the product's user manual, please visit [www.invacare.co.uk](http://www.invacare.co.uk)

	Various sizes	Height	Total product weight	Max. user weight <sup>1</sup>	Max. cleaning temperature	Warranty	
							Manufactured to European Standards: BS EN 1021-1 Cig, BS EN 1021-2 Match and BS 7176 (Crib 5)
<b>Matrx Flo-tech Solution</b>	A range of 80 sizes	78 mm	3.08 kg	140 kg	80°	3 years	


1. Guide maximum user weight limits are set assuming that the cushion selected is of appropriate width and depth for the intended user/s.

## Colours


Black Smooth



Black wipedown



Black spacer fabric



Baltic Blue

