



AQUATEC®

Transfer aids




Yes, you can.®



Aquatec *Trans*

- Extra large seat surface
- Safe: turntable is only unlocked by the weight of the user
- Seat surface turns and slides easily
- Soft, skin-friendly cover: easy-care, wipes clean



Aquatec *Disk*

- Extra large seat surface
- Seat surface turns very easily
- Soft, skin-friendly cover: easy-care, wipes clean
- Within seconds, the **Disk**, the shower stool **Dot** and an adapter are combined to a shower stool with a swivel seat



Aquatec *Step*

- Stable frame
- Feet with non-slip rubber caps
- Non-slip, removable cover, machine-washable

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

Order numbers	Aquatec <i>Trans</i>	Aquatec <i>Disk</i>	Aquatec <i>Step</i>
Turning and transfer aid, blue	4.03.002		
Turning and transfer aid, grey	4.03.008		
Rotating disk, blue		4.02.002	
Rotating disk, grey		4.02.008	
Step, blue			3.00.002
Dimensions			
Diameter	365 mm	365 mm	-
Length	516 mm	-	410 mm
Width	-	-	310 mm
Height	25 mm	25 mm	165 mm
More info			
Weight	2.2 kg	1.6 kg	1.5 kg
Load capacity	140 kg	140 kg	150 kg
Colours	Blue, grey	Blue, grey	Blue, white

Invacare Ltd
 Pencoed Technology Park - Pencoed
 Bridgend CF35 5AQ - United Kingdom
 Tel: +44 1656 776 222
 E-mail: uk@invacare.com
 Sales Order e-mail: ordersuk@invacare.com
 www.invacare.co.uk

Invacare Ireland Ltd
 Unit 5 - Seatown Business Campus Seatown
 Road - SWORDS County Dublin - Ireland
 Tel. +353 1 8107084
 Fax +353 1 8107085
 E-mail: Ireland@invacare.com
 www.invacare.ie



Yes, you can.®

