

Setting new standards in Bariatric care



With a safe working load up to 60 Stone (385 kgs), the Invacare **Octave** is a robust yet stylish profiling bed solution that assists in the care and comfort of bariatric clients. With a fully electrical, 4-sectioned profiling mattress support as standard, the Octave provides optimal comfort and support to the client, as well as an ideal working environment for the carer. The frameless construction of the mattress platform allows improved side access to the client, and the bed ends can be detached independently of the side rails for quick and easy access to both head and feet.

Features and Options







Support handle

For increased support whilst bed leaving.

Lifting pole

The handle can be adjusted in both height and depth, to offer support whilst re-positioning.

Extendable bedend

The bedend is extendable to offer mattress platform extension of 10 or 20 cm.

Great comfort and stylish design

The **Octave** is designed specifically with the 'plus sized' bariatric client in mind. The bed frame extends to the edges of the bed base for enhanced support and ease of transfer. An easy to use hand control allows the client to operate the bed's height and profiling function. The *Octave's* heavy duty motors enable a smooth and almost silent operation. The *Octave* is a substantial bed, however, its design includes all the wellknown qualities of Invacare's Scandinavian design.

Delivery and storage


To make delivery as easy as possible, the Octave is delivered on a transportation rack. Dismounting and assembly of the **Octave** couldn't be simpler, despite the size of the bed. Furthermore, the transportation rack makes it possible to economise on storage space and makes the **Octave** very easy to move between locations. Ideal for bed rental needs.


Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please visit www.invacare.co.uk.

















Wood: 199 ka

Invacare Octave

Invacare Octave

Max. user

2000 mm

Total height: 1720 - 2120mm

330-520 mm (4-sectioned)

0 - 30°

Steel 180 kg



60 stone 350 kg/

56 stone



Complies with IEC 60601-2-52.

The bed must not be used by patients under 12 years of age, or by patients with body size equivalent to an average 12 year old or smaller.





T: +44 1656 776222 uk@invacare.com www.invacare.co.uk www.invacare.ie

Invacare UK & Ireland Unit 4 Pencoed Technology Park Bridgend CF35 5AQ

