



Up to 220 kg



**Powerful, Comfortable and Reliable:
The *Invacare Comet^{HD}***

The *Invacare Comet^{HD}* is the latest addition to the *Invacare Comet* family, renowned for its safe, fast and enjoyable drive. This rugged outdoor vehicle is ideal for persons weighing between 150 kg to 220 kg (34 stone). Packed with safety features, this powerful scooter has a large contoured extra padded seat together with four individual adjustable suspension units giving an exceptional driving experience. The *Comet^{HD}* has been designed to allow for additional weight capacity and has passed thorough tests in the German independent TÜV test laboratory.



Yes, you can.®



Features

Durable design and power in one

Even for those with a weight of up to 220 kg (34 stone), safety is never a concern with the **Comet^{HD}**. It has been designed specifically to combine stability and safety with power and comfort.

That comfortable feeling

The XL contoured seat has been developed in conjunction with therapists to ensure optimum seating comfort. The seat base and height- and width-adjustable armrests have been reinforced to allow for additional stresses and can easily be used as supports when getting in and out the scooter. Seating comfort can be increased by using the pneumatically adjustable lumbar support.



Unique suspension system

Proven in the professional mountain biking industry, it allows for simple independent adjustment.



Reinforced seat column

Giving you the confidence of stability and safety.



Reinforced armrest

Strengthened key points allow for the additional stress of getting in and out of the **Comet^{HD}**.

Technical data

						
Invacare Comet^{HD} 660 mm	530 mm	440 - 515 mm	610 mm	660 mm	1475 mm	148 kg
					Invacare Limited Pencoed Technology Park - Pencoed - Bridgend CF35 5AQ - United Kingdom Tel: +44 1 656 776 222 Fax: +44 1656 776 220 E-mail: uk@invacare.com Sales Order E-mail: ordersuk@invacare.com www.invacare.co.uk	
Invacare Comet^{HD} 96 kg	220 kg (34 stone)	2 x 75 Ah	3260 mm	100 mm		
					Invacare Ireland Ltd Unit 5 - Seatown Business Campus Seatown Road - SWORDS - County Dublin - Ireland Tel. +353 1 8107084 Fax +353 1 8107085 Email: Ireland@invacare.com www.invacare.ie	
Invacare Comet^{HD} 135 mm	10°	650 W / 1800 W**	50* km	10 km/h		


* 10 km/h
** Peak performance

Colour of shrouds



Slate Grey

Upholstery



Black vinyl



Weight capacity
220 kg - 34 stone

© 2013 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Comet HD - UK - 07/2013 - 1567312



Yes, you can.®