



Safe, solid, secure: the *Invacare Leo*

The *Invacare Leo* is a scooter designed for individuals who value their independence and wish to get out and about unaided. Safety is a key feature of the **Leo**, which does not detract from its stylish and sporty looks. **Leo** offers individuals the freedom and confidence to enjoy their essential daily outings and leisure excursions. With high standards of design and safety engineering, the **Leo** sets a new market standard.

Trusted Solutions, Passionate People





For a safe, powerful drive

The *Invacare Leo* offers a smooth, safe drive in a scooter that is easy to control and manoeuvre. The *Leo's* intelligent lighting package includes indicators that automatically switch off and a rear brake light for safe driving. In addition to a visual warning for the battery level display, there is an audible warning (beep) when the battery needs recharging.

Compact

The *Leo* is a compact mini-scooter that can be easily dismantled for transportation during longer trips or weekend breaks. Easy manoeuvrability for getting around corners in the local shopping centre, museum or leisure complex becomes refreshingly easy. Access to the battery and electronics for maintenance or servicing is much easier with the *Leo* as only one part needs to be removed.

Reliable

Individuals can trust the *Leo* to take them wherever they need to go. From a trip to the shops or a ride to the park, the *Leo* is sturdy and reliable thanks to its exceptionally high quality electronic engineering and power. These elements are protected from corrosion and

water, (splash protection) which safeguard the vital workings of the scooter against dirt and water. A two-step disengaging lever protects the scooter from free-wheeling if inadvertently knocked out of drive mode. The lighting system is cleverly designed to keep functioning even when one light fails.


Features and Accessories



Finite tiller movement
Individuals can easily adjust the tiller to suit their needs.



Brake light
Rear brake lighting warns other people that the scooter is slowing down, thus preventing collisions.



Comfort seating
Swivel functionality and sliding rails ensure optimal seating position.



Two-step disengaging lever
This unique concept prevents the scooter from accidentally freewheeling.

Features and Accessories



Built in anti-splash guards
Protects electronics and the transaxle from dirt and water.



Disassembly into five parts
No tools are required to dismantle *Invacare Leo*



Crutch & cane holder
Keeps your crutch or cane stable during the drive.



Lockable front box
Secure box for personal valuables.



Seatbelt



Large mirror
Can be mounted left or right.




Heavy-duty storage cover
Keeps the scooter clean and dry.



Charging off-board
Simply remove batteries to charge with the off-board charger kit.









Seat suspension module
Adds extra comfort during driving.



Windscreen
Protects the individuals from the wind.

Technical data

						
4 wheels	470 mm	410 mm	460 - 560 mm	475 mm	590 mm	1220 mm

						
4 wheels	59 kg (excl. batteries) 83 kg (incl. batteries)	136 kg	2 x 36 Ah	2620 mm	60 mm	70 mm


				
4 wheels	10°	240 W	38 km	8 km/h

* Seat in maximum fifth rearward locking position from the front.

Colour of shrouds



Upholstery



Invacare Limited
Pencoed Technology Park -
Pencoed - Bridgend CF35 5AQ -
United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS -
County Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2022 Invacare International Sàrl All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Leo - UK - 04/2022 - 1503875



Follow us on Twitter:
[@Think_Mobility](https://twitter.com/Think_Mobility)



Like us on Facebook:
Invacare UK



Follow us on YouTube:
Invacare Ltd - UK



Yes, you can.®