



INVACARE® TOP END®
ELIMINATOR™ OSR
RACING WHEELCHAIR



The **Invacare® Top End® Eliminator™ OSR Racing Wheelchair** is designed for speed. The super-stiff roll-formed horizontal maintube, combined with a relaxed headtube angle for exceptional stability and aerodynamics, make it the choice of champion racers like Ernst van Dyk, Krige Schabort, Kurt Fearnley, Chantal Petitclerc and Edith Hunkeler. So, if you have your sights set on winning, the Eliminator OSR racing wheelchair will help you get to the finish line first.

Features

6061 T6 oval tube custom aluminum frame
Integral wrap-around fenders
4 cage designs: V, Open-V, U or I-cage
Kneeling, amputee or traditional upholstery
700C or 26" high performance rear wheels
Caliper hand brake
Aluminum side panels
Straight axle tube with rear wheel alignment
Aluminum fork
Precision track compensator system

20" front wheel with tubular tire
Rubber-coated 3/8" diameter aluminum rims
Aero bars
Powder coat paint
250 lb. weight capacity



1. I-cage
2. Open-V cage shown with cut-out axle option
3. V cage
4. U cage shown with knee plate
5. Oversized rake for maximum speed and cornering
6. Precision track compensator system



Options

Carbon fiber rear wheels
Corima™ carbon fiber front wheel
Alignment gauge
Wheel bag
Racing suit
Solid seat/kneel plate
Fixed footrest

Cushion
Cordless speedometer
Helmet
18" front wheel

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Top End

www.topendwheelchair.com
Sports and Recreation
Products
4501 63rd Circle North
Pinellas Park, Florida
33781-5914
(727) 522-8677
(800) 532-8677
Fax: (727) 522-1007

©2013 Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and ®. All trademarks are
owned by or licensed to
Invacare Corporation unless
otherwise noted. New York
City Marathon is a registered
trademark of New York Road
Runners Club, Inc. Corima
is a registered trademark
of Corima, S.A. Zipp is
a registered trademark
of Compositel, Inc.
Specifications are subject to
change without notification.
Form No. 01-192
Rev 02/13
130228

