

# Invacare® / Aquatec® | Toilet seat raisers



Invacare **Finesse**

- Toilet seat raiser made of soft, skin-friendly material
- Raises the toilet seat by 90 mm
- Fits on the toilet easily and safely, no assembly required
- Large hygiene recess at the front for easy personal hygiene
- Now available in Blue - Dementia friendly colour



Invacare **Maris**

- Toilet seat raiser
- Very easy to attach to the toilet
- Smooth surface for easy cleaning
- Large hygiene recess (250 x 190 mm)



Aquatec® **90**

- Toilet seat raiser
- Raises the toilet seat by 100 mm
- Attached safely to the toilet with two side clamps
- Large turning knobs for quick attachment and detachment
- Particularly high load capacity: 225 kg (35 stone)



**High weight limit  
225kg (35 stone)**

**3 year warranty**






**AQUATEC®**



**Yes, you can.®**

**Technical data**

Order number	Invacare <i>Finesse</i>	Invacare <i>Maris</i>	Aquatec® <i>90</i>
Toilet seat raiser	<b>1575059 (Single pack)</b> <b>1575061 (Pack of 5)</b>		
Toilet seat raiser with lid			<b>8.01.001</b>
Toilet seat raiser without lid		<b>1522262</b>	<b>8.00.001</b>
Dimensions			
Seat width	385 mm (A)	365 mm (A)	365 mm (A)
Seat depth	430 mm (B)	385 mm (B)	395 mm (B)
Seat height	90 mm (C)	–	–
Seat height front	–	125 mm (C)	105 mm (C)
Seat height rear	–	135 mm (D)	115 mm (D)
Total width	400 mm	380 mm	–
Total width (with lid)	–	–	380 mm (E)
Total width (without lid)	–	–	375 mm
Total depth	435 mm (F)	400 mm (F)	410 mm (F)
Total height	135 mm (G)	185 mm (G)	–
Total height (with lid)	–	–	135 mm
Total height (without lid)	–	–	125 mm (G)
Width of opening	225 mm (H)	190 mm (H)	210 mm (H)
Depth of opening	275 mm (I)	250 mm (I)	270 mm (I)
More info			
Weight	2.2kg	0,9 kg	–
Weight (with lid)	–	–	1,6 kg
Weight (without lid)	–	–	1,1 kg
Load capacity	125 kg (19.5 stone)	120 kg (19 stone)	225 kg! (35.5 stone)
Colour	Blue	White	White
Warranty	2 years	2 years	3 years

 Invacare **H304 - Finesse**

 Invacare **H301 - Maris**

 Aquatec **90**




**3 year warranty**

**Aquatec® 900**

- Securely attached toilet seat raiser
- Individual adjustment to the user at any time, no tools required:
  - three different seat heights: 2, 4 and 6"
  - three different angles, lower at the front
- Fold-way, ergonomic armrests for additional support when sitting down and standing up
- Available with or without armrests.
- Seat can be folded up in all setting positions – ease of cleaning



**Aquatec 90000**

- Toilet seat raiser
- Available in three different seat heights: 2, 4, 6"
- Swing-away, ergonomic armrests for additional support when sitting down and standing up
- Innovative clip system: The lid and toilet seat can be removed from the base; no tools are required


**Technical data**

Order numbers	Aquatec® 900	Aquatec 90000
Toilet seat raiser with armrests	<b>10128-10</b>	<b>8.31.101 (seat height 100 mm)</b> <b>8.31.111 (seat height 60 mm)</b> <b>8.31.121 (seat height 20 mm)</b>
Toilet seat raiser without armrests	<b>10129-10</b>	
Accessories: Adapter for forward-inclined assembly		<b>1470825</b>

Dimensions		
Seat width	360 mm (A)	365 mm (A)
Seat depth	425 mm (B)	390 mm (B)
Seat height 1	60 mm (C)	20 / 60 / 100 mm (C)
Seat height 2	100 mm (D)	–
Seat height 3	150 mm (E)	–
Total width with armrests	590 mm (F)	540 mm (D)
Total width without armrests	430 mm (G)	380 mm (E)
Total depth	490 mm (H)	475 mm (F)
Total height	–	–
Total height with armrests	380 mm (I)	230 / 270 / 310 mm (G)
Total height without armrests	180 mm (J)	30 / 70 / 110 mm (H)
Height, seat to armrests	200 mm (K)	200 mm (I)
Width between the armrests	510 mm (L)	450 mm (J)
Width of the splash guard	205 mm (M)	–
Depth of the splash guard	260 mm (N)	–
Width of opening	225 mm (O)	205 mm (K)
Depth of opening	285 mm (P)	275 mm (L)
Length, fit 1	420 mm (Q)	–
Length, fit 2	440 mm (R)	–

More info		
Weight	–	–
Weight with armrests	5.0 kg	4.4 kg (SH:100 mm) 4 kg (SH:60 mm) 3.6 kg (SH:20 mm)
Weight without armrests	2.7 kg	–
Load capacity	120 kg (19 stone)	150 kg (23 stone)
Colour	White	White
Warranty	3 years	3 years

Aquatec 900



Aquatec 90000

