Invacare | Alu Lite



A lighter way to get around

The *Invacare Alu Lite* is the basic-entry level wheelchair in the Invacare manual range. The lightweight aluminium construction of this transit wheelchair makes it easy to manoeuvre, even for an elderly carer, as well as providing a wheelchair that is easy to fold and transport.





Easy to transport and store

Being lightweight (aluminium chassis), compact and with a foldable backrest, the *Invacare Alu Lite* is very easy to handle, transport and store. Simply remove the footrests, fold the backrest and chassis and store away.

Pushing? No problem!

The overall manoeuvrability of a wheelchair is crucial - which is why the **Alu Lite** is designed to be very easy to push and manoeuvre, both indoors and outdoors. Due to the smaller overall dimensions and general compactness of this model, getting around in tight spaces is never a problem.



Folding backrestReduces overall size when folded for transportation or storage.



Swing away footrest
Footplates fold up as shown, and
footrests swing away, allowing
easier access in or out of the
chair.



Foldable cross frame Lightweight frame folds flat for easier portability.

Technical data



Invacare **Alu Lite** 405 / 455 mm



430 mm



490 mm



.....

450 mm

945 mm

Upholstery



230 mm

1030 mm



340 - 460 mm

760 mm

INVACARE® Limited
Pencoed Technology Park - Pencoed Bridgend CF35 5HZ
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
www.invacare.co.uk

Sales Order E-mail: ordersuk@invacare.com INVACARE® Ireland Ltd Unit 5 - Seatown Business Campus

Seatown Road - SWORDS - County

Dublin Tel. +353 1 8107084 Fax +353 1 8107085 Email: Ireland@invacare.com www.invacare.ie

© 2010 Invacare® International Sarl All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare® reserves the right to alter product specifications without prior consultation.

AluLite - UK - 11/2010 - 1507059





100 kg



520 / 570 mm

10.2 kg

Frame coloui

Invacare Alu Lite



Blue

Black nylon



Blac