

Invacare® | **Ben_{NG}**

NEW



Simple, practical and modular

The Invacare **Ben_{NG}** is the entry level chair in the Invacare manual range. This next generation of the well known **Ben 9 Plus** and **Duo** combines a reliable and proven construction with a modular concept that will allow you to reconfigure the propulsion type within a few minutes. Attendant and occupant propelled versions are available both in standard (127kg) or heavy duty versions for clients weighing up to 160 kg.



Yes, you can.®



Features and Options

Modularity

Invacare has designed the **Ben_{NG}** to make life easier. In addition to the well known key features that have given the **Ben 9** such a good reputation, the next generation proposes a new rear wheel support that allows you to change the propulsion type of your chair from a transit to a self propel, or vice versa. This increases adaptability to meet your individual needs.



Modular rear wheel position

Thanks to the removable rear wheel support you can transform your **Ben_{NG}** from a transit to a self propel, or vice versa.

Added value

The **Ben_{NG}**'s standard features include a swing away legrest and detachable armrests for easier transfer. The lightweight folding frame and backrest reduce the overall size when folded for transportation or storage. A wide range of options are available, from cushions to extended backrests and headrests for added support and comfort. Simply choose the options that suit your individual needs and tastes.



Compactness

The lightweight frame, backrest and legrests fold away easily to make the chair very compact for storage or transportation.



Heavy Duty version

A heavy duty version is available, ideal for people weighing up to 160 kg, in both self-propulsion and transit versions. Available in seat widths of 510, 560 and 610 mm with seat depths of 430 or 450 mm.

Technical data

Ben_{NG}	380 - 455 mm	405 / 430 mm	5°	485 mm	8°	420 mm	220 mm	330 - 430 mm
Ben_{NG} HD	510 - 610 mm	430 / 450 mm	5°	485 mm	8°	420 mm	220 mm	330 - 430 mm

Invacare Ben_{NG}	Seat width + 200 / 130* mm	1040 / 910* mm	780 / 670* mm	965 / 720 (folded) mm	11,5°
				965 / 720 (folded) mm	

Invacare Ben_{NG}	127 / 160** kg	19,3 / 17* / 22** kg

* Transit version

** Heavy Duty version

All product dimensions are in millimetres to apply EN 12183 standards. Our wheelchairs comply with ISO 7176-19 standards and with ISO 10542, and have additional features providing increased levels of occupant security and safety whilst travelling in a motor vehicle. Wheelchair users should however transfer to the vehicle seat and use the vehicle installed restraint system whenever it is feasible.

Invacare Limited
Pencoed Technology Park - Pencoed - Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail: ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2012 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

BenNG - UK - 10/2012 - 1551542

Frame colour



Black

Please note colours may vary slightly from those displayed above.

Upholstery



Black Nylon



Yes, you can.®