

alber

Simply mobile.

viamobil[®] eco

Push and brake aid



Pushing the easy way.



The new easy way to push.

► **Pushing and braking your wheelchair can now be as easy as this.** Uphill or downhill, on steep ramps, sloping drives or on longer journeys, viamobil eco, the electric pushing and braking aid for wheelchairs, helps carers with pushing at the times they need it most.

► **viamobil eco is also a safe partner when travelling downhill, with the automatic braking function providing a big advantage in safety.** The intelligent spring system ensures a comfortable and quiet ride, automatically compensating for bumps or dips in the road surface, making sure that the wheel always remains in contact with the ground.

These benefits not only relieve the strain on the carer, they enhance the personal freedom of carer and user alike and thereby improve the opportunity for social contact.



viamobil eco: with puncture-proof tyres, battery pack, charger and adaptable ErgoDrive handlebar control (all included in the price).

Just simple – the viamobil eco.

Wheelchair pushing can be as easy as this – a fact endorsed by TÜV, an independent organisation which tests quality and reliability. Under the guidance of TÜV, the user friendliness of the viamobil eco was exhaustively tested with the help of users and given the rating 'Simple to operate'.





ErgoDrive – Pushing the easy way.

ErgoDrive means easy, convenient operation. The ergonomic handlebar control makes pushing possible without effort. To move forwards or backwards, simply press the handlebar gently down or up. It's childishly simple and convenient.



Everything under control.

With the single-hand combination control switch, all functions are easy to see and are controlled with a handy, easy grip rotary switch that only needs one hand to operate. Optimal user convenience from the outset.



ErgoDrive: benefits at a glance.

- ▶ New patented handlebar control for guaranteed, fatigue-free and ergonomic wheelchair pushing.
- ▶ Single-hand combination switch for switching on and off, freewheeling and individual setting of speed.
- ▶ Intelligent pneumatic spring system that adapts to bumps and dips in the ground and allows simple lowering and raising of the drive wheels by foot-pedal



► **Technical data**

Load capacity	120 kg (weight of person)	
Speed	Infinitely adjustable 1-6 kph	
Range	10 km *	
Climbing capability	10 %	
Weight	Battery pack	6,0 kg
	Drive	10,5 kg
	Control unit	2,0 kg
Overall	18,5 kg	

* depending on load, terrain, ambient temperature to ISO 7176-4

Weight information V14

max. occupant's weight:	120 kg
max. overall weight:	170 kg
(incl. viamobil eco, wheelchair, wheelchair occupant)	

Ideal to transport.

In a matter of seconds, viamobil eco dismantles into three handy parts, without tools, and easily fits any car boot or luggage compartment.

Fits your wheelchair.

viamobil eco combines with most standard wheelchairs and, thanks to its well designed coupling system, attaches easily to the wheelchair without the need for tools.



INVACARE® Limited
 Pencoed Technology Park,
 Pencoed, Bridgend CF35 5HZ
 Phone +44 (0) 1656 776222
 Fax: +44 (0) 1656 776220
 E-mail: uk@invacare.com
 www.invacare.co.uk
 Sales Order E-mail:
 ordersuk@invacare.com

INVACARE® Ireland Ltd
 Unit 5
 Seatown Business Campus
 Seatown Road
 SWORDS, County Dublin
 Phone (353) 181 070 84
 Fax (353) 181 070 85
 Email: eire@invacare.com
 www.invacare.ie



Our quality claim = Your safety!
 Certified to the latest standards: ISO 9001:2000
 and ISO 13485:2003 for medical products.

www.alber.de