

# INVACARE® TOP END® Recreation Handcycles



Yes, you can.®

# INVACARE® TOP END® Recreation Handcycles




Safety, comfort and performance are the distinguishing characteristics of all our handcycles. The *Invacare Top End Exceleator* and *Invacare Top End XLT Series* are ideal recreational models. Even children are catered for with the *Invacare Top End Li'l Exceleator* and *XLT Jr.*

Featuring a comfortable recumbent design and 27 external gears which allow you to adjust the slightest changes encountered on the road (including variations of wind, inclines and descents), the *Invacare Top End XLT Pro* is the perfect touring and long-distance riding handcycle.


The *Invacare Top End XLT* and *XLT Jr.* models are designed with the recreational handcycling enthusiast in mind. Equipped with low maintenance 7 speed hub with reverse hands-on braking, the *Top End XLT* will easily handle most terrains, whether you're simply riding around town or taking a spin in the country.

For a great way to exercise, cross train or simply have fun, the *Top End Exceleator* will help sports enthusiasts 'excel' at all that they do. Designed to accommodate taller or larger frames, this upright and extremely manoeuvrable three-wheel handcycle is an 'excel-lent' choice!

## OPTIONS AND FEATURES



*Invacare Top End XLT Pro* handcycle with standard V crankset.



The *Invacare Top End XLT Jr.* handcycle is designed for children and small adults and features 20" cruiser wheels.



*Invacare Top End Exceleator* handcycle.



The *Top End Li'l Exceleator* handcycle is designed for children and has 20" cruiser wheels and tyres.



## TECHNICAL DATA



Max. 250 lbs  
(approx. 113 kg)

Invacare Limited  
Pencoed Technology Park - Pencoed  
Bridgend CF35 5AQ - United Kingdom  
Tel: +44 1 656 776 222  
Fax: +44 1656 776 220  
E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
Sales Order E-mail:  
[ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS  
County Dublin - Ireland  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

© 2012 Invacare International Sàrl  
All rights reserved.

All information quoted is believed  
to be correct at time of print.  
Invacare reserves the right to alter product  
specifications without prior consultation.

TopEnd Recreation - UK - 03/2012 - 1524493

Please refer to order form for all  
technical specifications.



Yes, you can.®