



# Invacare Reliant 350

## Product Data Sheet

**Product Code:** RPS350-1E

**Description:** Reliant 350 Stand Assist Lifter, with manual legs and detachable battery

All products are delivered with two batteries and off board charger and lead (please note the battery cannot be charged on the hoist).

### Key Features

- Steel powder coated frame
- SWL - 158kg (25 stone)
- Externally tested verified to ISO 10535
- Manually operated legs
- Padded, height adjustable knee pad
- Large stable footplate
- Multiple attachment points for sling
- Warranty 3 years

### Technical Specification

- 158kg (25 stone) weight limit
- Overall product weight – 48.9kg
- Castors - 76mm front/127mm rear (braked castors at rear)
- Maximum leg opening – 940mm
- Base width closed – 660mm (internal)
- Base width open - 940mm (internal)
- Lift height range – 1020 - 1670mm
- Base length – 900mm
- Number of cycles per full charge – 100-200 (dependent on client weight and range)
- Battery – 24V
- Emergency stop
- Emergency lower – Electrical and Manual

### Accessories:

- Stand Assist and Transfer Stand Assist Sling (Size - Small to Extra Large)  
(compatible with other manufacturers slings - for sling compatibility see bulletin PU661)



## Supporting Documentation

- Reliant 350 sales brochure
- Reliant 350 User manual
- Sling and Hoist Assessment Guidelines
- Sling brochure
- Sling application video ( go to [www.invacare.co.uk](http://www.invacare.co.uk), click on 'Invacare Videos' link at top of home page)

All above documentation can be found on [www.invacare.co.uk](http://www.invacare.co.uk)

Invacare's Business Development Managers can provide in field support and training on Hoists and Slings

---

### Manufacturers Contact Details

Invacare Ltd, Unit 4 Pencoed Technology Park, Bridgend, CF35 5AQ

**Tel:** 01656 776222 **E-mail:** [UK@invacare.com](mailto:UK@invacare.com)

**Web:** [www.invacare.co.uk](http://www.invacare.co.uk) / [www.thinkindependence.co.uk](http://www.thinkindependence.co.uk)



| safepatienthandling