



Invacare Sling Range

Product Data Sheet

High quality range of standard sling designs

Sling types available:

- Universal high
- Universal high plus
- Universal standard
- Universal low
- Dress toileting (High)
- Easy fit
- Comfort & Comfort toileting
- Amputee
- Stand assist
- Transfer Stand Assist
- Standing transfer vest



Key Features

- Range of fabric types – spacer, solid polyester and net
- SWL - 200kg (32 stone) Extra - large slings 250kg (39 stone)
- Sizes from extra small to extra large
- Machine washable to 95 degrees
- Can be tumble dried on low setting
- Warranty 1 year

Net fabric – dries quickly, ideal for bathing transfers

Polyester – smooth and easy to apply

Spacer – multi stretch, highly comfortable fabric that can be left in situ*




Net



Spacer



Polyester



Accessories:

Fleece padding – slides onto leg sections to provide additional protection for thighs
Leg insert – foam piece to stiffen leg section
Extension loops – provides additional strap length
Head support stiffener - additional head support

Supporting Documentation

- Sling sales brochure
- Sling user manual
- Sling and Hoist Assessment Guidelines
- Sling brochure
- Sling application video (go to www.invacare.co.uk, click on 'Invacare Videos' link at top of home page)

All above documentation can be found on www.invacare.co.uk

Invacare's Business Development Managers can provide in field support and training on Hoists and Slings

** A risk assessment must be carried before any sling is left in situ*

Manufacturers Contact Details

Invacare Ltd, Unit 4 Pencoed Technology Park, Bridgend, CF35 5AQ

Tel: 01656 776222 E-mail: UK@invacare.com

Web: www.invacare.co.uk / www.thinkindependence.co.uk



safepatienthandling