

# Invacare Octave Community Bed

## Product Data Sheet

---



### Product Range Available

- Octave

### Standard Features

- Robust bariatric community bed
- Frameless construction for easier manual handling access
- Comply to the Bed standard IEC 60601-2-52 **Third Edition**
- Four sectioned profiling mattress support
- Delivered on transportation rack
- Wooden side rails and side rail covers available
- Lifting pole as standard



**Delivered on transportation rack**

### Technical Specification

- Width outside: 1190mm
- Length outside: 2210mm
- Length inside: 2070mm
- Mattress support dimensions: 810mm (Head section) 340mm (Seat plate) x 330mm (thigh section) x 520mm (lower leg section)
- Backrest Angle: 0-57°
- Thigh Angle: 0-30°
- Full Legtest lift 0-30°
- Maximum user weight 350kg (56 stone)
- Safe Working Load 385kg (60 stone)

### **Optional accessories available including:**

- Support handles
- Side rail covers

### **Supporting Documentation**

- Octave range brochure PN: 1487891
- Bed Accessory brochure\*
- Beds & PAC leaflet PN: 1568402
- Bed rail ruler PN: 1556575

*\*Contact Marketing department*

---

### **Manufacturers Contact Details**

Invacare Ltd, Unit 4 Pencoed Technology Park, Bridgend, CF35 5AQ

**Tel:** 01656 776222 **E-mail:** [UK@invacare.com](mailto:UK@invacare.com)

**Web:** [www.invacare.co.uk](http://www.invacare.co.uk)