

Invacare Aquatec Ocean VIP Shower Chair Commode

Product Data Sheet



Product Range Available

- Ocean VIP
Product code: 1517060
- Ocean VIP XL
Product code: 1517061

All products are delivered with soft seat, insert for soft seat, pan and lid as standard

Cross compatibility on a wide range of accessories available on all products

Standard Features

- Modular shower chair commode
- Tension adjustable back
- Seat tilt can be individually adjusted (tilt angle 0°-35°)
- Stainless steel gas spring
- Extremely stable even at maximum tilt angle
- Stable, powder coated stainless steel frame – 100% corrosion free
- 4 Locking castors
- Angle adjustable, flip up and textured footplates with heel loops
- Height adjustable, swing away and detachable footrests
- Ergonomic swing away armrests, can be adjusted to two different height
- Height and depth adjustable headrest as standard

Technical Specification

- Adjustable seat to ground height 510 – 610mm (no tools required)
- Height, footrest to seat 345 – 465mm (no tools required)

- Armrests height 215 or 245mm
- Seat tilt individually adjusted 0° - 35°
- 3 Year Warranty
- 100% corrosion free frame
- Maximum user weight 150kg (23 stone)

Optional accessories available including:

- Soft seat and Visco soft seat
- Variable multi adjustable soft seat (can be retrofitted – no tools required)
- Solid plastic backrest – clips onto frame
- Hip and chest belts
- Armrest lock – lock in upright and down position
- Lateral supports
- Amputee and calf supports – no tools required

All accessories can be found on the product brochure

Supporting Documentation

- Aquatec Ocean range brochure PN: 1504395
- Ocean A5 quick guide PN: 1575042
- Ocean Poster PN: 1575974
- Ocean Order guide*
- Ocean packs are also available email UK@invacare.com to order your pack/packs
- Invacare's Business Development Managers can provide support with client assessments and provide training for OT's and carers.

**Contact marketing department*

Manufacturers Contact Details

Invacare Ltd, Unit 4 Pencoed Technology Park, Bridgend, CF35 5AQ

Tel: 01656 776222 **E-mail:** UK@invacare.com

Web: www.invacare.co.uk