

Invacare Aquatec Ocean Dual VIP Shower Chair Commode

Product Data Sheet

Product Range Available



- Ocean Dual VIP
Product code: 1517062

All products are delivered with soft seat, insert for soft seat, pan and lid as standard

Cross compatibility on a wide range of accessories available on all products

Standard Features

- Modular shower chair commode
- Tension adjustable back
- Stable, powder coated stainless steel frame – 100% corrosion free
- 4 Locking castors
- Angle adjustable, flip up and textured footplates with heel loops
- Height adjustable, swing away and detachable footrests
- Ergonomic swing away armrests, can be adjusted to two different height
- Seat tilt can be individually adjusted (tilt angle 0°-30°)
- Second operating lever with safety lock to lower backrest (Recline angle 10°-50°)
- Anti-tippers included as standard
- Only “off the shelf” shower chair with seat tilt and backrest recline on the market
- Height and depth adjustable headrest as standard

Technical Specification

- Adjustable seat to ground height 510 – 610mm (no tools required)
- Height, footrest to seat 345 – 465mm (no tools required)
- Armrests height 215 or 245mm

- Seat tilt individually adjusted 0° - 35°
- Backrest recline angle 10° - 50°
- 3 Year Warranty
- 100% corrosion free frame
- Maximum user weight 150kg (23 stone)

Optional accessories available including:

- Sot seat and Visco soft seat
- Variable multi adjustable soft seat (can be retrofitted – no tools required)
- Solid plastic backrest – clips onto frame
- Hip and chest belts
- Armrest lock – lock in upright and down position
- Lateral supports
- Amputee and calf supports – no tools required

All accessories can be found on the product brochure

Supporting Documentation

- Aquatec Ocean range brochure PN: 1504395
- Ocean A5 quick guide PN: 1575042
- Ocean Poster PN: 1575974
- Ocean Order guide*
- Ocean packs are also available email UK@invacare.com to order your pack/packs
- Invacare's Business Development Managers can provide support with client assessments and provide training for OT's and carers.

** Contact marketing department*

Manufacturers Contact Details

Invacare Ltd, Unit 4 Pencoed Technology Park, Bridgend, CF35 5AQ

Tel: 01656 776222 **E-mail:** UK@invacare.com

Web: www.invacare.co.uk

Trusted Solutions, Passionate People

