

Invacare Aquatec Ocean E-VIP Shower Chair Commode Product Data Sheet



Product Range Available

- Ocean E-VIP
Product code: 1518424

All products are delivered with soft seat, insert for soft seat, pan and lid as standard

Cross compatibility on a wide range of accessories available on all products

Standard Features

- Modular shower chair commode
- Tension adjustable back
- Stable, powder coated stainless steel frame – 100% corrosion free
- 4 Locking castors
- Angle adjustable, flip up and textured footplates with heel loops
- Height adjustable, swing away and detachable footrests
- Ergonomic swing away armrests, can be adjusted to two different height
- Seat height can be steplessly adjusted between 500 – 1050mm electrically
- Seat tilt can be electrically tilted forwards and backwards (tilt angle -5° to 35°)
- Clearly laid out hand control for simple adjustment of the seat height and seat tilt (4 buttons only)
- Sealed electronic components: ideal for use in shower
- Detachable battery pack can be recharged easily

Technical Specification

- Seat width 480mm
- Seat depth 450mm
- Armrests height 230 -270mm

- Height, footrest to seat 345 – 465mm (no tools required)
- Adjustable seat to ground height 500 - 1050mm (no tools required)
- 3 Year Warranty
- 100% corrosion free frame
- Maximum user weight 150kg (23 stone)

Optional accessories available including:

- Sot seat and Visco soft seat
- Hip and chest belts
- Armrest lock – lock in upright and down position
- Lateral supports
- Amputee and calf supports – no tools required

Manufacturers Contact Details

Invacare Ltd, Unit 4 Pencoed Technology Park, Bridgend, CF35 5AQ

Tel: 01656 776222 **E-mail:** UK@invacare.com

Web: www.invacare.co.uk