

Invacare® | *iPortal2*®

NEW



***iPortal 2*: media and mobility that keeps you connected**

iPortal2 is our most advanced wheelchair media-connectivity system yet. This pioneering solution allows wheelchair users the ability to connect with personal computers and portable computing devices from their power wheelchair joystick or special input device.

The new *iPortal2* module includes:

- A real-time display for a wheelchair's settings and battery status.
- Bluetooth® pairing to PC and iOS devices.
- Simple switching between devices.



Yes, you can.®

Dashboard App

The **iPortal2** Dashboard Application displays wheelchair information in real-time using Bluetooth® technology to connect with the wheelchair's operating system. It displays the speed range and seat adjustment indicators as well as a large battery status gauge. When connected to the iPhone 3GS (or later versions) users can also access a compass heading.



Features and Options



Accessibility

iPortal2's ground-breaking functionality makes a significant difference to users' lives - accessibility. It's a pioneering solution that gives wheelchair users the ability to control an iPhone, iPod Touch or iPad via their wheelchair joystick or special input device. Apple's Voiceover and Assistive Touch enables wheelchair users to control the Apple device. Surfing the net, making calls, playing music is now easy and accessible without the need to touch the Apple device.



Mouse Mover

With the Mouse Mover feature enabled with a simple Bluetooth pairing, the wheelchair's joystick can be used as a computer mouse pointer, allowing the user to navigate and interact with their PC or laptop.

Invacare International Sàrl
Route de Cité-Ouest 2 - 1196 Gland
Switzerland
Tel: +41 (0)22 354 60 10
Fax: +41 (0)22 354 60 11
hqeurope@invacare.com
www.invacare.eu.com

© 2013 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. iPhone and iPod Touch are trademarks of Apple Inc., registered in the U.S. and other countries. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth® SIG, Inc. and any use of such marks by Apple is under license. Apple is not a participant in or sponsor of this promotion. "Made for iPod," and "Made for iPhone," mean that an electronic accessory has been designed to connect specifically to iPod or iPhone, respectively, and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Please note that the use of this accessory with iPod or iPhone may affect wireless performance.

Compatibility

The **iPortal2** uses standard Bluetooth technology, whether built-in or via a dongle (Bluetooth 2.1 or higher) and is compatible with the iOS (iPod, iPhone, iPad), Mac, Windows and Linux operating systems as well as tablet computers. It is available to Dynamic Controls' Shark, DX & DX2 power wheelchairs controllers.



iPortal2 - EU - 11/2013



Yes, you can.®