

## On Track with EC Track

Modern, discreet and safe - EC Track offers free movement around the home.



**Move securely from one room to another**  
Two rooms can be connected allowing the hoist unit to glide smoothly between two rooms, for example bedroom to bathroom.



**Between hung traverse**  
EC Track has a discreet and stylish finish.



**Side hung**  
Side hung tracks can be mounted very close to the wall.

Invacare EC Track is part of the 'Safe Patient Handling' range.




safepatienthandling

Invacare Limited  
Pencoed Technology Park -  
Pencoed - Bridgend CF35 5AQ  
- United Kingdom  
Tel: +44 1656 776 222  
Fax: +44 1656 776 220  
E-mail: uk@invacare.com  
Sales Order E-mail:  
ordersuk@invacare.com  
www.invacare.co.uk

Invacare Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS -  
County Dublin - Ireland  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: Ireland@invacare.com  
www.invacare.ie

EC TRACK. UK:10/16 - @Invacare International Srl 2016. All rights reserved. Patents pending - ASPHYTE 334 860 285 R.C.S. TOURS



**EC Track system with traverse track under hung**  
 This configuration allows for more space between the track and the ceiling, which is an advantage if lighting fixtures or other elements are mounted on the ceiling.



**EC Track system with end hung wall mounted parallel tracks**  
 Where it is not possible to mount the parallel tracks onto the ceiling, the track or tracks can be directly mounted to the wall.

**EC Track brings two rooms together**  
 In order to connect two or more rooms with the same motor unit, a transit coupling can be mounted.



**Electronically operated for seamless transition**  
 The switch for the transit coupling is electronically operated. For optimal security the transit coupling will be locked before transferring from one room to the other.