



# NO SUCH THING AS UNBEATABLE CHALLENGES



## **STRONG**

The sturdy, yet flexible chassis, efficient shock absorption and advanced four-wheel drive form a solid base for smooth movements even on challenging surfaces and areas.

## **FLEXIBLE**

The Four X is comprehensively equipped with additional features such as seat tilt, balance adjustment, legrest and armrest adjustments to enable an individual to lead as active a lifestyle as possible. Upgrade to lift or stand support to get an even more tailored experience.

## STURDY

Every Chasswheel is manufactured using quality parts and assembled by highly trained experts. Four X comes with the industry's most comprehensive 10 year warranty.

# **GOING PLACES**

**FOUR** 

Chasswheel produce top of the range all terrain power chairs for individuals who choose an active lifestyle. Its innovative design features 4-wheel suspension and flexible chassis that makes it ideal for challenging terrain. With Chasswheel, going off the beaten track provides an exciting, all round driving experience.









packs

system

## **STANDARD FEATURES**

- Patented flexible chassis
- Tilt
- Gravity point adjustment
- LED indicators
- LED driving lights
- LED rear lights
- Headrest
- Seatbelt
- Adjustable armrest
- Adjustable legrest
- Recaro Specialist seats







## **TECHNICAL SPECIFICATIONS**

	FOUR X URBAN	FOUR X
Width:	620 mm	695 mm
Height (normal seating position w/out headrest):	1000 mm	700-850 mm
Length:	920 mm (1170mm footrest down)	1250 mm (footrest down)
Weight with battery:	135 Kg	140 Kg
	With lift 150 Kg	Four X DL Lift 170 Kg
		Four X DL SSS 170 Kg
Maximum speed:	10km/h	10 Km/h
Drive distance:	~40km	~40km
Batteries:	2 x 73Ah, gel	2 x 73Ah, gel
Motors:	4 x 250W	4 x 250W
Max user weight:	125 kg	125 Kg (160Kg with HD kit)
Seat width:	400-480 mm	400-480 mm
Seat depth:	390-510 mm	390-510 mm
Backrest height:	500-600 mm	500-600 mm
Ground clearance:	6" tyres 70mm	6" tyres 100mm
		8" tyres 120 mm
Slope climbing capacity:	20°	20°
Maximum obstacle climbing ability:	150 mm	150 mm
Seating height without cushions;	6" tyres - 500mm	6" tyres - Four X DL 440 mm
	With lift - 590 mm	Four X Lift - 530 mm
		Four X SSS - 530 mm
		8" tyres - Four X DL 465 mm
		Four X Lift - 555 mm
		Four X SSS - 555 mm
Gravity point adjustment:	200 mm	250 mm
Tilt angle:	30°	30°
Footrest angle:	80°	80°
Backrest angle:	84°-125°	84°-125°
Turning radius:	1.15m	1.15m
Electronics:	Dynamic DX2 REM420/550/	Four X, Four X DL, Four X Lift Dynamic DX2

FOLIR X LIRBAN

### **Invacare Ltd**

Four X SSS - R-Net

**FOLIR X** 

Unit 4 Pencoed Technology Park
Bridgend CF35 5AQ
T: 01656 776200 F: 01656 776201
e: uk@invacare.com www.invacare.co.uk



### Chasswheel brochure - 2015 - Pt no. 1589554

Invacare Ltd - Think Mobility

Think\_Mobility