



Compact and Versatile

The Invacare **Dragon** is a solid and cost effective power chair that provides excellent indoor and outdoor manoeuvrability. The new seat height and tilt feature can easily be adapted to suit the needs of the client and can be achieved without taking the seat apart.



Trusted Solutions, Passionate People





Features and Options

Mobility and Comfort

The robust frame is available with a choice of seating that makes the **Invacare Dragon** a reliable wheelchair for daily use. The compact dimensions provide good manoeuvrability whilst the spring suspension of the motors absorb shocks for increased comfort. Along with the longer range that is achieved from the larger 60Ah batteries. The powerful **Shark II** controller enhances driving with its fingertip precision and is programmable to suit individual needs.



Tilt and Seat Height

With this new feature seat height and tilt can be easily adapted to the clients needs.



Shark II Control

The version of the **Shark** remote offers the possibility to control lights with indicators and up to two electric adjustments.

Adjustability and Easy Service

Invacare Dragon features a range of seating, a manual settable seat tilt and 6 km/h motors. It is adjustable to suit the needs of clients and all parts of the power chair are accessible for ease of servicing.

Technical data

Standard Seat	440 - 490 mm ¹	410, 460, 510 mm	425 - 500 mm	480 / 540 mm	290 - 360 mm	300 - 430 mm	up to 22° ²
Comfort Seat	440 - 490 mm ¹	410, 460, 510 mm	425 - 500 mm	560 mm	290 - 360 mm	300 - 430 mm	up to 22° ²
Standard Seat	-10° - 30° ³	610 mm ⁴	980 / 1040 mm	1030 mm ⁵	850 mm	100 kg ⁸	140 kg
Comfort Seat	-10° - 30° ³	610 mm ⁴	1060 mm	1030 mm ⁵	850 mm	105 kg ⁸	140 kg
	2 x 12 V 60 Ah	1660 mm	60 mm ⁶	15% (8.5°)	2 x 180 W	30 km ⁷	6 km/h

Invacare Limited
Pencoed Technology Park -
Pencoed - Bridgend CF35 5AQ -
United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS -
County Dublin - Ireland
Tel: +353 1 8107084
Fax: +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2016 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Dragon - UK - 01/2016 - 1539041

Frame colour



Blue metallic

- 1/ seat cushion 430 mm
- 2/ varies with seat height
- 3/ in 6 steps, 5° - 30° with gas spring, option
- 4/ depending on seat adjustment

- 5/ incl. leg rest and antitippers
- 6/ 100 mm with kerb climber
- 7/ acc. to ISO 7176-4:2008
- 8/ total weight is depending on configuration

More detailed information, technical data, references for transport or description of certain options can be found in the wheelchair's user manual.



Follow us on Twitter:
[@Think_Mobility](https://twitter.com/Think_Mobility)



Follow us on YouTube:
Invacare Ltd - Think Mobility



Yes, you can.®