



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

MARCH 2021
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Action 3 JR

Standard Features

- Max user weight: 80kg (12.5 stone)
- Crash tested ISO 7176-Part 19
- Self Propulsion or Attendant Propulsion
- Choice of 3 seat widths 12" - 14" (305mm to 355mm)
- Choice of 2 seat depths (or depth adjustable upholstery)
- Choice of 3 seat heights 16" - 18" (410mm to 460mm)
- Black nylon upholstery
- Folding backrest set at 16" (400mm)
- Width and height adjustable backrest upholstery
- 80° detachable legrest hangers
- Flip up, angle and depth adjustable footplates
- Removable and height adjustable desk armrests
- Adjustable heel straps
- Choice of 7 colours (with 8 special colour options)
- Push or pull to lock brakes (configuration dependant)
- Transportation kit and buckle belt
- User tool kit

Standard Base price: £598



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description		
BASE PRODUCT*			
CNB0000	ACTION [®] JUNIOR		○
CNBSPEC	SPECIAL		POA
PROPULSION*			
CNB0100	SELF PROPULSION	(NC 1010, 1170)	£598
CNB0110	TRANSIT PROPULSION	(NC 1020, 1030, 1180, 1190, 1195, 1275, 1390, 1430, 1440)	£598
SEAT WIDTH*			
CNB0180	12" SW	(NC 0250, 0252, 0520, 0570, 1190, 1195)	○
CNB0190	13" SW	(NC 0240, 0251, 0530, 0580)	○
CNB0200	14" SW	(NC 0240, 0251, 0530, 0580)	◆

SEAT DEPTH*			
CNB0240		12" SD	○
CNB0250		14" SD	○
CNB0251		ADJ SD 10" - 12"	○
CNB0252		ADJ SD 12" - 14"	◆
SEAT HEIGHT*			
CNB0260		18" SH	○
CNB0270		17" SH	○
SEAT UPHOLSTERY*			
CNB0330		NYLON WITH HOOK AND LOOP STRIPS	◆
SEAT CUSHION*			
CNB33333		WITHOUT SEAT CUSHION	○
CNB0362		BASIC SEAT CUSHION (ESSENTIAL JUNIOR)	£47
BACKREST*			
CNB0370		FIXED 0° BACKREST 3 HEIGHTS 15" / 370MM, 16" / 400MM OR 17" / 430MM	○
CNB0380		FIXED 7° BACKREST 3 HEIGHTS 15" / 370MM, 16" / 400MM OR 17" / 430MM	○
CNB0381		FIXED 7° BACKREST 20"	○
CNB0390		RECLINING (RATCHET) 0/30° 5 HEIGHTS 13" - 17"	(NC 1190, 1195) £103
CNB0395		RECLINING (GAS) 0/30° 16"	(NC 1190, 1195) £170
CNB0400		FOLDING BACKREST 3 HEIGHTS 15" / 370MM, 16" / 400 OR 17" / 430MM	◆
CNB0401		FOLIDNG BACKREST 20"	○
BACKREST UPHOLSTERY*			
CNB0465		W/O BACK UPH	(NC 0855) ○
CNB0475		WIDTH AND HEIGHT ADJUSTABLE NYLON SLING TYPE	◆
CNB0480		CONTOUR (TENSION) ADJUSTABLE NYLON TYPE	(NC 0855) £100
LEGREST HANGER RH*			
CNB0505		90° STRAIGHT RH REMOVABLE HANGER UL STD 125 TO 260 MM	(NC 0600) £27
CNB0510		80° SHORT RH SWING IN/OUT HANGER 0° STANDARD POSITION; UL STD 220 TO 320 MM BLACK COLOUR ONLY	(NC 0562) ◆
CNB0512		80° LONG RH SWING IN/OUT HANGER 0° STANDARD POSITION; UL STD 300 TO 450 MM	(NC 0560) ○
CNB0520		ELEVATING LEGREST (0° TO 80°) RH SWING IN/SWING OUT HANGER (INCLUDING CALFPAD); UL STD 275 TO 375 MM	(NC 0600, 0611, 0880, 0886) £82
CNB0530		ELEVATING LEGREST (0° TO 80°) RH SWING IN/SWING OUT HANGER (WITHOUT CALFPAD, WITH CALFSTRAP); UL STD 235 TO 295 MM	(NC 0600, 0611, 0880, 0886) £82
LEGREST HANGER LH*			
CNB0555		90° STRAIGHT LH REMOVABLE HANGER UL STD 125 TO 260 MM	(NC 0620) £27
CNB0560		80° SHORT LH SWING IN/OUT HANGER 0° STANDARD POSITION; UL STD 220 TO 320 MM BLACK COLOUR ONLY	◆
CNB0562		80° LONG LH SWING IN/OUT HANGER 0° STANDARD POSITION; UL STD 300 TO 450 MM	○
CNB0570		ELEVATING LEGREST (0° TO 80°) LH SWING IN/SWING OUT HANGER (INCLUDING CALFPAD); UL ST 275 TO 375 MM	(NC 0611, 0620, 0880, 0886) £82
CNB0580		ELEVATING LEGREST (0° TO 80°) LH SWING IN/SWING OUT HANGER (WITHOUT CALFPAD, WITH CALFSTRAP); UL STD 235 TO 295 MM	(NC 0611, 0620, 0880, 0886) £82

FOOTREST RH*			
CNB0600		FLIP UP RH FOOTPLATE	(NC 0882, 0883, 0884, 0885, 0886) ○
CNB0605		FLIP UP DEPTH AND ANGLE ADJ RH	◆
CNB0611		ONE PIECE FOOTBOARD	(NC 0620, 0625, 0882, 0883, 0884, 0885) ○
FOOTREST LH*			
CNB0620		FLIP UP LH FOOTPLATE	(NC 0882, 0883, 0884, 0885, 0886) ○
CNB0625		FLIP UP DEPTH AND ANGLE ADJ LH	◆
ARMREST RH*			
CNB0640		REMOVABLE RH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED DESK, STD ARMPAD (60 MM) ADJ. 180 TO 250 MM	(NC 0745, 0746, 0750, 0755) ◆
CNB0645		REMOVABLE RH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED JUNIOR, STD ARMPAD (60 MM) ADJ. 150 TO 250 MM	(NC 0740, 0746, 0750, 0755) ○
CNB0646		REMOVABLE RH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED NARROW ARMPAD JUNIOR (50 MM) ADJ. 150 TO 250 MM	(NC 0745, 0750, 0755) ○
CNB0650		MUDGUARD RH HEIGHT ADJUSTABLE (155 TO 220 MM)	(NC 0740, 0745, 0746, 0755, 0910, 1201, 1471, 1472) £32
CNB0655		SKIRTGUARD RH JUNIOR	(NC 0740, 0745, 0746, 0750, 0910, 1471, 1472) ○
ARMREST LH*			
CNB0740		REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED DESK, STD ARMPAD (60 MM) ADJ. 180 TO 250 MM	◆
CNB0745		REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED JUNIOR, STD ARMPAD (60 MM) ADJ. 150 TO 250 MM	○
CNB0746		REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED NARROW ARMPAD JUNIOR (50 MM) ADJ. 150 TO 250 MM	○
CNB0750		MUDGUARD LH HEIGHT ADJUSTABLE (155 TO 220 MM)	(NC 0910, 1201, 1471, 1472) £32
CNB0755		SKIRTGUARD LH JUNIOR	(NC 0910, 1471, 1472) ○
OTHER SEATING OPTIONS			
CNB0855		BACKREST EXTENSION	(NC 0860, 0865) £67
CNB0860		HEADREST STANDARD MULTI - ADJUSTABLE	(NC 0865, 0902) £124
CNB0865		HEADREST ROUND MULTI - ADJUSTABLE	(NC 0902) £140
CNB0880		ADJUSTABLE CALF STRAP	£10
CNB0881		FRONT STABILISER KIT + 100 MM	○
CNB0882		ANKLE GRIPS BRACKET	(NC 0883, 0884, 0885) £61
CNB0883		ANKLE GRIPS "XS"	(NC 0884, 0885) £173
CNB0884		ANKLE GRIPS "S"	(NC 0885) £173
CNB0885		ANKLE GRIPS "M"	£173
CNB0886		FOOTPLATE ADJUSTMENT OPTIMIZER KIT	£19
CNB0902		WIDTH ADJUSTABLE BACK BRACE	£38
CNB0910		WITHOUT AUTOLOCK SAFETY DEVICE ON ARMRESTS N°1	○

COLOUR*			
CNB0952		PEARL WHITE	○
CNB0953		MATT BLACK	○
CNB0954		MATT RED	○
CNB0955		INTENSE YELLOW	○
CNB0956		POLAR BLUE	○
CNB0957		CANDY PINK	○
CNB0958		POP BLUE	○
CNB0991		SPECIAL COLOUR (ON DEMAND ONLY)	£156
SWO0001		SLATE GREY - 3 WEEK LEAD TIME	£156
SWO0002		ELECTRIC GREEN - 3 WEEK LEAD TIME	£156
SWO0003		HYPNOTIC PURPLE - 3 WEEK LEAD TIME	£156
SWO0004		ANTHRACITE BLACK - 3 WEEK LEAD TIME	£156
SWO0005		DARK RED - 3 WEEK LEAD TIME	£156
SWO0006		PEARL GREY - 3 WEEK LEAD TIME	£156
SWO0007		VIBRANT ORANGE - 3 WEEK LEAD TIME	£156
FRONT CASTORS*			
CNB0999		5" X 1" BLACK SOLID TYRE	◆
CNB1001		5" X 1" FLASHING LIGHTS SOLID TYRE	£43
CNB1007		5" X 1.3/4" BLACK SOLID TYRE	○
REAR WHEEL SIZE*			
CNB1010		12"	(NC 1045, 1130, 1140, 1141, 1180, 1253, 1254, 1255, 1256, 1380, 1390) ○
CNB1020		20"	(NC 1055, 1050, 1170) ◆
CNB1030		22"	(NC 1045, 1055, 1170) ○
REAR WHEELS*			
CNB1045		REINFORCED STANDARD SELF-PROPEL REAR WHEEL (COLORLESS HUB, CLEAR ANODIZED RIM AND STAINLESS STEEL SPOKES)	(NC 1280) ◆
CNB1050		DRUM BRAKE SELF-PROPEL 22 AND 24" or TRANSIT 12"	(NC 1270, 1275) ○
CNB1055		TRANSIT BLACK PLASTIC	(NC 1280) ○
TYRES*			
CNB1120		SOLID BLACK	◆
CNB1130		PNEUMATIC PROFILE BLACK	○
HANDRIM - A SELECTION MUST BE MADE FOR SELF-PROPULSION			
CNB1140		ALUMINIUM PRE-ANODIZED HANDRIM	○
CNB1141		HANDRIM COATED IN BLACK ANTI-SLIP PLASTIC	£69
REAR WHEEL AXLE*			
CNB1170		FIXED AXLE	◆
CNB1180		Q/R AXLE	○
REAR WHEEL POSITION*			
CNB1190		STANDARD REAR WHEELS POSITION	◆
CNB1195		MEDIUM REAR WHEEL POSITION	○
CNB1200		PASSIVE REAR WHEEL POSITION	○
CNB1201		DUAL AMPUTEE REAR WHEEL POSITION	(NC 1275, 1280) £63
PUMP			
CNB1230		PUMP AND SUPPORT	£9

CHASSIS OPTION			
CNB1253	<input type="checkbox"/>	SPOKE GUARDS CLEAR	£75
CNB1254	<input type="checkbox"/>	SPOKE GUARDS PIXEL	£92
CNB1255	<input type="checkbox"/>	SPOKE GUARDS SNOWFLAKES	£92
CNB1256	<input type="checkbox"/>	SPOKE GUARDS COLOR ME	£92
BRAKES*			
CNB1270	<input type="checkbox"/>	PUSH TO LOCK BRAKES	◆
CNB1275	<input type="checkbox"/>	PULL TO LOCK BRAKES	○
CNB1280	<input type="checkbox"/>	ATTENDANT HUB BRAKES	£136
TRANSPORTATION KIT*			
CNB1345	<input type="checkbox"/>	TRANSPORTATION BRACKETS	◆
SAFETY OPTIONS*			
CNB1360	<input type="checkbox"/>	POSTURE BELT	◆
CNB1361	<input type="checkbox"/>	PEDIATRIC BELT	£84
CNB1370	<input type="checkbox"/>	REAR ANTI-TIPPERS WITH TIPPING AID	(NC 1430, 1440) £32
CNB1380	<input type="checkbox"/>	REFLECTORS KIT	£6
CNB1390	<input type="checkbox"/>	BRAKE EXTENSION PAIR	£13
OTHER OPTIONS			
CNB1430	<input type="checkbox"/>	RH TIPPING AID	£13
CNB1440	<input type="checkbox"/>	LH TIPPING AID	£13
CNB1471	<input type="checkbox"/>	RH SWING AWAY ACRYLIC TRAY	£181
CNB1472	<input type="checkbox"/>	LH SWING AWAY ACRYLIC TRAY	£181

TOTAL PRICE	
--------------------	--

Non Conforming Options:

CNB0240	12" SD
CNB0250	14" SD
CNB0251	ADJ SD 10" - 12"
CNB0252	ADJ SD 12" - 14"
CNB0520	ELEVATING LEGREST (0° TO 80°) RH SWING IN/SWING OUT HANGER (INCLUDING CALFPAD); UL STD 275 TO 375 MM
CNB0530	ELEVATING LEGREST (0° TO 80°) RH SWING IN/SWING OUT HANGER (WITHOUT CALFPAD, WITH CALFSTRAP); UL STD 235 TO 295 MM
CNB0560	80° SHORT LH SWING IN/OUT HANGER 0° STANDARD POSITION; UL STD 220 TO 320 MM BLACK COLOUR ONLY
CNB0562	80° LONG LH SWING IN/OUT HANGER 0° STANDARD POSITION; UL STD 300 TO 450 MM
CNB0570	ELEVATING LEGREST (0° TO 80°) LH SWING IN/SWING OUT HANGER (INCLUDING CALFPAD); UL ST 275 TO 375 MM
CNB0580	ELEVATING LEGREST (0° TO 80°) LH SWING IN/SWING OUT HANGER (WITHOUT CALFPAD, WITH CALFSTRAP); UL STD 235 TO 295 MM
CNB0600	FLIP UP RH FOOTPLATE
CNB0611	ONE PIECE FOOTBOARD
CNB0620	FLIP UP LH FOOTPLATE
CNB0625	FLIP UP DEPTH AND ANGLE ADJ LH
CNB0740	REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED DESK, STD ARMPAD (60 MM) ADJ. 180 TO 250 MM
CNB0745	REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED JUNIOR, STD ARMPAD (60 MM) ADJ. 150 TO 250 MM
CNB0746	REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED NARROW ARMPAD JUNIOR (50 MM) ADJ. 150 TO 250 MM
CNB0750	MUDGUARD LH HEIGHT ADJUSTABLE (155 TO 220 MM)
CNB0755	SKIRTGUARD LH JUNIOR
CNB0855	BACKREST EXTENSION
CNB0860	HEADREST STANDARD MULTI - ADJUSTABLE
CNB0865	HEADREST ROUND MULTI - ADJUSTABLE
CNB0880	ADJUSTABLE CALF STRAP
CNB0882	ANKLE GRIPS BRACKET
CNB0883	ANKLE GRIPS "XS"
CNB0884	ANKLE GRIPS "S"
CNB0885	ANKLE GRIPS "M"
CNB0886	FOOTPLATE ADJUSTMENT OPTIMIZER KIT
CNB0902	WIDTH ADJUSTABLE BACK BRACE
CNB0910	WITHOUT AUTOLOCK SAFETY DEVICE ON ARMRESTS N°1
CNB1010	12"
CNB1020	20"
CNB1030	22"
CNB1045	REINFORCED STANDARD SELF-PROPEL REAR WHEEL (COLORLESS HUB, CLEAR ANODIZED RIM AND STAINLESS STEEL SPOKES)
CNB1050	DRUM BRAKE SELF-PROPEL 22 AND 24" or TRANSIT 12"
CNB1055	TRANSIT BLACK PLASTIC
CNB1130	PNEUMATIC PROFILE BLACK
CNB1140	ALUMINIUM PRE-ANODIZED HANDRIM
CNB1141	HANDRIM COATED IN BLACK ANTI-SLIP PLASTIC
CNB1170	FIXED AXLE
CNB1180	Q/R AXLE
CNB1190	STANDARD REAR WHEELS POSITION

CNB1195	MEDIUM REAR WHEEL POSITION
CNB1201	DUAL AMPUTEE REAR WHEEL POSITION
CNB1253	SPOKE GUARDS CLEAR
CNB1254	SPOKE GUARDS PIXEL
CNB1255	SPOKE GUARDS SNOWFLAKES
CNB1256	SPOKE GUARDS COLOR ME
CNB1270	PUSH TO LOCK BRAKES
CNB1275	PULL TO LOCK BRAKES
CNB1280	ATTENDANT HUB BRAKES
CNB1380	REFLECTORS KIT
CNB1390	BRAKE EXTENSION PAIR
CNB1430	RH TIPPING AID
CNB1440	LH TIPPING AID
CNB1471	RH SWING AWAY ACRYLIC TRAY
CNB1472	LH SWING AWAY ACRYLIC TRAY



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):