

# viamobil<sup>®</sup> eco

Invacare Ltd  
 Unit 4 Pencoed Technology Park  
 Pencoed, Bridgend  
 CF35 5AQ  
 T: 01656 776222  
 E: ordersuk@invacare.com  
 www.invacare.co.uk



Please fill in all specifications with each order.

Account code

Customer order No.

Company name

Address

<b>viamobil V14</b>	<b>1489118</b>	<b>£1034</b>
Included as standard:		
<ul style="list-style-type: none"> <li>• Push &amp; brake aid</li> <li>• Control unit with thumb control</li> <li>• Range: up to 12km</li> <li>• Removable battery 6kg</li> <li>• Max. user weight 120kg</li> <li>• Drive unit 10.5kg with automatic suspension</li> <li>• Anti-slip puncture proof rubber and pedal for raising/lowering the wheels</li> <li>• Suitable for slopes up to 10%</li> <li>• Total weight 18.5kg</li> </ul>		



## Wheelchair details

### Wheelchair manufacture

Model (e.g. Compact)

Seat width



cm

Seat height (e.g.Hemi)



cm

Self propel

Transit

## Accessories

Anti-tippers (pair)*	1489339	£177
Action 2NG conversion kit	1519963	£47
In-car battery charger	1520413	£288
Puncture proof tyres	1466421	£142

## Brackets required for installation onto wheelchair

Fitting on wheelchair at Invacare factory & delivery\* **£75**

Wheelchair order number

Collection & Fitting to 3rd party wheelchair & delivery \*\* **£75 + £96**

Fitting by trained dealer\*\* **£0**

\*Only when ordered with an Invacare wheelchair.

\*\*Mounting kits can only be sent to selected trained dealers.

## Other

Equipped with brake drums


Equipped with anti-tippers

Equipped with vehicle retention system

Rear wheel extension available

Available wheel position:

1=rear      2=centre      3=front



\*NB Always use anti-tippers (Alber or existing)

**Please note there may be some incompatibility with headrest options and rigidizer/tension adjustable backrest with the V14 Viamobil Eco**

Comments (mounting etc.)