viamobil eco

Invacare Ltd Unit 4 Pencoed Technology Park Pencoed, Bridgend **CF35 5AQ** T: 01656 776222 ordersuk@invacare.com www.invacare.co.uk







Please fill in all specifications with each order.

Account code

Customer order No.

Company name

Address

viamobil V14 1489118 £1034

Included as standard:

- Push & brake aid
- Control unit with thumb control
- Range: up to 12km
- Removable battery 6kg
- Max. user weight 120kg
- Drive unit 10.5kg with automatic suspension
- Anti-slip puncture proof rubber and pedal for raising/lowering the wheels
- Suitable for slopes up to 10%
- Total weight 18.5kg

Wheelchair details

Wheelchair manufacture

Model (e.g. Compact)

Seat width Seat height (e.g.Hemi)



cm

cm

Self propel

Transit

Accessories

Anti-tippers (pair)*	1489339	£165
Action 2NG conversion kit	1519963	£44
In-car battery charger	1520413	£269
Puncture proof tyres	1466421	£133

Please note there may be some incompatibility with headrest options and rigidizer/tension adjustable backrest with the V14 Viamobil Eco

Comments (mounting etc.)



Brackets required for installation onto wheelchair

Fitting on wheelchair at Invacare £70 factory & delivery*

Wheelchair order number

£70 + £90 Collection & Fitting to 3rd party wheelchair & delivery ** Fitting by trained dealer** £0

*Only when ordered with an Invacare wheelchair. **Mounting kits can only be sent to selected trained dealers.

Other

Equipped with brake drums

Equipped with anti-tippers

Equipped with vehicle retention

Rear wheel extension available

Available wheel position:

1=rear 2=centre 3=front

*NB Always use anti-tippers (Alber or existing)



