



Invacare Ltd
Unit 4 Pencoed Technology Park
Bridgend CF35 5AQ
T: 01656 776222
E: ordersuk@invacare.com
www.invacare.co.uk



Yes, you can.®

Please fill in all specifications with each order.

Account code

Customer order number

Company name

Address

e-fix E35

£4800

Included as standard:

22"

24"

- 2 e-fix drive wheels
- LED display unit
- Mobility Plus package
- Range 16km (10 miles)
- Max. user weight 120kg
- 6.0 Ah battery (2kg)
- Wheel weight 7.8kg
- Total weight 18.9kg

7.5 Ah battery (2kg)

£235

e-fix E36

£5565

Included as standard:

24"

- 2 e-fix drive wheels
- LED display unit
- Mobility Plus package
- Range 20km (12.4 miles)
- Max. user weight 160kg
- 7.5 Ah battery (2kg)
- Wheel weight 8kg
- Total weight 19.3kg

6Ah battery (2kg)

£235



Wheelchair details

Wheelchair manufacturer: (e.g. Kuschall)

Serial no:

Year:

Model (e.g. Compact):

Seat width:



cm

Seat height (e.g. Hemi):



cm

Armrest & Brake type:

Degree of wheel camber:

°

(Drum brakes must be removed)

Brackets required for installation onto wheelchair:

1457841 £94

Omit brackets £0

Fitting on wheelchair at Invacare factory & delivery* £70

Wheelchair order number:

Collection & fitting to third party:

wheelchair & delivery** £70 + £90

Fitting by trained dealer*** £0

Control Unit Mounting:

Right side £63

Left side £63

Right & Left (pluggable) £63

Bracket on: Frame mount Armrest mount

*Only when ordered with an Invacare wheelchair.

**Pending compatibility

***Mounting kits can only be sent to selected trained dealers.

Do you want to use both your standard and e-fix wheels in the future? If so, please order the retrofit stub axle.

Manual wheel axle retrofit for e-fix * / **

Retrofit, stub axle to 1/2" diameter £96

*If the manual wheels are used in combination with e-fix, these must be fitted with 1/2" stub axles, without drum brakes and 1" wide tyres.

**If drum brakes are installed, these must be removed as they are not compatible with e-fix. Standard lever brakes are used instead.



Yes, you can.®

Accessories

Anti-tippers with jack-up function (pair)*	1489214	£263
--	---------	------

Control Unit Accessories

Swivel arm	1467901	£334
Protection Bumper	1565836	£48

Joystick Options

Option special 40mm/1.57"	1565669	£25
Option soft 70mm/2.76"	1566196	£25
Option T - handle	1565670	£82

*NB Always use anti-tippers (Alber or existing)

Tray table incl. control unit bracket

Control unit centered SW 38-43cm	1566207	£235
Control unit centered SW 43-48cm	1566208	£235
Installation kit for control unit with existing tray table	1566206	£138

Extension for speed pre-selection

Set for speed pre-selection	1566197	£38
-----------------------------	---------	-----

Extension for coupling ring

Set for one drive wheel (single)	1566179	£96
Set for two drive wheels (pair)	1566178	£192

Spoke Covers (pair)

22" wheel, decor 'Alber'	1566201	£60
22" wheel, transparent	1566203	£60
24" wheel, decor 'Alber'	1566670	£60
24" wheel, transparent	1566666	£60

Attendant Control Bracket

Right	1457754	£117
Left	1457755	£117

Intuitive Attendant Control

Control unit intuitive attendant control	1566165	£428
Adapter cable intuitive attendant control	1566177	£88
Additional handle	1566166	£42
Opti-box	1566184	£329
Speed/Stop limit switch	1467698	£70

Battery pack

Spare battery 6.0 Ah	1565639	£517
Spare battery 7.5 Ah	1566153	£679
Transport bag for battery	1566168	£68
Attachment set for transport bag	1566169	£20

Battery holder fixation

Fixed assembly set for cross brace	1566164	£163
Fixed assembly set for cross tube	1566163	£140

Mobility Plus Package*

Allows the individual to connect to their smartphone via Bluetooth and access functions which include battery status, speed, range, drive mode. It also allows speed up to 8km/h, remote control function and on Android devices, the ability to monitor incoming phone calls and SMS via the e-fix display screen.

*iPhone and Android compatible


Other

External charger socket incl. clip	1566177	£91
------------------------------------	---------	-----

Brake extension	1467780	£29
-----------------	---------	-----

Rear wheel extension available

Available wheel position:



1=rear

2=centre

3=front

*NB Always use anti-tippers (Alber or existing)

Comments (mounting etc.)