

# twion®

Invacare Ltd  
Unit 4 Pencoed Technology Park  
Pencoed, Bridgend  
CF35 5AQ  
T: 01656 776222 F: 01656 776220  
E: ordersuk@invacare.com  
www.invacare.co.uk



Please fill in all specifications with each order.

Account code

Customer order number

Company name

Address

## twion M24

£4639

### Included as standard:

- 2 twion push rim 24" drive wheels
- Speed: 6km/h to 10km/h
- Weight: 6.3kg per wheel
- Built-in lithium-ion batteries in the wheel hubs
- Charger with magnetic connectors
- Mobility Plus Package

## twion T24

£4842

### Standard options same as above

- 1591624 White cover, Tyre Marathon, Curve L push rim
- 1591622 Black cover, Tyre Marathon, Curve L push rim
- 1592494 Black cover, Puncture-proof, Curve L push rim
- 1592495 White cover, Puncture-proof, Curve L push rim



### Brackets required for installation onto wheelchair

Bracket as standard	STD
Spare parts set	£213
Omit brackets	-£67
Fitting twion on wheelchair at Invacare factory*	£47
Wheelchair order number	
Fitting to third party wheelchair	£47
Fitting by trained dealer**	£0

\*Only when ordered with an Invacare wheelchair.

\*\*Mounting kits can only be sent to selected trained dealers.

### Wheelchair details

Wheelchair manufacturer (e.g. Kuschall)

Serial no

Year

Model (e.g. Compact)

Seat width

cm

Seat height (e.g. Hemi)

cm

Armrest & Brake type

Degree of wheel camber

°

(Drum brakes must be removed)



Do you want to use both your standard and twion wheels in the future? If so, please order the retrofit stub axle.

Manual wheel axle retrofit for twion \* / \*\*

Retrofit, stub axle to 1/2" diameter

£96



\*If the manual wheels are used in combination with twion, these must be fitted with 1/2" stub axles, without drum brakes and 1" wide tyres.

\*\*If drum brakes are installed, these must be removed as they are not compatible with e-fix. Standard lever brakes are used instead.



**Accessories**

Anti-tippers with jack-up function (pair)\* 1489214 £263

**Push rims - M24 ONLY**

Matt aluminium, coated black STD  
 Stainless steel, brushed (+400g each) £0



*\*NB Always use anti-tippers (Alber or existing)*

**Covers (pairs) - M24 ONLY**

White covers STD  
 Black covers £27



**Tyres (pairs) - M24 ONLY**

Schwalbe Marathon Plus Evolution STD  
 Solid infill black tyre 1/38 £50



**Extra options**

Transportation axle 1591362 £62  
 Loop axle 1591500 £117  
 Brake extension\* 1467780 £24



*\*maybe required depending on wheelchair type*

**Mobility Plus Package\***

Comprises of two wheel hub outer rings which allows the individual to connect via Bluetooth and the licence code information. Functions now include: Cruise control function, battery status, speed and distance, ECO / SPORT mode, adjust speed up to 10km/h, BlueDrive remote control mode, trip profile record.

*\*iPhone and Android compatible*


**Other**

**Comments (mounting etc.)**

External charger socket incl. clip 1566177 £91

Rear wheel extension available

Available wheel position:



1=rear  
 2=centre  
 3=front

*\*NB Always use anti-tippers (Alber or existing)*