



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Powerchairs Retail Prescription Form

NOVEMBER 2020
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

AVIVA RX20 Modulite

Standard Features

- Max user weight 136kg (21 Stone)
- Crash tested ISO 7176 -Part 19
- Modulite highly adjustable seating system
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 510mm stepless adjustment (16" - 20")
- Tension adjustable backrest
- Standard T-Armrests, detachable with short armpad
- Manual settable seat Tilt (0° - 15°)
- Manual settable backrest Recline (0° - 30° in 3.75° steps)
- 80° hangers with fixed footplate
- LiNX control REM110 with standard bracket
- Choice of 5 Shroud colours, 5 Rim insert colours
- 8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre
- 14" (355 mm) x 3" (75 mm), Solid inlay black drive wheels
- CTC Suspension Control - Traction - Comfort
- 6 km/h (4mph) 4 pole Durawatt motors and 52Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)



Standard Base price: £4,944

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
JBC0000	AVIVA RX 20 Base Model Max user weight up to 136 kg	£4,944
JBX1300	Special Extended lead time: + 10 days	<input type="radio"/>
CHASSIS WIDTH*		
JBL2810	Chassis Narrow (560 mm/22")	<input checked="" type="checkbox"/>
SEATING TYPE*		
JBC2100	Invacare® Modulite Seat Unit	<input checked="" type="checkbox"/>

SEAT FRAME*			
JBC2111	<input type="checkbox"/>	Telescopic Seat Frame and Plate with Side Rails	◆
JBC2112	<input type="checkbox"/>	Telescopic Seat Frame + Sling Seat with Side Rails	○
SEAT WIDTH*			
JBC2121	<input type="checkbox"/>	380 mm Seat Plate (NC 2174, 2175, 2176)	○
JBC2122	<input type="checkbox"/>	430 mm Seat Plate	○
JBC2123	<input type="checkbox"/>	480 mm Seat Plate (NC 2704, 2706)	○
JBC2124	<input type="checkbox"/>	530 mm Seat Plate (NC 2175, 2176, 2177, 2178, 2143, 2702, 2704, 2706, 2106)	○
Body Fit Kit			
JBC2106	<input type="checkbox"/>	Body fit Kit (NC 2192, 2194, 2171)	○
SEAT DEPTH*			
JBC2132	<input type="checkbox"/>	Seat Depth 410 mm (NC 2197, 2198)	○
JBC2134	<input type="checkbox"/>	Seat Depth 460 mm	○
JBC2136	<input type="checkbox"/>	Seat Depth 510 mm	○
SEAT CUSHION*			
JBC2142	<input type="checkbox"/>	Matrx Flo-tech Cushion Flo-tech Lite	£107
JBC2143	<input type="checkbox"/>	Matrx Flo-tech Contour Visco	£150
JBC2147	<input type="checkbox"/>	Matrx Libra Cushion	£354
JBC2144	<input type="checkbox"/>	Matrx Vi Cushion	£328
JBC2140	<input type="checkbox"/>	Without seat cushion	○
BACKREST SUPPORT*			
JBC2163	<input type="checkbox"/>	Backposts short (480 mm) (NC 2173)	○
JBC2162	<input type="checkbox"/>	Backposts tall (540 mm)	○
PUSH HANDLES FOR BACKPOSTS			
JBC2108	<input type="checkbox"/>	Backposts without push handles (NC 2219, 2220, 2227, 2228, 2245, 2291, 2292)	○
BACKREST UPHOLSTERY / CUSHION*			
JBC2171	<input type="checkbox"/>	Tension Adjustable Upholstery (NC 2230, 2247, 2106)	◆
JBC2179	<input type="checkbox"/>	Flex3 Tension Adjustable (NC 2192, 2193, 2194, 2199, 1446, 2229, 2155, 2230, 2247)	£140
JBC2172	<input type="checkbox"/>	Flex3 adjustable Backrest Plates without Backrest Cushion (NC 2128)	£228
JBC2173	<input type="checkbox"/>	Flex3 adjustable Backrest Plates with Laguna Backrest cushion (NC 2128)	£370
JBC2174	<input type="checkbox"/>	Matrx Posture Back, height 400 mm (up to 136 kg max user weight)	£508
JBC2175	<input type="checkbox"/>	Matrx Posture Back, height 500 mm (up to 136 kg max user weight)	£508
JBC2176	<input type="checkbox"/>	Matrx Elite Back, height 500 mm (up to 136 kg max user weight)	£508
JBC2177	<input type="checkbox"/>	Matrx Elite TR Back, height 500 mm	£508
JBC2178	<input type="checkbox"/>	Matrx Elite Deep Back, height 500 mm (up to 136 kg max user weight)	£590
JBC2170	<input type="checkbox"/>	Without Backrest Upholstery (NC 2323, 2324, 2325, 2326, 2327, 2230, 2247)	○
UPHOLSTERY TYPE*			
JBC2151	<input type="checkbox"/>	Black Fabric Upholstery	◆
JBC2152	<input type="checkbox"/>	Black multi stretch PU coated fabric Upholstery	○
JBC2155	<input type="checkbox"/>	Black Spacer Fabric Upholstery (NC 2171)	○
JBX3002	<input type="checkbox"/>	No Upholstry Colour required - No backrest upholstery (NC 2179, 2173, 2174, 2175, 2176, 2177, 2178)	○

ARMREST (Pair)*			
JBC2191		T-Arm detachable, Standard	◆
JBC2181		T-Arm detachable, Extended height	£37
JBC2198		Padded Armrest with padded side guard (WATER FALL TYPE) (NC 2195, 2196, 2197, 2199)	£100
JBC2183		Inrigged Armrests	○
JBC2193		Flip-up Armrest (NC 1445, 1446, 2179)	£59
JBC2192		Following Armrest (NC 1446, 2179, 2106)	£168
JBC2194		Following Armrest (NC 1446, 2179, 2106)	£168
ARMPAD (Pair)*			
JBC2195		Short Armpad (325 mm)	◆
JBC2196		Long Armpad (375 mm) (NC 2198)	○
JBC2197		Long Ergonomic Armpad (397 mm) (NC 2198, 2230, 2247)	○
JBC2199		Soft Armpad	£87
JBX3000		Please select with Option 2198 - Padded Armrest with padded side guard (WATER FALL TYPE) (NC 2181, 2183, 2193, 2192)	○
HIP SUPPORTS LH			
JBC2184		Hip Support LH, pad size 140 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£230
JBC2186		Hip Support LH, pad size 230 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£230
JBC2188		Hip Support LH, pad size 330 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£230
HIP SUPPORTS RH			
JBC2185		Hip Support RH, pad size 140 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£230
JBC2187		Hip Support RH, pad size 230 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£230
JBC2189		Hip Support RH, pad size 330 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£230
MOUNTING HARDWARE FOR HIP SUPPORTS			
JBC2164		Mounting hardware for hip supports (6,1" / 155 mm high)	○
JBC2165		Mounting hardware for hip supports (8,1" / 206 mm high)	○
HEADREST / NECKREST - Adequate headrest option must be selected*			
JBC1350		Neckrest depth adjustable (NC 2247)	£74
JBC1432		Headrest sideways adjustable (NC 2247)	£206
JBC1427		Headrest with cheek support (NC 2247)	£206
JBC2323		Elan Standard Pad 150 mm / 6" wide w/standard cover (NC 2247, 2170)	£329
JBC2324		Elan Standard Pad 255 mm / 10" wide w/standard cover (NC 2247, 2170)	£329
JBC2322		Elan Standard Pad 355 mm / 14" wide w/standard cover (NC 2247, 2170)	£329
JBC2325		Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (NC 2247, 2170)	£329
JBC2326		Elan Occipital Pad 305x200 mm / 12"x8" w/standard cover (NC 2247, 2170)	£329
JBC2327		Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard (NC 2247, 2170)	£446
JBC1496		Non Invacare Headrest but profile mount fitting Flex3 hard back (NC 2171, 2174, 2175, 2176, 2177, 2178, 1288, 1289)	○
JBC1495		Non Invacare Headrest	○
LATERAL SUPPORTS			
JBC1446		Lateral Support with extended arm, rigid pad (NC 2171, 2174, 2175, 2176, 2177, 2178, 2179, 2192, 2193, 2194)	£153
JBC1445		Lateral Support with extended arm, multi-adjustable pad (NC 2171, 2174, 2175, 2176, 2177, 2178, 2192, 2193, 2194)	£153
SEAT LIFTER / SEAT TILT*			
JBC1091		Manually settable Seat Tilt 0 - 15°	○
JBC1092		Powered "fixed pivot" Tilt 20° (NC 2201)	£368

SEAT TO GROUND HEIGHT*			
JBC2857		Seat-to-Floor Height 435 mm	○
JBC2858		Seat-to-Floor Height 460 mm	○
JBC2859		Seat-to-Floor Height 485 mm	○
BACKREST RECLINER*			
JBC0975		Manual Settable Recliner (in 3.75° steps) 90° - 120°	◆
JBC0977		Powered Recliner 30°	£460
CENTRE MOUNTED LEGREST*			
JBC0556		Centre mount legrest, angle settable Length: 280 - 385 mm (NC 2366, 0517, 1244, 1416)	£340
JBC0589		LNX Power centre mount legrest, length compensating Length: 340 - 410 mm (NC 2201, 1244, 1416)	£1,495
JBC2855		Without Legrests	○
FOOTPLATES / FOOT PLATFORM*			
JBC2365		Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")	£94
JBC0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")	○
JBC0516		Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5")	£123
JBC2366		Angle adjustable foot platform S (305 W x 215 mm D / 12"x8.5")	£246
JBC0517		Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")	£246
JBX3001		Please select if Centre mount footplate is not required. (NC 0589)	○
SIDE MOUNTED HANGERS LH			
JBC0555		Amputee Support LH	£177
SIDE MOUNTED HANGERS RH			
JBC0554		Amputee Support RH	£177
FOOTREST HANGER TYPE*			
JBC0526		Standard hanger 80° pair Lower leg length: 350 - 450 mm (NC 0522)	◆
JBC0503		Hanger Vari F pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with tools	£109
JBC0533		Hanger Vari A pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with handwheel and calf pad	£213
JBC0532		Hanger ADM pair Lower leg length: 290 - 460 mm Angle adjustment: 80 - 0° with lever and calf pad (NC 0522)	£202
JBC0531		Hanger ADE pair Lower leg length: 290 - 460 mm Powered Angle adjustment: 80 - 0° with calf pad (NC 2201, 0522)	£610
JBC0501		Without Hanger pair (NC 0504, 0513, 0522)	○
FOOTPLATE TYPE*			
JBC0512		Fixed Footplate pair (NC 0533, 0532, 0531)	◆
JBC0504		Angle Adjustable Footplate pair (NC 0526)	£29
JBC0513		Angle and depth adjustable Footplate pair	£108
JBC0522		Angle adjustable Footboard (NC 0532, 0531, 0501, 1244)	£140
JBX3003		No Side Mounted Footplate required (NC 0503, 0533, 0532, 0531)	○

CONTROLS*			
JBC2201		Remote LiNX DLX-REM110 = Drive only 1 Drive function with Drive only power module	(NC 1092, 0977, 0589, 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2258, 2219, 2245, 2270, 1908) ◆
JBC2203		Remote LiNX DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators	(NC 1437, 2231, 2232, 2229, 2233, 2230, 2247, 2257, 2240, 2270, 2258, 1908) ○
JBC2204		Remote LiNX DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights	(NC 2232, 2229, 2233, 2230, 2247, 2257, 2240, 2270, 2258) £138
JBC2205		Rehab Remote LiNX, touch-screen user interface DLX-REM400 Configurable jack sockets	(NC 2220) £311
JBC2206		Rehab Remote LiNX without joystick for alternative access users DLX-REM500 Requires a user input/alternative joystick: 2231, 2232, 2229, 2230, 2247, 2233	(NC 2220, 2270) £707
REM400 TOGGLE OPTION			
JBC2270		Assignable Toggle Switch Module f. REM400 (DLX-SW400TGL)	(NC 2203, 2204, 2206) £288
CONTROL BRACKETRY*			
JBC1207		Standard Bracket RH	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404) ○
JBC1208		Standard Bracket LH	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405) ○
JBC1707		Swing Away Bracket RH	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404) £121
JBC1708		Swing Away Bracket LH	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405) £121
JBC2213		Centre Mount, swing-away to RH, Nucleus	(NC 2283, 2287, 2207, 2265, 2275, 1415, 2217, 1404, 1405) £454
JBC2214		Centre Mount, swing-away to LH, Nucleus	(NC 2284, 2288, 2208, 2266, 2276, 1339, 2218, 1404, 1405) £454
JBC2227		Bracket RH, mounted for Attendant	(NC 1287, 2220, 2128, 2108, 2219, 1221, 1312) £58
JBC2228		Bracket LH, mounted for Attendant	(NC 1286, 2220, 2128, 2108, 2219, 1222, 1363) £58
JBX2239		Individual mounting option per request (Special) Requires additional information/specification	POA
DUAL CONTROL			
JBC2220		Dual Control LiNX (DLX-REM050 as secondary remote)	(NC 2273, 2108, 2205, 2206, 1221, 1222, 2227, 2228) £360
JBC2219		Compact Dual Control LiNX (DLX-ACU200)	(NC 2108, 1221, 1222, 2227, 2228) £360
INVACARE LINX IDC DUAL CONTROL			
JBC2245		Invacare LiNX IDC Dual control	(NC 2108, 2259, 2260) £684
DUAL CONTROL BRACKETRY			
JBC2259		Bracket RH, Attendant	(NC 2291, 2245, 1312) ○
JBC2260		Bracket LH, Attendant	(NC 2292, 2245, 1363) ○
INVACARE LINX IDC DUAL CONTROL BRACKETRY			
JBC1221		Bracket RH, detachable	(NC 2291, 2220, 2219, 2227, 1312) ○
JBC1222		Bracket LH, detachable	(NC 2292, 2220, 2219, 2228, 1363) ○
LINX COMPACT REMOTE			
JBC2231		Proportional Compact Remote DLX-CR400	(NC 2203, 2204) £655
JBC2232		Proportional Compact Remote Low Force DLX-CR-400LF	(NC 2203, 2204) £875
LINX COMPACT REMOTE BRACKET			
JBC2207		RH Mount	(NC 1707, 2213, 2217, 1404) ○
JBC2208		LH Mount	(NC 1708, 2214, 2218, 1405) ○
JBC2265		Centre Mount, swing-away to RH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!	(NC 1707, 2213, 2217, 1404, 1405, 1415) £454
JBC2266		Centre Mount, swing-away to LH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!	(NC 1708, 2214, 2218, 1404, 1405, 1339) £454
JBX2215		Individual mounting option per request (Special) Requires additional information/specification	POA

CHIN CONTROL			
JBC2229		Proportional Chin Joystick (ASL138) SELECTED HEADREST ALWAYS WITH STD MOUNTING HARDWARE	(NC 2201, 2203, 2204, 2230, 2247, 2233, 2257, 2240, 2179, 2128) £1,580
CHIN CONTROL BRACKETRY MANUAL			
JBC1289		Swing Away Bracket RH, side manual (ASL628) Requires Elan Headrest	(NC 1350, 1432, 1427, 2322, 1496) ○
JBC1288		Swing Away Bracket LH, side manual (ASL628) Requires Elan Headrest	(NC 1350, 1432, 1427, 2322, 1496) ○
CHIN CONTROL BRACKETRY POWERED			
JBC1287		Swing Away Bracket RH, powered	(NC 2259, 1221, 2291, 1312, 1415) £1,990
JBC1286		Swing Away Bracket LH, powered	(NC 2260, 1222, 2292, 1363, 1339) £1,990
LINX HEAD CONTROL			
JBC2230		Head Control with ultrasonic proximity sensors (ATOM, ASL104) with standard wings M3	(NC 2203, 2204, 2229, 2233, 2323, 2324, 2322, 2325, 2326, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427) £3,390
JBC2247		Head Control with ultrasonic proximity sensors (PROTON, ASL104P) with PROTON wings	(NC 2203, 2204, 2229, 2233, 2323, 2324, 2322, 2325, 2326, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427) £3,970
LINX FOUR -SWITCH PROXIMITY ARRAY			
JBC2257		Four-Switch Proximity Array (ASL106)	(NC 2203, 2204, 2229, 2230, 2247, 2233, 2240, 2128) £5,250
JBC2275		Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1707, 2213, 2217, 1404, 1405, 1415) £454
JBC2276		Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1708, 2214, 2218, 1404, 1405, 1339) £454
JBX2269		Individual mounting option per request (Special) Requires additional information/specification	POA
LINX WIRELESS ACCESSORIES FOR HEAD CONTROL AND PROXIMITY CONTROL			
JBC2249		Wireless Triple Switch Receiver (ASL557-3)	£311
JBC2250		Wireless Mouse Emulator (ASL558)	£745
LINX MICRO EXTREMITY CONTROL BRACKET			
JBC2233		MEC Proportional Joystick (Micro Extremity Control ASL130)	(NC 2203, 2204, 2229, 2230, 2247, 2257, 2128) £3,292
LINX HEAVY DUTY REMOTE			
JBC2240		Proportional Heavy Duty Remote Joystick w mode switch (ASL133)	(NC 2203, 2204, 2229, 2230, 2247, 2257, 2128) £1,835
LINX MICRO EXTREMITY AND HEAVY DUTY CONTROL BRACKET			
JBC2285		Tray Mount RH (100x150 mm)	(NC 1707, 2213, 2217, 1404, 1405, 1415) ○
JBC2286		Tray Mount LH (100x150 mm)	(NC 1708, 2214, 2218, 1404, 1405, 1339) ○
JBC2283		Centre Mount, swing-away to RH, Nucleus	(NC 1707, 2213, 2217, 1404, 1405, 1415) £454
JBC2284		Centre Mount, swing-away to LH, Nucleus	(NC 1708, 2214, 2218, 1404, 1405, 1339) £454
JBC2287		Centre Mount, swing-away to RH, Nucleus Incl. Tray ASL603, 100x150 mm	(NC 1707, 2213, 2217, 1404, 1405, 1415) £454
JBC2288		Centre Mount, swing-away to LH, Nucleus Incl. Tray ASL603, 100x150 mm	(NC 1708, 2214, 2218, 1404, 1405, 1339) £454
JBX2289		Individual mounting option per request (Special) Requires additional information/specification	POA
LINX 10-WAY SWITCH			
JBC2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access	£460
LINX 10-WAY SWITCH BRACKET			
JBC2217		Bracket RH, mounted for user	(NC 2207, 2287, 2283, 2275, 2213, 2265, 1404) ○
JBC2218		Bracket LH, mounted for user	(NC 2208, 2288, 2284, 2286, 2276, 2214, 2266, 1405) ○
JBC2291		Bracket RH, mounted for attendant	(NC 2259, 1221, 1287, 2108) £58
JBC2292		Bracket LH, mounted for attendant	(NC 2260, 1222, 1288, 2108) £58
LINX USB CHARGER PORT			
JBC2253		USB Charger Port DLX-USB02-A	(NC 2702, 2703, 2704, 2705, 2706) £150

LINX USB CHARGER PORT BRACKET			
JBC2281		RH Mount	○
JBC2282		LH Mount	○
REMOTE STOP SWITCH			
JBC2258		Remote Stop Switch (ASL504) (NC 2203, 2204)	£957
LINX IN200 MODULE			
JBC2273		IN200 Module (Jack port for external power button + control inputs via DB9)	£220
FURTHER LINX COMPONENTS / ACCESSORIES			
JBC2300		To be enclosed only	○
JBX2301		To be connected and programmed (Special)	POA
BUDDY BUTTONS			
JBC1685		Red Buddy Button	£80
JBC1686		Yellow Buddy Button	£80
JBC1687		Green Buddy Button	£80
JBC1690		Blue Buddy Button	£80
EGG SWITCH			
JBC2254		Egg Switch, black	£80
PIKO BUTTONS			
JBC1711		Red Piko Button	£144
JBC1712		Yellow Piko Button	£144
JBC1713		Green Piko Button	£144
JBC1714		Orange Piko Button	£144
JBC1715		Blue Piko Button	£144
JBC1716		Black Piko Button	£144
BODYPOINT JOYSTICK KNOBS			
JBC2721		3 inch U-Shaped, Joystick Handle, Bodypoint PC101A	£57
JBC2722		4 inch U-Shaped, Joystick Handle, Bodypoint PC102A	£57
JBC2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£57
JBC2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£57
5:1 SWITCH ADAPTOR			
JBC1681		5:1 Switch Adapter (9-pin sub D to 5 jack ports)	£110
SHROUD COLOUR*			
JBC2871		Mid range - Black/Anthracite Shroud	○
JBC2872		Mid range - Silver Shroud	○
JBC2873		Mid range - Bengal Red Shroud	○
JBC2874		Mid range - Mercury Blue Shroud	○
JBC2875		Mid range - Midnight Purple Shroud	○
RIM COLOUR*			
JBC2882		Mid range - Black Rim Inserts	○
JBC2883		Mid range - Silver Rim Inserts	○
JBC2884		Mid range - Bengal Red Rim Inserts	○
JBC2885		Mid range - Mercury Blue Rim Inserts	○
JBC2886		Mid range - Midnight Purple Rim Inserts	○
FORKS*			
JBC2632		Double sided fork Without castor rim insert	◆

FRONT WHEEL*			
JBC0828	<input type="checkbox"/>	8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre	◆
CASTOR LOCKS			
JBC1490	<input type="checkbox"/>	Castor Lock	£180
DRIVE WHEEL*			
JBC0704	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Pneumatic	◆
JBC0705	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	○
C.T.C SUSPENSION*			
JBC2870	<input type="checkbox"/>	CTC Suspension Control - Traction - Comfort	◆
SUSPENSION SETTING*			
JBC1910	<input type="checkbox"/>	Soft (suggested for up to 70 kg User Weight)	◆
JBC1912	<input type="checkbox"/>	Medium (suggested for up to 90 kg User Weight)	○
JBC1911	<input type="checkbox"/>	Firm (suggested for up to 110 kg User Weight)	○
JBC1915	<input type="checkbox"/>	Very firm (suggested for up to 136 kg User Weight)	○
MOTOR TYPE*			
JBC1294	<input type="checkbox"/>	Indoor setting, reduced speed approx. 3 km/h	○
JBC1282	<input type="checkbox"/>	6 km/h	◆
JBC1908	<input type="checkbox"/>	10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights	(NC 2203, 1437, 1254) £495
BATTERY TYPE*			
JBC2860	<input type="checkbox"/>	52 Ah AGM (Top-Gel)	◆
CHARGER*			
JBC1343	<input type="checkbox"/>	Opti Charge Charger 8A (IP44)	◆
JBC1313	<input type="checkbox"/>	Opti Charge Charger 10A (IP44)	£138
ANTI-TIPPERS*			
JBC1205	<input type="checkbox"/>	Anti-tipper	◆
LIGHTS			
JBC1437	<input type="checkbox"/>	LED Light system with Indicators 30 Lux headlights	(NC 2203, 1908) £437
JBC1254	<input type="checkbox"/>	Without Light system	○
KERB RIDER			
JBC1244	<input type="checkbox"/>	Kerb Rider With park position feature	(NC 0556, 0589, 1294) £193
TRANSPORTATION KIT*			
JBC1450	<input type="checkbox"/>	Preparation for DAHL-Docking station (Adequate headrest is compulsory) Includes a mounting bracket in the battery box Docking Station for car implementation has to be ordered separately with DAHL Engineering	£534
JBC1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory	◆
LAP BELT*			
JBC1329	<input type="checkbox"/>	Lap belt with Auto Lock Standard mount on backrest	○
JBC2113	<input type="checkbox"/>	Lap belt with Auto Lock Flexible mount in side rails	◆

POSITIONING BELT				
JBC2701		BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£70
JBC2702		BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£70
JBC2703		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£72
JBC2704		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£72
JBC2705		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£72
JBC2706		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£72
REAR VIEW MIRROR				
JBC1415		Rear view mirror RH	(NC 2217, 2285, 2213, 2265, 2275, 2283, 2287)	£40
JBC1339		Rear view mirror LH	(NC 2218, 2286, 2214, 2266, 2276, 2284, 2288)	£40
OTHER OPTIONS				
JBC1438		Luggage rack		£300
JBC2128		Bag Hooks on Modulite Seat Max load up to 6 kg = impact on max user weight	(NC 2172, 2173, 2227, 2228, 2229, 2233, 2240, 2257, 2291, 2292)	£70
JBC1312		Cane holder RH	(NC 2227, 2259, 1221)	£38
JBC1363		Cane holder LH	(NC 2228, 2260, 1222)	£38
JBC1416		Fender for Castors		£80
JBC1417		Fender for Drive wheels		£154
JBC1404		Tray, Swing Away RH	(NC 1707, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2207, 2213, 2214, 2217)	£152
JBC1405		Tray, Swing Away LH	(NC 1708, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2208, 2213, 2214, 2218)	£152

TOTAL PRICE	
--------------------	--

Non Conforming Options:

JBC2121	380 mm Seat Plate
JBC2122	430 mm Seat Plate
JBC2123	480 mm Seat Plate
JBC2124	530 mm Seat Plate
JBC2106	Body fit Kit
JBC2132	Seat Depth 410 mm
JBC2143	Matrx Flo-tech Contour Visco
JBC2108	Backposts without push handles
JBC2171	Tension Adjustable Upholstery
JBC2179	Flex3 Tension Adjustable
JBC2172	Flex3 adjustable Backrest Plates without Backrest Cushion
JBC2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion
JBC2174	Matrx Posture Back, height 400 mm (up to 136 kg max user weight)
JBC2175	Matrx Posture Back, height 500 mm (up to 136 kg max user weight)
JBC2176	Matrx Elite Back, height 500 mm (up to 136 kg max user weight)
JBC2177	Matrx Elite TR Back, height 500 mm
JBC2178	Matrx Elite Deep Back, height 500 mm (up to 136 kg max user weight)
JBC2170	Without Backrest Upholstery
JBC2155	Black Spacer Fabric Upholstery
JBC2181	T-Arm detachable, Extended height
JBC2198	Padded Armrest with padded side guard (WATER FALL TYPE)
JBC2183	Inrigged Armrests
JBC2193	Flip-up Armrest
JBC2192	Following Armrest
JBC2194	Following Armrest
JBC2195	Short Armpad (325 mm)
JBC2196	Long Armpad (375 mm)
JBC2197	Long Ergonomic Armpad (397 mm)
JBC2199	Soft Armpad
JBC2184	Hip Support LH, pad size 140 x 100 mm
JBC2186	Hip Support LH, pad size 230 x 100 mm
JBC2188	Hip Support LH, pad size 330 x 100 mm
JBC2185	Hip Support RH, pad size 140 x 100 mm
JBC2187	Hip Support RH, pad size 230 x 100 mm
JBC2189	Hip Support RH, pad size 330 x 100 mm
JBC1350	Neckrest depth adjustable
JBC1432	Headrest sideways adjustable
JBC1427	Headrest with cheek support
JBC2323	Elan Standard Pad 150 mm / 6" wide w/standard cover
JBC2324	Elan Standard Pad 255 mm / 10" wide w/standard cover
JBC2322	Elan Standard Pad 355 mm / 14" wide w/standard cover
JBC2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover

JBC2326	Elan Occipital Pad 305x200 mm / 12"x8" w/standard cover
JBC2327	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard
JBC1496	Non Invacare Headrest but profile mount fitting Flex3 hard back
JBC1446	Lateral Support with extended arm, rigid pad
JBC1445	Lateral Support with extended arm, multi-adjustable pad
JBC1092	Powered "fixed pivot" Tilt 20°
JBC0977	Powered Recliner 30°
JBC0556	Centre mount legrest, angle settable Length: 280 - 385 mm
JBC0589	LNX Power centre mount legrest, length compensating Length: 340 - 410 mm
JBC2366	Angle adjustable foot platform S (305 W x 215 mm D / 12"x8.5")
JBC0517	Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")
JBC0526	Standard hanger 80° pair Lower leg length: 350 - 450 mm
JBC0503	Hanger Vari F pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with tools
JBC0533	Hanger Vari A pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with handwheel and calf pad
JBC0532	Hanger ADM pair Lower leg length: 290 - 460 mm Angle adjustment: 80 - 0° with lever and calf pad
JBC0531	Hanger ADE pair Lower leg length: 290 - 460 mm Powered Angle adjustment: 80 - 0° with calf pad
JBC0501	Without Hanger pair
JBC0504	Angle Adjustable Footplate pair
JBC0513	Angle and depth adjustable Footplate pair
JBC0522	Angle adjustable Footboard
JBC2201	Remote LiNX DLX-REM110 = Drive only 1 Drive function with Drive only power module
JBC2203	Remote LiNX DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators
JBC2204	Remote LiNX DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights
JBC2205	Rehab Remote LiNX, touch-screen user interface DLX-REM400 Configurable jack sockets
JBC2206	Rehab Remote LiNX without joystick for alternative access users DLX-REM500 Requires a user input/alternative joystick: 2231, 2232, 2229, 2230, 2247, 2233
JBC2270	Assignable Toggle Switch Module f. REM400 (DLX-SW400TGL)
JBC1707	Swing Away Bracket RH
JBC1708	Swing Away Bracket LH
JBC2213	Centre Mount, swing-away to RH, Nucleus
JBC2214	Centre Mount, swing-away to LH, Nucleus
JBC2227	Bracket RH, mounted for Attendant
JBC2228	Bracket LH, mounted for Attendant
JBC2220	Dual Control LiNX (DLX-REM050 as secondary remote)
JBC2219	Compact Dual Control LiNX (DLX-ACU200)
JBC2245	Invacare LiNX IDC Dual control
JBC2259	Bracket RH, Attendant
JBC2260	Bracket LH, Attendant
JBC1221	Bracket RH, detachable
JBC1222	Bracket LH, detachable
JBC2231	Proportional Compact Remote DLX-CR400

JBC2232	Proportional Compact Remote Low Force DLX-CR-400LF
JBC2207	RH Mount
JBC2208	LH Mount
JBC2265	Centre Mount, swing-away to RH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!
JBC2266	Centre Mount, swing-away to LH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!
JBC2229	Proportional Chin Joystick (ASL138) SELECTED HEADREST ALWAYS WITH STD MOUNTING HARDWARE
JBC1289	Swing Away Bracket RH, side manual (ASL628) Requires Elan Headrest
JBC1288	Swing Away Bracket LH, side manual (ASL628) Requires Elan Headrest
JBC1287	Swing Away Bracket RH, powered
JBC1286	Swing Away Bracket LH, powered
JBC2230	Head Control with ultrasonic proximity sensors (ATOM, ASL104) with standard wings M3
JBC2247	Head Control with ultrasonic proximity sensors (PROTON, ASL104P) with PROTON wings
JBC2257	Four-Switch Proximity Array (ASL106)
JBC2275	Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JBC2276	Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JBC2233	MEC Proportional Joystick (Micro Extremity Control ASL130)
JBC2240	Proportional Heavy Duty Remote Joystick w mode switch (ASL133)
JBC2285	Tray Mount RH (100x150 mm)
JBC2286	Tray Mount LH (100x150 mm)
JBC2283	Centre Mount, swing-away to RH, Nucleus
JBC2284	Centre Mount, swing-away to LH, Nucleus
JBC2287	Centre Mount, swing-away to RH, Nucleus Incl. Tray ASL603, 100x150 mm
JBC2288	Centre Mount, swing-away to LH, Nucleus Incl. Tray ASL603, 100x150 mm
JBC2217	Bracket RH, mounted for user
JBC2218	Bracket LH, mounted for user
JBC2291	Bracket RH, mounted for attendant
JBC2292	Bracket LH, mounted for attendant
JBC2253	USB Charger Port DLX-USB02-A
JBC2258	Remote Stop Switch (ASL504)
JBC2273	IN200 Module (Jack port for external power button + control inputs via DB9)
JBC1294	Indoor setting, reduced speed approx. 3 km/h
JBC1908	10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBC1437	LED Light system with Indicators 30 Lux headlights
JBC1254	Without Light system
JBC1244	Kerb Rider With park position feature
JBC2701	BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm
JBC2702	BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm
JBC2703	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm
JBC2704	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm
JBC2705	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm

JBC2706	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm
JBC1415	Rear view mirror RH
JBC1339	Rear view mirror LH
JBC2128	Bag Hooks on Modulite Seat Max load up to 6 kg = impact on max user weight
JBC1312	Cane holder RH
JBC1363	Cane holder LH
JBC1416	Fender for Castors
JBC1404	Tray, Swing Away RH
JBC1405	Tray, Swing Away LH



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):