



rea

Confidence in posture solutions

rea | **Azalea[®] Tall**

No limits on comfort
for taller individuals



Yes, you can.

Rea Azalea Tall is specially designed to meet the needs of tall people who require a tilt-in-space wheelchair with a longer seat support. Adapted from the *Rea Azalea*, the *Rea Azalea Tall* boasts all the advantages of a reliable, tilting wheelchair and offers a unique weight-shifting mechanism. The *Rea Azalea Tall*: at the top of its class for comfort and stability, offers multi adjustable Flex 3 backrest to allow for individual adaptability. 5 frame colours and three upholstery options are available.



Real compatibility

The *Rea Azalea Tall* has all the features of the *Rea Azalea* such as tilt-in-space weight-shifting for a really stable tilt. Available as a transit or self propelling version, all accessories are cross compatible with the Azalea range.

Easy manoeuvrability

Attendants will appreciate the adjustable push bar, which offers good ergonomic grip and can be adopted to fit the height of the attendant, making the chair easier to manoeuvre.

Features and Options



Longer seat plate

An extra 50 mm of seat depth is designed to support longer legs.



Chassis length

An extra 50 mm in length enables greater forward stability.



Seat height

An extra 50 mm of height allows for longer legrests.


Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please visit www.invacare.co.uk


	Seat width	Seat depth	Seat height	Backrest height	Armrest height	Length of legrests	Seat angle	Backrest angle
<i>Rea Azalea Tall</i>	390-450/440-500/ 490-550 mm (390-490/440-540/ 490-590 mm)	480 - 550 mm	500 mm	Backrest plate: 560 - 790 mm Tension adjustable backrest: 540 - 715 mm	240 - 360 mm	330 - 500 mm	-1° to 23° (+1° to 25°)	0° - 30°
	Driving unit width	Total height	Total length inc. Footrests	Total weight	Max. User weight	Transport weight	Crash tested	Maximum slope for use of parking brake
<i>Rea Azalea Tall</i>	SW + 220 mm (SW + 230 mm)	1000 - 1300 mm	950 - 1070 mm	34 kg (SW 440 mm)	Max 135 kg	18,7 kg ¹	ISO 7176-19	7°

1. Transport weight increases by 4,2 kg when Power tilt & recline option is selected.

Frame colour



Upholstery



Invacare Limited
Pencoed Technology Park - Pencoed
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail: ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

