



**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

# Manual Retail Prescription Form

JANUARY 2022  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Küschall Champion 2.0

### Standard Features

- Max User Weight 120kg (19 Stone)
- Lightweight Folding Fully Hydroformed Aluminium Wheelchair
- Multiple adjustments, Front, Rear Seat Heights and Centre of Gravity
- Aluminium Front Frame with 75° or 90° Front Frame Angles
- 24" Self Propel Q/R Wheels with Grey Aluminium Handrims
- Schwalbe Black and Grey Pneumatic Tyres
- Choice of 7 Seat Widths 360mm to 480mm (14" - 19")
- Choice of 8 Seat Depths 340mm to 480mm (13" - 19")
- Choice of 5 Standard Frame Colours
- Range of Front Seat Heights 450mm to 540mm (17" -21")
- Range of Rear Seat Heights 390mm to 480mm (15" -19")
- Folding Backrest, Angle Adjustable and Soft Lock Feature
- Black Nylon Light Upholstery
- Angle and Height Adjustable Rigid Footrest, with Calf Strap
- Without Push Handles, Flip up Carbon Clothesguards
- Unfolding Help Cord. Black Comfort Front Castor
- Standard Brake Lever Bent
- Tool Kit



Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

Code	Description		
<b>BASE PRODUCT*</b>			
DHA0031	Küschall Champion 2.0	<input type="radio"/>	
DHA0000	MYKÜSCHALL Special Request - Please include quotation reference in special description at the end of the form where applicable	<input type="radio"/>	
<b>FRAME*</b>			
DHA1150	Frame Angle 75° Aluminium	(NC 1350)	£2,549
DHA1160	Frame Angle 90° Aluminium	(NC 1350)	£2,549
DHA1170	Frame Angle 75° Aluminium - Tapered Front Frame	(NC 1320, 1325, 7010)	£2,549
DHA1180	Frame Angle 90° Aluminium - Tapered Front Frame	(NC 1320, 1325, 7010)	£2,549

SEAT WIDTH*				
DHA1320		360 mm	(NC 7000)	<input type="radio"/>
DHA1325		380 mm	(NC 7000)	<input type="radio"/>
DHA1330		400 mm		<input type="radio"/>
DHA1335		420 mm		<input type="radio"/>
DHA1340		440 mm		<input type="radio"/>
DHA1345		460 mm		<input type="radio"/>
DHA1350		480 mm		<input type="radio"/>
SEAT DEPTH*				
DHA1405		340 mm	(NC 2630, 2720, 2730, 2750)	<input type="radio"/>
DHA1415		360 mm	(NC 2720, 2730)	<input type="radio"/>
DHA1425		380 mm	(NC 2720, 2730)	<input type="radio"/>
DHA1435		400 mm		<input type="radio"/>
DHA1445		420 mm		<input type="radio"/>
DHA1460		440 mm		<input type="radio"/>
DHA1470		460 mm		<input type="radio"/>
DHA1480		480 mm		<input type="radio"/>
BACK HEIGHT*				
DHA1600		Back Height for 3rd party backrest. Final height depends on the height of the positioning backrest, select if MatrX/Fixed backrest is required.	(NC 2200, 2220, 2230)	<input type="radio"/>
DHA1615		300 mm	(NC 2200, 2240, 2400, 2595)	<input type="radio"/>
DHA1620		315 mm	(NC 2200, 2240, 2400, 2595)	<input type="radio"/>
DHA1625		330 mm	(NC 2240, 2595)	<input type="radio"/>
DHA1630		345 mm	(NC 2240, 2595)	<input type="radio"/>
DHA1635		360 mm	(NC 2240, 2595)	<input type="radio"/>
DHA1640		375 mm	(NC 2240, 2595)	<input type="radio"/>
DHA1645		390 mm	(NC 2240, 2595)	<input type="radio"/>
DHA1650		405 mm	(NC 2240)	<input type="radio"/>
DHA1655		420 mm	(NC 2240)	<input type="radio"/>
DHA1660		435 mm	(NC 2240)	<input type="radio"/>
DHA1665		450 mm	(NC 2240)	<input type="radio"/>
DHA1670		465 mm	(NC 2240)	<input type="radio"/>
FRONT SEAT HEIGHT*				
DHA1740		450 mm	(NC 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960)	<input type="radio"/>
DHA1745		460 mm	(NC 1948, 1950, 1952, 1954, 1956, 1958, 1960)	<input type="radio"/>
DHA1750		470 mm	(NC 1950, 1952, 1954, 1956, 1958, 1960)	<input type="radio"/>
DHA1755		480 mm	(NC 1952, 1954, 1956, 1958, 1960)	<input type="radio"/>
DHA1760		490 mm	(NC 1954, 1956, 1958, 1960)	<input type="radio"/>
DHA1765		500 mm	(NC 1956, 1958, 1960)	<input type="radio"/>
DHA1770		510 mm	(NC 1958, 1960)	<input type="radio"/>
DHA1775		520 mm	(NC 1960)	<input type="radio"/>
DHA1780		530 mm	(NC 5000)	<input type="radio"/>
DHA1785		540 mm	(NC 5000, 5050, 5100, 5150, 5200, 5250)	<input type="radio"/>

**REAR SEAT HEIGHT. On 75° frame, Rear Seat Height can be 0mm - 70mm lower then front seat height and 90° can be 30mm to 100mm lower\***

DHA1815		390 mm	(NC 6420, 6430)	○
DHA1820		400 mm	(NC 6420, 6430)	○
DHA1825		410 mm	(NC 6430)	○
DHA1830		420 mm	(NC 6430)	○
DHA1835		430 mm		○
DHA1840		440 mm		○
DHA1845		450 mm	(NC 6400)	○
DHA1850		460 mm	(NC 6400)	○
DHA1855		470 mm	(NC 6400, 6410)	○
DHA1860		480 mm	(NC 6400, 6410)	○

**LOWER LEG LENGTH. Cushion height must be deducted\***

DHA1924		320 mm		○
DHA1926		330 mm		○
DHA1928		340 mm		○
DHA1930		350 mm		○
DHA1932		360 mm		○
DHA1934		370 mm		○
DHA1936		380 mm		○
DHA1938		390 mm		○
DHA1940		400 mm		○
DHA1942		410 mm		○
DHA1944		420 mm		○
DHA1946		430 mm		○
DHA1948		440 mm		○
DHA1950		450 mm		○
DHA1952		460 mm		○
DHA1954		470 mm		○
DHA1956		480 mm		○
DHA1958		490 mm		○
DHA1960		500 mm		○

**BACKREST\***

DHA2140		Angle Adustable, Forward Folding with Soft Lock. Backrest remains locked for transportation		◆
---------	--	---	--	---

**BACKREST UPHOLSTERY\***

DHA2200		Black Standard Light Upholstery, without Tension Adjustment	(NC 2430, 2595)	◆
DHA2220		Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener		£130
DHA2230		Tension Adjustable Backrest, with Black Aero, Breathable Spacer Upholstery Cover, Back Pocket/hook and loop fastener		£212
DHA2240		Without Backrest Upholstery Select if fitting a Matrix/Fixed backrest, no upholstery required	(NC 2430, 2595)	○

<b>PUSH HANDLES*</b>			
DHA2400		Without Push Handles	◆
DHA2410		Push Handles Fixed Short/Mini Grip	£126
DHA2420		Push Handles Fixed Long Grip	○
DHA2430		Push Handles Height Adjustable, Rear Set (NC 3480)	£241
DHA2440		Push Handles Height Adjustable, Integrated in backrest (NC 2595)	£208
DHA2450		Push Handles Foldable (NC 2595)	£188
<b>BACKREST ANGLE*</b>			
DHA2510		76° angle from the seat	○
DHA2530		80.5° angle from the seat	○
DHA2550		85° angle from the seat	◆
DHA2570		89.5° angle from the seat	○
DHA2590		94° angle from the seat	○
<b>BACKREST STABILIZATION BAR</b>			
DHA2595		Backrest Stabilization Bar, Swingaway and Removable	£156
<b>SEAT UPHOLSTERY*</b>			
DHA2600		Standard Seat Upholstery	◆
DHA2630		Light Seat Upholstery, mesh material to reduce weight	£66
<b>SEAT CUSHION</b>			
DHA2700		Standard Seat Cushion, 60 mm Double Foam with Visco Top Layer	£63
DHA2720		Matrx Libra Cushion, Black Smooth Cover - Very High Risk	£379
DHA2730		Matrx Vi Cushion, Black Smooth Cover - High Risk	£351
DCA2750		Vicair Active O2 lightweight cushion 6cm deep	£495
<b>SEAT BAG</b>			
DHA2800		Seat Bag, Click in and removable, includes the fixation clips	£70
<b>FOOTREST*</b>			
DHA3240		Rigid, Plastic Footplate, Set to front, Height and Angle Adjustable. Note: This does not flip up for transfer, but folds easy for transport and storage	◆
DHA3250		Rigid, Plastic Footplate, Set to rear, Height and Angle Adjustable. Note: This does not flip up for transfer, but folds easy for transport and storage	○
<b>CALF STRAP*</b>			
DHA3320		Calf Strap Mounted on Supporter	◆
DHA3330		Calf Strap Mounted on Frame	○
<b>CLOTHESGUARDS/MUDGUARDS*</b>			
DHA3420		Carbon Clothesguards, Flip Up	◆
DHA3430		Carbon Mudguards, Fixed (NC 3490)	£235
DHA3440		Carbon Mudguards, Removable (NC 3490)	£292
<b>ARMREST</b>			
DHA3480		Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included	£229
DHA3490		Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included	£101

FRAME COLOUR*			
DHA4000		Space Black (Standard)	○
DHA4010		Anthracite Grey (Standard)	○
DHA4020		Pure White (Standard) (NC 4220)	○
DHA4030		Carmin Red (Standard)	○
DHA4040		Azure Blue (Standard)	○
DHA4060		Flashy Yellow (Additional 10 day lead time)	£29
DHA4070		Atomic Violet (Additional 10 day lead time)	£29
DHA4100		Black Glossy (Additional 10 day lead time)	£29
DHA4130		Silver Grey (Additional 10 day lead time)	£29
DHA4140		Royal Blue (Additional 10 day lead time)	£29
DHA4150		Ruby Red (Additional 10 day lead time)	£29
DHA4170		Dynamite Orange (Additional 10 day lead time)	£29
DHA4180		Forest Green (Additional 10 day lead time)	£29
DHA4050		Personal frame colour. Provide standard RAL colour code of your choice	£295
FRAME DECAL TYPE*			
DHA4190		Bold Decal (Full coloured transfer on frame) (NC 4300, 4310, 1111)	£29
DHA4195		Tag (Small Kuschall brand transfer on front of frame) (NC 4200, 4210, 4220, 4230, 4240, 4250, 2222)	◆
DHA4197		Without Tag (Kuschall branding only) (NC 4200, 4210, 4220, 4230, 4240, 4250, 4300, 4310)	○
BOLD DECAL COLOUR*			
DHA4200		Black	○
DHA4210		Dark Grey	○
DHA4220		Light Grey	○
DHA4230		Red	○
DHA4240		Blue	○
DHA4250		Gold	○
DHA1111		No bold colour, select only with Option DHA4195 Tag or DHA4197 Without Tag (NC 4190)	○
TAG COLOUR*			
DHA4300		Black	○
DHA4310		Red	○
DHA2222		No tag colour, select only with Option DHA4190 Bold Decal or DHA4197 Without Tag (NC 4195)	○

FRONT CASTORS*			
DHA5000		75mm x 25mm (3" x 1") Supersport	○
DHA5050		100mm x 30mm (4" x 1.2") Everyday Black hub and tyre	○
DHA5100		100mm x 30mm (4" x 1.2") Comfort Black hub and tyre	○
DHA5150		100mm x 25mm (4" x 1") Starec Black hub and tyre	£79
DHA5200		100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs)	£88
DHA5250		100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs)	£88
DHA5300		125mm x 30mm (5" x 1.2") Everyday Black hub and tyre	○
DHA5350		125mm x 30mm (5" x 1.2") Comfort Black hub and tyre	◆
DHA5400		125mm x 25mm (5" x 1") Starec Black hub and tyre	£79
DHA5450		125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs)	£88
DHA5500		125mm x 36mm (5" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs)	£88
REAR WHEEL POSITION (CENTRE OF GRAVITY)*			
DHA6100		1st from front Most Active Position	○
DHA6110		2nd from front	○
DHA6120		3rd from front	○
DHA6130		4th from front	○
DHA6140		5th from front	◆
DHA6150		6th from front	○
DHA6160		7th from front	○
DHA6170		8th from front	○
DHA6180		9th from front	○
QUAD QUICK RELEASE AXLE			
DHA6300		Quad Quick Release/Tetra Clip (NC 0907, 0909, 6550)	£81
REAR WHEEL SIZE*			
DHA6400		22" (NC 6520, 6550, 6600, 6620, 6630, 6660, 6710, 6720, 6740, 6750, 6760, 6770)	○
DHA6410		24" (NC 6640)	◆
DHA6420		25" (NC 6510, 6600, 6640, 6660, 6700, 6710, 6720, 6760, 6770, 6900)	○
DHA6430		26" (NC 0907, 0909, 6520, 6550, 6600, 6630, 6640, 6660, 6720, 6760, 6770)	○
REAR WHEEL TYPE*			
DHA6510		22", 24", 26" Kuschall Starec Wheel, 24" & 26" Black Cross Spokes, Rim & Hub. 22" Silver Cross-Spokes, Colourless hub and Black rim	◆
DHA6520		24", 25" High Performance Wheel Black hub and rim (NC 6660)	£173
DHA6550		24", 25" Spinergy LXX, Black hub, rim and 18 spokes in Black (NC 6660)	£580
SWO0907		24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. Coloured spokes are available via MyKuschall (NC 6660)	£462
SWO0909		24", 25" Spinergy LX, Chrome hub and 12 spokes in Black. Coloured spokes are available via MyKuschall (NC 6660)	£698

<b>REAR WHEEL TYRE*</b>			
DHA6600		24" Schwalbe Rightrun Black/Grey Pneumatic tyre	◆
DHA6620		24", 25", 26" Primo Sentinel tyre. Pneumatic and puncture resistant	£95
DHA6630		24", 25" Schwalbe One V-Guard Puncture Resistant tyre (Cut-Resistant High-Tech Fibre)	£92
DHA6640		22" Profile Grey Pneumatic tyre 7 Bar (NC 7230)	○
DHA6660		24" Solid Profile Black Tyre (NC 7230)	£63
<b>HANDRIM*</b>			
DHA6700		22", 24", 26" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance	£212
DHA6710		24", 26" Steel, short tab with 1 hole, mounted at 3 cm distance	£59
DHA6720		24" Titanium, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance	£182
DHA6740		24" , 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip	£304
DHA6750		24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands	£304
DHA6760		24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion	£270
DHA6770		24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim	£338
DHA6780		22", 24", 25", 26" Light Aluminium Grey, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance	◆
DHA6790		22", 24", 25", 26" Light Aluminium Black, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3 cm distance	£40
<b>REAR WHEEL CAMBER*</b>			
DHA6820		0° Camber	◆
DHA6830		1° Camber	○
DHA6840		3° Camber (NC 7230)	○
<b>SPOKEGUARD</b>			
DHA6900		22", 24", 26" Transparent Spokeguards	£120
<b>FIXATION CLIP*</b>			
DHA7000		Fixation Strap (Hook and loop strap under the seat holds the frame together when folded for easy transportation, designed for clients with limited hand function)	○
DHA7010		Fixation Clip (Plastic clip on either side of the seat holds the frame together when folded for easy transportation)	◆
DHA3333		No fixation clip, select with Tapered frame options 1170, 1180	○
<b>UNFOLDING HELP*</b>			
DHA7100		Unfolding Help (Cord to help open the folding mechanism from the front of the chair)	◆

<b>BRAKES*</b>			
DHA7200		Pull Brakes	£74
DHA7210		Standard Push to Lock Brakes, Lever Bent	◆
DHA7220		Standard Push to Lock Brakes, Lever Straight	○
DHA7230		Performance Brakes, Push to Lock. Brake barrel flips away from tyre (NC 7300)	£74
DHA7240		Active/Scissor Brakes (NC 7300)	£74
DHA7250		Light/Compact Active, Aluminium Brakes (NC 7300)	£129
<b>BRAKE LEVER EXTENSION</b>			
DHA7300		Brake Lever Extension	£47
<b>SAFETY OPTIONS</b>			
DHA7500		Anti-Tipper Right	£116
DHA7510		Anti-Tipper Left	£116
DHA7520		Anti-Tipper Both sides	£232
DHA7800		Standard Pelvic belt	£70
<b>OTHER OPTIONS</b>			
DHA8010		Front Frame Protectors	£35
DHA8100		Tipper-Aid, Mounted Right Side Max. 2 adaptations per side (Adaptations = Anti tipper, Tipper aid and Cane holder)	£44
DHA8200		Cane Holder, Mounted Right Side Max. 2 adaptations per side (Adaptations = Anti tipper, Tipper aid and cane holder)	£71
DHA8500		Air Pump	£12
DHA8550		Tool Kit	◆
DHA8680		24" Mountain Bike Wheels with bag. Fitted with Specific Black Tyres, Grey Aluminium Handrims; Quick Release Axles Included	£304
DHA8700		Travel bag for Champion, Wheels can remain on the wheelchair	£241
<b>TOTAL PRICE</b>			



**Non Conforming Options:**

DHA1320	360 mm
DHA1325	380 mm
DHA1350	480 mm
DHA1946	430 mm
DHA1948	440 mm
DHA1950	450 mm
DHA1952	460 mm
DHA1954	470 mm
DHA1956	480 mm
DHA1958	490 mm
DHA1960	500 mm
DHA2200	Black Standard Light Upholstery, without Tension Adjustment
DHA2220	Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener
DHA2230	Tension Adjustable Backrest, with Black Aero, Breathable Spacer Upholstery Cover, Back Pocket/hook and loop fastener
DHA2240	Without Backrest Upholstery Select if fitting a Matrix/Fixed backrest, no upholstery required
DHA2400	Without Push Handles
DHA2430	Push Handles Height Adjustable, Rear Set
DHA2595	Backrest Stabilization Bar, Swingaway and Removable
DHA2630	Light Seat Upholstery, mesh material to reduce weight
DHA2720	Matrix Libra Cushion, Black Smooth Cover - Very High Risk
DHA2730	Matrix Vi Cushion, Black Smooth Cover - High Risk
DCA2750	Vicair Active O2 lightweight cushion 6cm deep
DHA3480	Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included
DHA3490	Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included
DHA4190	Bold Decal (Full coloured transfer on frame)
DHA4195	Tag (Small Kuschall brand transfer on front of frame)
DHA4200	Black
DHA4210	Dark Grey
DHA4220	Light Grey
DHA4230	Red
DHA4240	Blue
DHA4250	Gold
DHA1111	No bold colour, select only with Option DHA4195 Tag or DHA4197 Without Tag
DHA4300	Black
DHA4310	Red
DHA2222	No tag colour, select only with Option DHA4190 Bold Decal or DHA4197 Without Tag
DHA5000	75mm x 25mm (3" x 1") Supersport
DHA5050	100mm x 30mm (4" x 1.2") Everyday Black hub and tyre
DHA5100	100mm x 30mm (4" x 1.2") Comfort Black hub and tyre
DHA5150	100mm x 25mm (4" x 1") Starec Black hub and tyre
DHA5200	100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs)
DHA5250	100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs)
DHA6400	22"

DHA6410	24"
DHA6420	25"
DHA6430	26"
DHA6510	22", 24", 26" KÜschall Starec Wheel, 24" & 26" Black Cross Spokes, Rim & Hub. 22" Silver Cross-Spokes, Colourless hub and Black rim
DHA6520	24", 25" High Performance Wheel Black hub and rim
DHA6550	24", 25" Spinergy LXX, Black hub, rim and 18 spokes in Black
SWO0907	24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. Coloured spokes are available via MyKÜschall
SWO0909	24", 25" Spinergy LX, Chrome hub and 12 spokes in Black. Coloured spokes are available via MyKÜschall
DHA6600	24" Schwalbe Rightrun Black/Grey Pneumatic tyre
DHA6620	24", 25", 26" Primo Sentinel tyre. Pneumatic and puncture resistant
DHA6630	24", 25" Schwalbe One V-Guard Puncture Resistant tyre (Cut-Resistant High-Tech Fibre)
DHA6640	22" Profile Grey Pneumatic tyre 7 Bar
DHA6660	24" Solid Profile Black Tyre
DHA6700	22", 24", 26" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance
DHA6710	24", 26" Steel, short tab with 1 hole, mounted at 3 cm distance
DHA6720	24" Titanium, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance
DHA6740	24" , 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip
DHA6750	24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands
DHA6760	24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion
DHA6770	24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim
DHA6900	22", 24", 26" Transparent Spokeguards
DHA7000	Fixation Strap (Hook and loop strap under the seat holds the frame together when folded for easy transportation, designed for clients with limited hand function)
DHA7010	Fixation Clip (Plastic clip on either side of the seat holds the frame together when folded for easy transportation)
DHA7230	Performance Brakes, Push to Lock. Brake barrel flips away from tyre
DHA7300	Brake Lever Extension



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):