



Manual Retail Prescription Form

Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JULY 2021
Customer Ref:

| | |
|--------------------------|--------------------|
| Customer Ref* | Account No* |
| Contact No* | Order No* |
| Delivery Address* | Order Qty* |
| | Order date* |

Küschall Champion 2.0

Standard Features

- Max User Weight 120kg (19 Stone)
- Lightweight Folding Fully Hydroformed Aluminium Wheelchair
- Multiple adjustments, Front, Rear Seat Heights and Centre of Gravity
- Aluminium Front Frame with 75° or 90° Front Frame Angles
- 24" Self Propel Q/R Wheels with Grey Aluminium Handrims
- Schwalbe Black and Grey Pneumatic Tyres
- Choice of 7 Seat Widths 360mm to 480mm (14" - 19")
- Choice of 8 Seat Depths 340mm to 480mm (13" - 19")
- Choice of 5 Standard Frame Colours
- Range of Front Seat Heights 450mm to 540mm (17" -21")
- Range of Rear Seat Heights 390mm to 480mm (15" -19")
- Folding Backrest, Angle Adjustable and Soft Lock Feature
- Nylon Light Upholstery
- Angle and Height Adjustable Rigid Footrest, with Calf Strap
- Without Push Handles, Flip up Carbon Clothesguards
- Unfolding Help Cord
- Black Comfort Front Castor
- Standard Brake Lever Bent
- Tool Kit



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

| Code | Description | | |
|----------------------|--|-----------------------------|--------|
| BASE PRODUCT* | | | |
| DHA0031 | Küschall Champion 2.0 | <input type="radio"/> | |
| DHA0000 | MYKÜSCHALL Special Request - Please include quotation reference in special description at the end of the form where applicable | <input type="radio"/> | |
| FRAME* | | | |
| DHA1150 | Frame Angle 75° Aluminium | (NC 1350) | £2,320 |
| DHA1160 | Frame Angle 90° Aluminium | (NC 1350) | £2,320 |
| DHA1170 | Frame Angle 75° Aluminium - Tapered Front Frame | (NC 1320, 1325, 7000, 7010) | £2,320 |
| DHA1180 | Frame Angle 90° Aluminium - Tapered Front Frame | (NC 1320, 1325, 7000, 7010) | £2,320 |

| SEAT WIDTH* | | | | |
|--------------------|--|---|---|-----------------------|
| DHA1320 | | 360 mm | (NC 7000) | <input type="radio"/> |
| DHA1325 | | 380 mm | (NC 7000) | <input type="radio"/> |
| DHA1330 | | 400 mm | | <input type="radio"/> |
| DHA1335 | | 420 mm | | <input type="radio"/> |
| DHA1340 | | 440 mm | | <input type="radio"/> |
| DHA1345 | | 460 mm | | <input type="radio"/> |
| DHA1350 | | 480 mm | | <input type="radio"/> |
| SEAT DEPTH* | | | | |
| DHA1405 | | 340 mm | (NC 2630, 2720, 2730, 2750) | <input type="radio"/> |
| DHA1415 | | 360 mm | (NC 2720, 2730) | <input type="radio"/> |
| DHA1425 | | 380 mm | (NC 2720, 2730) | <input type="radio"/> |
| DHA1435 | | 400 mm | | <input type="radio"/> |
| DHA1445 | | 420 mm | | <input type="radio"/> |
| DHA1460 | | 440 mm | | <input type="radio"/> |
| DHA1470 | | 460 mm | | <input type="radio"/> |
| DHA1480 | | 480 mm | | <input type="radio"/> |
| BACK HEIGHT* | | | | |
| DHA1600 | | Back Height for 3rd party backrest. Final height depends on the height of the positioning backrest, select if Matrx/Fixed backrest is required. | (NC 2200, 2220, 2230) | <input type="radio"/> |
| DHA1615 | | 300 mm | (NC 2200, 2240, 2400, 2595) | <input type="radio"/> |
| DHA1620 | | 315 mm | (NC 2200, 2240, 2400, 2595) | <input type="radio"/> |
| DHA1625 | | 330 mm | (NC 2240, 2595) | <input type="radio"/> |
| DHA1630 | | 345 mm | (NC 2240, 2595) | <input type="radio"/> |
| DHA1635 | | 360 mm | (NC 2240, 2595) | <input type="radio"/> |
| DHA1640 | | 375 mm | (NC 2240, 2595) | <input type="radio"/> |
| DHA1645 | | 390 mm | (NC 2240, 2595) | <input type="radio"/> |
| DHA1650 | | 405 mm | (NC 2240) | <input type="radio"/> |
| DHA1655 | | 420 mm | (NC 2240) | <input type="radio"/> |
| DHA1660 | | 435 mm | (NC 2240) | <input type="radio"/> |
| DHA1665 | | 450 mm | (NC 2240) | <input type="radio"/> |
| DHA1670 | | 465 mm | (NC 2240) | <input type="radio"/> |
| FRONT SEAT HEIGHT* | | | | |
| DHA1740 | | 450 mm | (NC 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960) | <input type="radio"/> |
| DHA1745 | | 460 mm | (NC 1948, 1950, 1952, 1954, 1956, 1958, 1960) | <input type="radio"/> |
| DHA1750 | | 470 mm | (NC 1950, 1952, 1954, 1956, 1958, 1960) | <input type="radio"/> |
| DHA1755 | | 480 mm | (NC 1952, 1954, 1956, 1958, 1960) | <input type="radio"/> |
| DHA1760 | | 490 mm | (NC 1954, 1956, 1958, 1960) | <input type="radio"/> |
| DHA1765 | | 500 mm | (NC 1956, 1958, 1960) | <input type="radio"/> |
| DHA1770 | | 510 mm | (NC 1958, 1960) | <input type="radio"/> |
| DHA1775 | | 520 mm | (NC 1960) | <input type="radio"/> |
| DHA1780 | | 530 mm | (NC 5000) | <input type="radio"/> |
| DHA1785 | | 540 mm | (NC 5000, 5050, 5100, 5150, 5200, 5250) | <input type="radio"/> |

REAR SEAT HEIGHT. On 75° frame, Rear Seat Height can be 0mm - 70mm lower then front seat height and 90° can be 30mm to 100mm lower*

| | | | | |
|---------|--|--------|-----------------|---|
| DHA1815 | | 390 mm | (NC 6420, 6430) | ○ |
| DHA1820 | | 400 mm | (NC 6420, 6430) | ○ |
| DHA1825 | | 410 mm | (NC 6430) | ○ |
| DHA1830 | | 420 mm | (NC 6430) | ○ |
| DHA1835 | | 430 mm | | ○ |
| DHA1840 | | 440 mm | | ○ |
| DHA1845 | | 450 mm | (NC 6400) | ○ |
| DHA1850 | | 460 mm | (NC 6400) | ○ |
| DHA1855 | | 470 mm | (NC 6400, 6410) | ○ |
| DHA1860 | | 480 mm | (NC 6400, 6410) | ○ |

LOWER LEG LENGTH. Cushion height must be deducted*

| | | | | |
|---------|--|--------|--|---|
| DHA1924 | | 320 mm | | ○ |
| DHA1926 | | 330 mm | | ○ |
| DHA1928 | | 340 mm | | ○ |
| DHA1930 | | 350 mm | | ○ |
| DHA1932 | | 360 mm | | ○ |
| DHA1934 | | 370 mm | | ○ |
| DHA1936 | | 380 mm | | ○ |
| DHA1938 | | 390 mm | | ○ |
| DHA1940 | | 400 mm | | ○ |
| DHA1942 | | 410 mm | | ○ |
| DHA1944 | | 420 mm | | ○ |
| DHA1946 | | 430 mm | | ○ |
| DHA1948 | | 440 mm | | ○ |
| DHA1950 | | 450 mm | | ○ |
| DHA1952 | | 460 mm | | ○ |
| DHA1954 | | 470 mm | | ○ |
| DHA1956 | | 480 mm | | ○ |
| DHA1958 | | 490 mm | | ○ |
| DHA1960 | | 500 mm | | ○ |

BACKREST*

| | | | | |
|---------|--|---|--|---|
| DHA2140 | | Angle Adustable, Forward Folding with Soft Lock. Backrest remains locked for transportation | | ◆ |
|---------|--|---|--|---|

BACKREST UPHOLSTERY*

| | | | | |
|---------|--|--|-----------------|------|
| DHA2200 | | Black Standard Light Upholstery, without Tension Adjustment | (NC 2430, 2595) | ◆ |
| DHA2220 | | Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener | | £118 |
| DHA2230 | | Tension Adjustable Backrest, with Black Aero, Breathable Spacer Upholstery Cover, Back Pocket/hook and loop fastener | | £193 |
| DHA2240 | | Without Backrest Upholstery Select if fitting a Matrix/Fixed backrest, no upholstery required | (NC 2430, 2595) | ○ |

| PUSH HANDLES* | | | |
|-----------------------------------|--|---|------|
| DHA2400 | | Without Push Handles | ◆ |
| DHA2410 | | Push Handles Fixed Short/Mini Grip | £115 |
| DHA2420 | | Push Handles Fixed Long Grip | ○ |
| DHA2430 | | Push Handles Height Adjustable, Rear Set (NC 3480) | £219 |
| DHA2440 | | Push Handles Height Adjustable, Integrated in backrest (NC 2595) | £189 |
| DHA2450 | | Push Handles Foldable (NC 2595) | £171 |
| BACKREST ANGLE* | | | |
| DHA2510 | | 76° angle from the seat | ○ |
| DHA2530 | | 80.5° angle from the seat | ○ |
| DHA2550 | | 85° angle from the seat | ◆ |
| DHA2570 | | 89.5° angle from the seat | ○ |
| DHA2590 | | 94° angle from the seat | ○ |
| BACKREST STABILIZATION BAR | | | |
| DHA2595 | | Backrest Stabilization Bar, Swingaway and Removable | £142 |
| SEAT UPHOLSTERY* | | | |
| DHA2600 | | Standard Seat Upholstery | ◆ |
| DHA2630 | | Light Seat Upholstery, mesh material to reduce weight | £60 |
| SEAT CUSHION | | | |
| DHA2700 | | Standard Seat Cushion, 60 mm Double Foam with Visco Top Layer | £57 |
| DHA2720 | | Matrx Libra Cushion, Black Smooth Cover - Very High Risk | £338 |
| DHA2730 | | Matrx Vi Cushion, Black Smooth Cover - High Risk | £312 |
| DCA2750 | | Vicair Active O2 lightweight cushion 6cm deep | £495 |
| SEAT BAG | | | |
| DHA2800 | | Seat Bag, Click in and removable, includes the fixation clips | £64 |
| FOOTREST* | | | |
| DHA3240 | | Rigid, Plastic Footplate, Set to front, Height and Angle Adjustable. Note: This does not flip up for transfer, but folds easy for transport and storage | ◆ |
| DHA3250 | | Rigid, Plastic Footplate, Set to rear, Height and Angle Adjustable. Note: This does not flip up for transfer, but folds easy for transport and storage | ○ |
| CALF STRAP* | | | |
| DHA3320 | | Calf Strap Mounted on Supporter | ◆ |
| DHA3330 | | Calf Strap Mounted on Frame | ○ |
| CLOTHESGUARDS/MUDGUARDS* | | | |
| DHA3420 | | Carbon Clothesguards, Flip Up | ◆ |
| DHA3430 | | Carbon Mudguards, Fixed (NC 3490) | £214 |
| DHA3440 | | Carbon Mudguards, Removable (NC 3490) | £266 |
| ARMREST | | | |
| DHA3480 | | Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included | £208 |
| DHA3490 | | Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included | £92 |

| FRAME COLOUR* | | | |
|---------------|--|--|------|
| DHA4000 | | Space Black (Standard) | ○ |
| DHA4010 | | Anthracite Grey (Standard) | ○ |
| DHA4020 | | Pure White (Standard) (NC 4220) | ○ |
| DHA4030 | | Carmin Red (Standard) | ○ |
| DHA4040 | | Azure Blue (Standard) | ○ |
| DHA4060 | | Flashy Yellow (Additional 10 day lead time) | £26 |
| DHA4070 | | Atomic Violet (Additional 10 day lead time) | £26 |
| DHA4100 | | Black Glossy (Additional 10 day lead time) | £26 |
| DHA4130 | | Silver Grey (Additional 10 day lead time) | £26 |
| DHA4140 | | Royal Blue (Additional 10 day lead time) | £26 |
| DHA4150 | | Ruby Red (Additional 10 day lead time) | £26 |
| DHA4170 | | Dynamite Orange (Additional 10 day lead time) | £26 |
| DHA4180 | | Forest Green (Additional 10 day lead time) | £26 |
| DHA4050 | | Personal frame colour. Provide standard RAL colour code of your choice | £268 |

| FRAME DECAL TYPE* | | | |
|-------------------|--|---|-----|
| DHA4190 | | Bold Decal (Full coloured transfer on frame) (NC 4300, 4310, 1111) | £26 |
| DHA4195 | | Tag (Small Kuschall brand transfer on front of frame) (NC 4200, 4210, 4220, 4230, 4240, 4250, 2222) | ◆ |
| DHA4197 | | Without Tag (Kuschall branding only) (NC 4200, 4210, 4220, 4230, 4240, 4250, 4300, 4310) | ○ |

| BOLD DECAL COLOUR* | | | |
|--------------------|--|--|---|
| DHA4200 | | Black | ○ |
| DHA4210 | | Dark Grey | ○ |
| DHA4220 | | Light Grey | ○ |
| DHA4230 | | Red | ○ |
| DHA4240 | | Blue | ○ |
| DHA4250 | | Gold | ○ |
| DHA1111 | | No bold colour, select only with Option DHA4195 Tag or DHA4197 Without Tag (NC 4190) | ○ |

| TAG COLOUR* | | | |
|-------------|--|--|---|
| DHA4300 | | Black | ○ |
| DHA4310 | | Red | ○ |
| DHA2222 | | No tag colour, select only with Option DHA4190 Bold Decal or DHA4197 Without Tag (NC 4195) | ○ |

| FRONT CASTORS* | | | |
|---|--|---|------|
| DHA5000 | | 75mm x 25mm (3" x 1") Supersport | ○ |
| DHA5050 | | 100mm x 30mm (4" x 1.2") Everyday Black hub and tyre | ○ |
| DHA5100 | | 100mm x 30mm (4" x 1.2") Comfort Black hub and tyre | ○ |
| DHA5150 | | 100mm x 25mm (4" x 1") Starec Black hub and tyre | £72 |
| DHA5200 | | 100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs) | £80 |
| DHA5250 | | 100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs) | £80 |
| DHA5300 | | 125mm x 30mm (5" x 1.2") Everyday Black hub and tyre | ○ |
| DHA5350 | | 125mm x 30mm (5" x 1.2") Comfort Black hub and tyre | ◆ |
| DHA5400 | | 125mm x 25mm (5" x 1") Starec Black hub and tyre | £72 |
| DHA5450 | | 125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs) | £80 |
| DHA5500 | | 125mm x 36mm (5" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs) | £80 |
| REAR WHEEL POSITION (CENTRE OF GRAVITY)* | | | |
| DHA6100 | | 1st from front Most Active Position | ○ |
| DHA6110 | | 2nd from front | ○ |
| DHA6120 | | 3rd from front | ○ |
| DHA6130 | | 4th from front | ○ |
| DHA6140 | | 5th from front | ◆ |
| DHA6150 | | 6th from front | ○ |
| DHA6160 | | 7th from front | ○ |
| DHA6170 | | 8th from front | ○ |
| DHA6180 | | 9th from front | ○ |
| QUAD QUICK RELEASE AXLE | | | |
| DHA6300 | | Quad Quick Release/Tetra Clip (NC 0907, 0909, 6550) | £74 |
| REAR WHEEL SIZE* | | | |
| DHA6400 | | 22" (NC 6520, 6550, 6620, 6630, 6660, 6710, 6720, 6740, 6750, 6760, 6770) | ○ |
| DHA6410 | | 24" (NC 6510, 6640, 6660, 6700, 6710, 6720, 6760, 6770, 6900) | ◆ |
| DHA6420 | | 25" (NC 0907, 0909, 6520, 6550, 6630, 6640, 6660, 6720, 6760, 6770) | ○ |
| DHA6430 | | 26" (NC 0907, 0909, 6520, 6550, 6630, 6640, 6660, 6720, 6760, 6770) | ○ |
| REAR WHEEL TYPE* | | | |
| DHA6510 | | 22", 24", 26" Kuschall Starec Wheel, 24" & 26" Black Cross Spokes, Rim & Hub. 22" Silver Cross-Spokes, Colourless hub and Black rim | ◆ |
| DHA6520 | | 24", 25" High Performance Wheel Black hub and rim (NC 6660) | £157 |
| DHA6550 | | 24", 25" Spinergy L XK, Black hub, rim and spokes in Black (NC 6660) | £528 |
| SWO0907 | | 24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. Coloured spokes are available via MyKuschall (NC 6660) | £420 |
| SWO0909 | | 24", 25" Spinergy LX, Chrome hub and 12 spokes in Black. Coloured spokes are available via MyKuschall (NC 6660) | £635 |

| REAR WHEEL TYRE* | | | |
|---------------------------|--|--|------|
| DHA6600 | | 22", 24", 25", 26" Schwalbe Rightrun Black/Grey Pneumatic tyre | ◆ |
| DHA6620 | | 24", 25", 26" Schwalbe Marathon Plus Smart Guard (Puncture Protected tyre with 5mm layer of flexible rubber) | £86 |
| DHA6630 | | 24", 25" Schwalbe One V-Guard Puncture Resistant tyre (Cut-Resistant High-Tech Fibre) | £84 |
| DHA6640 | | 22", 24" Profile Grey Pneumatic tyre 7 Bar (NC 7230) | ○ |
| DHA6660 | | 24" Solid Profile Black Tyre (NC 7230) | £57 |
| HANDRIM* | | | |
| DHA6700 | | 24", 26" Supergripp, mounted at 3 cm distance | £193 |
| DHA6710 | | 24", 26" Steel, mounted at 3 cm distance | £54 |
| DHA6720 | | 24" Titanium, long tab with 2 holes for 3 cm or 4 cm mountings, mounted at 3cm distance | £166 |
| DHA6740 | | 24", 25" Carbolife Gekko mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip | £277 |
| DHA6750 | | 24", 25" Carbolife Gekko S, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands | £277 |
| DHA6760 | | 24" Carbolife Curve L, mounted at 3 cm distance. Unique shape for low-effort propulsion | £246 |
| DHA6770 | | 24" Carbolife Quadro, mounted at 3 cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim | £308 |
| DHA6780 | | 24", 25", 26" Light Aluminium Grey, long tab with 2 holes for 3 cm or 4cm mountings, mounted at 3cm distance | ◆ |
| DHA6790 | | 24", 25", 26" Light Aluminium Black, mounted at 3 cm distance | £36 |
| REAR WHEEL CAMBER* | | | |
| DHA6820 | | 0° Camber | ◆ |
| DHA6830 | | 1° Camber | ○ |
| DHA6840 | | 3° Camber (NC 7230) | ○ |
| SPOKEGUARD | | | |
| DHA6900 | | 22", 24", 26" Transparent Spokeguards | £109 |
| FIXATION CLIP* | | | |
| DHA7000 | | Fixation Strap (Hook and loop strap under the seat holds the frame together when folded for easy transportation, designed for clients with limited hand function) | ○ |
| DHA7010 | | Fixation Clip (Plastic clip on either side of the seat holds the frame together when folded for easy transportation) | ◆ |
| DHA3333 | | No fixation clip, select with Tapered frame options 1170, 1180 | ○ |
| UNFOLDING HELP* | | | |
| DHA7100 | | Unfolding Help (Cord to help open the folding mechanism from the front of the chair) | ◆ |
| BRAKES* | | | |
| DHA7200 | | Pull Brakes | £67 |
| DHA7210 | | Standard Push to Lock Brakes, Lever Bent | ◆ |
| DHA7220 | | Standard Push to Lock Brakes, Lever Straight | ○ |
| DHA7230 | | Performance Brakes, Push to Lock. Brake barrel flips away from tyre (NC 7300) | £67 |
| DHA7240 | | Active/Scissor Brakes (NC 7300) | £67 |
| DHA7250 | | Light/Compact Active, Aluminium Brakes (NC 7300) | £180 |

| BRAKE LEVER EXTENSION | | |
|-----------------------|--------------------------|---|
| DHA7300 | <input type="checkbox"/> | Brake Lever Extension £43 |
| SAFETY OPTIONS | | |
| DHA7500 | <input type="checkbox"/> | Anti-Tipper Right £106 |
| DHA7510 | <input type="checkbox"/> | Anti-Tipper Left £106 |
| DHA7520 | <input type="checkbox"/> | Anti-Tipper Both sides £211 |
| DHA7800 | <input type="checkbox"/> | Standard Pelvic belt £64 |
| OTHER OPTIONS | | |
| DHA8010 | <input type="checkbox"/> | Front Frame Protectors £32 |
| DHA8100 | <input type="checkbox"/> | Tipper-Aid, Mounted Right Side Max. 2 adaptations per side (Adaptations = Anti tipper, Tipper aid and Cane holder) £40 |
| DHA8200 | <input type="checkbox"/> | Cane Holder, Mounted Right Side Max. 2 adaptations per side (Adaptations = Anti tipper, Tipper aid and cane holder) £65 |
| DHA8500 | <input type="checkbox"/> | Air Pump £11 |
| DHA8550 | <input type="checkbox"/> | Tool Kit ◆ |
| DHA8680 | <input type="checkbox"/> | 24" Mountain Bike Wheels with bag. Fitted with Specific Black Tyres, Grey Aluminium Handrims; Quick Release Axles Included £277 |
| DHA8700 | <input type="checkbox"/> | Travel bag for Champion, Wheels can remain on the wheelchair £219 |

| | |
|--------------------|--|
| TOTAL PRICE | |
|--------------------|--|

Non Conforming Options:

| | |
|---------|--|
| DHA1320 | 360 mm |
| DHA1325 | 380 mm |
| DHA1350 | 480 mm |
| DHA1946 | 430 mm |
| DHA1948 | 440 mm |
| DHA1950 | 450 mm |
| DHA1952 | 460 mm |
| DHA1954 | 470 mm |
| DHA1956 | 480 mm |
| DHA1958 | 490 mm |
| DHA1960 | 500 mm |
| DHA2200 | Black Standard Light Upholstery, without Tension Adjustment |
| DHA2220 | Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener |
| DHA2230 | Tension Adjustable Backrest, with Black Aero, Breathable Spacer Upholstery Cover, Back Pocket/hook and loop fastener |
| DHA2240 | Without Backrest Upholstery Select if fitting a Matrix/Fixed backrest, no upholstery required |
| DHA2400 | Without Push Handles |
| DHA2430 | Push Handles Height Adjustable, Rear Set |
| DHA2595 | Backrest Stabilization Bar, Swingaway and Removable |
| DHA2630 | Light Seat Upholstery, mesh material to reduce weight |
| DHA2720 | Matrix Libra Cushion, Black Smooth Cover - Very High Risk |
| DHA2730 | Matrix Vi Cushion, Black Smooth Cover - High Risk |
| DCA2750 | Vicair Active O2 lightweight cushion 6cm deep |
| DHA3480 | Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included |
| DHA3490 | Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included |
| DHA4190 | Bold Decal (Full coloured transfer on frame) |
| DHA4195 | Tag (Small Kuschall brand transfer on front of frame) |
| DHA4200 | Black |
| DHA4210 | Dark Grey |
| DHA4220 | Light Grey |
| DHA4230 | Red |
| DHA4240 | Blue |
| DHA4250 | Gold |
| DHA1111 | No bold colour, select only with Option DHA4195 Tag or DHA4197 Without Tag |
| DHA4300 | Black |
| DHA4310 | Red |
| DHA2222 | No tag colour, select only with Option DHA4190 Bold Decal or DHA4197 Without Tag |
| DHA5000 | 75mm x 25mm (3" x 1") Supersport |
| DHA5050 | 100mm x 30mm (4" x 1.2") Everyday Black hub and tyre |
| DHA5100 | 100mm x 30mm (4" x 1.2") Comfort Black hub and tyre |
| DHA5150 | 100mm x 25mm (4" x 1") Starec Black hub and tyre |
| DHA5200 | 100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs) |
| DHA5250 | 100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs) |
| DHA6400 | 22" |

| | |
|---------|--|
| DHA6410 | 24" |
| DHA6420 | 25" |
| DHA6430 | 26" |
| DHA6510 | 22", 24", 26" KÜschall Starec Wheel, 24" & 26" Black Cross Spokes, Rim & Hub. 22" Silver Cross-Spokes, Colourless hub and Black rim |
| DHA6520 | 24", 25" High Performance Wheel Black hub and rim |
| DHA6550 | 24", 25" Spinergy LXX, Black hub, rim and spokes in Black |
| SWO0907 | 24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. Coloured spokes are available via MyKÜschall |
| SWO0909 | 24", 25" Spinergy LX, Chrome hub and 12 spokes in Black. Coloured spokes are available via MyKÜschall |
| DHA6620 | 24", 25", 26" Schwalbe Marathon Plus Smart Guard (Puncture Protected tyre with 5mm layer of flexible rubber) |
| DHA6630 | 24", 25" Schwalbe One V-Guard Puncture Resistant tyre (Cut-Resistant High-Tech Fibre) |
| DHA6640 | 22", 24" Profile Grey Pneumatic tyre 7 Bar |
| DHA6660 | 24" Solid Profile Black Tyre |
| DHA6700 | 24", 26" Supergripp, mounted at 3 cm distance |
| DHA6710 | 24", 26" Steel, mounted at 3 cm distance |
| DHA6720 | 24" Titanium, long tab with 2 holes for 3 cm or 4 cm mountings, mounted at 3cm distance |
| DHA6740 | 24", 25" Carbolife Gekko mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip |
| DHA6750 | 24", 25" Carbolife Gekko S, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands |
| DHA6760 | 24" Carbolife Curve L, mounted at 3 cm distance. Unique shape for low-effort propulsion |
| DHA6770 | 24" Carbolife Quadro, mounted at 3 cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim |
| DHA6900 | 22", 24", 26" Transparent Spokeguards |
| DHA7000 | Fixation Strap (Hook and loop strap under the seat holds the frame together when folded for easy transportation, designed for clients with limited hand function) |
| DHA7010 | Fixation Clip (Plastic clip on either side of the seat holds the frame together when folded for easy transportation) |
| DHA7230 | Performance Brakes, Push to Lock. Brake barrel flips away from tyre |
| DHA7300 | Brake Lever Extension |



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):