



**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

JANUARY 2022  
Customer Ref:

|                          |                    |
|--------------------------|--------------------|
| <b>Customer Ref*</b>     | <b>Account No*</b> |
| <b>Contact No*</b>       | <b>Order No*</b>   |
| <b>Delivery Address*</b> | <b>Order Qty*</b>  |
|                          | <b>Order date*</b> |

## Küschall Champion SK 2.0

### Standard Features

- Max User Weight 120kg (19 Stone)
- Active, Lightweight Folding Fully Hydroformed Aluminium Wheelchair
- Multiple adjustments, Front and Rear Seat Heights and Centre of Gravity
- Aluminium Folding Front Frame with 75° or 90° Front Frame Angles
- 24" Self Propel Q/R Wheels with Grey Aluminium Handrims and Schwalbe Black and Grey Pneumatic Tyres
- Choice of 6 Seat Widths 360mm to 460mm (14" - 18")
- Choice of 8 Seat Depths 340mm to 480mm (13" - 19")
- Choice of 5 Standard Frame Colours
- Range of Front Seat Heights 450mm to 540mm (17" -21")
- Range of Rear Seat Heights 390mm to 480mm (15" -19")
- Folding Backrest, Angle Adjustable and Soft Lock Feature with Nylon Light Upholstery
- Angle and Height Adjustable Rigid Footrest, with Calf Strap
- Without Push Handles
- Flip up Carbon Clothesguards
- Unfolding Help Cord
- Black Comfort Front Castor
- Standard Brake Lever Bent
- Tool Kit



**Standard Base price: £2,977**

Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

| Code                 | Description   |                       |
|----------------------|---|-----------------------|
| <b>BASE PRODUCT*</b> |   |                       |
| DHA0031              | <input type="checkbox"/> Küschall Champion SK 2.0   | <input type="radio"/> |
| DHA0000              | <input type="checkbox"/> MYKÜSCHALL Special Request | <input type="radio"/> |
| <b>FRAME*</b>        |   |                       |
| DHA1150              | <input type="checkbox"/> Frame Angle 75° Aluminium  | <input type="radio"/> |
| DHA1160              | <input type="checkbox"/> Frame Angle 90° Aluminium  | <input type="radio"/> |

| <b>SK FOLDING FRONT*</b>  |  |   |   |
|---------------------------|--|---|---|
| DHA1250                   |  | SK Folding Front Frame  | £2,977  |
| <b>SEAT WIDTH*</b>        |  |   |   |
| DHA1320                   |  | 360 mm  | (NC 7000) ○   |
| DHA1325                   |  | 380 mm  | (NC 7000) ○   |
| DHA1330                   |  | 400 mm  | ○   |
| DHA1335                   |  | 420 mm  | ○   |
| DHA1340                   |  | 440 mm  | ○   |
| DHA1345                   |  | 460 mm  | ○   |
| <b>SEAT DEPTH*</b>        |  |   |   |
| DHA1405                   |  | 340 mm  | (NC 2630, 2720, 2730, 2750) ○                         |
| DHA1415                   |  | 360 mm  | (NC 2720, 2730) ○                                     |
| DHA1425                   |  | 380 mm  | (NC 2720, 2730) ○                                     |
| DHA1435                   |  | 400 mm  | ○   |
| DHA1445                   |  | 420 mm  | ○   |
| DHA1460                   |  | 440 mm  | ○   |
| DHA1470                   |  | 460 mm  | ○   |
| DHA1480                   |  | 480 mm  | ○   |
| <b>BACK HEIGHT*</b>       |  |   |   |
| DHA1600                   |  | Back Height for 3rd party backrest. Final height depends on the height of the positioning backrest, select if Matrx/Fixed backrest is required. | (NC 2200, 2220, 2230) ◆                               |
| DHA1615                   |  | 300 mm  | (NC 2200, 2240, 2400, 2595) ○                         |
| DHA1620                   |  | 315 mm  | (NC 2200, 2240, 2400, 2595) ○                         |
| DHA1625                   |  | 330 mm  | (NC 2240, 2595) ○                                     |
| DHA1630                   |  | 345 mm  | (NC 2240, 2595) ○                                     |
| DHA1635                   |  | 360 mm  | (NC 2240, 2595) ○                                     |
| DHA1640                   |  | 375 mm  | (NC 2240, 2595) ○                                     |
| DHA1645                   |  | 390 mm  | (NC 2240, 2595) ○                                     |
| DHA1650                   |  | 405 mm  | (NC 2240) ○   |
| DHA1655                   |  | 420 mm  | (NC 2240) ○   |
| DHA1660                   |  | 435 mm  | (NC 2240) ○   |
| DHA1665                   |  | 450 mm  | (NC 2240) ○   |
| DHA1670                   |  | 465 mm  | (NC 2240) ○   |
| <b>FRONT SEAT HEIGHT*</b> |  |   |   |
| DHA1740                   |  | 450 mm  | (NC 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960) ○ |
| DHA1745                   |  | 460 mm  | (NC 1948, 1950, 1952, 1954, 1956, 1958, 1960) ○       |
| DHA1750                   |  | 470 mm  | (NC 1950, 1952, 1954, 1956, 1958, 1960) ○             |
| DHA1755                   |  | 480 mm  | (NC 1952, 1954, 1956, 1958, 1960) ○                   |
| DHA1760                   |  | 490 mm  | (NC 1954, 1956, 1958, 1960) ○                         |
| DHA1765                   |  | 500 mm  | (NC 1956, 1958, 1960) ○                               |
| DHA1770                   |  | 510 mm  | (NC 1958, 1960) ○                                     |
| DHA1775                   |  | 520 mm  | (NC 1960) ○   |
| DHA1780                   |  | 530 mm  | (NC 5000) ○   |
| DHA1785                   |  | 540 mm  | (NC 5000, 5050, 5100, 5150, 5200, 5250) ○             |

**REAR SEAT HEIGHT. On 75° frame, Rear Seat Height can be 0mm - 70mm lower then front seat height and 90° can be 30mm to 100mm lower\***

|         |  |        |                 |   |
|---------|--|--------|-----------------|---|
| DHA1815 |  | 390 mm | (NC 6420, 6430) | ○ |
| DHA1820 |  | 400 mm | (NC 6420, 6430) | ○ |
| DHA1825 |  | 410 mm | (NC 6430)       | ○ |
| DHA1830 |  | 420 mm | (NC 6430)       | ○ |
| DHA1835 |  | 430 mm |                 | ○ |
| DHA1840 |  | 440 mm |                 | ○ |
| DHA1845 |  | 450 mm | (NC 6400)       | ○ |
| DHA1850 |  | 460 mm | (NC 6400)       | ○ |
| DHA1855 |  | 470 mm | (NC 6400, 6410) | ○ |
| DHA1860 |  | 480 mm | (NC 6400, 6410) | ○ |

**LOWER LEG LENGTH. Cushion height must be deducted\***

|         |  |        |  |   |
|---------|--|--------|--|---|
| DHA1936 |  | 380 mm |  | ○ |
| DHA1938 |  | 390 mm |  | ○ |
| DHA1940 |  | 400 mm |  | ○ |
| DHA1942 |  | 410 mm |  | ○ |
| DHA1944 |  | 420 mm |  | ○ |
| DHA1946 |  | 430 mm |  | ○ |
| DHA1948 |  | 440 mm |  | ○ |
| DHA1950 |  | 450 mm |  | ○ |
| DHA1952 |  | 460 mm |  | ○ |
| DHA1954 |  | 470 mm |  | ○ |
| DHA1956 |  | 480 mm |  | ○ |
| DHA1958 |  | 490 mm |  | ○ |
| DHA1960 |  | 500 mm |  | ○ |

**BACKREST\***

|         |  |   |  |   |
|---------|--|---|--|---|
| DHA2140 |  | Angle Adustable, Forward Folding with Soft Lock. Backrest remains locked for transportation |  | ◆ |
|---------|--|---|--|---|

**BACKREST UPHOLSTERY\***

|         |  |  |                 |      |
|---------|--|--|-----------------|------|
| DHA2200 |  | Black Standard Light Upholstery, without Tension Adjustment  | (NC 2430, 2595) | ◆    |
| DHA2220 |  | Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener |                 | £130 |
| DHA2230 |  | Tension Adjustable Backrest, with Black Aero, Breathable Spacer Upholstery Cover, Back Pocket/hook and loop fastener           |                 | £212 |
| DHA2240 |  | Without Backrest Upholstery Select if fitting a Matrix/Fixed backrest, no upholstery required                                  | (NC 2430, 2595) | ○    |

**PUSH HANDLES\***

|         |  |  |           |      |
|---------|--|--|-----------|------|
| DHA2400 |  | Without Push Handles                                   |           | ◆    |
| DHA2410 |  | Push Handles Fixed Short/Mini Grip                     |           | £126 |
| DHA2420 |  | Push Handles Fixed Long Grip                           |           | ○    |
| DHA2430 |  | Push Handles Height Adjustable, Rear Set               | (NC 3480) | £241 |
| DHA2440 |  | Push Handles Height Adjustable, Integrated in backrest | (NC 2595) | £208 |
| DHA2450 |  | Push Handles Foldable                                  | (NC 2595) | £188 |

| <b>BACKREST ANGLE*</b>            |  |   |      |
|-----------------------------------|--|---|------|
| DHA2510                           |  | 76° angle from the seat   | ○    |
| DHA2530                           |  | 80.5° angle from the seat   | ○    |
| DHA2550                           |  | 85° angle from the seat   | ◆    |
| DHA2570                           |  | 89.5° angle from the seat   | ○    |
| DHA2590                           |  | 94° angle from the seat   | ○    |
| <b>BACKREST STABILIZATION BAR</b> |  |   |      |
| DHA2595                           |  | Backrest Stabilization Bar, Swingaway and Removable   | £156 |
| <b>SEAT UPHOLSTERY*</b>           |  |   |      |
| DHA2600                           |  | Standard Seat Upholstery  | ◆    |
| DHA2630                           |  | Light Seat Upholstery, mesh material to reduce weight   | £66  |
| <b>SEAT CUSHION</b>               |  |   |      |
| DHA2700                           |  | Standard Seat Cushion, 60 mm Double Foam with Visco Top Layer   | £63  |
| DHA2720                           |  | Matrx Libra Cushion, Black Smooth Cover - Very High Risk  | £338 |
| DHA2730                           |  | Matrx Vi Cushion, Black Smooth Cover - High Risk  | £312 |
| DCA2750                           |  | Vicair Active O2 lightweight cushion 6cm deep   | £495 |
| <b>SEAT BAG</b>                   |  |   |      |
| DHA2800                           |  | Seat Bag, Click in and removable, includes the fixation clips   | £70  |
| <b>FOOTREST*</b>                  |  |   |      |
| DHA3240                           |  | Rigid, Plastic Footplate, Set to front, Height and Angle Adjustable. Note: This does not flip up for transfer, but folds easy for transport and storage | ◆    |
| DHA3250                           |  | Rigid, Plastic Footplate, Set to rear, Height and Angle Adjustable. Note: This does not flip up for transfer, but folds easy for transport and storage  | ○    |
| <b>CALF STRAP*</b>                |  |   |      |
| DHA3320                           |  | Calf Strap Mounted on Supporter   | ◆    |
| DHA3330                           |  | Calf Strap Mounted on Frame   | ○    |
| <b>CLOTHESGUARDS/MUDGUARDS*</b>   |  |   |      |
| DHA3420                           |  | Carbon Clothesguards, Flip Back   | ◆    |
| DHA3430                           |  | Carbon Mudguards, Fixed (NC 3490)   | £235 |
| DHA3440                           |  | Carbon Mudguards, Removable (NC 3490)   | £292 |
| <b>ARMREST</b>                    |  |   |      |
| DHA3480                           |  | Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included                                     | £229 |
| DHA3490                           |  | Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included   | £101 |

| FRAME COLOUR* |  |  |      |
|---------------|--|--|------|
| DHA4000       |  | Space Black (Standard)   | ○    |
| DHA4010       |  | Anthracite Grey (Standard)   | ○    |
| DHA4020       |  | Pure White (Standard)  | ○    |
| DHA4030       |  | Carmine Red (Standard)   | ○    |
| DHA4040       |  | Azure Blue (Standard)  | ○    |
| DHA4060       |  | Flashy Yellow (Additional 10 day lead time)                            | £29  |
| DHA4070       |  | Atomic Violet (Additional 10 day lead time)                            | £29  |
| DHA4100       |  | Black Glossy (Additional 10 day lead time)                             | £29  |
| DHA4130       |  | Silver Grey (Additional 10 day lead time)                              | £29  |
| DHA4140       |  | Royal Blue (Additional 10 day lead time)                               | £29  |
| DHA4150       |  | Ruby Red (Additional 10 day lead time)                                 | £29  |
| DHA4170       |  | Dynamite Orange (Additional 10 day lead time)                          | £29  |
| DHA4180       |  | Forest Green (Additional 10 day lead time)                             | £29  |
| DHA4050       |  | Personal frame colour. Provide standard RAL colour code of your choice | £295 |

| FRAME DECAL TYPE* |  |   |                   |
|-------------------|--|---|-------------------|
| DHA4195           |  | Tag (Small Kuschall brand transfer on front of frame) | (NC 2222) ◆       |
| DHA4197           |  | Without Tag (Kuschall branding only)                  | (NC 4300, 4310) ○ |

| TAG COLOUR* |  |  |             |
|-------------|--|--|-------------|
| DHA4300     |  | Black  | ○           |
| DHA4310     |  | Red  | ○           |
| DHA2222     |  | No tag colour, select only with Option DHA4197 Without Tag | (NC 4195) ○ |

| FRONT CASTORS* |  |  |     |
|----------------|--|--|-----|
| DHA5000        |  | 75mm x 25mm (3" x 1") Supersport   | ○   |
| DHA5050        |  | 100mm x 30mm (4" x 1.2") Everyday Black hub and tyre                           | ○   |
| DHA5100        |  | 100mm x 30mm (4" x 1.2") Comfort Black hub and tyre                            | ○   |
| DHA5150        |  | 100mm x 25mm (4" x 1") Starec Black hub and tyre                               | £79 |
| DHA5200        |  | 100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs)     | £88 |
| DHA5250        |  | 100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs) | £88 |
| DHA5300        |  | 125mm x 30mm (5" x 1.2") Everyday Black hub and tyre                           | ○   |
| DHA5350        |  | 125mm x 30mm (5" x 1.2") Comfort Black hub and tyre                            | ◆   |
| DHA5400        |  | 125mm x 25mm (5" x 1") Starec Black hub and tyre                               | £79 |
| DHA5450        |  | 125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs)     | £88 |
| DHA5500        |  | 125mm x 36mm (5" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs) | £88 |

| <b>REAR WHEEL POSITION (CENTRE OF GRAVITY)*</b> |  |   |      |
|---|--|---|------|
| DHA6100   |  | 1st from front Most Active Position   | ○    |
| DHA6110   |  | 2nd from front  | ○    |
| DHA6120   |  | 3rd from front  | ○    |
| DHA6130   |  | 4th from front  | ○    |
| DHA6140   |  | 5th from front  | ◆    |
| DHA6150   |  | 6th from front  | ○    |
| DHA6160   |  | 7th from front  | ○    |
| DHA6170   |  | 8th from front  | ○    |
| DHA6180   |  | 9th from front  | ○    |
| <b>QUAD QUICK RELEASE AXLE</b>                  |  |   |      |
| DHA6300   |  | Quad Quick Release/Tetra Clip (NC 6550)   | £81  |
| <b>REAR WHEEL SIZE*</b>                         |  |   |      |
| DHA6400   |  | 22" (NC 6520, 6550, 6600, 6620, 6630, 6660, 6710, 6720, 6760, 6770)   | ○    |
| DHA6410   |  | 24" (NC 6640)   | ◆    |
| DHA6420   |  | 25" (NC 6510, 6600, 6640, 6700, 6710, 6720, 6760, 6770, 6900)   | ○    |
| DHA6430   |  | 26" (NC 6520, 6550, 6600, 6630, 6640, 6720, 6740, 6750, 6760, 6770)   | ○    |
| <b>REAR WHEEL TYPE*</b>                         |  |   |      |
| DHA6510   |  | 22", 24", 26" Kuschall Starec Wheel, 24" & 26" Black Cross Spokes, Rim & Hub. 22" Silver Cross-Spokes, Colourless hub and Black rim | ◆    |
| DHA6520   |  | 24", 25" High Performance Wheel Black hub and rim (NC 6660)   | £173 |
| DHA6550   |  | 24", 25" Spinergy LXX, Black hub, rim and 18 spokes in Black (NC 6660)  | £580 |
| SWO0907   |  | 24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. Coloured spokes are available via MyKuschall (NC 6660)                   | £462 |
| SWO0909   |  | 24", 25" Spinergy LX, Chrome hub and 12 spokes in Black. Coloured spokes are available via MyKuschall (NC 6660)                     | £698 |
| <b>REAR WHEEL TYRE*</b>                         |  |   |      |
| DHA6600   |  | 24" Schwalbe Rightrun Black/Grey Pneumatic tyre   | ◆    |
| DHA6620   |  | 24", 25", 26" Primo Sentinel tyre. Pneumatic and puncture resistant   | £95  |
| DHA6630   |  | 24", 25" Schwalbe One V-Guard Puncture Resistant tyre (Cut-Resistant High-Tech Fibre)   | £92  |
| DHA6640   |  | 22" Profile Grey Pneumatic tyre 7 Bar (NC 7230)   | ○    |
| DHA6660   |  | 24" Solid Profile Black Tyre (NC 7230)  | £63  |

| <b>HANDRIM*</b>              |  |   |      |
|------------------------------|--|---|------|
| DHA6700                      |  | 22", 24", 26" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance   | £212 |
| DHA6710                      |  | 24", 26" Steel, short tab with 1 hole, mounted at 3 cm distance   | £59  |
| DHA6720                      |  | 24" Titanium, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance  | £182 |
| DHA6740                      |  | 24" , 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip                                    | £304 |
| DHA6750                      |  | 24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands | £304 |
| DHA6760                      |  | 24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion  | £270 |
| DHA6770                      |  | 24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim  | £338 |
| DHA6780                      |  | 22", 24", 25", 26" Light Aluminium Grey, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance   | ◆    |
| DHA6790                      |  | 22", 24", 25", 26" Light Aluminium Black, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3 cm distance   | £40  |
| <b>REAR WHEEL CAMBER*</b>    |  |   |      |
| DHA6820                      |  | 0° Camber   | ◆    |
| DHA6830                      |  | 1° Camber   | ○    |
| DHA6840                      |  | 3° Camber (NC 7230)   | ○    |
| <b>SPOKEGUARD</b>            |  |   |      |
| DHA6900                      |  | 22", 24", 26" Transparent   | £120 |
| <b>FIXATION CLIP*</b>        |  |   |      |
| DHA7000                      |  | Fixation Strap (Hook and loop strap under the seat holds the frame together when folded for easy transportation, designed for clients with limited hand function)   | ○    |
| DHA7010                      |  | Fixation Clip (Plastic clip on either side of the seat holds the frame together when folded for easy transportation)  | ◆    |
| DHA3333                      |  | No fixation clip, select with Tapered frame options 1170, 1180  | ○    |
| <b>UNFOLDING HELP*</b>       |  |   |      |
| DHA7100                      |  | Unfolding Help (Cord to help open the folding mechanism from the front of the chair)  | ◆    |
| <b>BRAKES*</b>               |  |   |      |
| DHA7200                      |  | Pull Brakes   | £74  |
| DHA7210                      |  | Standard Push to Lock Brakes, Lever Bent  | ◆    |
| DHA7220                      |  | Standard Push to Lock Brakes, Lever Straight  | ○    |
| DHA7230                      |  | Performance Brakes, Push to Lock. Brake barrel flips away from tyre (NC 7300)   | £74  |
| DHA7240                      |  | Active/Scissor Brakes (NC 7300)   | £74  |
| DHA7250                      |  | Light/Compact Active, Aluminium Brakes (NC 7300)  | £129 |
| <b>BRAKE LEVER EXTENSION</b> |  |   |      |
| DHA7300                      |  | Brake Lever Extension   | £47  |

| SAFETY OPTIONS |  |  |      |
|----------------|--|--|------|
| DHA7500        |  | Anti-Tipper Right  | £116 |
| DHA7510        |  | Anti-Tipper Left   | £116 |
| DHA7520        |  | Anti-Tipper Both sides   | £232 |
| DHA7800        |  | Standard Pelvic belt   | £70  |
| OTHER OPTIONS  |  |  |      |
| DHA8010        |  | Front Frame Protectors; 2 Piece. Removable with hook and loop  | £35  |
| DHA8100        |  | Tipper-Aid, Mounted Right Side Max. 2 adaptations per side (Adaptations = Anti tipper, Tipper aid and Cane holder)         | £44  |
| DHA8200        |  | Cane Holder, Mounted Right Side Max. 2 adaptations per side (Adaptations = Anti tipper, Tipper aid and cane holder)        | £71  |
| DHA8500        |  | Air Pump   | £12  |
| DHA8550        |  | Tool Kit   | ◆    |
| DHA8680        |  | 24" Mountain Bike Wheels with bag. Fitted with Specific Black Tyres, Grey Aluminium Handrims; Quick Release Axles Included | £304 |
| DHA8750        |  | Travel bag for Champion SK,  | £354 |

|                    |  |
|--------------------|--|
| <b>TOTAL PRICE</b> |  |
|--------------------|--|



**Non Conforming Options:**

|         |   |
|---------|---|
| DHA1946 | 430 mm  |
| DHA1948 | 440 mm  |
| DHA1950 | 450 mm  |
| DHA1952 | 460 mm  |
| DHA1954 | 470 mm  |
| DHA1956 | 480 mm  |
| DHA1958 | 490 mm  |
| DHA1960 | 500 mm  |
| DHA2200 | Black Standard Light Upholstery, without Tension Adjustment   |
| DHA2220 | Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener      |
| DHA2230 | Tension Adjustable Backrest, with Black Aero, Breathable Spacer Upholstery Cover, Back Pocket/hook and loop fastener                |
| DHA2240 | Without Backrest Upholstery Select if fitting a Matrix/Fixed backrest, no upholstery required                                       |
| DHA2400 | Without Push Handles  |
| DHA2430 | Push Handles Height Adjustable, Rear Set  |
| DHA2595 | Backrest Stabilization Bar, Swingaway and Removable   |
| DHA2630 | Light Seat Upholstery, mesh material to reduce weight   |
| DHA2720 | Matrix Libra Cushion, Black Smooth Cover - Very High Risk   |
| DHA2730 | Matrix Vi Cushion, Black Smooth Cover - High Risk   |
| DCA2750 | Vicair Active O2 lightweight cushion 6cm deep   |
| DHA3480 | Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included                 |
| DHA3490 | Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included   |
| DHA4195 | Tag (Small Kuschall brand transfer on front of frame)   |
| DHA4300 | Black   |
| DHA4310 | Red   |
| DHA2222 | No tag colour, select only with Option DHA4197 Without Tag  |
| DHA5000 | 75mm x 25mm (3" x 1") Supersport  |
| DHA5050 | 100mm x 30mm (4" x 1.2") Everyday Black hub and tyre  |
| DHA5100 | 100mm x 30mm (4" x 1.2") Comfort Black hub and tyre   |
| DHA5150 | 100mm x 25mm (4" x 1") Starec Black hub and tyre  |
| DHA5200 | 100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs)  |
| DHA5250 | 100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs)  |
| DHA6400 | 22"   |
| DHA6410 | 24"   |
| DHA6420 | 25"   |
| DHA6430 | 26"   |
| DHA6510 | 22", 24", 26" Kuschall Starec Wheel, 24" & 26" Black Cross Spokes, Rim & Hub. 22" Silver Cross-Spokes, Colourless hub and Black rim |
| DHA6520 | 24", 25" High Performance Wheel Black hub and rim   |
| DHA6550 | 24", 25" Spinergy LXX, Black hub, rim and 18 spokes in Black  |
| DHA6600 | 24" Schwalbe Rightrun Black/Grey Pneumatic tyre   |
| DHA6620 | 24", 25", 26" Primo Sentinel tyre. Pneumatic and puncture resistant   |
| DHA6630 | 24", 25" Schwalbe One V-Guard Puncture Resistant tyre (Cut-Resistant High-Tech Fibre)   |
| DHA6640 | 22" Profile Grey Pneumatic tyre 7 Bar   |

|         |   |
|---------|---|
| DHA6660 | 24" Solid Profile Black Tyre  |
| DHA6700 | 22", 24", 26" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance   |
| DHA6710 | 24", 26" Steel, short tab with 1 hole, mounted at 3 cm distance   |
| DHA6720 | 24" Titanium, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance  |
| DHA6740 | 24" , 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip                                    |
| DHA6750 | 24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands |
| DHA6760 | 24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion  |
| DHA6770 | 24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim  |
| DHA6900 | 22", 24", 26" Transparent   |
| DHA7000 | Fixation Strap (Hook and loop strap under the seat holds the frame together when folded for easy transportation, designed for clients with limited hand function)   |
| DHA7230 | Performance Brakes, Push to Lock. Brake barrel flips away from tyre   |
| DHA7300 | Brake Lever Extension   |



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):