

Essential Plus

Quality Foam Mattress



A cost effective solution for comfort
and pressure redistribution



The *Invacare Essential Plus* is an entry level pressure ulcer prevention mattress designed for the community healthcare environment. Practical and durable, this mattress offers comfort, support and pressure redistribution at a cost effective price.

Effective Pressure Ulcer Prevention

- ▶ The core of the *Invacare Essential Plus* is constructed from a series of 'keyhole cut' cells, which move independently to evenly distribute weight, promoting greater comfort and facilitating pressure ulcer prevention.

Cover is resistant to fluid ingress

- ▶ The outer cover is a multi stretch, water resistant, polyurethane fabric providing excellent recovery and durability, which helps to reduce shear and friction forces. The cover also has a white inner liner allowing easy inspection for strikethrough or other cover damage.



Zip around three sides

- ▶ Allows easy inspection of the foam core and easy removal for laundering.

Turn and rotate

- ▶ Double sided castellated foam allows the mattress to be turned or rotated for increased product longevity.

Design

- ▶ The *Essential Plus* mattress has an innovative and modern cover design with the simple wave giving the mattress a fresh and current feel.

Technical data ▼

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.



Size range ▼



Max. user weight ▼



Max. cleaning temperature ▼

Nominal density range (kg/m³) ▼

Article no. ▼

Essential Plus
(1950 mm length)

1950 x 880 x 140 mm

140 kg

95°

40 - 42

1576200

Essential Plus
(2000 mm length)

2000 x 855 x 140 mm

110 kg

95°

40 - 42

1576201

Mattresses is manufactured to:
EN 597-1 Cig, BS EN 597-2 Match and BS 7177: Crib Source 5

1. *Essential nursing care is pivotal in pressure ulcer prevention. This mattress will positively contribute to the outcome of a pressure ulcer prevention care plan. Education, clinical judgement and action based planning based on vulnerability are fundamental factors in the prevention of pressure ulcers. A range of assessment scales can be used as a formal method of assessing risk from pressure ulcer development, and should be used in conjunction with an informal assessment (informed nursing judgement).*

Invacare Limited
Pencoed Technology Park - Pencoed
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS -
County Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie



Yes, you can.®

