



AQUATEC®

Ocean Ergo

an advanced range of modular shower chairs



Premium personal care solutions

Ergonomic, versatile and easy to manoeuvre



Yes, you can.®

Even greater versatility and ergonomic enhancements

Our latest Aquatec Ocean Ergo range has specifically been designed to make showering and toileting even easier.

From a standard attendant propelled shower chair, to self-propel or a tilt-in-space, Ocean Ergo offers a solution for almost everyone.

In collaboration with leading therapists, our design goals were focused on the ergonomics and the lifestyle and wellbeing of patients.

Every model in the Ocean Ergo range is optimised with new designs to achieve higher levels of safety, comfort and support.

Whether personal care is independent or with the help of a carer, each one of our highly modular chairs offers a dignified solution.

Tilt-in-space models providing additional supportive seating solutions

Ocean VIP Ergo and Ocean Dual VIP Ergo have a -5° to 40° seat tilt angle with an integrated weight shift mechanism that ensures optimal stability on a small footprint.

The Ocean Dual VIP Ergo is multi-functional benefiting from both tilt and recline, allowing the seat and backrest angle to open up between 5° and 35° , in a simple smooth adjustment thanks to the lightweight frame.

5° - 35°



NEW features

Improved armrest system

- > Newly designed armrests provide **increased support and comfort for the user**. The new and optional armrest pads offer individuals a softer surface area and resting position. Featuring an easy to operate optional integrated safety lock system.



Reductions in seat height

- > Additional **seat to ground height adjustment** has been introduced to promote safe and independent transfers. The frame can be lowered further accommodating people of a smaller stature, available on the Ocean Ergo, Ocean 24" Ergo and Ocean Ergo XL.



Tilt stop function

- > An easy to apply optional function which **prevents the frame tilting anteriorly past 0°**. This ensures the frame remains horizontal to the ground and is an aid for those who are hoisted in/out of the chair.



Ergonomic headrest holder

- > With multiple axis adjustment the ergonomic headrest holder offers **superior postural support for the head/neck region**. Adjustments are quick and easy by simply rotating the integrated hand-screws to achieve the correct position for the individual. Ideal for assessment purposes and set up by practitioners with no tools required.



Innovative ergonomic seat

- > **Insight-led design** encourages a more upright seating position for enhanced independence and comfort.
- > Designed around the **imprints** people leave behind when seated.
- > A **5° seat tilt** lifts the knees into a natural squat position for toileting.
- > Transfer supports at the front to allow the user to position themselves into a **comfortable position**.



Highly configurable

- > Available with an extensive range of positioning, toileting and seating accessories that can be used to support a wide range of conditions and postural needs.



Extremely stable

- > **Enhanced stability stainless steel frame** for additional safety and reassurance.



Easy to handle

- > A small, compact footprint and an **optimised wheel base** makes it easy to manoeuvre in small spaces.
- > Easy disassembly for transportation and storage.

Introducing our latest **Ocean Ergo** range

Discover the full range below including optional accessories for specific requirements.



Ocean Ergo

An entry level attendant propelled model



Ocean Ergo XL

An entry level attendant propelled model with a heavier maximum user weight



Ocean 24" Ergo

An entry level self-propelling model



Ocean 24" Ergo XL

An entry level self-propelling model with a heavier maximum user weight



Ocean VIP Ergo

A tilt-in-space model for additional seating and positioning support



Ocean Dual VIP Ergo

A tilt-in-space model with recline for additional seating and positioning support

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.



Seat width ▼



Seat depth ▼



Total width ▼



Total depth ▼



Width between armrests ▼

Ocean Ergo	480 mm	480 mm	565 mm	905 mm	475 mm
Ocean Ergo XL	480 mm	480 mm	645 mm	905 mm	555 mm
Ocean 24" Ergo	480 mm	480 mm	670 mm	1050 mm	475 mm
Ocean 24" Ergo XL	480 mm	480 mm	670 mm	1050 mm	555 mm
Ocean VIP Ergo	480 mm	480 mm	565 mm	1000 mm	475 mm
Ocean Dual VIP Ergo	480 mm	480 mm	565 mm	1000 mm	475 mm



Seat height ▼



Tilt angle seat ▼



Tilt angle back ▼



Weight ▼



Load capacity ▼

Ocean Ergo	440 - 600 mm	-	-	17 kg	150 kg
Ocean Ergo XL	440 - 600 mm	-	-	18 kg	180 kg
Ocean 24" Ergo	440 - 600 mm	-	-	23 kg	150 kg
Ocean 24" Ergo XL	475 - 600 mm	-	-	24 kg	180 kg
Ocean VIP Ergo	475 - 600 mm	-5° - 40°	-	23.5 kg	150 kg
Ocean Dual VIP Ergo	475 - 600 mm	-5° - 40°	5° - 35°	28.3 kg	150 kg



White

Frame colour

T: +44 1656 776222
 uk@invacare.com
 www.invacare.co.uk
 www.invacare.ie

Invacare UK & Ireland
 Unit 4 Pencoed Technology
 Park
 Bridgend CF35 5AQ



Yes, you can.®