



Yes, you can:

Powerchairs Retail Prescription Form

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

MARCH 2021
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

FOX Modulite Teen

Standard Features

- Max User Weight 127Kg (20 Stone)
- Modulite highly adjustable seating system
- Available in seat widths from 305mm - 380mm (12" - 15")
- Seat depths 3050mm - 560mm stepless adjustment (12" - 20")
- Tension adjustable backrest
- Standard T-armrest, detachable with short armpads
- Manual settable seat Tilt (from 0° to 8°)
- Manual settable backrest Recline (0° - 30°, in 3.75° steps)
- Standard 80° hangers with fixed flip-up footplates
- LiNX control REM110 with standard bracket
- Choice of 5 colour options
- Solid Grey castors 200mm x 50mm (8" x 2"), Silver Rim
- Grey drive wheels, puncture proof solid Inlay, (12.1/2" x 2 1/4"), Silver Rim
- 6 km/h (4mph) motors with 40Ah Gel battery
- 8A (200 - 250V) charger
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)
- Wheel lock / parking brake

Standard Base price: £3,034



* Configurations may vary from the image shown

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
AXB2100	FOX TEEN	£3,034
CHAIR TYPE*		
AXZ0000	Fox Teen Basic Model	○
AXZ0000	Special Fox Teen Basic Model	POA
FITTING*		
SWO0025	Invacare ONE to ONE fitting service	◆
SEATING TYPE*		
AXB2100	Modulite seating unit	◆

SEAT FRAME AND SEAT PLATE*			
AXB2111		Telescopic seat frame	◆
SEAT WIDTH*			
AXB0001		305mm (12")	○
AXB0002		330mm (13")	○
AXB2121		380mm (15")	○
SEAT DEPTH*			
AXB0003		305mm (12")	○
AXB0004		330mm (13")	○
AXB0005		355mm (14")	○
AXB0006		380mm (15")	○
AXB2132		410mm (16")	(NC 2196) ○
AXB2134		460mm (18")	○
AXB2136		510mm (20")	○
SEAT CUSHION*			
AXB2142		Matrx Flo-tech Lite	(NC 0001, 0002) £107
AXB2147		Matrx Libra	(NC 2155, 0001, 0002) £372
1626906		PS cushion with Kidabra print 310mm x 310mm, 12" x 12" (Max user weight 68kg)	£183
1626907		PS cushion with Kidabra print 310mm x 360mm, 12" x 14" (Max user weight 68kg)	£183
1626908		PS cushion with Kidabra print 310mm x 410mm, 12" x 16" (Max user weight 68kg)	£183
1627456		PS cushion with Kidabra print 330mm x 310mm, 13" x 12" (Max user weight 68kg)	£183
1627457		PS cushion with Kidabra print 330mm x 360mm, 13" x 14" (Max user weight 68kg)	£183
1627458		PS cushion with Kidabra print 330mm x 360mm, 13" x 14" (Max user weight 68kg)	£183
1626909		PS cushion with Kidabra print 360mm x 310mm, 14" x 12" (Max user weight 68kg)	£183
1626910		PS cushion with Kidabra print 360mm x 360mm, 14" x 14" (Max user weight 68kg)	£183
1626911		PS cushion with Kidabra print 360mm x 410mm, 14" x 16" (Max user weight 68kg)	£183
AXB2140		Without seat cushion	◆
BACKREST SUPPORT*			
AXB2163		Short 480mm	(NC 2172, 2173) ◆
AXB2162		Tall 540mm	○
PUSH HANDLES FOR BACKPOSTS*			
AXB2107		With push handles	◆
AXB2108		Without push handles	(NC 1307, 1495, 2219, 2220, 2227, 2228) £34

BACKREST UPHOLSTERY / CUSHION*			
AXB2171		Tension adjustable backrest	(NC 2152, 1496, 1445, 1446, 2155) ◆
AXB2172		Flex3 backrest without cushion	(NC 2163, 2128) £213
AXB2173		Flex3 backrest with Laguna Backrest	(NC 2163, 2128) £346
SP1562143		Mistral Flex3 backrest shell	(NC 0001, 0002) £129
SP1431668-TR26		Mistral backrest cushion Black multi stretch PU coated fabric	(NC 0001, 0002) £109
SP1431668-TR35		Mistral backrest cushion Black Plush	(NC 0001, 0002) £112
EMP		Extended mounting pins	£45
MINI1207		Mini Back with Kidabra print 310mm x 180mm, 12" x 7" (Max user weight 68kg) (extending mounting Pins required)	£485
MINI1210		Mini Back with Kidabra print 310mm x 260mm, 12" x 10" (Max user weight 68kg) (extending mounting Pins required)	£485
MINI1212		Mini Back with Kidabra print 310mm x 310mm, 12" x 12" (Max user weight 68kg) (extending mounting Pins required)	£485
MINI1214		Mini Back with Kidabra print 310mm x 360mm, 12" x 14" (Max user weight 68kg) (extending mounting Pins required)	£485
MINI1307		Mini Back with Kidabra print 330mm x 180mm, 13" x 7" (Max user weight 68kg) (extending mounting Pins required)	£485
MINI1310		Mini Back with Kidabra print 330mm x 260mm, 13" x 10" (Max user weight 68kg) (extending mounting Pins required)	£485
MINI1312		Mini Back with Kidabra print 330mm x 310mm, 13" x 12" (Max user weight 68kg) (extending mounting Pins required)	£485
MINI1314		Mini Back with Kidabra print 330mm x 360mm, 13" x 14" (Max user weight 68kg) (extending mounting Pins required)	£485
MINI1407		Mini Back with Kidabra print 360mm x 180mm, 14" x 7" (Max user weight 68kg) (extending mounting Pins required)	£485
MINI1410		Mini Back with Kidabra print 360mm x 260mm, 14" x 10" (Max user weight 68kg) (extending mounting Pins required)	£485
MINI1412		Mini Back with Kidabra print 360mm x 310mm, 14" x 12" (Max user weight 68kg) (extending mounting Pins required)	£485
MINI1414		Mini Back with Kidabra print 360mm x 360mm, 14" x 14" (Max user weight 68kg) (extending mounting Pins required)	£485
AXB2170		Without backrest upholstery	(NC 1445, 1446, 1496) ○
BACKREST SUPPORT BRACKET			
SP1571542		Angle adjustable backrest support bracket RH	(NC 0977) £36
SP1571553		Angle adjustable backrest support bracket LH	(NC 0977) £36
UPHOLSTERY TYPE*			
AXB2151		Black fabric	◆
AXB2152		Black multi stretch PU coated fabric	(NC 2171) ○
AXB2155		Black Spacer fabric	(NC 2147, 2171) ○
AXZ1002		Without	○
ARMREST*			
AXB2191		T-Armrest standard	◆
AXB2181		T-Armrest extended height	£34
AXB2183		In-rigged armrests	○
SP1582159		Inrigged armrests for 12" seat width	£30

ARMPADS*			
AXB2195		Short 325mm	◆
AXB2196		Long 375mm	(NC 2132) ○
HEADREST / NECKREST - Adequate headrest option must be selected*			
1603255		Elan standard pad small	£113
1603254		Elan standard pad medium	£113
1603291		Elan Occipital pad	£113
1629822		Elan 4 point headrest pad small	£210
AXB1432		Headrest	£193
AXB1496		All bracketry as per option 1432	(NC 2170, 2171) ○
AXB1495		Push bar 22mm	(NC 2108) ○
ELAN HEADREST MOUNTING HARDWARE			
1603250		Matrx Elan 6-point hardware	(NC 1432, 1496, 1495) £167
1603251		Matrx Elan 6-point Mini	(NC 1432, 1496, 1495) £167
1603249		Matrx standard multi-axis hardware	(NC 1432, 1496, 1495) £150
1603252		Elan adapter for Flex3 back	(NC 1432, 1496, 1495) £10
LATERAL SUPPORT (PAIR)			
AXB1446		Rigid pad	(NC 2170, 2171) £143
AXB1445		Multi-adjustable pad	(NC 2170, 2171) £143
LATERAL SUPPORTS			
1627466		Matrx Mini Lateral Support Swing-away	£124
EPAD-XS		Matrx Mini Lateral Support-Extra small	£60
1627465		Chest harness interface kit	£53
SEAT HEIGHT SETTING*			
AXB0312		Low Seat to Floor height 480mm	○
AXB0344		High Seat to Floor height 510mm	○
SEAT TILT*			
AXB1002		Manual settable	◆
BACKREST RECLINER*			
AXB0975		Manual settable	◆
AXB0977		Powered recliner (0° - 30°)	£431
FOOTREST HANGER TYPE*			
AXB0526		Standard 80° hangers	◆
1606173-0601		Black short 90°	(NC 0512, 0513) £63
1606172-0601		Black short 80°	(NC 0512, 0513) £73
AXB0501		Without hanger	(NC 0512, 0513) ○
LEGREST SUPPORT CARRIER			
SP1531195		LH	(NC 0526) £31
SP1531196		RH	(NC 0526) £31
AMPUTEE SUPPORTS			
AXB0554		Amputee Support RH	£166
AXB0555		Amputee Support LH	£166

FOOTPLATE TYPE*			
AXB0512		Fixed footplate	(NC 0501) ◆
AXB0513		Angle and depth adjustable footplate	(NC 0501) £76
1428660		Angle and depth adjustable footplate Action style LH	(NC 0526) £28
1428654		Angle and depth adjustable footplate Action style RH	(NC 0526) £28
1530620		Fixed footplate Action style LH	(NC 0526) £15
1530614		Fixed footplate Action style RH	(NC 0526) £15
1492150		Up mounted footplate Action style LH	(NC 0526) £60
1492149		Up mounted footplate Action style RH	(NC 0526) £60
ARZ1003		Without footplates	(NC 0526) ○
CONTROLS*			
AXB2201		REM110	(NC 2219) ◆
AXB2203		REM211	○
CONTROL BRACKETRY*			
AXB1207		Standard bracket RH	(NC 1404) ○
AXB1208		Standard bracket LH	(NC 1405) ○
AXB1707		Swing away RH	(NC 1404) £113
AXB1708		Swing away LH	(NC 1405) £113
AXB2227		RH bracket mounted for Attendant	(NC 2108, 2219, 2220, 2128) £54
AXB2228		LH bracket mounted for Attendant	(NC 2108, 2219, 2220, 2128) £54
DUAL CONTROL			
AXB2220		Dual control REM050	(NC 2108) £338
AXB2219		Compact (ACU200)	(NC 2108) £338
DUAL CONTROL BRACKET			
AXB1207		Standard bracket RH	(NC 2227, 2228) ○
AXB1208		Standard bracket LH	(NC 2227, 2228) ○
BODYPOINT JOYSTICK KNOBS			
AXB2721		3" U-shaped, joystick handle	£62
AXB2722		4" U-shaped, joystick handle	£62
AXB2723		Domed rubber, joystick handle	£62
AXB2724		Mushroom grooved, joystick handle	£62
WHEELBASE SETTING*			
AXB1226		Recommended Factory Setting Wheel setting positions can be changed as desired after performing the appropriate risk assessment	◆
SHROUD COLOUR*			
AXB0201		Pearl Night Blue	○
AXB0233		Metallic Silver	○
AXB0228		Silky Matt Black	○
AXB0204		Transparent Cherry Red	○
AXB0222		White	○
FRONT CASTORS*			
AXB0802		Castor Grey 200mm x 50mm (8" x 2") PU tyre	○
AXB0803		Castor Grey 200mm x 50mm (8" x 2") Solid Inlay	◆
DRIVE WHEELS*			
AXB0709		Drive Wheel grey, Puncture Proof Solid Inlay (12 1/2 * 2 1/4")	○
AXB0710		Drive Wheel grey, Puncture Proof Solid Inlay (12 1/2 * 2 1/4")	◆

MOTOR TYPE*			
AXB1294		3 km/h (2mph)	○
AXB1282		6 km/h (4mph)	◆
BATTERY*			
AXB1255		40 Ah Gel	◆
AXB2148		50 Ah AGM	£32
CHARGER*			
AXB1343		OptiCharge 8A (200 - 250V)	◆
AXB1313		OptiCharge 10A (200 - 250V)	£130
WHEEL LOCK / PARKING BREAK*			
AXB1206		Manual wheel locks	◆
KERB RIDER			
AXB1244		Kerbrider	£181
TRANSPORTATION KIT*			
AXB1494		4-Point tie-down	◆
FIXATION BELT*			
AXB1329		Lap belt with auto-lock	◆
HB216-S38-B2		2-point belt, dual pull, side release Small	£59
HB216-M46-B2		2-point belt, dual pull, side release, Medium	£59
HB219-S38-B2		2-point belt, rear pull, side release, Small	£70
HB19-M46-B2		2-point belt, rear pull, side release, Medium	£70
HB235-S32-B2		2-pointbelt, centre pull, 320mm pads	£70
HB235-S38-B2		2-point belt, centre pull, 380mm pads	£70
EB235S		Evoflex Small pelvic stabilizer, push button	£75
EB225M		Evoflex Medium pelvic stabilizer, rehab latch	£75
EB225L		Evoflex Large pelvic stabilizer, rehab latch	£75
BODYPOINT HARNESS			
SH290XXS-B6		X-Extra small, standard PivotFit	£102
SH290XS-B6		Extra small, standard PivotFit	£102
SH290S-B6		Small, standard PivotFit	£102
SH290M-B6		Medium, standard PivotFit	£102
SH295XXS-B6		X-Extra small, Dynamic PivotFit	£107
SH295XS-B6		Extra small, Dynamic PivotFit	£107
SH295S-B6		Small, Dynamic PivotFit	£107
SH295M-B6		Medium, Dynamic PivotFit	£107
1627465		Chest harness interface kit	£42

OTHER OPTIONS

AXB2128		Bag hook	(NC 2172, 2173, 2227, 2228)	£66
AXB1312		RH		£36
AXB1363		LH		£36
AXB1404		Swing away RH		£188
AXB1405		Swing away LH		£188
AXB1490		Castor lock levers		£27

TOTAL PRICE

Non Conforming Options:

AXB0001	305mm (12")
AXB0002	330mm (13")
AXB2132	410mm (16")
AXB2147	Matrx Libra
AXB2163	Short 480mm
AXB2108	Without push handles
AXB2171	Tension adjustable backrest
AXB2172	Flex3 backrest without cushion
AXB2173	Flex3 backrest with Laguna Backrest
MINI1307	Mini Back with Kidabra print 330mm x 180mm, 13" x 7" (Max user weight 68kg) (extending mounting Pins required)
AXB2170	Without backrest upholstery
AXB2152	Black multi stretch PU coated fabric
AXB2155	Black Spacer fabric
AXB2196	Long 375mm
AXB1432	Headrest
AXB1496	All bracketry as per option 1432
AXB1495	Push bar 22mm
AXB1446	Rigid pad
AXB1445	Multi-adjustable pad
AXB0977	Powered recliner (0° - 30°)
AXB0526	Standard 80° hangers
AXB0501	Without hanger
AXB0512	Fixed footplate
AXB0513	Angle and depth adjustable footplate
AXB2227	RH bracket mounted for Attendant
AXB2228	LH bracket mounted for Attendant
AXB2220	Dual control REM050
AXB2219	Compact (ACU200)
AXB2128	Bag hook
AXB1404	Swing away RH
AXB1405	Swing away LH



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):