



**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

JULY 2022  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

**Storm 4 Modulite Teen**

**Standard Features**

- Max user weight 150Kg (23.5 stone)
- Modulite highly adjustable seating system
- Available in Seat Widths from 305mm - 380mm (12" - 15")
- Seat Depths 305mm - 560mm Stepless Adjustment (12" - 22")
- Tension adjustable
- Standard T-Armrests, detachable with Short armpad
- Manual settable seat Tilt
- Manual Settable backrest Recline (0° - 30° in 3.75° steps)
- VARI F hangers with Fixed Footplate
- LiNX Control REM110 with standard bracket
- Choice of 6 Colours
- Puncture Proof Castors
- Puncture Proof Drive Wheels
- 73.5 Ah Batteries
- 6 km/h (4mph) Motors
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)
- Wheel lock / parking brake
- Anti-Tippers
- Suspension
- 8A (200 - 250V) Charger

**Standard Base price: £6,340**



Q'STRAIT  
QLK-150  
Compatible

\* Configurations may vary from the image shown

Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>CHAIR TYPE*</b>		
ATZ0000	Storm 4 Modulite	£6,340
ATZ0000	Special Storm 4 Modulite (Please provide details in special section)	POA
<b>FITTING*</b>		
SWO0025	Invacare ONE to ONE fitting service	◆

<b>SEATING OPTIONS*</b>			
ATB2100		Modulite Seating Unit (Flex3 seating unit)	◆
ATB2111		Telescopic Seat Frame with metal seat plate (with stepless adjustment on armrests, legrests and backrest)	◆
<b>SEAT WIDTH*</b>			
ATB0001		305mm (12")	○
ATB0002		330mm (13")	○
ATB2121		380mm (15")	○
<b>SEAT DEPTH*</b>			
ATB0003		305mm (12")	○
ATB0004		330mm (13")	○
ATB0005		355mm (14")	○
ATB0006		380mm (15")	○
ATB2132		410mm (16") (NC 2196)	○
ATB2134		460mm (18")	○
ATB2136		510mm (20")	○
<b>SEAT CUSHION*</b>			
ATB2142		Matrx Flo-tech Lite (NC 0001, 0002)	£107
ATB2147		Matrx Libra Cushion 90mm thick (NC 2155, 0001, 0002)	£372
1626906		Matrx Mini Series Cushion : PS cushion with Kidabra print 310mm x 310mm, 12" x 12" (Max user weight 68kg)	£182
1626907		Matrx Mini Series Cushion : PS cushion with Kidabra print 310mm x 360mm, 12" x 14" (Max user weight 68kg)	£182
1626908		Matrx Mini Series Cushion : PS cushion with Kidabra print 310mm x 410mm, 12" x 16" (Max user weight 68kg)	£182
1627456		Matrx Mini Series Cushion : PS cushion with Kidabra print 330mm x 310mm, 13" x 12" (Max user weight 68kg)	£182
1627457		Matrx Mini Series Cushion : PS cushion with Kidabra print 330mm x 360mm, 13" x 14" (Max user weight 68kg)	£182
1627458		Matrx Mini Series Cushion : PS cushion with Kidabra print 330mm x 360mm, 13" x 14" (Max user weight 68kg)	£182
1626909		Matrx Mini Series Cushion : PS cushion with Kidabra print 360mm x 310mm, 14" x 12" (Max user weight 68kg)	£182
1626910		Matrx Mini Series Cushion : PS cushion with Kidabra print 360mm x 360mm, 14" x 14" (Max user weight 68kg)	£182
1626911		Matrx Mini Series Cushion : PS cushion with Kidabra print 360mm x 410mm, 14" x 16" (Max user weight 68kg)	£182
ATB2140		Without seat cushion	○
<b>BACKREST SUPPORT*</b>			
ATB2162		Backposts 540mm (Tall) must be selected for Flex3 Laguna backrest	◆
ATB2163		Backposts 480mm (Short) (NC 2172, 2173)	○
<b>PUSH HANDLES FOR BACKPOSTS*</b>			
ATB2107		Backposts WITH push handles	◆
ATB0108		Backposts WITHOUT push handles, Push handles will be cut down on the vertical section (NC 1495, 2220, 2219, 2227, 2228, 2245)	£34

<b>BACKREST UPHOLSTERY / CUSHION*</b>			
ATB2171		Basic Tension Adjustable Backrest (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear, Black fabric only	(NC 2152, 1496, 1445, 1446, 2229, 2155, 0001, 0002) ◆
ATB2172		Flex3 adjustable Backrest Plates (without Backrest Cushion)	(NC 2163, 2128, 0001, 0002, 2128) £213
ATB2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion (maximum seat depth achievable 510mm (20"))	(NC 2163, 2128, 0001, 0002, 2128) £346
SP1562143		Mistral Flex3 Backrest shell - shorter & narrower than the Laguna Flex4	(NC 0001, 0002) £129
SP1431668-TR26		Mistral Backrest Cushion, Black multi stretch PU coated fabric	(NC 0001, 0002) £108
SP1431668-TR35		Mistral Backrest Cushion, Black Plush	(NC 0001, 0002) £112
EMP		Extended Mounting Pins	£45
MINI1207		Matrx Mini Series Backrest : Mini Back with Kidabra print 310mm x 180mm, 12" x 7" (Max user weight 68kg) (extending mounting Pins required)	£485
MINI1210		Matrx Mini Series Backrest : Mini Back with Kidabra print 310mm x 260mm, 12" x 10" (Max user weight 68kg) (extending mounting Pins required)	£485
MINI1212		Matrx Mini Series Backrest : Mini Back with Kidabra print 310mm x 310mm, 12" x 12" (Max user weight 68kg) (extending mounting Pins required)	£485
MINI1214		Matrx Mini Series Backrest : Mini Back with Kidabra print 310mm x 360mm, 12" x 14" (Max user weight 68kg) (extending mounting Pins required)	£485
MINI1307		Matrx Mini Series Backrest : Mini Back with Kidabra print 330mm x 180mm, 13" x 7" (Max user weight 68kg) (extending mounting Pins required)	£485
MINI1310		Matrx Mini Series Backrest : Mini Back with Kidabra print 330mm x 260mm, 13" x 10" (Max user weight 68kg) (extending mounting Pins required)	£485
MINI1312		Matrx Mini Series Backrest : Mini Back with Kidabra print 330mm x 310mm, 13" x 12" (Max user weight 68kg) (extending mounting Pins required)	£485
MINI1314		Matrx Mini Series Backrest : Mini Back with Kidabra print 330mm x 360mm, 13" x 14" (Max user weight 68kg) (extending mounting Pins required)	£485
MINI1407		Matrx Mini Series Backrest : Mini Back with Kidabra print 360mm x 180mm, 14" x 7" (Max user weight 68kg) (extending mounting Pins required)	£485
MINI1410		Matrx Mini Series Backrest : Mini Back with Kidabra print 360mm x 260mm, 14" x 10" (Max user weight 68kg) (extending mounting Pins required)	£485
MINI1412		Matrx Mini Series Backrest : Mini Back with Kidabra print 360mm x 310mm, 14" x 12" (Max user weight 68kg) (extending mounting Pins required)	£485
MINI1414		Matrx Mini Series Backrest : Mini Back with Kidabra print 360mm x 360mm, 14" x 14" (Max user weight 68kg) (extending mounting Pins required)	£485
ATB2170		Without Backrest Upholstery	(NC 1496, 1445, 1446, 2151, 2152, 2155) ○
<b>BACKREST SUPPORT BRACKET</b>			
SP1571542		Angle adjustable backrest support bracket RH	(NC 0977) £36
SP1571553		Angle adjustable backrest support bracket LH	(NC 0977) £36
<b>UPHOLSTERY TYPE*</b>			
ATB2151		Black fabric	(NC 2170) ◆
ATB2152		Black multi stretch PU coated fabric	(NC 2171, 2170) ○
ATB2155		Black Spacer fabric	(NC 2147, 2171, 2170) ○
ATB2170		Without Backrest Upholstery (for option 2170)	○

<b>ARMREST*</b>			
ATB2191		T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm	◆
ATB2181		T-Armrest Extended height, Range of adjustment from seat plate 295mm - 365mm.	£34
ATB2183		In-rigged Armrests (reduces seat width by 25mm each side) (NC 2184, 2185, 2186, 2187)	○
SP1582159		Inrigged Armrests for 12" seat width, Standard inrigged arm option <b>MUST BE SELECTED</b> with additional brackets (must be selected with 12" seat width)	£30
ATB2195		Short Armpad (325mm)	◆
ATB2196		Long Armpad (375mm) (NC 2132)	○
<b>HIP SUPPORTS LH (Not available with 410mm seat depth and Manual Centre legrests)</b>			
ATB2184		SMALL 140mm x 100mm, for flip up and following armrests only (NC 2183)	£215
ATB2186		MEDIUM 230mm x 100mm, for flip up and following armrests only (NC 2183)	£215
<b>HIP SUPPORTS RH (Not available with 410mm seat depth and Manual Centre legrests)</b>			
ATB2185		SMALL 140mm x 100mm, for flip up and following armrests only (NC 2183)	£215
ATB2187		MEDIUM 230mm x 100mm, for flip up and following armrests only (NC 2183)	£215
<b>MOUNTING HARDWARE FOR HIP SUPPORTS</b>			
ATB2165		Mounting hardware for hip supports (8.1" / 206mm high)	○
<b>HEADREST / NECKREST - Adequate headrest option must be selected*</b>			
1603255		Elan Series Headrest : Elan Standard Pad - Small (160mm x 90mm) Standard cover	£113
1603254		Elan Series Headrest : Elan Standard Pad - Medium (260mm x 130mm) Standard cover	£113
1603291		Elan Series Headrest : Elan Occipital Pad - (2300mm x 180mm) Standard cover	£113
1629822		Elan Series Headrest : Elan Standard Pad - Small (160mm x 90mm) Standard cover	£210
ATB1432		Standard Headrest Options : Headrest, height, depth and sideways adjustable (NC 2230, 2247)	£193
ATB1496		Standard Headrest Options : All bracketry as option 1432 without headrest cushion (NC 2171)	○
ATB1495		Standard Headrest Options : Push bar dia 22mm in conjunction with telescopic seat frame / backposts, no Headrest - no mounting hardware	○
<b>ELAN HEADREST MOUNTING HARDWARE*</b>			
1603250		Matrx Elan 6-point Hardware Infinite adjustments (NC 1432, 1496, 1495)	£162
1603251		Matrx Elan 6-point Mini (paediatric) 100mm shorter post (NC 1432, 1496, 1495)	£162
1603249		Matrx standard multi-axis hardware Limited adjustability (NC 1432, 1496, 1495)	£150
1603252		Elan adapter for Flex3 back (NC 1432, 1496, 1495)	£10
<b>LATERAL SUPPORTS</b>			
ATB1446		Rigid pad (Azalea Type, multi stretch PU coated fabric) (NC 2170, 2171)	£143
ATB1445		Multi-adjustable pad (Azalea type, multi stretch PU coated fabric) (NC 2170, 2171)	£143
1627466		Matrx Mini Lateral Support : Swing-away Lateral Hardware (each)	£124
EPAD-XS		Matrx Mini Lateral Support : Extra small Lateral pad (each)	£60
1627465		Matrx Mini Lateral Support : Chest harness interface kit	£53

<b>SEAT HEIGHT / SEAT RISER*</b>			
ATB1023		Fixed Seat Height 450mm (17.75")	○
ATB1025		Fixed Seat Height (19") (NC 2201)	○
ATB1030		Electrical Seat Riser (400mm - 650mm) (16" - 25.5") (NC 2201)	£1,037
ATB1040		Electrical Sear Riser (400mm - 690mm) includes adapter (17.25" - 27") (NC 2201)	£1,101
<b>SEAT TILT*</b>			
ATB1002		Manual Seat Tilt by spindle	◆
ATB1014		Electrical Tilt 0° to 25°	£514
<b>BACKREST RECLINER*</b>			
ATB0975		Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	◆
ATB0977		Powered Recliner (0° - 30°) Foldable (NC 2201)	£431
<b>FOOTREST HANGER TYPE*</b>			
ATB0503		VARI F Angle adjustable with tools 0° - 70°, length 290mm - 460mm (NC 0504)	○
1606173-0601		Black short 90° Legrest Hanger x 2 - Action style (NC 0512, 0504, 0513)	£62
1606172-0601		Black short 80° Legrest Hanger x 2 - Action style (NC 0512, 0504, 0513)	£60
ATB0501		Without Hanger (NC 0512, 0504, 0513)	○
<b>LEGREST SUPPORT CARRIER</b>			
SP1531195		Legrest Support Carrier LH (required for Action style Legrests above)	£31
SP1531196		Legrest Support Carrier RH (required for Action style Legrests above)	£31
<b>FOOTPLATE TYPE*</b>			
ATB0512		Fixed Footplate (NC 0501)	○
ATB0504		Angle adjustable Footplate (NC 0501)	£63
ATB0513		Angle and depth adjustable Footplate (NC 0501)	£101
1428660		Footplate Options for Action Style Hangers : Angle and depth adjustable Footplate LH including stem	£28
1428654		Footplate Options for Action Style Hangers : Angle and depth adjustable Footplate RH including stem	£28
1530620		Footplate Options for Action Style Hangers : Fixed Footplate LH including stem	£15
1530614		Footplate Options for Action Style Hangers : Fixed Footplate RH including stem	£15
1492150		Footplate Options for Action Style Hangers : Upmounted Footplate LH	£61
1492149		Footplate Options for Action Style Hangers : Upmounted Footplate RH	£61
ATBXXXX		Without footplates for option 0501 without hanger	○
<b>CONTROLS*</b>			
ATB2201		LiNX REM110 Drive only (NC 0977, 1437, 2209, 2210, 2270, 2219, 2245)	◆
ATB2203		LiNX REM211 Drive & Actuators (NC 1437, 2209, 2210, 2270)	○
ATB2204		LiNX REM216 Drive, Actuators & Lights (NC 2209, 2210, 2230, 2247, 2257, 2233, 2240, 2270)	£130
ATB2205		LiNX REM400, touch-screen Interface, built in Bluetooth HID (NC 2220, 2209, 2210)	£290
ATB2206		LiNX REM500, without joystick for alternative access, built in Bluetooth HID (Secondary User input required) (NC 2220, 2227, 2228, 2270)	£662

<b>CONTROL BRACKETRY*</b>				
ATB1207		Standard Bracket RH	(NC 1404, 2206)	○
ATB1208		Standard Bracket LH	(NC 1405, 2206)	○
ATB1707		Swing Away RH Bracket	(NC 1404, 2206)	£113
ATB1708		Swing Away LH Bracket	(NC 1405, 2206)	£113
ATB2213		Centre mount swing away to RH Nucleus Bracket	(NC 1404, 1405, 1415)	£425
ATB2214		Centre mount swing away to LH Nucleus Bracket	(NC 1404, 1405, 1339)	£425
ATB2227		Bracket RH mounted for Attendant	(NC 2219, 2220, 2245, 1287, 2128)	£54
ATB2228		Bracket LH mounted for Attendant	(NC 2219, 2220, 2245, 1286, 2128)	£54
ATZ2239		Individual mounting option per (special) - Requires additional information / specification		POA
ATB2209		LiNX REM500 Swing away RH multi-adjustable Bracket	(NC 2201, 2203, 2204, 2205, 1404)	○
ATB2210		LiNX REM500 Swing away LH multi-adjustable Bracket	(NC 2201, 2203, 2204, 2205, 1405)	○
<b>REM400 TOGGLE OPTION</b>				
ATB2270		Assignable Toggle Switch Module for REM400, (as default set to left ON / OFF, right speed dial)	(NC 2201, 2203, 2204)	£270
<b>DUAL CONTROL</b>				
ATB2220		LiNX DLX-REM050 Dual Control (Please select LiNX user control and bracket) (Please select dual control bracket below)	(NC 2201, 2205, 2206, 1221, 1222, 2227, 2228, 2245, 2128)	£338
ATB2219		LiNX Compact Dual Control (DLX-ACU200) (Please select Dual control bracket below)	(NC 2201, 1221, 1222, 2227, 2228, 2245, 2128)	£338
<b>DUAL CONTROL BRACKETRY</b>				
ATB1207		Dual Control Standard Bracket RH	(NC 1221, 1222, 2227, 2228, 2245)	○
ATB1208		Dual Control Standard Bracket LH	(NC 1221, 1222, 2227, 2228, 2245)	○
<b>INVACARE LiNX IDC DUAL CONTROL</b>				
ATB2245		Invacare IDC Dual control (Please select IDC dual control bracket below)	(NC 2219, 2220, 2201, 2128)	£641
<b>INVACARE LiNX IDC DUAL CONTROL BRACKETRY</b>				
ATB1221		IDC Dual control Detachable Bracket RH	(NC 2220, 2219)	○
ATB1222		IDC Dual control Detachable Bracket LH	(NC 2220, 2219)	○
<b>LINX COMPACT REMOTE</b>				
ATB2231		Proportional Compact Remote DLX-CR400 (Requires Remote REM400 or REM500) (Please select LiNX Proportional Bracket)	(NC 2201, 2203, 2204)	£705
ATB2232		Proportional Compact Remote Low Force DLX-CR400LF (Requires Primary Remote REM400 or REM500) (Please select LiNX Proportional Bracketry)	(NC 2201, 2203, 2204)	£943
<b>LINX COMPACT REMOTE BRACKET</b>				
ATB2207		Proportional Compact RH Mount Bracket	(NC 2201, 2203, 2204)	○
ATB2208		Proportional Compact LH Mount Bracket	(NC 2201, 2203, 2204)	○
ATB2213		Proportional Compact Centre Mount swing away to RH (Nucleus)	(NC 2201, 2203, 2204)	£425
ATB2214		LiNX Proportional Compact Centre Mount swing away to LH (Nucleus)	(NC 2201, 2203, 2204)	£425
ATZ2215		Individual mounting option per request (Special) Requires additional information / specification	(NC 2201, 2203, 2204)	POA
<b>CHIN CONTROL</b>				
ATB2229		Proportional Chin Joystick (ASL138) (Requires Primary Remote REM400 or REM500) (Please select chin control bracket below)	(NC 2201, 2203, 2204, 1287, 1286, 2171, 1496, 1495, 2128)	£1,702
<b>CHIN CONTROL BRACKETRY</b>				
ATB1289		Chin Control Swing Away Bracket RH Manual (requires ELAN Headrest)	(NC 1432, 1496, 1495)	○
ATB1288		Chin Control Swing Away Bracket RH Manual (requires ELAN Headrest)	(NC 1432, 1496, 1495)	○
ATB1287		Chin Control Swing Away Bracket RH Powered	(NC 2227, 1404, 1415, 1312)	£2,144
ATB1286		Chin Control Swing Away Bracket LH Powered	(NC 2228, 1405, 1339, 1363)	£2,144

<b>LINX HEAD CONTROL</b>			
ATB2230		Head Control with ultrasonic proximity sensors (ATOM , ASL104) with standard wings M3 (Requires Primary Remote REM400 or REM500)	(NC 2201, 2203, 2204) £3,652
ATB2247		Head Control with ultrasonic proximity sensors (PROTON , ASL104P) with PROTON wings (Requires Primary Remote REM400 or REM500)	(NC 2201, 2203, 2204) £4,277
<b>LINX FOUR-SWITCH PROXIMITY ARRAY (with ATOM technology)</b>			
ATZ2257		Four-switch Proximity Array (ASL106)	£5,655
<b>LINX FOUR-SWITCH PROXIMITY ARRAY (with ATOM technology) MOUNT</b>			
ATZ2275		Centre Mount, swing-away to RH, Nucleus (includes Eclipse Tray ASL619, 260mm x 180mm)	(NC 2201, 2203, 2204, 2220, 2229, 2230, 2247, 2233, 2240, 2128) £425
ATZ2276		Centre Mount, swing-away to LH, Nucleus (includes Eclipse Tray ASL619, 260mm x 180mm)	£425
ATZ2215		Individual mounting option per request (Special) Requires additional information / specification	POA
<b>LINX WIRELESS ACCESSORIES FOR ATOM</b>			
ATB2249		Wireless Triple Switch Receiver (ASL557-3) (Sub D plug)	£290
ATB2250		Wireless Mouse Emulator (ASL558)	£419
<b>LINX MICRO EXTREMITY CONTROL</b>			
ATB2233		MEC Proportional Joystick (Micro Extremity Control) (Please select MEC Bracketry) (Requires Primary Remote REM400 or REM500)	(NC 2201, 2203, 2204) £3,546
<b>LINX MICRO EXTREMITY CONTROL BRACKET</b>			
ATB2235		MEC Tray Mount RH (ASL603, 100mm x 150mm)	(NC 2201, 2203, 2204) ○
ATB2236		MEC Tray Mount LH (ASL603, 100mm x 150mm)	(NC 2201, 2203, 2204) ○
ATB2213		Centre Mount, swing away to RH, Nucleus	(NC 2201, 2203, 2204) £425
ATB2214		Centre Mount, swing away to LH, Nucleus	(NC 2201, 2203, 2204) £425
ATB2241		Centre Mount, swing away to RH, Nucleus (Includes Tray ASL603, 100mm x 150mm)	(NC 2201, 2203, 2204) £425
ATB2242		Centre Mount, swing away to LH, Nucleus (Includes Tray ASL603, 100mm x 150mm)	(NC 2201, 2203, 2204) £425
ATZ2215		Individual mounting option per request (Special) Requires additional information / specification	(NC 2201, 2203, 2204) POA
<b>LINX HEAVY DUTY REMOTE</b>			
ATZ2240		Proportional Heavy Duty Remote Joystick with mode switch (ASL133) (Please select Heavy Duty Bracketry) (Requires Primary Remote REM400 or REM500)	(NC 2201, 2203, 2204) £1,977
<b>LINX HEAVY DUTY REMOTE BRACKET</b>			
ATZ2235		Tray Mount RH, (ASL603, 100mm x 150mm)	(NC 2201, 2203, 2204) ○
ATZ2236		Tray Mount LH, (ASL603, 100mm x 150mm)	(NC 2201, 2203, 2204) ○
ATZ2213		Centre Mount, swing away to RH, Nucleus	(NC 2201, 2203, 2204) £425
ATZ2214		Centre Mount, swing away to LH, Nucleus	(NC 2201, 2203, 2204) £425
ATZ2241		Centre Mount, swing away to RH, Nucleus (Includes Tray ASL603, 100mm x 150mm)	(NC 2201, 2203, 2204) £425
ATZ2242		Centre Mount, swing away to LH, Nucleus (Includes Tray ASL603, 100mm x 150mm)	(NC 2201, 2203, 2204) £425
ATZ2215		Individual mounting option per request (Special) Requires additional information / specification	(NC 2201, 2203, 2204) POA
<b>LINX 10-WAY SWITCH</b>			
ATB2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access (included loom to connect + additional jack for mode switch) (Please select 10 way switch bracket below)	£431

<b>LINX 10-WAY SWITCH BRACKET</b>				
ATB2217		10-Way Switch Bracket RH, Mounted for User	(NC 2209, 2213, 2235, 2241)	○
ATB2218		10-Way Switch Bracket LH, Mounted for User	(NC 2210, 2214, 2236, 2242)	○
ATB2227		10-Way Switch Bracket RH, Mounted for Attendant	(NC 1287, 2128)	£54
ATB2228		10-Way Switch Bracket LH, Mounted for Attendant	(NC 1286, 2128)	£54
<b>LINX ACCESSORY</b>				
ATB2252		Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch)		£86
<b>LINX USB CHARGER PORT</b>				
ATB2253		USB Charger Port DLX-USB02-A		£136
<b>LINX USB CHARGER PORT BRACKET</b>				
ATB2207		USB Charger Port RH Mount		○
ATB2208		USB Charger Port LH Mount		○
<b>FURTHER LINX COMPONENTS / ACCESSORIES</b>				
ATB2300		To be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected)		○
ATZ2301		To be connected and programmed (Special)		POA
<b>BUDDY BUTTONS</b>				
ATB1685		Red Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)		£75
ATB1686		Yellow Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)		£75
ATB1687		Green Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)		£75
ATB1690		Blue Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)		£75
<b>EGG SWITCH</b>				
ATB2254		Egg Switch, black (to be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected)		£75
<b>PIKO BUTTONS</b>				
ATB1711		Red Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)		£135
ATB1712		Yellow Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)		£135
ATB1713		Green Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)		£135
ATB1714		Orange Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)		£135
ATB1715		Blue Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)		£135
ATB1716		Black Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)		£135
<b>TOGGLES / BUTTONS</b>				
ATZ2261		Stereo Toggle - direct control for Tilt only		£230
ATZ2262		Stereo Button - direct control for Tilt only		£219
ATZ2263		Four-Way Toggle - direct control for Tilt & Recline		£286
ATZ2264		Four-Way Button - direct control for Tilt & Recline		£306
ATZ2207		Toggle / Button RH Mount		○
ATZ2208		Toggle / Button LH Mount		○
ATZ2225		Without Mounting		○



<b>BODYPOINT JOYSTICK KNOBS</b>			
ATB2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A	£62
ATB2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A	£62
ATB2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£62
ATB2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£62
<b>REMOTE STOP SWITCH</b>			
ATZ2258		Remote Stop Switch (ASL504)	£897
<b>SHROUD COLOUR*</b>			
ATB0202		Blue	○
ATB0228		Black	○
ATB0220		Champagne	○
ATB0204		Transparent Cherry Red	○
ATB0222		White	○
ATB0216		Olive green	○
ATB0299		Special shroud colour (RAL code required)	POA
<b>FRONT WHEEL*</b>			
ATB0851		3.00 - 6 Grey Pneumatic	○
ATB0853		3.00 - 6 Grey Puncture Proof (sealed)	○
ATB0855		3.00 - 6 Grey Puncture Proof (solid inlay)	£92
<b>REAR WHEEL*</b>			
ATB0771		3.00 - 8 Grey Pneumatic	○
ATB0773		3.00 - 8 Grey Puncture Proof (sealed)	○
ATB0750		3.00 - 8 Grey Puncture Proof (solid inlay)	£92
<b>MOTOR TYPE*</b>			
ATB1277		6km/h (4mph)	◆
ATB1908		10km/h (6mph) (includes LED Lights system with indicators and rear wheel fender)	(NC 1437) £442
<b>BATTERY*</b>			
ATB1364		71 Ah Gel	◆
<b>SUSPENSION*</b>			
ATB1240		Suspension low (up to 80kg user weight)	○
ATB1241		Suspension middle (from 70kg up to 100kg user weight)	○
ATB1242		Suspension hard (from 90kg up to 150kg user weight)	○
<b>G TRACK</b>			
ATB1215		G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	£703
<b>CHARGER*</b>			
ATB1343		OptiCharge Charger 8A (200 - 250V)	◆
ATB1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	£130
<b>LIGHTS*</b>			
ATB1437		LED light system with Indicators	(NC 2201, 2203, 1908) £410
ATB1254		Without Light system	(NC 1908) ◆
<b>ANTI TIPPERS*</b>			
ATB1205		Anti-tipper (standard version)	○

KERB RIDER		
ATB1244		Central Kerb Rider (with parking position) capable of climbing kerbs up to 100mm (4") £181
TRANSPORTATION KIT*		
ATB1450		Preparation for DAHL-Docking station, (Adequate headrest is compulsory, see options below) (NC 1244) £349
ATB1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory) ◆

**FIXATION BELT\***

ATB1329		Lap belt (adjustable both sides) fitted to base of backrest	○
ATB2113		Lap belt (adjustable both sides) adjustable mounting on seat rail extrusion	◆
HB216-S38-B2		Bodypoint Belts : HB216 2-point belt Small x 380mm Pads, cinch mount	£56
HB216-M46-B2		Bodypoint Belts : HB216 2-point belt Medium x 460mm Pads, cinch mount	£56
HB219-S38-B2		Bodypoint Belts : HB219 2-point belt Small x 380mm Pads, cinch mount	£65
HB19-M46-B2		Bodypoint Belts : HB219 2-point belt Medium x 460mm Pads, cinch mount	£65
HB235-S32-B2		Bodypoint Belts : HB235 2-point belt Extra small x, 320mm Pads, cinch mount	£68
HB235-S38-B2		Bodypoint Belts : HB235 2-point belt Small x 380mm Pads, cinch mount	£68
EB235S		Bodypoint Harnesses : Evoflex Small pelvic stabilizer, push button	£80
EB225M		Bodypoint Harnesses : Evoflex Medium pelvic stabilizer, Latch	£80
EB225L		Bodypoint Harnesses : Evoflex Large pelvic stabilizer, Latch	£80
SH250XXS-B2		Bodypoint Harnesses : PivotFit XX-Small, front pull, cinch mount, dynamic trunk support	£97
SH250XS-B2		Bodypoint Harnesses : PivotFit X-Small, front pull, cinch mount, dynamic trunk support	£97
SH250S-B2		Bodypoint Harnesses : PivotFit Small, front pull, cinch mount, dynamic trunk support	£97
SH250M-B2		Bodypoint Harnesses : PivotFit Medium, front pull, cinch mount, dynamic trunk support	£97
SH250XXS-C2-B2		Bodypoint Harnesses : PivotFit XX-Small, shoulder harness, Plain top strap, Cinch mount	£102
SH250XS-C2-B2		Bodypoint Harnesses : PivotFit Small, shoulder harness, Plain top strap, Cinch mount	£102
SH250M-C2-B2		Bodypoint Harnesses : PivotFit Medium, shoulder harness, Plain top strap, Cinch mount	£102
SH260XXS-B2		Bodypoint Harnesses : PivotFit XX-small, rear pull, cinch mount, standard	£102
SH260XS-B2		Bodypoint Harnesses : PivotFit X-small, rear pull, cinch mount, standard	£102
SH260S-B2		Bodypoint Harnesses : PivotFit Small, rear pull, cinch mount, standard	£102
SH260M-B2		Bodypoint Harnesses : PivotFit Medium, rear pull, cinch mount, standard	£102
SH270XXS-B2		Bodypoint Harnesses : PivotFit XX-Small, Front-Pull cinch mount, dynamic trunk support	£102
SH270XS-B2		Bodypoint Harnesses : PivotFit X-Small, Front-Pull cinch mount, dynamic trunk support	£102
SH270S-B2		Bodypoint Harnesses : PivotFit Small, Front-Pull, cinch mount, dynamic trunk support	£102
SH270M-B2		Bodypoint Harnesses : PivotFit Medium, Front-Pull cinch mount, dynamic trunk support	£102
SH270XXXS-C2-B2		Bodypoint Harnesses : PivotFit XXS, Front pull, Plain top strap, cinch mount	£107
SH270XS-C2-B2		Bodypoint Harnesses : PivotFit XS, Front pull, Plain top strap, cinch mount	£107
SH270S-C2-B2		Bodypoint Harnesses : PivotFit Small, Front pull, Plain top strap, cinch mount	£107
SH270M-C2-B2		Bodypoint Harnesses : PivotFit Medium, Front pull, Plain top strap, cinch mount	£107
SH280XXS-B2		Bodypoint Harnesses : PivotFit XX-Small, rear pull, cinch mount, dynamic	£107
SH280XS-B2		Bodypoint Harnesses : PivotFit X-Small, rear pull, cinch mount, dynamic	£107
SH280S-B2		Bodypoint Harnesses : PivotFit Small, rear pull, cinch	£107

<b>REAR VIEW MIRROR</b>			
ATB1415	<input type="checkbox"/>	Rear view Mirror RH - Must be selected with 10 km/h Motor	£36
ATB1339	<input type="checkbox"/>	Rear view Mirror LH - Must be selected with 10 km/h Motor	£36
<b>BAG HOOK</b>			
ATB2128	<input type="checkbox"/>	Bag Hook (Max load up to 6 kg = impact on max user weight) (NC 2172, 2173, 2227, 2228, 2220, 2219, 2229, 2245)	£66
<b>LUGGAGE RACK</b>			
ATB1394	<input type="checkbox"/>	Luggage Rack	£281
<b>CANE HOLDER</b>			
ATB1312	<input type="checkbox"/>	Cane holder RH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote) (NC 2227, 1287)	£36
ATB1363	<input type="checkbox"/>	Cane holder LH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote) (NC 2228, 1286)	£36
<b>TRAY</b>			
ATB1404	<input type="checkbox"/>	Tray Swing Away RH (NC 2213, 2214, 2209, 1207, 1707, 1287)	£300
ATB1405	<input type="checkbox"/>	Tray Swing Away LH (NC 2213, 2214, 2210, 1208, 1707, 1286)	£300
<b>FENDER</b>			
ATB1416	<input type="checkbox"/>	Fender for front Wheels	£103
ATB1417	<input type="checkbox"/>	Fender for rear wheels (mandatory for Lighting system)	£103

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

ATB0001	305mm (12")
ATB0002	330mm (13")
ATB2132	410mm (16")
ATB2147	Matrx Libra Cushion 90mm thick
ATB2163	Backposts 480mm (Short)
ATB2171	Basic Tension Adjustable Backrest (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear, Black fabric only
ATB2172	Flex3 adjustable Backrest Plates (without Backrest Cushion)
ATB2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion (maximum seat depth achievable 510mm (20"))
MINI1207	Matrx Mini Series Backrest : Mini Back with Kidabra print 310mm x 180mm, 12" x 7" (Max user weight 68kg) (extending mounting Pins required)
MINI1312	Matrx Mini Series Backrest : Mini Back with Kidabra print 330mm x 310mm, 13" x 12" (Max user weight 68kg) (extending mounting Pins required)
ATB2170	Without Backrest Upholstery
ATB2151	Black fabric
ATB2152	Black multi stretch PU coated fabric
ATB2155	Black Spacer fabric
ATB2170	Without Backrest Upholstery (for option 2170)
ATB2183	In-rigged Armrests (reduces seat width by 25mm each side)
ATB2196	Long Armpad (375mm)
ATB2184	SMALL 140mm x 100mm, for flip up and following armrests only
ATB2186	MEDIUM 230mm x 100mm, for flip up and following armrests only
ATB2185	SMALL 140mm x 100mm, for flip up and following armrests only
ATB2187	MEDIUM 230mm x 100mm, for flip up and following armrests only
ATB1432	Standard Headrest Options : Headrest, height, depth and sideways adjustable
ATB1496	Standard Headrest Options : All bracketry as option 1432 without headrest cushion
ATB1495	Standard Headrest Options : Push bar dia 22mm in conjunction with telescopic seat frame / backposts, no Headrest - no mounting hardware
ATB1446	Rigid pad (Azalea Type, multi stretch PU coated fabric)
ATB1445	Multi-adjustable pad (Azalea type, multi stretch PU coated fabric)
ATB0977	Powered Recliner (0° - 30°) Foldable
ATB0501	Without Hanger
ATB0512	Fixed Footplate
ATB0504	Angle adjustable Footplate
ATB0513	Angle and depth adjustable Footplate
ATB2201	LiNX REM110 Drive only
ATB2203	LiNX REM211 Drive & Actuators
ATB2204	LiNX REM216 Drive, Actuators & Lights
ATB2205	LiNX REM400, touch-screen Interface, built in Bluetooth HID
ATB2206	LiNX REM500, without joystick for alternative access, built in Bluetooth HID (Secondary User input required)
ATB1207	Standard Bracket RH
ATB1208	Standard Bracket LH
ATB1707	Swing Away RH Bracket
ATB2213	Centre mount swing away to RH Nucleus Bracket
ATB2214	Centre mount swing away to LH Nucleus Bracket
ATB2227	Bracket RH mounted for Attendant

ATB2228	Bracket LH mounted for Attendant
ATB2209	LiNX REM500 Swing away RH multi-adjustable Bracket
ATB2210	LiNX REM500 Swing away LH multi-adjustable Bracket
ATB2270	Assignable Toggle Switch Module for REM400, (as default set to left ON / OFF, right speed dial)
ATB2220	LiNX DLX-REM050 Dual Control (Please select LiNX user control and bracket) (Please select dual control bracket below)
ATB2219	LiNX Compact Dual Control (DLX-ACU200) (Please select Dual control bracket below)
ATB1207	Dual Control Standard Bracket RH
ATB1208	Dual Control Standard Bracket LH
ATB2245	Invacare IDC Dual control (Please select IDC dual control bracket below)
ATB1221	IDC Dual control Detachable Bracket RH
ATB1222	IDC Dual control Detachable Bracket LH
ATB2213	Proportional Compact Centre Mount swing away to RH (Nucleus)
ATB2214	LiNX Proportional Compact Centre Mount swing away to LH (Nucleus)
ATB2229	Proportional Chin Joystick (ASL138) (Requires Primary Remote REM400 or REM500) (Please select chin control bracket below)
ATB1287	Chin Control Swing Away Bracket RH Powered
ATB1286	Chin Control Swing Away Bracket LH Powered
ATB2230	Head Control with ultrasonic proximity sensors (ATOM , ASL104) with standard wings M3 (Requires Primary Remote REM400 or REM500)
ATB2247	Head Control with ultrasonic proximity sensors (PROTON , ASL104P) with PROTON wings (Requires Primary Remote REM400 or REM500)
ATZ2257	Four-switch Proximity Array (ASL106)
ATB2233	MEC Proportional Joystick (Micro Extremity Control) (Please select MEC Bracketry) (Requires Primary Remote REM400 or REM500)
ATB2235	MEC Tray Mount RH (ASL603, 100mm x 150mm)
ATB2236	MEC Tray Mount LH (ASL603, 100mm x 150mm)
ATB2213	Centre Mount, swing away to RH, Nucleus
ATB2214	Centre Mount, swing away to LH, Nucleus
ATB2241	Centre Mount, swing away to RH, Nucleus (Includes Tray ASL603, 100mm x 150mm)
ATB2242	Centre Mount, swing away to LH, Nucleus (Includes Tray ASL603, 100mm x 150mm)
ATZ2240	Proportional Heavy Duty Remote Joystick with mode switch (ASL133) (Please select Heavy Duty Bracketry) (Requires Primary Remote REM400 or REM500)
ATZ2235	Tray Mount RH, (ASL603, 100mm x 150mm)
ATZ2236	Tray Mount LH, (ASL603, 100mm x 150mm)
ATZ2213	Centre Mount, swing away to RH, Nucleus
ATZ2214	Centre Mount, swing away to LH, Nucleus
ATZ2241	Centre Mount, swing away to RH, Nucleus (Includes Tray ASL603, 100mm x 150mm)
ATZ2242	Centre Mount, swing away to LH, Nucleus (Includes Tray ASL603, 100mm x 150mm)
ATB2227	10-Way Switch Bracket RH, Mounted for Attendant
ATB2228	10-Way Switch Bracket LH, Mounted for Attendant
ATB1908	10km/h (6mph) (includes LED Lights system with indicators and rear wheel fender)
ATB1437	LED light system with Indicators
ATB1244	Central Kerb Rider (with parking position) capable of climbing kerbs up to 100mm (4")
ATB1415	Rear view Mirror RH - Must be selected with 10 km/h Motor
ATB1339	Rear view Mirror LH - Must be selected with 10 km/h Motor
ATB2128	Bag Hook (Max load up to 6 kg = impact on max user weight)
ATB1312	Cane holder RH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote)

---

ATB1363	Cane holder LH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote)
ATB1404	Tray Swing Away RH
ATB1405	Tray Swing Away LH

---



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):