



**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

JANUARY 2022  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Action 3 Junior Small

### Standard Features

- Max user weight 80kg (12.5 stone)
- Crash tested ISO 7176-19
- Self Propulsion or Attendant Propulsion
- Choice of 4 seat widths 8" - 11" (205mm - 280mm)
- Choice of 3 seat depths (depth adjustable upholstery)
- Choice fo 3 seat heights 16" - 18" (410mm - 460mm)
- Black nylon upholstery
- Folding backrest set at 15" (370mm)
- Width and height adjustable backrest upholstery
- 80° Detachable swing away legrest hangers
- Flip up, angle and depth adjustable footplates
- Removable and height adjustable desk armrests
- Adjustable heel straps
- Choice of 7 colours (with 8 special colour options)
- +100mm outriggered front castors
- Push or pull to lock brakes (configuration dependant)
- Transportation kit and buckle belt (adjustable both sides)
- User tool kit



**Standard Base price: £616**

Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description		
<b>BASE PRODUCT*</b>			
CNC0001	ACTION <sup>3</sup> JUNIOR		○
CNCSPEC	SPECIAL		POA
<b>PROPULSION*</b>			
CNC0100	SELF PROPULSION	(NC 1055, 1170)	£616
CNC0110	TRANSIT PROPULSION	(NC 0280, 1020, 1030, 1045, 1140, 1141, 1180, 1190, 1195, 1275, 1390, 1430, 1440)	£616

SEAT WIDTH*				
CNC0140		8" SW	(NC 0465)	○
CNC0141		9" SW	(NC 0465)	○
CNC0142		10" SW	(NC 0465)	◆
CNC0143		11" SW	(NC 0465)	○
SEAT DEPTH*				
CNC0160		ADJ SD 8" - 10"	(NC 1030)	○
CNC0161		ADJ SD 10" - 12"		◆
CNC0162		ADJ SD 12" - 14"		○
SEAT HEIGHT*				
CNC0260		18" SH		○
CNC0270		17" SH		○
CNC0280		16" SH	(NC 1430, 1440)	○
SEAT UPHOLSTERY*				
CNC0335		NYLON WITH HOOK AND LOOP STRIPS - WIDTH AND DEPTH ADJUSTABLE		◆
SEAT CUSHION*				
CNC33333		WITHOUT SEAT CUSHION		○
CNC0362		BASIC SEAT CUSHION (ESSENTIAL JUNIOR)		£48
BACKREST*				
CNC0370		FIXED 0° BACKREST 3 HEIGHTS 15" / 370MM, 16" / 400MM OR 17" / 430MM		○
CNC0380		FIXED 7° BACKREST 3 HEIGHTS 15" / 370MM, 16" / 400MM OR 17" / 430MM		○
CNC0381		FIXED 7° BACKREST 20"		○
CNC0390		RECLINING (RATCHET) 0/30° 5 HEIGHTS 13" - 17"	(NC 1050, 1190, 1195)	£106
CNC0395		RECLINING (GAS) 0/30° 16"	(NC 1050, 1190, 1195)	£175
CNC0400		FOLDING BACKREST 3 HEIGHTS 15" / 370MM, 16" / 400 OR 17" / 430MM		◆
CNC0401		FOLIDNG BACKREST 20"		○
CNC0402		FOLDING BACKREST 2 HEIGHTS 15" / 20" REAR SET PUSH HANDLES	(NC 0465, 0855, 0902, 1426)	£133
CNC0420		FOLDING BACKREST 16" HEIGHT ADJUSTABLE PUSH HANDLES	(NC 0855)	£72
CNC0425		FIXED 7° BACKREST 16" HEIGHT ADJUSTABLE PUSH HANDLES	(NC 0855)	£72
BACKREST UPHOLSTERY*				
CNC0465		W/O BACK UPH - DELIVERED WITH SIMPLY SLEEVE	(NC 0855)	○
CNC0475		WIDTH AND HEIGHT ADJUSTABLE NYLON SLING TYPE		◆
LEGREST HANGER RH*				
CNC0505		90° STRAIGHT RH REMOVABLE HANGER UL STD 125 TO 260 MM		£28
CNC0510		80° SHORT RH SWING IN/OUT HANGER 0° STANDARD POSITION; UL STD 220 TO 320 MM BLACK COLOUR ONLY	(NC 0562)	◆
CNC0512		80° LONG RH SWING IN/OUT HANGER 0° STANDARD POSITION; UL STD 300 TO 450 MM	(NC 0560)	○

<b>LEGREST HANGER LH*</b>			
CNC0555		90° STRAIGHT LH REMOVABLE HANGER UL STD 125 TO 260 MM	£28
CNC0560		80° SHORT LH SWING IN/OUT HANGER 0° STANDARD POSITION; UL STD 220 TO 320 MM BLACK COLOUR ONLY	◆
CNC0562		80° LONG LH SWING IN/OUT HANGER 0° STANDARD POSITION; UL STD 300 TO 450 MM	○
<b>FOOTREST RH*</b>			
CNC0605		FLIP UP DEPTH AND ANGLE ADJ RH	◆
CNC0611		ONE PIECE FOOTBOARD (NC 0625)	○
<b>FOOTREST LH*</b>			
CNC0625		FLIP UP DEPTH AND ANGLE ADJ LH	◆
CNC66666		NO FOOTPLATE REQUIRED (PLEASE SELECT IN COMBINATION WITH OPTION 0611)	○
<b>ARMREST RH*</b>			
CNC0640		REMOVABLE RH (N°1) DESK, STD ARMPAD (60 MM) ADJ. 180 TO 250 MM AUTOLOCK SAFETY DEVICE INCLUDED (NC 0745, 0746, 0750, 0755)	◆
CNC0645		REMOVABLE RH (N°1) JUNIOR, STD ARMPAD (60 MM) ADJ. 150 TO 250 MM AUTOLOCK SAFETY DEVICE INCLUDED (NC 0740, 0746, 0750, 0755)	○
CNC0646		REMOVABLE RH (N°1) NARROW ARMPAD JUNIOR (50 MM) ADJ. 150 TO 250 MM AUTOLOCK SAFETY DEVICE INCLUDED (NC 0740, 0745, 0750, 0755)	○
CNC0650		MUDGUARD RH HEIGHT ADJUSTABLE (155 TO 220 MM) (NC 0740, 0745, 0746, 0755, 0910, 1201, 1471, 1472)	£33
CNC0655		SKIRTGUARD RH JUNIOR (NC 0740, 0745, 0746, 0750, 0910, 1471, 1472)	○
<b>ARMREST LH*</b>			
CNC0740		REMOVABLE LH (N°1) DESK, STD ARMPAD (60 MM) ADJ. 180 TO 250 MM AUTOLOCK SAFETY DEVICE INCLUDED	◆
CNC0745		REMOVABLE LH (N°1) JUNIOR, STD ARMPAD (60 MM) ADJ. 150 TO 250 MM AUTOLOCK SAFETY DEVICE INCLUDED	○
CNC0746		REMOVABLE LH (N°1) NARROW ARMPAD JUNIOR (50 MM) ADJ. 150 TO 250 MM AUTOLOCK SAFETY DEVICE INCLUDED	○
CNC0750		MUDGUARD LH HEIGHT ADJUSTABLE (155 TO 220 MM) (NC 0910, 1201, 1471, 1472)	£33
CNC0755		SKIRTGUARD LH JUNIOR (NC 0910, 1471, 1472)	○
<b>OTHER SEATING OPTIONS</b>			
CNC0855		BACKREST EXTENSION (NC 0860, 0865)	£69
CNC0860		HEADREST STANDARD MULTI - ADJUSTABLE (NC 0865, 0902)	£128
CNC0865		HEADREST ROUND MULTI - ADJUSTABLE (NC 0902)	£144
CNC0880		ADJUSTABLE CALF STRAP	£10
CNC0875		BODYPOINT ANKLE HUGGERS (PAIR) SIZE S - UNMOUNTED	£109
CNC0886		FOOTPLATE ADJUSTMENT OPTIMIZER KIT	£20
CNC0902		WIDTH ADJUSTABLE BACK BRACE	£39
CNC0910		WITHOUT AUTOLOCK SAFETY DEVICE ON ARMRESTS N°1	○

COLOUR*			
CNC0952		PEARL WHITE	○
CNC0953		MATT BLACK	○
CNC0954		MATT RED	○
CNC0955		INTENSE YELLOW	○
CNC0956		POLAR BLUE	○
CNC0957		CANDY PINK	○
CNC0958		POP BLUE	○
CNC0991		SPECIAL COLOUR (ON DEMAND ONLY)	£161
SWO0001		SLATE GREY - 3 WEEK LEAD TIME	£161
SWO0002		ELECTRIC GREEN - 3 WEEK LEAD TIME	£161
SWO0003		HYPNOTIC PURPLE - 3 WEEK LEAD TIME	£161
SWO0004		ANTHRACITE BLACK - 3 WEEK LEAD TIME	£161
SWO0005		DARK RED - 3 WEEK LEAD TIME	£161
SWO0006		PEARL GREY - 3 WEEK LEAD TIME	£161
SWO0007		VIBRANT ORANGE - 3 WEEK LEAD TIME	£161
FRONT CASTORS*			
CNC0999		5" X 1" BLACK SOLID TYRE	◆
CNC1001		5" X 1" FLASHING LIGHTS SOLID TYRE	£44
CNC1007		5" X 1.3/4" BLACK SOLID TYRE	○
REAR WHEEL SIZE*			
CNC1010		12"	(NC 1040, 1045, 1130, 1140, 1141, 1180, 1253, 1254, 1255, 1256, 1380, 1390) ○
CNC1020		20"	(NC 1040, 1055, 1050, 1170) ○
CNC1030		22"	(NC 1045, 1055, 1170) ○
REAR WHEELS*			
CNC1040		STANDARD SELF-PROPEL REAR WHEEL (BLACK HUB, COLORLESS RIM AND ZINC PLATTED SPOKES)	(NC 1280) ○
CNC1045		REINFORCED STANDARD SELF-PROPEL REAR WHEEL (COLORLESS HUB, CLEAR ANODIZED RIM AND STAINLESS STEEL SPOKES)	(NC 1280) ○
CNC1050		DRUM BRAKE SELF-PROPEL 22" or TRANSIT 12"	(NC 1270, 1275) ○
CNC1055		TRANSIT COMPOSITE (BLACK PLASTIC HUB AND RIM)	(NC 1280) ○
TYRES*			
CNC1120		SOLID BLACK	○
CNC1130		PNEUMATIC PROFILE BLACK	○
HANDRIM - A SELECTION MUST BE MADE FOR SELF PROPULSION			
CNC1140		ALUMINIUM PRE-ANODIZED HANDRIM	○
CNC1141		HANDRIM COATED IN BLACK ANTI-SLIP PLASTIC	£71
REAR WHEEL AXLE*			
CNC1170		FIXED AXLE	○
CNC1180		Q/R AXLE	○
REAR WHEEL POSITION*			
CNC1190		STANDARD REAR WHEELS POSITION	○
CNC1195		MEDIUM REAR WHEEL POSITION	○
CNC1200		PASSIVE REAR WHEEL POSITION	○
CNC1201		DUAL AMPUTEE REAR WHEEL POSITION	(NC 1275, 1280) £65

<b>PUMP</b>			
CNC1220	<input type="checkbox"/>	W/O PUMP	○
CNC1230	<input type="checkbox"/>	PUMP AND SUPPORT	£9
<b>CHASSIS OPTION</b>			
CNC1253	<input type="checkbox"/>	SPOKE GUARDS CLEAR	£77
CNC1254	<input type="checkbox"/>	SPOKE GUARDS PIXEL	£95
CNC1255	<input type="checkbox"/>	SPOKE GUARDS SNOWFLAKES	£95
CNC1256	<input type="checkbox"/>	SPOKE GUARDS COLOR ME	£95
<b>BRAKES*</b>			
CNC1270	<input type="checkbox"/>	PUSH TO LOCK BRAKES	◆
CNC1275	<input type="checkbox"/>	PULL TO LOCK BRAKES	○
CNC1280	<input type="checkbox"/>	ATTENDANT HUB BRAKES (NC 1426)	£140
<b>TRANSPORTATION KIT*</b>			
CNC1345	<input type="checkbox"/>	TRANSPORTATION BRACKETS	◆
<b>SAFETY OPTIONS*</b>			
CNC1360	<input type="checkbox"/>	POSTURE BELT	◆
CNC1361	<input type="checkbox"/>	PEDIATRIC BELT (Flotech dual pull)	£87
CNC1370	<input type="checkbox"/>	REAR ANTI-TIPPERS WITH TIPPING AID (NC 1430, 1440)	£33
CNC1380	<input type="checkbox"/>	REFLECTORS KIT	£6
CNC1390	<input type="checkbox"/>	BRAKE EXTENSION PAIR	£13
<b>OTHER OPTIONS</b>			
CNC1426	<input type="checkbox"/>	STROLLER PUSH HANDLES	£85
CNC1430	<input type="checkbox"/>	RH TIPPING AID	£13
CNC1440	<input type="checkbox"/>	LH TIPPING AID	£13
CNC1471	<input type="checkbox"/>	RH SWING AWAY ACRYLIC TRAY	£187
CNC1472	<input type="checkbox"/>	LH SWING AWAY ACRYLIC TRAY	£187

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

CNC0280	16" SH
CNC0465	W/O BACK UPH - DELIVERED WITH SIMPLY SLEEVE
CNC0560	80° SHORT LH SWING IN/OUT HANGER 0° STANDARD POSITION; UL STD 220 TO 320 MM BLACK COLOUR ONLY
CNC0562	80° LONG LH SWING IN/OUT HANGER 0° STANDARD POSITION; UL STD 300 TO 450 MM
CNC0625	FLIP UP DEPTH AND ANGLE ADJ LH
CNC0740	REMOVABLE LH (N°1) DESK, STD ARMPAD (60 MM) ADJ. 180 TO 250 MM AUTOLOCK SAFETY DEVICE INCLUDED
CNC0745	REMOVABLE LH (N°1) JUNIOR, STD ARMPAD (60 MM) ADJ. 150 TO 250 MM AUTOLOCK SAFETY DEVICE INCLUDED
CNC0746	REMOVABLE LH (N°1) NARROW ARMPAD JUNIOR (50 MM) ADJ. 150 TO 250 MM AUTOLOCK SAFETY DEVICE INCLUDED
CNC0750	MUDGUARD LH HEIGHT ADJUSTABLE (155 TO 220 MM)
CNC0755	SKIRTGUARD LH JUNIOR
CNC0855	BACKREST EXTENSION
CNC0860	HEADREST STANDARD MULTI - ADJUSTABLE
CNC0865	HEADREST ROUND MULTI - ADJUSTABLE
CNC0902	WIDTH ADJUSTABLE BACK BRACE
CNC0910	WITHOUT AUTOLOCK SAFETY DEVICE ON ARMRESTS N°1
CNC1020	20"
CNC1030	22"
CNC1040	STANDARD SELF-PROPEL REAR WHEEL (BLACK HUB, COLORLESS RIM AND ZINC PLATTED SPOKES)
CNC1045	REINFORCED STANDARD SELF-PROPEL REAR WHEEL (COLORLESS HUB, CLEAR ANODIZED RIM AND STAINLESS STEEL SPOKES)
CNC1050	DRUM BRAKE SELF-PROPEL 22" or TRANSIT 12"
CNC1055	TRANSIT COMPOSITE (BLACK PLASTIC HUB AND RIM)
CNC1130	PNEUMATIC PROFILE BLACK
CNC1140	ALUMINIUM PRE-ANODIZED HANDRIM
CNC1141	HANDRIM COATED IN BLACK ANTI-SLIP PLASTIC
CNC1170	FIXED AXLE
CNC1180	Q/R AXLE
CNC1190	STANDARD REAR WHEELS POSITION
CNC1195	MEDIUM REAR WHEEL POSITION
CNC1201	DUAL AMPUTEE REAR WHEEL POSITION
CNC1253	SPOKE GUARDS CLEAR
CNC1254	SPOKE GUARDS PIXEL
CNC1255	SPOKE GUARDS SNOWFLAKES
CNC1256	SPOKE GUARDS COLOR ME
CNC1270	PUSH TO LOCK BRAKES
CNC1275	PULL TO LOCK BRAKES
CNC1280	ATTENDANT HUB BRAKES
CNC1380	REFLECTORS KIT
CNC1390	BRAKE EXTENSION PAIR
CNC1426	STROLLER PUSH HANDLES
CNC1430	RH TIPPING AID
CNC1440	LH TIPPING AID
CNC1471	RH SWING AWAY ACRYLIC TRAY





For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):