



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Manual Retail Prescription Form

NOVEMBER 2020
 Customer Ref:

| | |
|--------------------------|--------------------|
| Customer Ref* | Account No* |
| Contact No* | Order No* |
| Delivery Address* | Order Qty* |
| | Order date* |

Action 3NG Comfort

Standard Features

- Max user weight 125kg (19.7 stone)
- Crash tested ISO 7176-Part 19
- Self Propulsion or Attendant Propulsion
- Choice of 4 seat widths (16" - 20") 405mm to 505mm
- Choice of 2 seat depths (16" or 18") 400mm or 450mm
- Choice of 3 seat height (16" - 18") 410mm to 460mm
- Black nylon seat upholstery
- Reclining (Ratchet) backrest 0° - 30° with Laguna shell 20° (510mm)
- Comfort Black multi stretch PU coated fabric Passad 2 backrest
- Elevating legrest hangers 0° to 75° including calf pad
- Flip up, angle and depth adjustable footplates
- Removable, height and depth adjustable comfort armrests
- Black multi stretch PU coated fabric trunk supports and headrest
- Choice of 8 colours (with 7 special colour options)
- Push or Pull to lock brakes (configuration dependant)
- Comfort transportation kit and buckle belt (adjustable both sides)
- Adjustable heel straps
- User tool kit



Standard Base price: £969

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

| Code | Description | |
|----------------------|--------------------|--|
| BASE PRODUCT* | | |
| CEG0005 | ACTION 3NG COMFORT | ○ |
| CEGSPEC | SPECIAL | POA |
| PROPULSION* | | |
| CEG0100 | SELF PROPULSION | (NC 1062, 1074) £969 |
| CEG0110 | TRANSIT PROPULSION | (NC 0290, 0300, 1068, 1150, 1180, 1260, 1275, 1390) £969 |

| SEAT WIDTH* | | | |
|-----------------------|--|---|------|
| CEG0190 | | 16" / 405 MM SEAT WIDTH | ○ |
| CEG0200 | | 17" / 430 MM SEAT WIDTH | ◆ |
| CEG0210 | | 18" / 455 MM SEAT WIDTH | ○ |
| CEG0230 | | 20" / 505 MM SEAT WIDTH | ○ |
| SEAT DEPTH* | | | |
| CEG0240 | | 16" / 400 MM SEAT DEPTH | ○ |
| CEG0250 | | 18" / 450 MM SEAT DEPTH | ◆ |
| SEAT HEIGHT* | | | |
| CEG0280 | | 18" / 460 MM SEAT HEIGHT | ◆ |
| CEG0290 | | 17" / 435 MM SEAT HEIGHT (NC 1010, 1062) | ○ |
| CEG0300 | | 16" / 410 MM SEAT HEIGHT (NC 1010, 1062, 1370, 1440, 1441) | ○ |
| SEAT UPHOLSTERY* | | | |
| CEG0340 | | NYLON WITH HOOK AND LOOP STRIPS | ◆ |
| CEG0372 | | CUSHION COMFORT MISTRAL 2 LAYER (INCLUDING SEAT UPHOLSTERY 3 VELCROS) | ○ |
| BACKREST* | | | |
| CEG0441 | | RECLINING (RATCHET) 0/30° LAGUNA SHELL 20" | ◆ |
| BACKREST UPHOLSTERY* | | | |
| CEG0490 | | COMFORT BACKREST LAGUNA BLACK MULTI STRETCH PU COATED 20" / 510 MM HEIGHT | ○ |
| CEG0491 | | COMFORT BACKREST PASSAD 2 BLACK MULTI STRETCH PU COATED 20" / 510 MM HEIGHT | ◆ |
| LEGREST HANGER RH* | | | |
| CEG0520 | | ELEVATING LEGREST RH SWING IN/SWING OUT (NC 66666) | ◆ |
| CEG0545 | | AMPUTEES LEGREST RH (NC 0610, 1111) | £115 |
| LEGREST HANGER LH* | | | |
| CEG0570 | | ELEVATING LEGREST LH SWING IN/SWING OUT (NC 55555) | ◆ |
| CEG0595 | | AMPUTEES LEGREST LH (NC 0630, 2222) | £115 |
| FOOTREST RH* | | | |
| CEG0610 | | FLIP UP DEPTH AND ANGLE ADJ RH | ◆ |
| SWO1111 | | FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY) (NC 0545) | ◆ |
| CEE66666 | | NO FOOTREST | ○ |
| FOOTREST LH* | | | |
| CEG0630 | | FLIP UP DEPTH AND ANGLE ADJ LH | ◆ |
| SWO2222 | | FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY) (NC 0595) | ◆ |
| CEE55555 | | NO FOOTREST | ○ |
| ARMREST RH* | | | |
| CEG0691 | | REMOVABLE RH (N°4) COMFORT ARM | ◆ |
| CEG0732 | | REMOVABLE RH (N°4) HEMI ARM (NC 1481) | £51 |
| ARMREST LH* | | | |
| CEG0791 | | REMOVABLE LH (N°4) COMFORT ARM | ◆ |
| CEG0832 | | REMOVABLE LH N°4 HEMI ARM (NC 1481) | £51 |
| OTHER SEATING OPTIONS | | | |
| CEG0881 | | TRUNK SUPPORT COMFORT BLACK MULTI STRETCH PU COATED FABRIC (PAIR) ONLY FOR LAGUNA SHELL | ◆ |
| CEG0906 | | HEADREST BLACK MULTI STRETCH PU COATED FABRIC | ◆ |

| COLOUR* | | | |
|-----------------------------|--|--|------|
| CEG0961 | | SLATE GREY (RAL7024 FINE TEXTURE) | ○ |
| CEG0964 | | PEARL WHITE | ○ |
| CEG0965 | | ELECTRIC GREEN | ○ |
| CEG0966 | | MATT BLACK | ○ |
| CEG0967 | | BLUE POP | ○ |
| CEG0968 | | HYPNOTIC PURPLE | ○ |
| CEG0969 | | RED MATT | ○ |
| CEG0980 | | ANTHRACITE BLACK | ○ |
| CEG0991 | | SPECIAL COLOUR ON DEMAND ONLY | £149 |
| SWO3333 | | DARK RED - SPECIAL | £149 |
| SWO4444 | | PEARL GREY - SPECIAL | £149 |
| SWO5555 | | INTENSE YELLOW | £149 |
| SWO6666 | | VIBRANT ORANGE | £149 |
| SWO7777 | | POLAR BLUE | £149 |
| SWO8888 | | CANDY PINK | £149 |
| CASTOR HOUSING* | | | |
| CEG0995 | | STANDARD | ◆ |
| CEG0996 | | LIGHT WEIGHT | £26 |
| FRONT CASTORS* | | | |
| CEG1010 | | 8" X 1.1/4" BLACK SOLID TYRE | ◆ |
| CEG1050 | | 6" X 1.1/4" BLACK SOLID TYRE | ○ |
| REAR WHEEL SIZE* | | | |
| CEG1062 | | 12" (NC 1071, 1072, 1084, 1088, 1140, 1150, 1180, 1260, 1390) | ○ |
| CEG1068 | | 24" | ◆ |
| REAR WHEELS* | | | |
| CEG1071 | | STANDARD (BASIC) (BLACK HUB, COLORLESS RIM AND ZINC PLATTED SPOKES) (NC 1280) | ◆ |
| CEG1072 | | SELF-PROPEL REINFORCED REAR WHEEL (NC 1280) | £20 |
| CEG1073 | | DRUM BRAKE SELF-PROPEL 22 AND 24" or TRANSIT 12" (NC 1270, 1275) | ○ |
| CEG1074 | | TRANSIT BLACK (NC 1088, 1260, 1280) | ○ |
| TYRES* | | | |
| CEG1082 | | SOLID BLACK | ◆ |
| CEG1084 | | PNEUMATIC PROFILE BLACK | £18 |
| CEG1088 | | PNEUMATIC LOW PROFILE BLACK | ○ |
| HANDRIM | | | |
| CEG1140 | | ALUMINIUM PRE-ANODIZED HANDRIM | ◆ |
| CEG1150 | | HANDRIM COATED IN BLACK ANTI-SLIP PLASTIC | £29 |
| REAR WHEEL AXLE* | | | |
| CEG1170 | | FIXED AXLE | ○ |
| CEG1180 | | Q/R AXLE | ◆ |
| REAR WHEEL POSITION* | | | |
| CEG1200 | | PASSIVE REAR WHEELS POSITION | ◆ |
| CEG1201 | | DUAL AMPUTEE REAR WHEEL POSITION (NC 1275, 1280) | £64 |
| PUMP | | | |
| CEG1230 | | PUMP AND SUPPORT | £6 |

| CHASSIS OPTION | | | |
|-----------------------------|--------------------------|--|------|
| CEG1241 | <input type="checkbox"/> | FRONT WHEELS EXTENSION | £19 |
| CEG1260 | <input type="checkbox"/> | SPOKE GUARDS PAIR 24" BLK | £68 |
| BRAKES* | | | |
| CEG1270 | <input type="checkbox"/> | PUSH TO LOCK | ◆ |
| CEG1275 | <input type="checkbox"/> | PULL TO LOCK BRAKES | ○ |
| CEG1280 | <input type="checkbox"/> | ATTENDANT HUB BRAKES | £204 |
| COMFORT TRANSPORTATION KIT* | | | |
| CEG1346 | <input type="checkbox"/> | COMFORT TRANSPORTATION KIT | ◆ |
| SAFETY OPTIONS | | | |
| CEG1361 | <input type="checkbox"/> | BODYPOINT 2-PT POSTURE BELT M (NC 1362) | £71 |
| CEG1362 | <input type="checkbox"/> | BODYPOINT 2-PT POSTURE BELT L | £71 |
| CEG1370 | <input type="checkbox"/> | REAR ANTI-TIPPERS PAIR (NC 1375, 1440, 1441) | ○ |
| CEG1375 | <input type="checkbox"/> | REAR ANTI-TIPPERS PAIR WITH TIPPING AID (NC 1440, 1441) | £62 |
| CEG1390 | <input type="checkbox"/> | BRAKE EXTENSION PAIR | £12 |
| OTHER OPTIONS | | | |
| CEG1430 | <input type="checkbox"/> | CANE HOLDER | £20 |
| CEG1440 | <input type="checkbox"/> | TIPPING AID RH | £12 |
| CEG1441 | <input type="checkbox"/> | TIPPING AID LH | £12 |
| CEG1481 | <input type="checkbox"/> | COMFORT ACRYLIC TRAY | £90 |
| SWO10001 | <input type="checkbox"/> | ABDUCTED LEGRESTS (PR) 0° TO 30° - INTERNAL REFERENCE DE259618 | £95 |

| | |
|--------------------|--|
| TOTAL PRICE | |
|--------------------|--|

Non Conforming Options:

| | |
|----------|---|
| CEG0290 | 17" / 435 MM SEAT HEIGHT |
| CEG0300 | 16" / 410 MM SEAT HEIGHT |
| CEG0545 | AMPUTEE LEGREST RH |
| CEG0595 | AMPUTEE LEGREST LH |
| CEG0610 | FLIP UP DEPTH AND ANGLE ADJ RH |
| SWO1111 | FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY) |
| CEE66666 | NO FOOTREST |
| CEG0630 | FLIP UP DEPTH AND ANGLE ADJ LH |
| SWO2222 | FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY) |
| CEE55555 | NO FOOTREST |
| CEG1010 | 8" X 1.1/4" BLACK SOLID TYRE |
| CEG1062 | 12" |
| CEG1068 | 24" |
| CEG1071 | STANDARD (BASIC) (BLACK HUB, COLORLESS RIM AND ZINC PLATTED SPOKES) |
| CEG1072 | SELF-PROPEL REINFORCED REAR WHEEL |
| CEG1074 | TRANSIT BLACK |
| CEG1084 | PNEUMATIC PROFILE BLACK |
| CEG1088 | PNEUMATIC LOW PROFILE BLACK |
| CEG1140 | ALUMINIUM PRE-ANODIZED HANDRIM |
| CEG1150 | HANDRIM COATED IN BLACK ANTI-SLIP PLASTIC |
| CEG1180 | Q/R AXLE |
| CEG1260 | SPOKE GUARDS PAIR 24" BLK |
| CEG1270 | PUSH TO LOCK |
| CEG1275 | PULL TO LOCK BRAKES |
| CEG1280 | ATTENDANT HUB BRAKES |
| CEG1362 | BODYPOINT 2-PT POSTURE BELT L |
| CEG1370 | REAR ANTI-TIPPERS PAIR |
| CEG1375 | REAR ANTI-TIPPERS PAIR WITH TIPPING AID |
| CEG1390 | BRAKE EXTENSION PAIR |
| CEG1440 | TIPPING AID RH |
| CEG1441 | TIPPING AID LH |
| CEG1481 | COMFORT ACRYLIC TRAY |



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):