



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

NOVEMBER 2020
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Action 3NG Self Propel Compact

Standard Features

- Max user weight 125kg (19.7 stone)
- Crash tested ISO 7176-Part 19
- Self Propulsion with 20" Quick release rear wheels
- Choice of 3 seat widths 380mm to 430mm (15" - 17")
- Choice of 2 seat depths 400mm or 450mm (16" or 18")
- Choice of 5 seat heights 360mm to 460mm (14" - 18")
- Black nylon upholstery
- Folding 7° 430mm backrest
- Lumbar support padded backrest with shaped cover
- 90° Swing in / Swing out legrest hangers
- Flip up, Angle and Depth Adjustable Footplates
- Swing away and detachable desk armrests
- Choice of 8 standard colours (with 7 special colour options)
- Push or pull to lock brakes (configuration dependant)
- Transportation kit and buckle belt (adjustable both sides)
- Adjustable heel straps
- User tool kit

Standard Base price: £618



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
CEN0012	ACTION 3NG COMPACT	○
CENSPEC	SPECIAL	POA
PROPULSION*		
CEN0100	SELF PROPULSION	£618

SEAT WIDTH*			
CEN44444		13" / 330 MM SEAT WIDTH - SELECT SEAT WIDTH OPTION CEE0180 AND IN-RIGGED ARMREST OPTION CEE0952 (lapbelt adjustment 420mm - 1520mm)	£36
CEN0180		15" / 380 MM SEAT WIDTH	○
CEN0190		16" / 405 MM SEAT WIDTH	◆
CEN0200		17" / 430 MM SEAT WIDTH	○
SEAT DEPTH*			
CEN0240		16" / 400 MM SEAT DEPTH (NC 0879)	◆
CEN0250		18" / 450 MM SEAT DEPTH (NC 0878)	○
SEAT HEIGHT*			
CEN0280		18" / 460 MM SEAT HEIGHT	○
CEN0290		17" / 435 MM SEAT HEIGHT	◆
CEN0300		16" / 410 MM SEAT HEIGHT (NC 1370, 1440, 1441)	○
CEN0310		15" / 385 MM SEAT HEIGHT (NC 0511, 0561, 0660, 0670, 0760, 0770, 1370, 1440, 1441)	○
CEN0315		14" / 360 MM SEAT HEIGHT, 0° SEAT ANGLE (AVAILABLE WITH LOW FORK, 5" FRONT WHEEL AND 20" REAR WHEELS) (NC 0511, 0561, 0611, 0660, 0670, 0760, 0770, 1066, 1375, 1440, 1441)	○
SEAT UPHOLSTERY*			
CEN0330		NYLON W/O HOOK AND LOOP STRIPS	◆
CEN0340		NYLON WITH HOOK AND LOOP STRIPS	○
BACKREST*			
CEN0380		FIXED 0° BACKREST 2 HEIGHTS 16" / 400 MM or 17" / 430 MM (DELIVERED 17" / 430 MM) FIXED PUSH HANDLES (NC 0471, 0477, 0945)	○
CEN0390		FIXED 0° BACKREST 20" / 510 MM HEIGHT FIXED PUSH HANDLES (NC 0471, 0477, 0945)	○
CEN0395		FIXED 7° BACKREST HEIGHT ADJUST 15" / 380 MM TO 18" / 460 MM (DELIVERED 15" / 380 MM) FIXED PUSH HANDLES (follow upholstery height) (NC 0474, 0477, 0478, 0481, 0905, 0906, 0907, 0940)	£54
CEN0396		FIXED 7° BACKREST HEIGHT ADJUST 15" / 380 TO 18" / 460 MM (DELIVERED 15" / 380 MM) HEIGHT ADJUSTABLE PUSH HANDLES (independantly from upholstery height) (NC 0474, 0477, 0478, 0481, 0905, 0906, 0907, 0940)	£71
CEN0410		FIXED 7° BACKREST 2 HEIGHTS 16" / 400 or 17" / 430 MM (DELIVERED 17" / 430 MM) FIXED PUSH HANDLES (NC 0471, 0477, 0945)	◆
CEN0420		FIXED 7° BACKREST 20" / 510 MM HEIGHT FIXED PUSH HANDLES (NC 0471, 0477, 0945)	○
CEN0425		FIXED 7° BACKREST 17" / 430 MM HEIGHT HEIGHT ADJUSTABLE PUSHHANDLE (6" / 150 MM) (NC 0471, 0477, 0945)	£38
CEN0428		FIXED 7° BACKREST ACTIVE 14" / 350 MM HEIGHT (ADJUSTABLE TO 16" / 406 MM) FIXED PUSH HANDLES (SHORT) (NC 0471, 0477, 0481, 0905, 0906, 0907, 0940)	○
CEN0432		ANGLE ADJUSTABLE AND FOLD DOWN (15° to + 15°) FIXED 7°, HEIGHT 17" / 425MM WITH FIXED LONG PUSHHANDLES HEIGHT 17" / 425 MM 18" / 450 MM SET AT 425MM / 17" (2 positions, 1" / 25mm steps) (NC 0465, 0471, 0474, 0478, 0481, 0882, 0952)	£77
CEN0445		RECLINING (GAS SPRING) 0/30° BACKREST INCL. BACKRACE AND PASSIVE REAR WHEEL POS. 3 HEIGHTS 15" to 17" / 390/410/430 MM (DELIVERED 17" / 430 MM) ANTITIPPERS REQUIRED (NC 0471, 0477, 0940, 0945, 1190)	£206
CEN0450		RECLINING (GAS SPRING) 0/30° BACKREST INCL. BACKRACE AND PASSIVE REAR WHEEL POS. 3 HEIGHTS 18" to 20" / 470/490/510 MM (DELIVERED 20" / 510 MM) ANTITIPPERS REQUIRED (NC 0471, 0477, 0940, 0945, 1190)	£206
CEN0455		FOLDING 7° BACKREST 17" / 430 MM HEIGHT FIXED PUSH HANDLES (NC 0471, 0477, 0945)	£55
CEN0460		FOLDING 7° BACKREST 20" / 510 MM HEIGHT FIXED PUSH HANDLES (NC 0471, 0477, 0481, 0945)	£55

BACKREST UPHOLSTERY*				
CEN0465		W/O BACK UPH	(NC 0882)	○
CEN0471		STANDARD NYLON SLING (JEMIMA)	(NC 0882, 0905, 0906, 0907)	○
CEN0474		LUMBAR SUPPORT BI MATERIAL PADDED COVER		○
CEN0477		CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER		£49
CEN0478		CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE	(NC 0940, 0945)	£55
CEN0481		CONTOUR ADJUSTABLE (7+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE	(NC 0940, 0945)	£106
LEGREST HANGER RH*				
CEN0511		90° RH SWING IN/SWING OUT	(NC 1111, 66666)	◆
CEN0520		ELEVATING LEGREST RH SWING IN/SWING OUT	(NC 0611, 0930, 66666)	£102
CEN0545		AMPUTEE LEGREST RH	(NC 0610, 0611, 0930, 1111)	£118
LEGREST HANGER LH*				
CEN0561		90° LH SWING IN/SWING OUT	(NC 2222, 55555)	◆
CEN0570		ELEVATING LEGREST LH SWING IN/SWING OUT	(NC 0611, 0930, 55555)	£102
CEN0595		AMPUTEE LEGREST LH	(NC 0611, 0630, 0930, 2222)	£118
FOOTREST RH*				
CEN0610		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.		○
SWO1111		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)		○
CEN0611		FLIP UP ONE PIECE	(NC 0630, 55555)	£69
CEN66666		NO FOOTREST		○
FOOTREST LH*				
CEN0630		FLIP UP DEPTH AND ANGLE ADJUSTABLE LH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.		○
SWO2222		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)		○
CEN55555		NO FOOTREST		○
ARMREST RH*				
CEN0640		REMOVABLE RH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED DESK ADJ. 200 TO 280 MM HEIGHT	(NC 0760, 0770, 0795, 0796, 1480, 1483)	○
CEN0650		REMOVABLE RH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED FULL LENGTH ADJ. 200 TO 280 MM HEIGHT	(NC 0760, 0770, 0795, 0796)	○
CEN0660		SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT	(NC 0740, 0750, 0795, 0796, 0952, 0955, 1480, 1483)	◆
CEN0670		SWING AWAY DETACHABLE RH (N°2) FULL LENGTH 240 MM HEIGHT	(NC 0740, 0750, 0795, 0796, 0952, 0955)	○
CEN0695		REMOVABLE RH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ.	(NC 0740, 0750, 0760, 0770, 0952, 0955, 1480, 1483)	£53
CEN0696		REMOVABLE RH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ.	(NC 0740, 0750, 0760, 0770, 0952, 0955, 1483)	£53

ARMREST LH*			
CEN0740		REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED DESK ADJ. 200 TO 280 MM HEIGHT (NC 1480, 1484)	○
CEN0750		REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED FULL LENGTH ADJ. 200 TO 280 MM HEIGHT	○
CEN0760		SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT (NC 0952, 0955, 1480, 1484)	◆
CEN0770		SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT (NC 0952, 0955)	○
CEN0795		REMOVABLE LH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ. (NC 0952, 0955, 1480, 1484)	£53
CEN0796		REMOVABLE LH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. (NC 0952, 0955, 1484)	£53
SEAT CUSHION			
CEN0849		BASIC 50 MM FOAM CUSHION - NYLON AND ANTI-SLIP COVER (NC 0952)	£21
SEAT CUSHION DEPTH			
CEN0878		FOR 400 MM SEAT DEPTH	○
CEN0879		FOR 450 MM SEAT DEPTH	○
OTHER SEATING OPTIONS			
CEN0882		LATERAL POSITIONING PADS (PAIR)	£54
CEN0905		NECKREST BLACK MULTI STRETCH PU COATED FABRIC (NC 0906, 0907, 0940, 0945)	£117
CEN0906		HEADREST BLACK MULTI STRETCH PU COATED FABRIC (NC 0907, 0940, 0945)	£140
CEN0907		HEADREST WITH CHEEK SUPPORT BLACK MULTI STRETCH PU COATED FABRIC (NC 0940, 0945)	£162
CEN0930		CALF STRAP	○
CEN0940		BACKBRACE (NC 0945)	£31
CEN0945		FOLDABLE BACKBRACE	£34
CEN0952		IN-RIGGED ARMREST (NC 1480, 1483, 1484)	£37
CEN0955		W/O AUTOLOCK SAFETY DEVICE FOR ARMREST N°1 ONLY	○
COLOUR*			
CEN0961		SLATE GREY (RAL7024 FINE TEXTURE)	○
CEN0964		PEARL WHITE	○
CEN0965		ELECTRIC GREEN	○
CEN0966		MATT BLACK	○
CEN0967		BLUE POP	○
CEN0968		HYPNOTIC PURPLE	○
CEN0969		RED MATT	○
CEN0980		ANTHRACITE BLACK	○
CASTOR HOUSING*			
CEN0995		STANDARD	◆
CEN0996		LIGHT WEIGHT	£26
FRONT CASTORS*			
CEN1054		5" X 1" BLACK SOLID TYRE	◆
REAR WHEEL SIZE*			
CEN1064		20" (NC 1071, 1076, 1088)	◆
CEN1066		22"	○

REAR WHEELS*			
CEN1071	<input type="checkbox"/>	STANDARD (BASIC) (BLACK HUB, COLORLESS RIM AND ZINC PLATED SPOKES)	◆
CEN1072	<input type="checkbox"/>	SELF-PROPEL REINFORCED REAR WHEEL	£20
CEN1076	<input type="checkbox"/>	LIGHTWEIGHT (NC 1230, 1261)	£52
TYRES*			
CEN1082	<input type="checkbox"/>	SOLID BLACK	◆
CEN1084	<input type="checkbox"/>	PNEUMATIC PROFILE BLACK	£32
CEN1088	<input type="checkbox"/>	PNEUMATIC LOW PROFILE BLACK	○
HANDRIM*			
CEN1140	<input type="checkbox"/>	ALUMINIUM PRE-ANODIZED HANDRIM	◆
REAR WHEEL AXLE*			
CEN1180	<input type="checkbox"/>	Q/R AXLE	◆
REAR WHEEL POSITION*			
CEN1190	<input type="checkbox"/>	STANDARD REAR WHEELS POSITION	◆
CEN1200	<input type="checkbox"/>	PASSIVE REAR WHEELS POSITION	○
PUMP			
CEN1230	<input type="checkbox"/>	PUMP AND SUPPORT	£12
CEN1235	<input type="checkbox"/>	PUMP - HIGH PRESSURE	£36
CHASSIS OPTION			
CEN1241	<input type="checkbox"/>	FRONT WHEELS EXTENSION	£19
CEN1261	<input type="checkbox"/>	SPOKE GUARDS PAIR 560 AND 515 MM (22" AND 20") CLR	£86
BRAKES*			
CEN1270	<input type="checkbox"/>	PUSH TO LOCK	◆
CEN1275	<input type="checkbox"/>	PULL TO LOCK BRAKES	○
TRANSPORTATION KIT*			
CEN1345	<input type="checkbox"/>	TRANSPORTATION KIT	◆
SAFETY OPTIONS			
CEN1361	<input type="checkbox"/>	BODYPOINT 2-PT POSTURE BELT M (NC 1362)	£71
CEN1362	<input type="checkbox"/>	BODYPOINT 2-PT POSTURE BELT L	£71
CEN1370	<input type="checkbox"/>	REAR ANTI-TIPPERS PAIR (NC 1375, 1440, 1441)	£32
CEN1375	<input type="checkbox"/>	REAR ANTI-TIPPERS PAIR WITH TIPPING AID (NC 1440, 1441)	£64
CEN1380	<input type="checkbox"/>	REFLECTORS KIT	£6
CEN1390	<input type="checkbox"/>	BRAKE EXTENSION PAIR	£12
OTHER OPTIONS			
CEN1440	<input type="checkbox"/>	TIPPING AID RH	£12
CEN1441	<input type="checkbox"/>	TIPPING AID LH	£12
CEN1480	<input type="checkbox"/>	ACRYLIC TRAY (NC 1483, 1484, 1485)	£102
CEN1483	<input type="checkbox"/>	HALF TABLE TRAY RH (NC 1484)	£78
CEN1484	<input type="checkbox"/>	HALF TABLE TRAY LH	£78
CEN1485	<input type="checkbox"/>	COMFORT PADDING FOR HALF TABLE TRAY	£39
SWO10002	<input type="checkbox"/>	ABDUCTED LEGRESTS (PR) 0° TO 30° - INTERNAL REFERENCE DE259618	£95

TOTAL PRICE	
--------------------	--

Non Conforming Options:

CEN0465	W/O BACK UPH
CEN0471	STANDARD NYLON SLING (JEMIMA)
CEN0474	LUMBAR SUPPORT BI MATERIAL PADDED COVER
CEN0477	CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER
CEN0478	CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE
CEN0481	CONTOUR ADJUSTABLE (7+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE
CEN0511	90° RH SWING IN/SWING OUT
CEN0561	90° LH SWING IN/SWING OUT
CEN0610	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.
SWO1111	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CEN0611	FLIP UP ONE PIECE
CEN66666	NO FOOTREST
CEN0630	FLIP UP DEPTH AND ANGLE ADJUSTABLE LH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.
SWO2222	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CEN55555	NO FOOTREST
CEN0660	SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT
CEN0670	SWING AWAY DETACHABLE RH (N°2) FULL LENGTH 240 MM HEIGHT
CEN0740	REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED DESK ADJ. 200 TO 280 MM HEIGHT
CEN0750	REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED FULL LENGTH ADJ. 200 TO 280 MM HEIGHT
CEN0760	SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT
CEN0770	SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT
CEN0795	REMOVABLE LH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ.
CEN0796	REMOVABLE LH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ.
CEN0878	FOR 400 MM SEAT DEPTH
CEN0879	FOR 450 MM SEAT DEPTH
CEN0882	LATERAL POSITIONING PADS (PAIR)
CEN0905	NECKREST BLACK MULTI STRETCH PU COATED FABRIC
CEN0906	HEADREST BLACK MULTI STRETCH PU COATED FABRIC
CEN0907	HEADREST WITH CHEEK SUPPORT BLACK MULTI STRETCH PU COATED FABRIC
CEN0930	CALF STRAP
CEN0940	BACKBRACE
CEN0945	FOLDABLE BACKBRACE
CEN0952	IN-RIGGED ARMREST
CEN0955	W/O AUTOLOCK SAFETY DEVICE FOR ARMREST N°1 ONLY
CEN1066	22"
CEN1071	STANDARD (BASIC) (BLACK HUB, COLORLESS RIM AND ZINC PLATTED SPOKES)
CEN1076	LIGHTWEIGHT
CEN1088	PNEUMATIC LOW PROFILE BLACK
CEN1190	STANDARD REAR WHEELS POSITION

CEN1230	PUMP AND SUPPORT
CEN1261	SPOKE GUARDS PAIR 560 AND 515 MM (22" AND 20") CLR
CEN1362	BODYPOINT 2-PT POSTURE BELT L
CEN1370	REAR ANTI-TIPPERS PAIR
CEN1375	REAR ANTI-TIPPERS PAIR WITH TIPPING AID
CEN1440	TIPPING AID RH
CEN1441	TIPPING AID LH
CEN1480	ACRYLIC TRAY
CEN1483	HALF TABLE TRAY RH
CEN1484	HALF TABLE TRAY LH
CEN1485	COMFORT PADDING FOR HALF TABLE TRAY



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):