



# Manual Retail Prescription Form

**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

NOVEMBER 2020  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Action 3NG Transit

### Standard Features

- Max user weight 125kg (19.7 stone)
- Crash tested ISO 7176-Part 19
- Attendant Propulsion
- Choice of 6 seat widths 380mm to 505mm (15" to 20")
- Choice of 2 seat depths 400mm or 450mm (16" or 18")
- Choice of 3 seat heights 460mm to 510mm (18" to 20")
- Black nylon upholstery
- Folding 7° 430mm backrest
- Lumbar support padded backrest with shaped cover
- 80° Swing in / swing out legrest hangers
- Flip up footplates
- Swing away and detachable desk armrests
- Choice of 8 colours (with 7 special colours)
- Push to lock brakes
- Transportation kit and buckle belt (adjustable both sides)
- Adjustable heel straps
- User tool kit



**Standard Base price: £567**

Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE PRODUCT*</b>		
CEF0004	ACTION 3NG TRANSIT (NC 1074)	<input type="radio"/>
CEFSPEC	SPECIAL	POA
<b>PROPULSION*</b>		
CEF0110	TRANSIT PROPULSION	£567

SEAT WIDTH*			
CEF44444		13" / 330 MM SEAT WIDTH - SELECT SEAT WIDTH OPTION CEE0180 AND IN-RIGGED ARMREST OPTION CEE0952 (lapbelt adjustment 420mm - 1520mm)	£36
CEF0180		15" / 380 MM SEAT WIDTH	○
CEF0190		16" / 405 MM SEAT WIDTH	○
CEF0200		17" / 430 MM SEAT WIDTH	◆
CEF0210		18" / 455 MM SEAT WIDTH	○
CEF0220		19" / 480 MM SEAT WIDTH	○
CEF0230		20" / 505 MM SEAT WIDTH	○
SEAT DEPTH*			
CEF0240		16" / 400 MM SEAT DEPTH (NC 0879)	○
CEF0250		18" / 450 MM SEAT DEPTH (NC 0878)	◆
CEF0252		16" TO 18" / 400 TO 450 MM ADJUSTABLE SEAT DEPTH (NC 0330, 0350)	£18
SEAT HEIGHT*			
CEF0260		20" / 510 MM SEAT HEIGHT	○
CEF0270		19" / 485 MM SEAT HEIGHT	○
CEF0280		18" / 460 MM SEAT HEIGHT	○
SEAT UPHOLSTERY*			
CEF0330		NYLON W/O HOOK AND LOOP STRIPS	◆
CEF0340		NYLON WITH HOOK AND LOOP STRIPS	○
CEF0350		PADDED NYLON	£26
BACKREST - Anti-tippers should be selected with reclining backrests 0435, 0440, 0445 and 0450*			
CEF0380		FIXED 0° BACKREST 2 HEIGHTS 16" / 400 MM or 17" / 430 MM (NC 0471, 0477, 0945)	○
CEF0390		FIXED 0° BACKREST 20" / 510 MM HEIGHT (NC 0471, 0477, 0945)	○
CEF0395		FIXED 7° BACKREST HT ADJ 15" / 380 MM TO 18" / 460 MM FIX PUSH HANDLES (NC 0474, 0477, 0478, 0481, 0905, 0906, 0907, 0940)	£59
CEF0396		FIXED 7° BACKREST HEIGHT ADJUST 15" / 380 TO 18" / 460 MM HT ADJ PUSH HANDLES (NC 0474, 0477, 0478, 0481, 0905, 0906, 0907, 0940)	£71
CEF0410		FIXED 7° BACKREST 2 HEIGHTS 16" / 400 or 17" / 430 MM FIX PUSH HANDLES (NC 0471, 0477, 0945)	◆
CEF0420		FIXED 7° BACKREST 20" / 510 MM HEIGHT FIX PUSH HANDLES (NC 0471, 0477, 0945)	○
CEF0425		FIXED 7° BACKREST 17" / 430 MM HEIGHT HT ADJ PUSHHANDLE (NC 0471, 0477, 0945)	£38
CEF0432		ANGLE ADJUSTABLE AND FOLD DOWN (15° to + 15°) 17" (NC 0465, 0471, 0474, 0478, 0481, 0882, 0952)	£77
CEF0450		RECLINING (GAS SPRING) 0/30° 3 HEIGHTS 18" 19" 20" (NC 0471, 0477, 0940, 0945)	£206
CEF0455		FOLDING 7° BACKREST 17" (NC 0471, 0477, 0945)	£55
CEF0460		FOLDING 7° BACKREST 20" (NC 0471, 0477, 0481, 0945)	£55
BACKREST UPHOLSTERY*			
CEF0465		W/O BACK UPH (NC 0882)	○
CEF0471		STANDARD NYLON SLING (JEMIMA) (NC 0882, 0905, 0906, 0907)	○
CEF0474		LUMBAR SUPPORT BI MATERIAL PADDED COVER	○
CEF0477		CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER	£48
CEF0478		CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE (NC 0940, 0945)	£55
CEF0481		CONTOUR ADJUSTABLE (7+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE (NC 0940, 0945)	£106

LEGREST HANGER RH*				
CEF0510		80° RH SWING IN/SWING OUT	(NC 1111, 66666)	◆
CEF0520		ELEVATING LEGREST RH SWING IN/SWING OUT	(NC 0611, 0930, 0935, 0936, 66666)	£102
CEF0545		AMPUTEES LEGREST RH	(NC 0600, 0610, 0611, 0930, 0936, 1111)	£118
LEGREST HANGER LH*				
CEF0560		80° LH SWING IN/SWING OUT	(NC 2222)	◆
CEF0595		AMPUTEES LEGREST LH	(NC 0611, 0620, 0630, 0930, 0936, 2222)	£118
FOOTREST RH*				
CEF0600		FLIP UP RH FOOTPLATE		◆
CEF0610		FLIP UP DEPTH AND ANGLE AJUSTABLE RH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.		£38
SWO1111		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)		£38
CEF0611		FLIP UP ONE PIECE	(NC 0620, 0630, 0935, 0936)	£69
CEF66666		NO FOOTREST		○
FOOTREST LH*				
CEF0620		FLIP UP LH FOOTPLATE		◆
CEF0630		FLIP UP DEPTH AND ANGLE AJUSTABLE LH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.		£38
SWO2222		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)		£38
CEF55555		NO FOOTREST		○
ARMREST RH*				
CEF0640		REMOVABLE RH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED DESK ADJ. 200 TO 280 MM HEIGHT	(NC 0760, 0770, 0780, 0790, 0795, 0796, 1480, 1483)	○
CEF0650		REMOVABLE RH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED FULL LENGTH ADJ. 200 TO 280 MM HEIGHT	(NC 0760, 0770, 0780, 0790, 0795, 0796)	○
CEF0660		SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT	(NC 0740, 0750, 0780, 0790, 0795, 0796, 0952, 0955, 1480, 1483)	◆
CEF0670		SWING AWAY DETACHABLE RH (N°2) FULL LENGTH 240 MM HEIGHT	(NC 0740, 0750, 0780, 0790, 0795, 0796, 0952, 0955)	○
CEF0680		SWING AWAY DETACHABLE RH (N°3) DESK ADJ. IN HEIGHT (170-295 MM)	(NC 0740, 0750, 0760, 0770, 0795, 0796, 0952, 0955, 1480, 1483)	£49
CEF0690		SWING AWAY DETACHABLE RH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH	(NC 0740, 0750, 0760, 0770, 0795, 0796, 0952, 0955, 1483)	£49
CEF0695		REMOVABLE RH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ.	(NC 0740, 0750, 0760, 0770, 0780, 0790, 0952, 0955, 1480, 1483)	£53
CEF0696		REMOVABLE RH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ.	(NC 0740, 0750, 0760, 0770, 0780, 0790, 0952, 0955, 1483)	£53
CEF0732		REMOVABLE RH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION	(NC 0952, 1480, 1483)	£52

<b>ARMREST LH*</b>			
CEF0740		REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED DESK ADJ. 200 TO 280 MM HEIGHT (NC 1480, 1484)	○
CEF0750		REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED FULL LENGTH ADJ. 200 TO 280 MM HEIGHT	○
CEF0760		SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT (NC 0952, 0955, 1480, 1484)	◆
CEF0770		SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT (NC 0952, 0955)	○
CEF0780		SWING AWAY DETACHABLE LH (N°3) DESK ADJ. IN HEIGHT (170-295 MM) (NC 0952, 0955, 1480, 1484)	£49
CEF0790		SWING AWAY DETACHABLE LH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH (NC 0952, 0955, 1484)	£49
CEF0795		REMOVABLE LH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ. (NC 0952, 0955, 1480, 1484)	£53
CEF0796		REMOVABLE LH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. (NC 0952, 0955, 1484)	£53
CEF0832		REMOVABLE LH N°4 ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION ARMPAD HEMI MULTI POSITION (NC 0952, 1480, 1484)	£52
<b>SEAT CUSHION</b>			
CEF0849		BASIC 50 MM FOAM CUSHION - NYLON AND ANTI-SLIP COVER (NC 0952)	£21
<b>SEAT CUSHION DEPTH</b>			
CEF0878		FOR 400 MM SEAT DEPTH	○
CEF0879		FOR 450 MM SEAT DEPTH	○
<b>OTHER SEATING OPTIONS</b>			
CEF0882		LATERAL POSITIONING PADS (PAIR)	£54
CEF0905		NECKREST BLACK MULTI STRETCH PU COATED FABRIC (NC 0906, 0907, 0940, 0945)	£117
CEF0906		HEADREST BLACK MULTI STRETCH PU COATED FABRIC (NC 0907, 0940, 0945)	£140
CEF0907		HEADREST WITH CHEEK SUPPORT BLACK MULTI STRETCH PU COATED FABRIC (NC 0940, 0945)	£162
SWO3333		MULT ADJUSTABLE LATERALS	£252
CEF0930		CALF STRAP	◆
CEF0935		LONGER FOOTREST TUBE	£27
CEF0936		SHORTER FOOTREST SUPPORT	£37
CEF0940		BACKBRACE (NC 0945)	£31
CEF0945		FOLDABLE BACKBRACE	£34
CEF0952		IN-RIGGED ARMREST (NC 1480, 1483, 1484)	£37
CEF0955		W/O AUTOLOCK SAFETY DEVICE FOR ARMREST N°1 ONLY	○

COLOUR*			
CEF0961		SLATE GREY (RAL7024 FINE TEXTURE)	○
CEF0964		PEARL WHITE	○
CEF0965		ELECTRIC GREEN	○
CEF0966		MATT BLACK	○
CEF0967		BLUE POP	○
CEF0968		HYPNOTIC PURPLE	○
CEF0969		RED MATT	○
CEF0980		ANTHRACITE BLACK	○
CEF0991		SPECIAL COLOUR ON DEMAND ONLY	£149
SWO4444		DARK RED - SPECIAL	£149
SWO5555		PEARL GREY - SPECIAL	£149
SWO6666		INTENSE YELLOW	£149
SWO7777		VIBRANT ORANGE	£149
SWO8888		POLAR BLUE	£149
SWO9999		CANDY PINK	£149
CASTOR HOUSING*			
CEF0995		STANDARD	◆
CEF0996		LIGHT WEIGHT	£26
FRONT CASTORS*			
CEF1010		8" X 1.1/4" BLACK SOLID TYRE	◆
CEF1030		8" X 2" BLACK SOLID TYRE	£25
REAR WHEEL SIZE*			
CEF1062		12"	◆
REAR WHEELS*			
CEF1073		DRUM BRAKE - TRANSIT 12" (NC 1270)	○
CEF1074		TRANSIT COMPOSITE (BLACK PLASTIC HUB AND RIM) (NC 1280)	◆
TYRES*			
CEF1082		SOLID BLACK	◆
REAR WHEEL AXLE*			
CEF1170		FIXED AXLE	◆
REAR WHEEL POSITION*			
CEF1200		PASSIVE REAR WHEELS POSITION	◆
CEF1201		DUAL AMPUTEE REAR WHEEL POSITION (NC 1280)	£64
CHASSIS OPTION			
CEF1241		FRONT WHEELS EXTENSION	£19
BRAKES*			
CEF1270		PUSH TO LOCK	◆
CEF1280		ATTENDANT HUB BRAKES	£204
TRANSPORTATION KIT*			
CEF1345		TRANSPORTATION KIT	◆
SAFETY OPTIONS			
CEF1361		BODYPOINT 2-PT POSTURE BELT M (NC 1362)	£71
CEF1362		BODYPOINT 2-PT POSTURE BELT L	£71
CEF1370		REAR ANTI-TIPPERS PAIR (NC 1375, 1440, 1441)	£32
CEF1375		REAR ANTI-TIPPERS PAIR WITH TIPPING AID (NC 1440, 1441)	£64

OTHER OPTIONS			
CEF1430		CANE HOLDER	£20
CEF1440		TIPPING AID RH	£12
CEF1441		TIPPING AID LH	£12
CEF1480		ACRYLIC TRAY (NC 1483, 1484, 1485)	£102
CEF1483		HALF TABLE TRAY RH (NC 1484)	£79
CEF1484		HALF TABLE TRAY LH	£79
CEF1485		COMFORT PADDING FOR HALF TABLE TRAY	£39
SWO10001		ABDUCTED LEGRESTS (PR) 0° TO 30° - INTERNAL REFERENCE DE259618	£95

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

CEF0330	NYLON W/O HOOK AND LOOP STRIPS
CEF0350	PADDED NYLON
CEF0465	W/O BACK UPH
CEF0471	STANDARD NYLON SLING (JEMIMA)
CEF0474	LUMBAR SUPPORT BI MATERIAL PADDED COVER
CEF0477	CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER
CEF0478	CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE
CEF0481	CONTOUR ADJUSTABLE (7+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE
CEF0600	FLIP UP RH FOOTPLATE
CEF0610	FLIP UP DEPTH AND ANGLE AJUSTABLE RH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.
SWO1111	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CEF0611	FLIP UP ONE PIECE
CEF66666	NO FOOTREST
CEF0620	FLIP UP LH FOOTPLATE
CEF0630	FLIP UP DEPTH AND ANGLE AJUSTABLE LH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.
SWO2222	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CEF0740	REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED DESK ADJ. 200 TO 280 MM HEIGHT
CEF0750	REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED FULL LENGTH ADJ. 200 TO 280 MM HEIGHT
CEF0760	SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT
CEF0770	SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT
CEF0780	SWING AWAY DETACHABLE LH (N°3) DESK ADJ. IN HEIGHT (170-295 MM)
CEF0790	SWING AWAY DETACHABLE LH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH
CEF0795	REMOVABLE LH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ.
CEF0796	REMOVABLE LH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ.
CEF0878	FOR 400 MM SEAT DEPTH
CEF0879	FOR 450 MM SEAT DEPTH
CEF0882	LATERAL POSITIONING PADS (PAIR)
CEF0905	NECKREST BLACK MULTI STRETCH PU COATED FABRIC
CEF0906	HEADREST BLACK MULTI STRETCH PU COATED FABRIC
CEF0907	HEADREST WITH CHEEK SUPPORT BLACK MULTI STRETCH PU COATED FABRIC
CEF0930	CALF STRAP
CEF0935	LONGER FOOTREST TUBE
CEF0936	SHORTER FOOTREST SUPPORT
CEF0940	BACKBRACE
CEF0945	FOLDABLE BACKBRACE
CEF0952	IN-RIGGED ARMREST
CEF0955	W/O AUTOLOCK SAFETY DEVICE FOR ARMREST N°1 ONLY
CEF1074	TRANSIT COMPOSITE (BLACK PLASTIC HUB AND RIM)
CEF1270	PUSH TO LOCK
CEF1280	ATTENDANT HUB BRAKES

CEF1362	BODYPOINT 2-PT POSTURE BELT L
CEF1375	REAR ANTI-TIPPERS PAIR WITH TIPPING AID
CEF1440	TIPPING AID RH
CEF1441	TIPPING AID LH
CEF1480	ACRYLIC TRAY
CEF1483	HALF TABLE TRAY RH
CEF1484	HALF TABLE TRAY LH
CEF1485	COMFORT PADDING FOR HALF TABLE TRAY





For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):