



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

May 2026
Customer Ref:

| | |
|--------------------------|--------------------|
| Customer Ref* | Account No* |
| Contact No* | Order No* |
| Delivery Address* | Order Qty* |
| | Order date* |

Action 4NG Self Propel

Standard Features

- Max user weight: 125kg (19.7 stone)
- Crash tested ISO 7176 - Part 19
- Self Propulsion
- Choice of 5 seat widths (15" - 19") 380mm - 480mm
- Choice of 3 depths (16" - 20") 400mm - 505mm
- Choice of 5 seat heights (16" - 20") 400mm - 510mm
- Black nylon seat upholstery
- Fixed 7° 380mm backrest
- Standard padded nylon sling type backrest upholstery
- 80° swing in / swing out legrest hangers
- Flip up footplates
- Swing away detachable desk armrests
- Choice of 8 colours (with 7 special colours)
- Push or pull to lock brakes (configuration dependant)
- Transportation kit and lap belt (adjustable both sides)
- Adjustable heel straps
- User tool kit

Standard Base price: £700



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

| Code | Description | |
|----------------------|--------------------------------------|-----------------------|
| BASE PRODUCT* | | |
| CKE0000 | <input type="checkbox"/> ACTION 4NG | <input type="radio"/> |
| CKESPEC | <input type="checkbox"/> SPECIAL | POA |
| PROPULSION* | | |
| CKE0100 | <input type="checkbox"/> SELF PROPEL | £700 |

| SEAT WIDTH* | | | | |
|----------------------|--|--|---|------|
| CKF0119 | | 14" (355mm) SW | (NC 0180, 0190, 0195, 0290, 0295, 0360, 0364, 0367, 0368, 0375, 0400, 0420, 0450, 0470, 0480, 0491, 0500, 0659, 0661, 0681, 0691, 0693, 0703, 0656, 1270, 1271, 1272) | ○ |
| CKE0120 | | 15" (380mm) SW | (NC 0656) | ◆ |
| CKE0130 | | 16" (405mm) SW | | ◆ |
| CKE0140 | | 17" (430mm) SW | | ◆ |
| CKE0150 | | 18" (455mm) SW | | ◆ |
| CKE0160 | | 19" (480mm) SW | | ◆ |
| SEAT DEPTH* | | | | |
| CKE0180 | | 16" (400mm) SD | (NC 0702, 0703) | ◆ |
| CKE0190 | | 18" (450mm) SD | (NC 0701, 0703) | ◆ |
| CKE0195 | | 20" (500mm) SD | (NC 0295, 0659, 0661, 0681, 0691, 0701, 0702) | ◆ |
| CKE0197 | | 14" - 16" ADJ SD | (NC 0295) | £21 |
| CKE0200 | | 16" - 18" ADJ SD | (NC 0295, 0703) | £21 |
| SEAT HEIGHT* | | | | |
| CKE0210 | | 20" (510mm) SH | (NC 0556, 0616, 0890, 0895, 0946) | ◆ |
| CKE0220 | | 19" (485mm) SH | | ◆ |
| CKE0230 | | 18" (460mm) SH | | ◆ |
| CKE0240 | | 17" (430mm) SH | (NC 0850, 0870) | ◆ |
| CKE0250 | | 16" (410mm) SH | (NC 0556, 0616, 0770, 0850, 0870, 1130, 1196, 1245, 1246) | ◆ |
| SEAT UPHOLSTERY* | | | | |
| CKE0290 | | NYLON WITH HOOK AND LOOP STRIPS | | ◆ |
| CKE0295 | | PADDED NYLON | | £22 |
| CKE0298 | | LOOP BAND TENSION ADJUSTABLE IN DEPTH AND SHAPE | | ○ |
| BACKREST* | | | | |
| CKE0300 | | FIXED 7° BACKREST HT ADJ 15" - 18" FIXED P/HANDLES | (NC 0364, 0367, 0368, 0375, 0537, 0538, 0539, 0540, 0597, 0598, 0599, 0600, 0716, 0717, 0718, 0780) | ◆ |
| CKE0310 | | FIXED 7° BACKREST HT ADJ 15" - 18" HT ADJ P/HANDLES | (NC 0367, 0368, 0375, 0537, 0538, 0539, 0540, 0597, 0598, 0599, 0600, 0716, 0717, 0718, 0780) | £80 |
| CKE0311 | | FIXED 7° BACKREST 2 HEIGHTS 16" OR 17" | (NC 0360, 0367, 0782) | ◆ |
| CKE0312 | | FIXED 7° BACKREST 20" | (NC 0360, 0367, 0782) | ◆ |
| CKE0314 | | FIXED 7° BACKREST 17" ADJ P/HANDLES 6" | (NC 0360, 0367, 0782) | £43 |
| CKE0318 | | ANGLE ADJUSTABLE (-15° TO +15°) | (NC 0350, 0360, 0364, 0368, 0375, 0716, 0717, 0718, 0721) | £86 |
| CKE0320 | | RECLINING (RATCHET) 0/30° HEIGHTS 15"/16"/17" | (NC 0360, 0367, 0556, 0616, 0780, 0782, 1055, 1060) | £175 |
| CKE0334 | | RECLINING (GAS) 0/30° HEIGHTS 18"/19"/20" | (NC 0360, 0367, 0556, 0616, 0780, 0782, 1055, 1060) | £233 |
| CKE0345 | | RECLINING/FOLDING (GAS) 0/30° HEIGHTS 15"/16"/17" | (NC 0360, 0367, 0556, 0616, 0780, 0782, 1055, 1060) | £233 |
| CKE0341 | | FOLDING 7° BACKREST 17" | (NC 0360, 0367, 0782) | ◆ |
| CKE0342 | | FOLDING 7° BACKREST 20" | (NC 0360, 0367, 0375, 0782) | ◆ |
| BACKREST UPHOLSTERY* | | | | |
| CKE0350 | | WITHOUT BACK UPHOLSTERY | (NC 0721) | ◆ |
| CKE0360 | | STANDARD JEMIMA | (NC 0721) | ◆ |
| CKE0364 | | LUMBAR +1 STRAP | | ◆ |
| CKE0367 | | CONTOUR 3+1 STRAPS | | £55 |
| CKE0368 | | CONTOUR 3+1 STRAPS | (NC 0780, 0782) | £62 |
| CKE0375 | | CONTOUR 7+1 STRAPS | (NC 0780, 0782) | £120 |

| RH LEGREST HANGER* | | | | |
|---------------------------|--|--|---|------|
| CKE0390 | | 80° RH | (NC 1111, 66666) | ◆ |
| CKE0400 | | ELEVATING LEGREST RH | (NC 0491, 0760, 0770, 0771, 66666) | £116 |
| CKE0420 | | AMPUTEE LEGREST RH | (NC 0480, 0490, 0491, 0760, 0771, 1111) | £133 |
| LH LEGREST HANGER* | | | | |
| CKE0440 | | 80° LH | (NC 2222) | ◆ |
| CKE0450 | | ELEVATING LEGREST LH | (NC 0491, 0760, 0770, 0771, 55555) | £116 |
| CKE0470 | | AMPUTEE LEGREST LH | (NC 0491, 0500, 0510, 0760, 0771, 2222) | £133 |
| RH FOOTREST* | | | | |
| CKE0480 | | FLIP UP RH | | ◆ |
| CKE0490 | | DEPTH AND ANGLE ADJUSTABLE RH | | £43 |
| SWO1111 | | FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY) | | £43 |
| CKE0491 | | ONE PIECE | (NC 0500, 0510, 0770, 0771) | £78 |
| CKE66666 | | NO FOOTREST | | ○ |
| LH FOOTREST* | | | | |
| CKE0500 | | FLIP UP LH | | ◆ |
| CKE0510 | | DEPTH AND ANGLE ADJUSTABLE LH | | £43 |
| SWO2222 | | FLIP UP DEPTH AND ANGLE ADJUSTABLE LH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY) | | £43 |
| CKE55555 | | NO FOOTREST | | ○ |
| RH ARMREST* | | | | |
| CKE0520 | | REMOVABLE RH "RIO" (N°1) DESK ADJ. 200 TO 280 MM HEIGHT + SKIRTGUARD DEPTH AND HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED | (NC 0593, 0596, 0597, 0598, 0599, 0600, 0616, 1270, 1271) | ◆ |
| CKE0530 | | REMOVABLE RH "RIO"(N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT + SKIRTGUARD DEPTH AND HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED | (NC 0593, 0596, 0597, 0598, 0599, 0600, 0616) | ◆ |
| CKE0533 | | REMOVABLE RH (N°4) DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED | (NC 0580, 0590, 0597, 0598, 0599, 0600, 0616, 1270, 1271) | £33 |
| CKE0536 | | REMOVABLE RH (N°4) FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED | (NC 0580, 0590, 0597, 0598, 0599, 0600, 0616, 1271) | £33 |
| CKE0537 | | SWING AWAY DETACHABLE RH (N°2) DESK | (NC 0580, 0590, 0593, 0596, 0599, 0600, 0616, 1055, 1270, 1271) | ◆ |
| CKE0538 | | SWING AWAY DETACHABLE RH (N°2) FULL LENGTH | (NC 0580, 0590, 0593, 0596, 0599, 0600, 0616, 1055) | ◆ |
| CKE0539 | | SWING AWAY DETACHABLE RH (N°3) DESK ADJ. IN HEIGHT (170-295 MM) | (NC 0580, 0590, 0593, 0596, 0597, 0598, 0616, 1055, 1270, 1271) | £56 |
| CKE0540 | | SWING AWAY DETACHABLE RH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH | (NC 0580, 0590, 0593, 0596, 0597, 0598, 0616, 1055, 1271) | £60 |
| CKE0555 | | REMOVABLE RH (N°4) ARMPAD HEMI MULTI STRETCH PU COATED FABRIC MULTI POSITION AUTOLOCK SAFETY DEVICE INCLUDED | (NC 0616, 1270, 1271) | £60 |
| CKE0556 | | REMOVABLE MUDGUARD RH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) HEIGHT ADJUSTABLE (155 TO 220 MM) | (NC 0580, 0590, 0593, 0596, 0597, 0598, 0599, 0600, 0615, 0946, 1055, 1075, 1270, 1271) | £35 |

| LH ARMREST* | | | |
|--|--|---|------|
| CKE0580 | | REMOVABLE LH "RIO" (N°1) DESK ADJ. 200 TO 280 MM HEIGHT + SKIRTGUARD DEPTH AND HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED (NC 1270, 1272) | ◆ |
| CKE0590 | | REMOVABLE LH "RIO"(N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT + SKIRTGUARD DEPTH AND HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED | ◆ |
| CKE0593 | | REMOVABLE LH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ. (NC 1270, 1272) | £33 |
| CKE0596 | | REMOVABLE LH (N°4) FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED (NC 1272) | £33 |
| CKE0597 | | SWING AWAY DETACHABLE LH (N°2) DESK (NC 1055, 1270, 1272) | ◆ |
| CKE0598 | | SWING AWAY DETACHABLE LH (N°2) FULL LENGTH (NC 1055) | ◆ |
| CKE0599 | | SWING AWAY DETACHABLE LH (N°3) DESK ADJ. IN HEIGHT (170-295 MM) (NC 1055, 1270, 1272) | £56 |
| CKE0600 | | SWING AWAY DETACHABLE LH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH (NC 1055, 1272) | £60 |
| CKE0615 | | REMOVABLE LH (N°4) ARMPAD HEMI MULTI STRETCH PU COATED FABRIC MULTI POSITION AUTOLOCK SAFETY DEVICE INCLUDED (NC 1270, 1272) | £60 |
| CKE0616 | | REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) HEIGHT ADJUSTABLE (155 TO 220 MM) (NC 0946, 1270, 1272) | £35 |
| SEAT CUSHION | | | |
| CKE0654 | | BASIC 50MM - NYLON AND ANT-SLIP COVER (NC 0706, 0707, 0708, 0709, 0706, 0707, 0708, 0709) | £24 |
| CKE0659 | | FLO-TECH LITE CUSHION - Low Risk (NC 0703, 0707) | £121 |
| CKE0661 | | FLO-TECH IMAGE CUSHION - High Risk (NC 0703, 0707) | £169 |
| CKE0681 | | FLO-TECH CONTOUR VISCO NG CUSHION - HIGH RISK (NC 0703, 0706, 0707, 0708, 0709) | £288 |
| CKE0691 | | FLO-TECH SOLUTION XTRA CUSHION HIGH RISK (NC 0703) | £398 |
| CKE0693 | | MATRX LIBRA CUSHION - Very High Risk | £410 |
| SEAT CUSHION DEPTH | | | |
| CKE0701 | | 16" SD | ◆ |
| CKE0702 | | 18" SD | ◆ |
| CKE0703 | | 20" SD (NC 0656) | ◆ |
| 1st POSITIONING ACCESSORIES MATRX LIBRA & FLOTECH | | | |
| CKE0706 | | TAPERED WEDGES (PAIR) | £30 |
| CKE0707 | | ABDUCTOR | £23 |
| CKE0708 | | LATERAL WEDGE | £20 |
| CKE0709 | | ANTERIOR SEAT WEDGE | £27 |
| 2nd POSITIONING ACCESSORIES MATRX LIBRA & FLOTECH | | | |
| CKE0706 | | TAPERED WEDGES (PAIR) | £30 |
| CKE0707 | | ABDUCTOR | £23 |
| CKE0708 | | LATERAL WEDGE | £20 |
| CKE0709 | | ANTERIOR SEAT WEDGE | £27 |

| OTHER SEATING OPTIONS | | | |
|-----------------------|--|---|------|
| CKE0656 | | SEAT RIGIDIZER | £39 |
| CKE0716 | | NECKREST INCLUDING BACKBRACE (NC 0717, 0718, 0780, 0782) | £131 |
| CKE0717 | | HEADREST INCLUDING BACKBRACE (NC 0718, 0780, 0782) | £159 |
| CKE0718 | | HEADREST WITH CHEEK SUPPORT (NC 0780, 0782) | £183 |
| CKE0721 | | LATERAL OSITIONING PADS PAIR | £61 |
| SWO3333 | | MULT ADJUSTABLE LATERALS | £284 |
| CKE0760 | | CALF STRAP | ◆ |
| CKE0770 | | LONGER FOOTREST TUBE (NC 0771) | £30 |
| CKE0771 | | SHORTER FOOTREST TUE | £42 |
| CKE0780 | | BACKBRACE (NC 0782) | £35 |
| CKE0782 | | FOLDABLE BACKBRACE | £37 |
| COLOUR* | | | |
| CKE0801 | | SLATE GREY | ◆ |
| CKE0804 | | PEARL WHITE | ◆ |
| CKE0807 | | MATT BLACK | ◆ |
| CKE0811 | | RED MATT | ◆ |
| CKE0815 | | WATER BLUE | ◆ |
| CKE0816 | | BLUE | ◆ |
| CKE0817 | | DARK PURPLE | ◆ |
| CKE0820 | | ANTHRACITE BLACK | ◆ |
| CKE0821 | | SPECIAL DEMAND | £169 |
| SWO4444 | | DARK RED - SPECIAL | £169 |
| SWO5555 | | PEARL GREY - SPECIAL | £169 |
| SWO6666 | | INTENSE YELLOW | £169 |
| SWO7777 | | VIBRANT ORANGE | £169 |
| SWO8888 | | POLAR BLUE | £169 |
| SWO9999 | | CANDY PINK | £169 |
| FRONT CASTORS* | | | |
| CKE0850 | | 8" X 1.1/4" BLACK SOLID | ◆ |
| CKE0870 | | 8" X 2" BLACK SOLID | £28 |
| CKE0890 | | 6" X 1.1/4" BLACK SOLID | ○ |
| CKE0895 | | 6" X 1.1/4" BLACK SOLID MEDIUM HARD | £28 |
| REAR WHEEL SIZE* | | | |
| CKE0946 | | 22" (NC 0982, 1110) | ○ |
| CKE0948 | | 24" (NC 1111) | ◆ |
| REAR WHEEL* | | | |
| CKE0962 | | STANDARD (NC 0982, 1130) | ◆ |
| CKE0963 | | DRUM BRAKE (NC 0982, 1120, 1125) | ○ |
| CKE0966 | | LIGHTWEIGHT (NC 1075, 1090, 1110, 1111, 1130) | £59 |
| TYRES* | | | |
| CKE0972 | | BLACK SOLID 1"3/8 | ◆ |
| CKE0974 | | PNEUMATIC BLACK 1"3/8 | £21 |
| CKE0978 | | PNEUMATIC BLACK 1" | ○ |
| CKE0982 | | PNEUMATIC AND PUNCTURE RESISTANT MAXXIS APLHA PROTECT TYRE 1" / 25 MM | £39 |

| HANDRIM* | | | |
|-----------------------------|--------------------------|---|-------------------------------|
| CKE1020 | <input type="checkbox"/> | ALUMINIUM H/R | ◆ |
| CKE1030 | <input type="checkbox"/> | H/R BLACK ANTI-SLIP | £34 |
| REAR WHEEL AXLE* | | | |
| CKE1050 | <input type="checkbox"/> | Q/R AXLE | ◆ |
| REAR WHEEL POSITION* | | | |
| CKE1055 | <input type="checkbox"/> | ACTIVE REAR WHEELS POSITION 52,5 MM FRONTWARD/BACKREST CANE AXLE | (NC 1130) ◆ |
| CKE1060 | <input type="checkbox"/> | STANDARD REAR WHEELS POSITION 35 MM FRONTWARD/BACKREST CANE AXLE | ◆ |
| CKE1070 | <input type="checkbox"/> | PASSIVE REAR WHEELS POSITION 35 MM BACKWARD/BACKREST CANE AXLE | ◆ |
| CKE1075 | <input type="checkbox"/> | DUAL AMPUTEE POSITION 52,5 MM BACKWARD/BACKREST CANE AXLE,ADJUSTABLE TO 75 MM BY THE DEALER | (NC 1125, 1130) £72 |
| PUMP | | | |
| CKE1090 | <input type="checkbox"/> | PUMP | £7 |
| CHASSIS OPTION | | | |
| CKE1101 | <input type="checkbox"/> | FRONT WHL EXTENSION | £22 |
| CKE1110 | <input type="checkbox"/> | SPOKE GUARDS PAIR (24") BLACK | £78 |
| CKE1111 | <input type="checkbox"/> | SPOKE GUARDS PAIR (22"/20") CLEAR | £98 |
| BRAKES* | | | |
| CKE1120 | <input type="checkbox"/> | PUSH TO LOCK BRAKES | ◆ |
| CKE1125 | <input type="checkbox"/> | PULL TO LOCK BRAKES | ◆ |
| CKE1130 | <input type="checkbox"/> | ATTENDANT HUB BRAKES (INCLUDING PUSH TO LOCK BRAKES) | £231 |
| TRANSPORTATION KIT* | | | |
| CKE1175 | <input type="checkbox"/> | TRANSPORTATION BRACKETS | ◆ |
| SAFETY OPTIONS* | | | |
| CKE1190 | <input type="checkbox"/> | LAP BELT | (NC 1189, 1192, 1193, 1194) ◆ |
| CKE1189 | <input type="checkbox"/> | BODYPOINT 2 POINT CENTER-PULL PUSH- BUTTON MEDIUM | (NC 1192, 1193, 1194) £83 |
| CKE1192 | <input type="checkbox"/> | BODYPOINT 2 POINT CENTER-PULL PUSH- BUTTON LARGE | (NC 1193, 1194) £83 |
| CKE1193 | <input type="checkbox"/> | BODYPOINT 4 POINT CENTER-PULL PUSH- BUTTON MEDIUM | (NC 1194) £98 |
| CKE1194 | <input type="checkbox"/> | BODYPOINT 4 POINT CENTER-PULL PUSH- BUTTON LARGE | £98 |
| CKE1196 | <input type="checkbox"/> | ANTI-TIPPER (PAIR) | (NC 1200, 1245, 1246) £36 |
| CKE1200 | <input type="checkbox"/> | REAR ANTI-TIPPERS PAIR WITH TIPPING AID | (NC 1245, 1246) £72 |
| CKE1210 | <input type="checkbox"/> | REFLECTORS KIT | £8 |
| CKE1220 | <input type="checkbox"/> | BRAKE EXTENSION | £14 |

| OTHER OPTIONS | | | |
|---------------|--|---|--------|
| CKE1240 | | CANE HOLDER | £23 |
| CKE1245 | | TIPPING AID RH | £14 |
| CKE1246 | | TIPPING AID LH | £14 |
| CKE1270 | | ACRYLIC TRAY (NC 1271, 1272, 1273) | £116 |
| CKE1271 | | HALF TRAY RH (NC 1272) | £89 |
| CKE1272 | | HALF TRAY LH | £89 |
| CKE1273 | | COMFORT PADDING FOR HALF TRAY | £44 |
| SWO10002 | | ABDUCTED LEGRESTS (PR) 0° TO 30° - INTERNAL REFERENCE DE259618 | £107 |
| 1595002 | | ViaGo Push aid. Includes as standard: Drive Unit, Charger, Control unit and Brackets (Note: 10 minute fitting time) | £1,995 |

| | |
|--------------------|--|
| TOTAL PRICE | |
|--------------------|--|

Non Conforming Options:

| | |
|----------|--|
| CKE0180 | 16" (400mm) SD |
| CKE0190 | 18" (450mm) SD |
| CKE0195 | 20" (500mm) SD |
| CKE0290 | NYLON WITH HOOK AND LOOP STRIPS |
| CKE0295 | PADDED NYLON |
| CKE0350 | WITHOUT BACK UPHOLSTERY |
| CKE0360 | STANDARD JEMIMA |
| CKE0364 | LUMBAR +1 STRAP |
| CKE0367 | CONTOUR 3+1 STRAPS |
| CKE0368 | CONTOUR 3+1 STRAPS |
| CKE0375 | CONTOUR 7+1 STRAPS |
| CKE0400 | ELEVATING LEGREST RH |
| CKE0420 | AMPUTEE LEGREST RH |
| CKE0450 | ELEVATING LEGREST LH |
| CKE0470 | AMPUTEE LEGREST LH |
| CKE0480 | FLIP UP RH |
| CKE0490 | DEPTH AND ANGLE ADJUSTABLE RH |
| SWO1111 | FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY) |
| CKE0491 | ONE PIECE |
| CKE66666 | NO FOOTREST |
| CKE0500 | FLIP UP LH |
| CKE0510 | DEPTH AND ANGLE ADJUSTABLE LH |
| SWO2222 | FLIP UP DEPTH AND ANGLE ADJUSTABLE LH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY) |
| CKE55555 | NO FOOTREST |
| CKE0537 | SWING AWAY DETACHABLE RH (N°2) DESK |
| CKE0538 | SWING AWAY DETACHABLE RH (N°2) FULL LENGTH |
| CKE0539 | SWING AWAY DETACHABLE RH (N°3) DESK ADJ. IN HEIGHT (170-295 MM) |
| CKE0540 | SWING AWAY DETACHABLE RH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH |
| CKE0556 | REMOVABLE MUDGUARD RH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) HEIGHT ADJUSTABLE (155 TO 220 MM) |
| CKE0580 | REMOVABLE LH "RIO" (N°1) DESK ADJ. 200 TO 280 MM HEIGHT + SKIRTGUARD DEPTH AND HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED |
| CKE0590 | REMOVABLE LH "RIO"(N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT + SKIRTGUARD DEPTH AND HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED |
| CKE0593 | REMOVABLE LH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ. |
| CKE0596 | REMOVABLE LH (N°4) FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED |
| CKE0597 | SWING AWAY DETACHABLE LH (N°2) DESK |
| CKE0598 | SWING AWAY DETACHABLE LH (N°2) FULL LENGTH |
| CKE0599 | SWING AWAY DETACHABLE LH (N°3) DESK ADJ. IN HEIGHT (170-295 MM) |
| CKE0600 | SWING AWAY DETACHABLE LH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH |
| CKE0615 | REMOVABLE LH (N°4) ARMPAD HEMI MULTI STRETCH PU COATED FABRIC MULTI POSITION AUTOLOCK SAFETY DEVICE INCLUDED |
| CKE0616 | REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) HEIGHT ADJUSTABLE (155 TO 220 MM) |
| CKE0659 | FLO-TECH LITE CUSHION - Low Risk |

| | |
|---------|---|
| CKE0661 | FLO-TECH IMAGE CUSHION - High Risk |
| CKE0681 | FLO-TECH CONTOUR VISCO NG CUSHION - HIGH RISK |
| CKE0691 | FLO-TECH SOLUTION XTRA CUSHION HIGH RISK |
| CKE0693 | MATRX LIBRA CUSHION - Very High Risk |
| CKE0701 | 16" SD |
| CKE0702 | 18" SD |
| CKE0703 | 20" SD |
| CKE0706 | TAPERED WEDGES (PAIR) |
| CKE0707 | ABDUCTOR |
| CKE0708 | LATERAL WEDGE |
| CKE0709 | ANTERIOR SEAT WEDGE |
| CKE0706 | TAPERED WEDGES (PAIR) |
| CKE0707 | ABDUCTOR |
| CKE0708 | LATERAL WEDGE |
| CKE0709 | ANTERIOR SEAT WEDGE |
| CKE0656 | SEAT RIGIDIZER |
| CKE0716 | NECKREST INCLUDING BACKBRACE |
| CKE0717 | HEADREST INCLUDING BACKBRACE |
| CKE0718 | HEADREST WITH CHEEK SUPPORT |
| CKE0721 | LATERAL OSITIONING PADS PAIR |
| CKE0760 | CALF STRAP |
| CKE0770 | LONGER FOOTREST TUBE |
| CKE0771 | SHORTER FOOTREST TUE |
| CKE0780 | BACKBRACE |
| CKE0782 | FOLDABLE BACKBRACE |
| CKE0850 | 8" X 1.1/4" BLACK SOLID |
| CKE0870 | 8" X 2" BLACK SOLID |
| CKE0890 | 6" X 1.1/4" BLACK SOLID |
| CKE0895 | 6" X 1.1/4" BLACK SOLID MEDIUM HARD |
| CKE0946 | 22" |
| CKE0982 | PNEUMATIC AND PUNCTURE RESISTANT MAXXIS APLHA PROTECT TYRE 1" / 25 MM |
| CKE1055 | ACTIVE REAR WHEELS POSITION 52,5 MM FRONTWARD/BACKREST CANE AXLE |
| CKE1060 | STANDARD REAR WHEELS POSITION 35 MM FRONTWARD/BACKREST CANE AXLE |
| CKE1075 | DUAL AMPUTEE POSITION 52,5 MM BACKWARD/BACKREST CANE AXLE,ADJUSTABLE TO 75 MM BY THE DEALER |
| CKE1090 | PUMP |
| CKE1110 | SPOKE GUARDS PAIR (24") BLACK |
| CKE1111 | SPOKE GUARDS PAIR (22"/20") CLEAR |
| CKE1120 | PUSH TO LOCK BRAKES |
| CKE1125 | PULL TO LOCK BRAKES |
| CKE1130 | ATTENDANT HUB BRAKES (INCLUDING PUSH TO LOCK BRAKES) |
| CKE1189 | BODYPOINT 2 POINT CENTER-PULL PUSH-BUTTON MEDIUM |
| CKE1192 | BODYPOINT 2 POINT CENTER-PULL PUSH-BUTTON LARGE |
| CKE1193 | BODYPOINT 4 POINT CENTER-PULL PUSH-BUTTON MEDIUM |
| CKE1194 | BODYPOINT 4 POINT CENTER-PULL PUSH-BUTTON LARGE |

| | |
|---------|---|
| CKE1196 | ANTI-TIPPER (PAIR) |
| CKE1200 | REAR ANTI-TIPPERS PAIR WITH TIPPING AID |
| CKE1245 | TIPPING AID RH |
| CKE1246 | TIPPING AID LH |
| CKE1270 | ACRYLIC TRAY |
| CKE1271 | HALF TRAY RH |
| CKE1272 | HALF TRAY LH |
| CKE1273 | COMFORT PADDING FOR HALF TRAY |



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):