



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

MARCH 2021
Customer Ref:

Customer Ref*, Contact No*, Delivery Address*, Account No*, Order No*, Order Qty*, Order date*

Fox Standard

Standard Features

- Max User Weight 127Kg (20 Stone)
Canvas sling seat with Black upholstery and backrest height 490mm
Seat widths from 350mm - 500mm
Seat depths from 370mm - 460mm
Seat to floor heights 480mm or 510mm
Detachable T-armrest with standard armpads
Manually settable seat Tilt (Set at 4°, 0° - 8° achievable)
Settable backrest Recline (85°- 115°)
Standard 80° legrest hangers with fixed flip-up footplates
LiNX control REM110 with standard bracket
Choice of 5 colour options
Solid Grey castors 200mm x 50mm (8" x 2"), Silver Rim
Grey drive wheels, puncture proof solid inlay, (12 1/2" x 2 1/4"), Silver Rim (aluminium)
6km/h (4mph) motors with 34Ah AGM battery
8A (200 - 250V) charger
Bolted on anti tipper
Wheel lock / parking brake
4 Point tie down system
Lapbelt with autolock (adjustable both sides)
Turning radius 1600mm - 1630mm depending on configuration



Standard Base price: £2,514

Please make your selection with a cross in the box [] Free of charge [O] standard feature [◆] POA special request

Customer weight = KG

Table with columns Code and Description. Rows include BASE PRODUCT* (AXA1140 FOX £2,514), CHAIR TYPE* (AXA0000 Fox, AXZ0000 Special Fox POA), and SEATING TYPE* (AXA1140 Fox seat unit).

SEAT WIDTH*			
AXA0368	<input type="checkbox"/>	350mm (14")	(NC 0606) ○
AXA0361	<input type="checkbox"/>	400mm (16")	○
AXA0345	<input type="checkbox"/>	450mm (18")	○
AXA0350	<input type="checkbox"/>	500mm (20")	(NC 1187, 1184, 2174, 2175, 2176) ○
SEAT DEPTH*			
AXA0324	<input type="checkbox"/>	370mm	○
AXA0308	<input type="checkbox"/>	400mm	○
AXA0327	<input type="checkbox"/>	430mm	○
AXA0319	<input type="checkbox"/>	460mm	○
SEAT TYPE*			
AXA1187	<input type="checkbox"/>	Sling seat	(NC 0350, 2147) ◆
AXA1188	<input type="checkbox"/>	Seat board	○
SEAT CUSHION*			
AXA1409	<input type="checkbox"/>	Without seat cushion	◆
AXA1336	<input type="checkbox"/>	Seat Cushion 50mm thick	£40
AXA2147	<input type="checkbox"/>	Matrx Libra	(NC 1187) £371
BACKREST UPHOLSTERY / CUSHION*			
AXA1184	<input type="checkbox"/>	Sling upholstery	(NC 0350) ◆
AXA1178	<input type="checkbox"/>	Tension adjustable backrest	£46
AXA2174	<input type="checkbox"/>	Matrx PB 400mm	(NC 0350) £471
AXA2175	<input type="checkbox"/>	Matrx PB 500mm	(NC 0350) £471
AXA2176	<input type="checkbox"/>	Matrx Elite	(NC 0350) £471
AXA2177	<input type="checkbox"/>	Matrx Elite TR	£471
AXA2178	<input type="checkbox"/>	Matrx Elite Deep	£524
BACKREST HEIGHT for Sling Upholstery and Tension Adjustable			
AXA0449	<input type="checkbox"/>	490mm	◆
AXA0457	<input type="checkbox"/>	570mm	£46
BACKREST SUPPORT FOR MATRX BACKS			
AXA0413	<input type="checkbox"/>	Backrest frame short (470mm)	○
AXA0414	<input type="checkbox"/>	Backrest frame long (540mm)	○
ARMREST*			
AXA0606	<input type="checkbox"/>	T-armrest	(NC 0368, 0645) ◆
AXA0677	<input type="checkbox"/>	In-rigged armrests	(NC 0650) £68
ARMPAD*			
AXA0650	<input type="checkbox"/>	Long Ergonomic (397mm)	(NC 0677) ○
AXA0647	<input type="checkbox"/>	Standard (325mm)	◆
AXA0649	<input type="checkbox"/>	Long (375mm)	○
AXA0645	<input type="checkbox"/>	Short (220mm)	(NC 0606) ○
HEADREST / NECKREST (multi stretch PU coated fabric)			
AXA1350	<input type="checkbox"/>	Neckrest	£145
AXA1432	<input type="checkbox"/>	Headrest	£192
SEAT HEIGHT SETTING*			
AXA0312	<input type="checkbox"/>	Low Seat To floor height 480mm	○
AXA0344	<input type="checkbox"/>	High Seat To floor height 510mm	○
SEAT TILT*			
AXA1002	<input type="checkbox"/>	Manually settable	◆

BACKREST RECLINER*			
AXA0961	<input type="checkbox"/>	Settable in steps (85° /95° /105° /115°)	◆
FOOTREST HANGER TYPE*			
AXA0526	<input type="checkbox"/>	Standard 80° hangers	◆
AXA0529	<input type="checkbox"/>	Manual elevating legrests	£181
AMPUTEE SUPPORT			
AXA0554	<input type="checkbox"/>	RH	£166
AXA0555	<input type="checkbox"/>	LH	£166
FOOTPLATE TYPE*			
AXA0512	<input type="checkbox"/>	Fixed	◆
AXA0513	<input type="checkbox"/>	Angle and depth adjustable	£76
CONTROLS*			
AXA2201	<input type="checkbox"/>	REM110	◆
CONTROL BRACKETS*			
AXA1207	<input type="checkbox"/>	Standard bracket RH (NC 1404)	○
AXA1208	<input type="checkbox"/>	Standard bracket LH (NC 1405)	○
AXA1707	<input type="checkbox"/>	Swing away RH (NC 1404)	£113
AXA1708	<input type="checkbox"/>	Swing away LH (NC 1405)	£113
AXA2227	<input type="checkbox"/>	RH bracket mounted for Attendant (NC 2219, 2220)	£54
AXA2228	<input type="checkbox"/>	LH bracket mounted for Attendant (NC 2219, 2220)	£54
DUAL CONTROLS			
AXA2220	<input type="checkbox"/>	Dual control REM050	£336
AXA2219	<input type="checkbox"/>	Compact (ACU200)	£336
DUAL CONTROL BRACKET			
AXA1207	<input type="checkbox"/>	Standard bracket RH	○
AXA1208	<input type="checkbox"/>	Standard bracket LH	○
BODYPOINT JOYSTICK KNOBS			
AXA2721	<input type="checkbox"/>	3" U-shaped, joystick handle	£62
AXA2722	<input type="checkbox"/>	4" U-shaped, joystick handle	£62
AXA2723	<input type="checkbox"/>	Domed rubber, joystick handle	£62
AXA2724	<input type="checkbox"/>	Mushroom grooved, joystick handle	£62
WHEELBASE SETTING*			
AXA1226	<input type="checkbox"/>	Recommended factory setting	◆
SHROUD COLOUR*			
AXA0201	<input type="checkbox"/>	Pearl Night Blue	○
AXA0233	<input type="checkbox"/>	Metallic Silver	○
AXA0228	<input type="checkbox"/>	Silky Matt Black	○
AXA0204	<input type="checkbox"/>	Transparent Cherry Red	○
AXA0222	<input type="checkbox"/>	White	○
FRONT CASTORS*			
AXA0802	<input type="checkbox"/>	Castor Grey 200mm x 50mm (8" x 2") PU tyre	○
AXA0803	<input type="checkbox"/>	Castor Grey 200mm x 50mm (8" x 2") Solid Inlay	◆
DRIVE WHEELS*			
AXA0709	<input type="checkbox"/>	Drive wheel Grey, Puncture Proof Solid Inlay, (12 1/2 * 2 1/4")	○
AXA0710	<input type="checkbox"/>	Drive wheel Grey, Puncture Proof Solid Inlay, (12 1/2 * 2 1/4")	◆

MOTOR TYPE*			
AXA1294	<input type="checkbox"/>	3 km/h (2mph)	○
AXA1282	<input type="checkbox"/>	6 km/h (4mph)	◆
BATTERY*			
AXA1486	<input type="checkbox"/>	34 Ah AGM	◆
AXA1255	<input type="checkbox"/>	40 Ah Gel	£103
CHARGER*			
AXA1343	<input type="checkbox"/>	OptiCharge 8A (200 - 250V)	◆
AXA1313	<input type="checkbox"/>	OptiCharge 10A (200 - 250V)	£130
PARKING BREAK*			
AXA1206	<input type="checkbox"/>	Manual wheel locks	◆
KERB RIDER			
AXA1244	<input type="checkbox"/>	Kerbrider	£181
TRANSPORTATION KIT*			
AXA1494	<input type="checkbox"/>	4-Point tie-down	◆
LAP BELT*			
AXA1329	<input type="checkbox"/>	Lap belt with auto-lock	◆
OTHER OPTIONS			
AXA1312	<input type="checkbox"/>	RH	£36
AXA1363	<input type="checkbox"/>	LH	£36
AXA1404	<input type="checkbox"/>	Tray, Swing Away RH	£188
AXA1405	<input type="checkbox"/>	Tray, Swing Away LH	£188
AXA1454	<input type="checkbox"/>	Oxygen bottle holder	£92
AXA1490	<input type="checkbox"/>	Castor lock levers	£28

TOTAL PRICE	
--------------------	--

Non Conforming Options:

AXA0368	350mm (14")
AXA0350	500mm (20")
AXA1187	Sling seat
AXA2147	Matrx Libra
AXA1184	Sling upholstery
AXA2174	Matrx PB 400mm
AXA2175	Matrx PB 500mm
AXA2176	Matrx Elite
AXA0606	T-armrest
AXA0677	In-rigged armrests
AXA0650	Long Ergonomic (397mm)
AXA0645	Short (220mm)
AXA2220	Dual control REM050
AXA2219	Compact (ACU200)
AXA1404	Tray, Swing Away RH
AXA1405	Tray, Swing Away LH



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions
Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):