



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Power Retail Prescription Form

01 DECEMBER 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Storm 4 MAX

Standard Features

- Max User Weight 220KG (34.5 Stone)
- Max Seating
- Seat widths from 490mm - 610mm
- Seat depths from 560mm - 630mm
- Black Multi Stretch PU Coated Wipe Down Fabric TR26
- Seat heights from 480mm - 510mm
- Manual Settable Tilt with Manual Settable Recline
- Reinforced Armrests
- Reinforced Legrests with angle adjustable flip-back footboard
- Headrest, sideways and depth adjustable
- LiNX DLX-REM110 with standard bracket
- Choice of 3 colours
- Puncture Proof Castors
- Puncture Proof Drive Wheels
- 6 km/h (4mph) Motors
- 8A (200 - 250V) Charger
- 4 Point tie down system
- 73.5 Ah (Gel) batteries
- Bolted on anti tipers

Standard Base price: £7,789



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
AWA0000	Storm 4 MAX	£7,789
AWZ0000	Special Storm 4 MAX (Extended Lead Time). Please provide details in comments section.	POA
SEATING TYPE*		
AWA1145	MAX SEAT UNIT	<input checked="" type="checkbox"/>
BACK CUSHION*		
AWA1106	Tension adjustable 630mm - 620mm high with SW55 / 610mm 570mm High with SW490mm	<input checked="" type="checkbox"/>

SEAT WIDTH*			
AWA0378	<input type="checkbox"/>	Seat Width 550mm	<input type="radio"/>
AWA0380	<input type="checkbox"/>	Seat Width 610mm (NC 0503)	<input type="radio"/>
SEAT CUSHION*			
AWA1459	<input type="checkbox"/>	Seat cushion Max (550mm x 635mm)	<input type="radio"/>
AWA1460	<input type="checkbox"/>	Seat cushion Max (610mm x 635mm)	<input type="radio"/>
SEAT DEPTH*			
AWA0396	<input type="checkbox"/>	Seat Depth 560mm	<input type="radio"/>
AWA0397	<input type="checkbox"/>	Seat Depth 595mm	<input type="radio"/>
ADAPTER TO REDUCE SEAT DEPTH			
AWA0399	<input type="checkbox"/>	Adapter to shorten the seat depth by 130mm	£49
UPHOLSTERY COLOUR*			
AWA0278	<input type="checkbox"/>	Black Multi Stretch PU Coated Wipe Down Fabric TR26	◆
ARMRESTS INCLUDING REINFORCED DEPTH ADJUSTABLE BRACKET*			
AWA0635	<input type="checkbox"/>	Fixed armrest 250mm - 340mm height adjustable removable	<input type="radio"/>
ARMPADS*			
AWA0647	<input type="checkbox"/>	Standard Armpad 325mm	<input type="radio"/>
AWA0649	<input type="checkbox"/>	Armpad extended plus 50mm	<input type="radio"/>
PUSH BAR*			
AWA1307	<input type="checkbox"/>	Push bar - height adjustable	<input type="radio"/>
HEADREST / NECKREST (Adequate headrest option must be selected)*			
AWA1432	<input type="checkbox"/>	Headrest sideways and depth adjustable	◆
AWA1427	<input type="checkbox"/>	Headrest with cheek support sideways adjustable	£195
AWA1495	<input type="checkbox"/>	Equivalent non Invacare Headrest supporting the requirements of ISO7176-19:2001. Provided independently from this order	<input type="radio"/>
SEAT HEIGHT*			
AWA1025	<input type="checkbox"/>	Fixed seat height 480mm	<input type="radio"/>
AWA1006	<input type="checkbox"/>	Fixed seat height 510mm	<input type="radio"/>
MANUAL / POWERED SEATING OPTIONS*			
AWA1002	<input type="checkbox"/>	Manual Seat Tilt by spindle	<input type="radio"/>
AWA1014	<input type="checkbox"/>	Electrical Tilt 0° to 15°	£514
BACKREST RECLINER*			
AWA0958	<input type="checkbox"/>	Settable recliner by spindle (0° - 25°)	<input type="radio"/>
AWA0947	<input type="checkbox"/>	Electrical recliner (with high pivot, shearing compensation)	£454
LEG REST TYPE*			
AWA0564	<input type="checkbox"/>	Side Mounted: Reinforced hanger fixed angle (not detachable without tools) (NC 0522)	<input type="radio"/>
AWA0503	<input type="checkbox"/>	Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools (NC 0599)	£102
CALF PADS			
AWA0583	<input type="checkbox"/>	Calf Pads adjustable	<input type="radio"/>
FOOTBOARD*			
AWA0599	<input type="checkbox"/>	Angle adjustable flip back heavy duty footboard (NC 0503)	◆
AWA0522	<input type="checkbox"/>	Angle adjustable footboard (for max user weight 200kg) (NC 0564)	<input type="radio"/>

LINX CONTROLS*				
AWA2201		DLX-REM110 (Drive Only)	(NC 1437)	◆
AWA2203		DLX-REM211 (Drive & Actuators)	(NC 1437)	○
AWA2204		DLX-REM216 (Drive, Actuators & Lights)		£130
CONTROL BRACKETRY*				
AWA1207		Standard Bracket RH		○
AWA1208		Standard Bracket LH		○
AWA1707		Swing Away RH Bracket		£113
AWA1708		Swing Away LH Bracket		£113
AWA1709		Matrx Resolve Swing in/out Bracket RH		£219
AWA1710		Matrx Resolve Swing in/out Bracket LH		£219
AWA2227		Attendant Mounted RH	(NC 2219, 2220)	£54
AWA2228		Attendant Mounted LH	(NC 2219, 2220)	£54
DUAL CONTROLS				
AWA2220		LiNX DLX-REM050 (Dual Control used as secondary remote)	(NC 2227, 2228, 2245, 1221, 1222)	£338
AWA2219		LiNX DLX-ACU200 (Compact Dual Control)	(NC 2227, 2228, 2245, 2201, 1221, 1222)	£338
AWA2245		LiNX DLX-IDC (Intuitive Dual Control)	(NC 2227, 2228, 2220, 2219, 2201, 1207, 1208)	£641
DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)				
AWA1207		Dual Control Standard Bracket RH	(NC 2227, 2228, 2245)	○
AWA1208		Dual Control Standard Bracket LH		○
AWA1221		IDC Dual Control Detachable Bracket RH	(NC 2227, 2228, 2220, 2219)	○
AWA1222		IDC Dual Control Detachable Bracket LH	(NC 2227, 2228, 2220, 2219)	○
LINX 10-WAY SWITCH				
AWA2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access		£431
LINX 10-WAY SWITCH BRACKET				
AWA2217		Bracket RH, mounted for user		○
AWA2218		Bracket LH, mounted for user		○
AWA2227		Bracket RH, mounted for attendant		£54
AWA2228		Bracket LH, mounted for attendant		£54
LINX ACCESSORY CABLE				
AWA2252		Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch)		£86
REMOTE STOP SWITCH				
AWZ2258		Remote Stop Switch (ASL504)		£897
FURTHER LINX COMPONENTS / ACCESSORIES				
AWA2302		Not Necessary	(NC 1686, 1687, 1690, 2254, 1711, 1712, 1713, 1714, 1715, 1716)	◆
AWA2300		To be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected)		○
BUDDY BUTTONS				
AWA1685		Red Buddy Button		£75
AWA1686		Yellow Buddy Button		£75
AWA1687		Green Buddy Button		£75
AWA1690		Blue Buddy Button		£75
EGG SWITCH				
AWA2254		Egg Switch, black		£75

PIKO BUTTONS			
AWA1711		Red Piko Button	£135
AWA1712		Yellow Piko Button	£135
AWA1713		Green Piko Button	£135
AWA1714		Orange Piko Button	£135
AWA1715		Blue Piko Button	£135
AWA1716		Black Piko Button	£135
TOGGLES / BUTTONS			
AWZ2261		Stereo Toggle - direct control for Tilt only	£230
AWZ2262		Stereo Button - direct control for Tilt only	£219
AWZ2263		Four-Way Toggle - direct control for Tilt & Recline	£286
AWZ2264		Four-Way Button - direct control for Tilt & Recline	£306
AWZ2207		Toggle / Button RH Mount	○
AWZ2208		Toggle / Button LH Mount	○
AWZ2225		Without Mounting	○
BODYPOINT JOYSTICK KNOBS			
AWA2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A	£62
AWA2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A	£62
AWA2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£62
AWA2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£62
SHROUD COLOUR*			
AWA0202		Blue	○
AWA0228		Black	○
AWA0204		Transparent Cherry Red	○
AWZ0299		Special Shroud Colour (RAL Code must be supplied)	POA
CASTORS*			
AWA0840		10" Grey Pneumatic	○
AWA0858		10" Grey Puncture Proof Sealed	○
AWA0842		10" Grey Puncture Proof Solid Inlay	◆
DRIVE WHEELS*			
AWA0771		14" Grey Pneumatic	○
AWA0773		14" Grey Puncture Proof Sealed	○
AWA0750		14" Grey Puncture Proof Solid Inlay	◆
MOTOR TYPE*			
AWA1277		7 km/h	◆
BATTERY*			
AWA1364		73.5 AH Gel	◆
SUSPENSION*			
AWA1203		Heavy Duty Suspension	◆
G TRAC			
AWA1215		G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	£703

CHARGER*			
AWA1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	◆
AWA1313	<input type="checkbox"/>	OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	£130
LIGHTS*			
AWA1254	<input type="checkbox"/>	No Lights	◆
AWA1437	<input type="checkbox"/>	LED Light system with indicators (30 LUX)	£410
TRANSPORTATION KIT*			
AWA1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory.	◆
FIXATION BELT*			
AWA1455	<input type="checkbox"/>	Fixation Belt Extra long with Auto Lock	◆
REAR VIEW MIRROR			
AWA1415	<input type="checkbox"/>	Rear view Mirror RH - Must be selected with 10 km/h Motor (NC 2217)	£36
AWA1339	<input type="checkbox"/>	Rear view Mirror LH - Must be selected with 10 km/h Motor (NC 2218)	£36
CANE HOLDER			
AWA1312	<input type="checkbox"/>	Cane holder RH	£36
AWA1363	<input type="checkbox"/>	Cane holder LH	£36
WHEEL FENDERS			
AWA1417	<input type="checkbox"/>	Fender for Drive Wheels	£104

			TOTAL PRICE	
--	--	--	--------------------	--

Non Conforming Options:

AWA0564	Side Mounted: Reinforced hanger fixed angle (not detachable without tools)
---------	--

AWA0503	Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools
---------	--

AWA0599	Angle adjustable flip back heavy duty footboard
---------	---

AWA0522	Angle adjustable footboard (for max user weight 200kg)
---------	--

AWA2201	DLX-REM110 (Drive Only)
---------	-------------------------

AWA1207	Standard Bracket RH
---------	---------------------

AWA1208	Standard Bracket LH
---------	---------------------

AWA2227	Attendant Mounted RH
---------	----------------------

AWA2228	Attendant Mounted LH
---------	----------------------

AWA2220	LiNX DLX-REM050 (Dual Control used as secondary remote)
---------	---

AWA2219	LiNX DLX-ACU200 (Compact Dual Control)
---------	--

AWA2245	LiNX DLX-IDC (Intuitive Dual Control)
---------	---------------------------------------

AWA1207	Dual Control Standard Bracket RH
---------	----------------------------------

AWA1208	Dual Control Standard Bracket LH
---------	----------------------------------

AWA1221	IDC Dual Control Detachable Bracket RH
---------	--

AWA1222	IDC Dual Control Detachable Bracket LH
---------	--

AWA2217	Bracket RH, mounted for user
---------	------------------------------

AWA2218	Bracket LH, mounted for user
---------	------------------------------

AWA2227	Bracket RH, mounted for attendant
---------	-----------------------------------

AWA2228	Bracket LH, mounted for attendant
---------	-----------------------------------

AWA1686	Yellow Buddy Button
---------	---------------------

AWA1687	Green Buddy Button
---------	--------------------

AWA1690	Blue Buddy Button
---------	-------------------

AWA2254	Egg Switch, black
---------	-------------------

AWA1711	Red Piko Button
---------	-----------------

AWA1712	Yellow Piko Button
---------	--------------------

AWA1713	Green Piko Button
---------	-------------------

AWA1714	Orange Piko Button
---------	--------------------

AWA1715	Blue Piko Button
---------	------------------

AWA1716	Black Piko Button
---------	-------------------

AWA1437	LED Light system with indicators (30 LUX)
---------	---



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):