



**Yes, you can:**

Tel: 01656 776222  
Email: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
Online Spares available at: [www.invacare.co.uk](http://www.invacare.co.uk)

# Power Retail Prescription Form

01 DECEMBER 2025  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## TDX SP2 MAX

### Standard Features

- Max User Weight 180kg (28 stone)
- Crash tested ISO7176 - 19
- Seat widths from 490mm - 610mm
- Seat depths from 560mm - 630mm
- Tension Adjustable backrest
- Armrest with reinforced depth adjustable bracket with short armpads
- Manual settable Seat Tilt
- Manual settable Backrest Recline
- VARI F hangers with fixed footplates
- LiNX DLX-REM110 with standard bracket
- Selection of 5 shroud and rim colours
- Centre (Driving Wheels) 355mm x 76mm (14" x 3") Black Tyres
- Suspension and Sure Step Stability Technology
- 6 km/h (4mph) 4 Pole Motors
- 8A (200 - 250V) Charger
- 73.5 Ah Batteries
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: £5,382**



**Q'STRAIT**  
**QLK-150**  
Compatible

Please make your selection with a cross in the box    ☐ Free of charge    ☒ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE MODEL*</b>		
ARH0000	<input type="checkbox"/> TDX SP2 MAX	£5,382
ARZ0001	<input type="checkbox"/> Special TDX SP2 MAX (Extended Lead Time). Please provide details in comments section.	POA
<b>USER WEIGHT*</b>		
ARH1924	<input type="checkbox"/> User weight up to 180 kg	<input checked="" type="checkbox"/>
<b>SEATING TYPE*</b>		
ARH1145	<input type="checkbox"/> MAX Seating System	<input checked="" type="checkbox"/>

SEAT WIDTH*			
ARH0378	<input type="checkbox"/>	Seat Width of Seat Frame 550 mm (+100 mm adjustment in arms)	<input type="radio"/>
ARH0380	<input type="checkbox"/>	Seat Width of Seat Frame 610 mm (+100 mm adjustment in arms) (NC 0533)	<input type="radio"/>
SEAT DEPTH*			
ARH0396	<input type="checkbox"/>	560 mm	<input type="radio"/>
ARH0397	<input type="checkbox"/>	595 mm	<input type="radio"/>
ARH0398	<input type="checkbox"/>	630 mm	<input type="radio"/>
SEAT CUSHION*			
ARH1484	<input type="checkbox"/>	MAX Seat Cushion	<input type="radio"/>
ARH1114	<input type="checkbox"/>	Without Seat Cushion (10 mm seat board only)	<input type="radio"/>
ADAPTOR TO REDUCE THE SEAT DEPTH BY 130MM			
ARH0399	<input type="checkbox"/>	Adaptor to shorten the seat depth by 130 mm. (Required if seat width 490 mm is selected in combination with seat depth 595 mm or 630 mm)	<input type="radio"/>
BACKREST UPHOLSTERY*			
ARH1106	<input type="checkbox"/>	Tension Adjustable	<input checked="" type="radio"/>
UPHOLSTERY COLOUR*			
ARH0278	<input type="checkbox"/>	Black Multi Stretch PU Coated Wipe Down Fabric TR26	<input checked="" type="radio"/>
ARMRESTS*			
ARH0635	<input type="checkbox"/>	Fixed Armrest, 250mm - 340mm tool adjustable in height, detachable. (Includes reinforced depth adjustable bracket).	<input checked="" type="radio"/>
ARMPADS*			
ARH0647	<input type="checkbox"/>	Short armpad (325mm)	<input checked="" type="radio"/>
ARH0649	<input type="checkbox"/>	Long armpad (375mm)	<input type="radio"/>
PUSH BAR*			
ARH1307	<input type="checkbox"/>	Push bar	<input checked="" type="radio"/>
HEADREST / NECKREST (Adequate headrest option must be selected)*			
ARH1432	<input type="checkbox"/>	Headrest sideways adjustable	<input checked="" type="radio"/>
ARH1427	<input type="checkbox"/>	Headrest with cheek support	£199
ARH1495	<input type="checkbox"/>	Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts	<input type="radio"/>
MANUAL / POWERED SEATING OPTIONS*			
ARH1091	<input type="checkbox"/>	Manually settable Seat Tilt (max user weight 180kg)	<input checked="" type="radio"/>
ARH1095	<input type="checkbox"/>	Powered Heavy Duty Basic Tilt (0° - 20°). Select either 470mm or 495mm Seat to Ground Height. (max user weight 180kg)	£384
SEAT TO GROUND HEIGHT*			
ARH1024	<input type="checkbox"/>	Fixed seat height 470 mm (18.5")	<input type="radio"/>
ARH1046	<input type="checkbox"/>	Fixed seat height 495 mm (19.5")	<input type="radio"/>
BACKREST RECLINER*			
ARH0958	<input type="checkbox"/>	Manually settable recliner by spindle (0° - 25°)	<input checked="" type="radio"/>
ARH0947	<input type="checkbox"/>	Electrical recliner (0° - 25°)	£444

LEG REST TYPE*				
ARH0556		Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm	(NC 0504, 0513, 0599)	£328
ARH0589		Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm	(NC 0504, 0513, 0599)	£1,442
ARH0564		Side Mounted: Reinforced hangers, fixed angle (not detachable without tools). Includes calf pads.		○
ARH0503		Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools		◆
ARH0533		Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad	(NC 0559)	£206
ARH0531		Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad	(NC 0513, 0599, 0559)	£588
FOOTPLATE TYPE*				
ARH0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")	(NC 0503, 0533, 0531)	○
ARH0512		Fixed Footplate pair	(NC 0533, 0556, 0589)	○
ARH0504		Angle Adjustable Footplate pair	(NC 0589)	£65
ARH0513		Angle and depth adjustable Footplate pair	(NC 0589)	£104
ARH0599		Angle adjustable flip-back, heavy-duty footboard (Recommended for SW 550mm)	(NC 0503, 0531, 0533, 0556, 0589)	£294
LiNX CONTROLS*				
ARH2201		DLX-REM110 (Drive Only)	(NC 1437, 1095, 0947, 0589, 0531, 2219, 2245)	◆
ARH2203		DLX-REM211 (Drive & Actuators)	(NC 1437)	○
ARH2204		DLX-REM216 (Drive, Actuators & Lights)		£134
CONTROL BRACKETRY*				
ARH1207		Standard Bracket RH		○
ARH1208		Standard Bracket LH		○
ARH1707		Swing Away RH Bracket		£116
ARH1708		Swing Away LH Bracket		£116
ARH1709		Matrx Resolve Swing in/out Bracket RH		£219
ARH1710		Matrx Resolve Swing in/out Bracket LH		£219
ARH2227		Bracket RH mounted for Attendant		£56
ARH2228		Bracket LH mounted for Attendant		£56
DUAL CONTROLS				
ARH2220		LiNX DLX-REM050 (Dual Control used as secondary remote)	(NC 1221, 1222)	£348
ARH2219		LiNX DLX-ACU200 (Compact Dual Control)	(NC 1222, 2201)	£348
ARH2245		LiNX DLX-IDC (Intuitive Dual Control)	(NC 2219, 2227, 2228, 1207, 1208)	£660
DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)				
ARH1207		Dual Control Standard Bracket RH		○
ARH1208		Dual Control Standard Bracket LH		○
ARH1221		IDC Dual Control Detachable Bracket RH	(NC 2219, 2227, 2228)	○
ARH1222		IDC Dual Control Detachable Bracket LH	(NC 2219, 2227, 2228)	○
LiNX 10-WAY SWITCH				
ARH2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access		£444

LiNX 10-WAY SWITCH BRACKETRY (Only applicable if 10-Way Switch is selected)			
ARH2217		Bracket RH, mounted for user	○
ARH2218		Bracket LH, mounted for user	○
ARH2227		Bracket RH, mounted for attendant	£56
ARH2228		Bracket LH, mounted for attendant	£56
LiNX ACCESSORY CABLE			
ARH2252		Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch)	£89
FURTHER LiNX COMPONENTS / ACCESSORIES			
ARH2302		Not Necessary	(NC 1686, 1687, 1690, 2254, 1711, 1712, 1713, 1714, 1715, 1716) ◆
ARH2300		To be enclosed only	○
ARZ2301		To be connected and programmed (Special)	POA
BUDDY BUTTONS			
ARH1685		Red Buddy Button	£77
ARH1686		Yellow Buddy Button	£77
ARH1687		Green Buddy Button	£77
ARH1690		Blue Buddy Button	£77
EGG SWITCH			
ARH2254		Egg Switch, black	£77
PIKO BUTTONS			
ARH1711		Red Piko Button	£139
ARH1712		Yellow Piko Button	£139
ARH1713		Green Piko Button	£139
ARH1714		Orange Piko Button	£139
ARH1715		Blue Piko Button	£139
ARH1716		Black Piko Button	£139
TOGGLES / BUTTONS			
ARZ2261		Stereo Toggle - direct control for Tilt only	£237
ARZ2262		Stereo Button - direct control for Tilt only	£226
ARZ2263		Four-Way Toggle - direct control for Tilt & Recline	£295
ARZ2264		Four-Way Button - direct control for Tilt & Recline	£315
ARZ2207		Toggle / Button RH Mount	○
ARZ2208		Toggle / Button LH Mount	○
ARZ2225		Without Mounting	○
BODYPOINT JOYSTICK KNOBS			
ARH2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A	£64
ARH2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A	£64
ARH2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£64
ARH2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£64
REMOTE STOP SWITCH			
ARZ2258		Remote Stop Switch (ASL504)	£924

SHROUD COLOUR*			
ARH2640		Obsidian Black Frame Shrouds	○
ARH2641		Ruby Red Frame Shrouds	○
ARH2642		Purple Passion Frame Shrouds	○
ARH2643		Peacock Blue Frame Shrouds	○
ARH2644		Arctic Blue Frame Shrouds	○
RIM COLOUR*			
ARH2650		Obsidian Black Rim Inserts	○
ARH2651		Ruby Red Rim Inserts	○
ARH2652		Purple Passion Rim Inserts	○
ARH2653		Peacock Blue Rim Inserts	○
ARH2654		Arctic Blue Rim Inserts	○
FORKS*			
ARH2631		Single sided fork (New) includes Black castor	£310
ARH2632		Double sided fork	○
CENTRE WHEELS (DRIVE WHEELS)*			
ARH0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	○
ARH0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam	◆
ARH0706		DW 14" (355 mm) x 3" (75 mm), Puncture protected, gel-sealed	○
MOTOR TYPE*			
ARH1294		3 km/h (2mph)	○
ARH1282		6 km/h (4mph) (4 pole)	◆
BATTERY*			
ARH1260		73.5 Ah	◆
G TRAC			
ARH1215		G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	£724
CHARGER*			
ARH1343		OptiCharge Charger 8A (200 - 250V)	◆
ARH1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	£134
LIGHTS*			
ARH1437		LED Light System with Indicators (30 LUX)	£422
ARH1254		Without light system	◆
TRANSPORTATION KIT			
ARH1450		TEMPORARILY UNAVAILABLE ON TDX MODELS UNTIL FURTHER NOTICE - Preparation for DAHL-Docking station. (NC 0001)	£514
ARH1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory.	◆
FIXATION BELT*			
ARH1329		MAX Lapbelt (adjustable both sides) fitted to base of backrest	◆
REAR VIEW MIRROR			
ARH1415		Rear-view Mirror RH	£37
ARH1339		Rear-view Mirror LH	£37

CANE HOLDER

ARH1312		Cane holder RH	£37
ARH1363		Cane holder LH	£37

TOTAL PRICE

**Non Conforming Options:**

ARZ0001	Special TDX SP2 MAX (Extended Lead Time). Please provide details in comments section.
ARH1095	Powered Heavy Duty Basic Tilt (0° - 20°). Select either 470mm or 495mm Seat to Ground Height. (max user weight 180kg)
ARH0947	Electrical recliner (0° - 25°)
ARH0556	Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm
ARH0589	Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm
ARH0503	Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools
ARH0533	Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad
ARH0531	Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad
ARH0559	Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")
ARH0504	Angle Adjustable Footplate pair
ARH0513	Angle and depth adjustable Footplate pair
ARH0599	Angle adjustable flip-back, heavy-duty footboard (Recommended for SW 550mm)
ARH2201	DLX-REM110 (Drive Only)
ARH1207	Standard Bracket RH
ARH1208	Standard Bracket LH
ARH2227	Bracket RH mounted for Attendant
ARH2228	Bracket LH mounted for Attendant
ARH2219	LiNX DLX-ACU200 (Compact Dual Control)
ARH2245	LiNX DLX-IDC (Intuitive Dual Control)
ARH1207	Dual Control Standard Bracket RH
ARH1208	Dual Control Standard Bracket LH
ARH1221	IDC Dual Control Detachable Bracket RH
ARH1222	IDC Dual Control Detachable Bracket LH
ARH2227	Bracket RH, mounted for attendant
ARH2228	Bracket LH, mounted for attendant
ARH1686	Yellow Buddy Button
ARH1687	Green Buddy Button
ARH1690	Blue Buddy Button
ARH2254	Egg Switch, black
ARH1711	Red Piko Button
ARH1712	Yellow Piko Button
ARH1713	Green Piko Button
ARH1714	Orange Piko Button
ARH1715	Blue Piko Button
ARH1716	Black Piko Button
ARH1437	LED Light System with Indicators (30 LUX)



For a more complex customer requirement,  
please email: [orders@invacare.com.au](mailto:orders@invacare.com.au)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):