



# Manual Retail Prescription Form

**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

MARCH 2021  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## XLT Max

### Standard Features

- MAX USER WEIGHT 180KG/28 STONE
- CRASH TESTED ISO7176-PART 19
- LIGHTWEIGHT TITANIUM RIGID FRAME WITH RIGID FOOTREST TUBE
- LOCKDOWN FOLDING BACKREST, ANGLE ADJUSTABLE
- TENSION ADJUSTABLE UPHOLSTERY
- CHOICE OF 3 SEAT WIDTHS 505MM - 605MM (20" -24")
- 3 ADJUSTABLE SEAT DEPTHS 460MM - 600MM (18" - 24")
- REMOVABLE CLOTHESGUARDS
- CHOICE OF 6 FRAME COLOURS
- 24" QUICK RELEASE SELF PROPEL REAR WHEELS



**Standard Base price: £2,437**

Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE PRODUCT*</b>		
BTA0010	INVACARE XLT	£2,437
BTA8000	SPECIAL REQUEST - PLEASE INCLUDE QUOTATION REFERENCE IN SPECIAL DESCRIPTION AT THE END OF THE FORM WHERE APPLICABLE	POA
<b>PROPULSION TYPE*</b>		
BTA0020	STANDARD WHEEL HUB	◆
BTA0030	ASSISTANT DRUM BRAKES (NC 0570, 0580, 0630, 0970, 1050)	£276
<b>FRAME*</b>		
BTA0040	ACTIVE FRAME, KNEE ANGLE 75°	◆

SEAT WIDTH*				
BTA0050		505 MM		○
BTA0060		555 MM		○
BTA0070		605 MM		○
SEAT DEPTH*				
BTA0080		500 MM ADJUSTABLE FROM SEAT DEPTH 460 TO 500 MM		◆
BTA0090		550 MM ADJUSTABLE FROM SEAT DEPTH 510 TO 550 MM		○
BTA0100		600 MM ADJUSTABLE FROM SEAT DEPTH 560 TO 600 MM	(NC 0650)	○
FRONT SEAT HEIGHT*				
BTA0110		440 MM	(NC 0200, 0260, 0270, 0410, 0420, 0430, 0440, 0895, 0900, 0920)	○
BTA0120		450 MM	(NC 0200, 0270, 0410, 0420, 0430, 0440, 0895, 0900, 0920)	○
BTA0130		460 MM	(NC 0200, 0270, 0410, 0420, 0430, 0440, 0880, 0920)	○
BTA0140		470 MM	(NC 0200, 0420, 0430, 0440, 0920)	○
BTA0150		480 MM	(NC 0430, 0440)	◆
BTA0160		490 MM	(NC 0430, 0440)	○
BTA0170		500 MM	(NC 0440)	○
BTA0180		510 MM	(NC 0440)	○
BTA0190		520 MM	(NC 0200, 0210, 0880)	○
REAR SEAT HEIGHT*				
BTA0200		390 MM		○
BTA0210		400 MM		○
BTA0220		410 MM		○
BTA0230		420 MM		○
BTA0240		430 MM		○
BTA0250		440 MM		◆
BTA0260		450 MM		○
BTA0270		460 MM		○
LOWER LEG LENGTH*				
BTA0280		270 MM ONLY WITH FOOTREST HIGH MOUNTED	(NC 0670, 0680)	○
BTA0290		280 MM ONLY WITH FOOTREST HIGH MOUNTED	(NC 0670, 0680)	○
BTA0300		290 MM ONLY WITH FOOTREST HIGH MOUNTED	(NC 0670, 0680)	○
BTA0310		300 MM ONLY WITH FOOTREST HIGH MOUNTED	(NC 0670, 0680)	○
BTA0320		310 MM ONLY WITH FOOTREST HIGH MOUNTED	(NC 0670, 0680)	○
BTA0330		320 MM ONLY WITH FOOTREST HIGH MOUNTED	(NC 0670, 0680)	○
BTA0340		330 MM ONLY WITH FOOTREST HIGH MOUNTED	(NC 0670, 0680)	○
BTA0350		340 MM ONLY WITH FOOTREST HIGH MOUNTED	(NC 0670, 0680)	○
BTA0360		350 MM ONLY WITH FOOTREST HIGH MOUNTED	(NC 0670, 0680)	○
BTA0370		360 MM ONLY WITH FOOTREST HIGH MOUNTED	(NC 0670, 0680)	○
BTA0380		370 MM ONLY WITH FOOTREST HIGH MOUNTED	(NC 0670, 0680)	○
BTA0390		380 MM	(NC 0700)	○
BTA0400		400 MM	(NC 0700)	○
BTA0410		420 MM	(NC 0700)	○
BTA0420		440 MM	(NC 0700)	◆
BTA0430		460 MM	(NC 0700)	○
BTA0440		480 MM	(NC 0700)	○

<b>BACKREST TYPE*</b>			
BTA0450		VELCRO TENSION ADJUSTABLE, HEIGHT ADJUSTABLE, INFINITELY VARIABLE ANGLE ADJUSTABLE AND FOLDABLE BACKREST	◆
<b>BACKREST HEIGHT*</b>			
BTA0510		400 MM HEIGHT ADJUSTABLE FROM 400 TO 440 MM BY 20 MM STEPS	○
BTA0520		420 MM HEIGHT ADJUSTABLE FROM 400 TO 440 MM BY 20 MM STEPS	○
BTA0530		440 MM HEIGHT ADJUSTABLE FROM 400 TO 440 MM BY 20 MM STEPS	○
BTA0540		450 MM HEIGHT ADJUSTABLE FROM 450 TO 490 MM BY 20 MM STEPS	○
BTA0550		470 MM HEIGHT ADJUSTABLE FROM 450 TO 490 MM BY 20 MM STEPS	○
BTA0560		490 MM HEIGHT ADJUSTABLE FROM 450 TO 490 MM BY 20 MM STEPS	○
<b>PUSH HANDLES*</b>			
BTA0570		FIXED PUSH HANDLES	◆
BTA0580		HEIGHT ADJUSTABLE PUSH HANDLES, EXTERNAL, SHORT GRIP	£53
BTA0590		HEIGHT ADJUSTABLE PUSH HANDLES, EXTERNAL, LONG GRIP FOR DRUM BRAKES	£53
BTA0600		HEIGHT ADJUSTABLE PUSH HANDLES, EXTERNAL, INCLUDING NECKREST, LONG GRIP FOR DRUM BRAKES	£88
BTA0610		HEIGHT ADJUSTABLE PUSH HANDLES, EXTERNAL, INCLUDING HEADREST, LONG GRIP FOR DRUM BRAKES	£111
BTA0630		WITHOUT PUSH HANDLES	○
<b>SEAT TYPE*</b>			
BTA0640		VELCRO TENSION- AND DEPTH ADJUSTABLE	◆
<b>SEAT CUSHION*</b>			
BTA0650		FLO-TECH CONTOUR	£140
BTA0660		WITHOUT SEAT CUSHION	○
<b>FOOTREST*</b>			
BTA0670		FOOT BOW RIGID	◆
BTA0680		FOOT BOW RIGID WITH COVER	£13
BTA0700		FOOTREST, HIGH MOUNTED, FOR KNEE TO HEEL LENGTH 270 TO 370 MM	○
<b>CLOTHESGUARDS/ARMRESTS LH*</b>			
BTA0710		CLOTHESGUARD REMOVABLE, LH	◆
BTA0720		MUDGUARD, REMOVABLE, LH	£45
BTA0730		ARMREST HEIGHT ADJUSTABLE, SHORT PAD (280 MM), LH	£38
<b>CLOTHESGUARDS/ARMRESTS RH*</b>			
BTA0760		CLOTHESGUARD REMOVABLE, RH	◆
BTA0770		MUDGUARD, REMOVABLE, RH	£45
BTA0780		ARMREST HEIGHT ADJUSTABLE, SHORT PAD (280 MM), RH	£38

FRAME COLOUR*			
BTA0810	<input type="checkbox"/>	ANTHRACITE BLACK	<input type="radio"/>
BTA0820	<input type="checkbox"/>	GREY	<input type="radio"/>
BTA0830	<input type="checkbox"/>	BLUE	<input type="radio"/>
BTA0840	<input type="checkbox"/>	RED	<input type="radio"/>
BTA0850	<input type="checkbox"/>	ORANGE	<input type="radio"/>
BTA0860	<input type="checkbox"/>	BLASTED TITANIUM	<input type="radio"/>
CASTOR ATTACHMENT*			
BTA0870	<input type="checkbox"/>	HOUSING WITH ANGLE ADJUSTMENT	<input checked="" type="checkbox"/>
FRONT CASTORS*			
BTA0880	<input type="checkbox"/>	140X40 SOFT (5,5")	<input type="radio"/>
BTA0895	<input type="checkbox"/>	150X30 MEDIUM HARD (6")	<input type="radio"/>
BTA0900	<input type="checkbox"/>	150X30 SOFT (6")	<input type="radio"/>
BTA0920	<input type="checkbox"/>	180X40 SOFT (7")	<input type="radio"/>
CASTOR FORK*			
BTA0940	<input type="checkbox"/>	FIXED STEPS	<input checked="" type="checkbox"/>
REAR WHEEL*			
BTA0950	<input type="checkbox"/>	24" STANDARD	<input checked="" type="checkbox"/>
TYRE TYPE*			
BTA0960	<input type="checkbox"/>	TYRES SOLID	<input checked="" type="checkbox"/>
BTA0970	<input type="checkbox"/>	TYRES BLACK MARATHON EVOLUTION	£75
HANDRIM TYPE*			
BTA0980	<input type="checkbox"/>	HANDRIM ALUMINIUM	<input checked="" type="checkbox"/>
BTA0990	<input type="checkbox"/>	HANDRIM STAINLESS	<input type="radio"/>
BTA1000	<input type="checkbox"/>	HANDRIM PLASTIC COATED	£26
BTA1010	<input type="checkbox"/>	HANDRIM FOAM COATED	£26
BTA1020	<input type="checkbox"/>	HANDRIM BLACK GRIP	£26
REAR WHEEL FIXATION*			
BTA1030	<input type="checkbox"/>	STANDARD REAR WHEEL FIXATION	<input checked="" type="checkbox"/>
CAMBER ANGLE*			
BTA1040	<input type="checkbox"/>	0° CAMBER	<input checked="" type="checkbox"/>
BTA1050	<input type="checkbox"/>	2° CAMBER	<input type="radio"/>
SPOKEGUARD			
BTA1060	<input type="checkbox"/>	SPOKEGUARD FROST	£45
BRAKE*			
BTA1070	<input type="checkbox"/>	STANDARD BRAKE	<input checked="" type="checkbox"/>
SAFTEY OPTIONS*			
BTA1080	<input type="checkbox"/>	PELVIC BELT	<input type="radio"/>
BTA1090	<input type="checkbox"/>	BRAKE LEVER EXTENSION, FOLD AWAY, LH	£13
BTA1100	<input type="checkbox"/>	BRAKE LEVER EXTENSION, FOLD AWAY, RH	£13
BTA1110	<input type="checkbox"/>	AUTO LOCK FOR ARMREST	£9
BTA1120	<input type="checkbox"/>	ANTI-TIP DEVICES, LH FOR SAFETY 2 ANTI-TIPS ARE RECOMMENDED	£34
BTA1130	<input type="checkbox"/>	ANTI-TIP DEVICES, RH FOR SAFETY 2 ANTI-TIPS ARE RECOMMENDED	£34

**OTHER OPTIONS**

BTA1150		PUMP	£13
BTA1160		CANE HOLDER, LH	£18
BTA1170		CANE HOLDER, RH	£18
BTA1200		TOOL KIT	◆
BTA1210		FRAME PROTECTOR	£26

**TOTAL PRICE**

**Non Conforming Options:**

BTA0200	390 MM
BTA0210	400 MM
BTA0260	450 MM
BTA0270	460 MM
BTA0410	420 MM
BTA0420	440 MM
BTA0430	460 MM
BTA0440	480 MM
BTA0570	FIXED PUSH HANDLES
BTA0580	HEIGHT ADJUSTABLE PUSH HANDLES, EXTERNAL, SHORT GRIP
BTA0630	WITHOUT PUSH HANDLES
BTA0650	FLO-TECH CONTOUR
BTA0670	FOOT BOW RIGID
BTA0680	FOOT BOW RIGID WITH COVER
BTA0700	FOOTREST, HIGH MOUNTED, FOR KNEE TO HEEL LENGTH 270 TO 370 MM
BTA0880	140X40 SOFT (5,5")
BTA0895	150X30 MEDIUM HARD (6")
BTA0900	150X30 SOFT (6")
BTA0920	180X40 SOFT (7")
BTA0970	TYRES BLACK MARATHON EVOLUTION
BTA1050	2° CAMBER



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):