

Contact: \_\_\_\_\_  
 Service/PCT: \_\_\_\_\_  
 City: \_\_\_\_\_  
 Phone Number: \_\_\_\_\_  
 Marked For: \_\_\_\_\_  
 Purchase Order # \_\_\_\_\_

Client reference \_\_\_\_\_

Assessment Date \_\_\_\_\_



- Special request required
- Quote
- Order
- E-mail Address \_\_\_\_\_
- Fax Number \_\_\_\_\_

**Wheelchair Specifications:**

Model: \_\_\_\_\_  
 Width: \_\_\_\_\_  
 Depth: \_\_\_\_\_

**Seating Specifications:**


A - Back Height \_\_\_\_\_  
 B - Lateral Height \_\_\_\_\_  
 C - Lateral Width \_\_\_\_\_  
 D - Hip Width \_\_\_\_\_  
 E - Set Up Seat Depth \_\_\_\_\_  
 G - "E" + 3" for Growth \_\_\_\_\_  
 H - Seat Width \_\_\_\_\_  
 I - Back Width \_\_\_\_\_



- 1 \_\_\_\_\_
- 2 \_\_\_\_\_
- 3 \_\_\_\_\_
- 4 \_\_\_\_\_
- 5 \_\_\_\_\_

**Supine**

- 6 \_\_\_\_\_
- 7 \_\_\_\_\_
- HL \_\_\_\_\_
- HR \_\_\_\_\_





Freedom Seating  
March 2021

**Upholstery and Monograms:**

Client \_\_\_\_\_

Total system covered in:  Black Wipedown  Black Smooth

Contact Surface:  Black Wipedown  Black Smooth Non-contact Surface:  Naugahyde  Lycra

Naugahyde Color Choices:

- |  |  |  |   |
|--|--|--|---|
| <input type="checkbox"/> K-16 Nutmeg (Red Sparkle) | <input type="checkbox"/> K-09 Bright Yellow  | <input type="checkbox"/> F-168 Deep Blue | <input type="checkbox"/> K-10 Bright Green    |
| <input type="checkbox"/> F-117 Orange Sizzle       | <input type="checkbox"/> K-17 Sparkle Purple | <input type="checkbox"/> F-28 Bright Red | <input type="checkbox"/> F-112 Purple Passion |
| <input type="checkbox"/> F-118 Lime Punch          | <input type="checkbox"/> F-119 Rosalee       | <input type="checkbox"/> F-75 Space Blue | <input type="checkbox"/> F-124 Shimmer        |

Monogram Name: \_\_\_\_\_

Font:  Script OR  **BLOCK**


Name Colour \_\_\_\_\_

**Step 1. Choose Hardware for Planar / Flat Seat**  
Note: Price includes 1/2" wood & 1-1/2" 4# foam.

- |                                     |   |         |
|-------------------------------------|---|---------|
| <input type="checkbox"/> EE1        | Seat - Econo-eze 2 point hardware (Add Back Hardware)                 | £364.70 |
| <input type="checkbox"/> Fastclamp1 | Seat - Fast Clamp 2 point hardware (Add back hardware)                | £307.20 |
|                                     | <input type="checkbox"/> 1" posts <input type="checkbox"/> 7/8" posts |         |
| <input type="checkbox"/> 101-7500   | Seat -Econo-eze 4 point hardware                                      | £398.10 |

**Step 2. Choose Modifications to Planar / Flat Seat**

- |   |             |   |             |
|---|-------------|---|-------------|
| <input type="checkbox"/> Mod-1K Anti-Thrust Seat Foam | Add £125.40 | <input type="checkbox"/> Mod-1Y Contour Abd/Anti Thrust/Add | Add £214.20 |
|---|-------------|---|-------------|



(see page 2)

- |   |             |   |             |
|---|-------------|---|-------------|
| <input type="checkbox"/> Mod-1V Contour Abd/Anti Thrust     | Add £167.20 | <input type="checkbox"/> Mod-1W Contour Add/Anti Thrust | Add £167.20 |
| <input type="checkbox"/> Mod-1U Removable Seat Cover        | Add £61.70  |   |             |
| <input type="checkbox"/> Mod - 1I Incontinent Liner—Plastic | £18.80      |   |             |

**2 Point Pelvic Positioning Belt**

Add £35.50

- SB-3331-S Centre pull/pls Push Button - Small  
 SB-3331-M Centre pull/pls Push Button - Medium  
 SB-3331-L Centre pull/pls Push Button - Large

**4 Point Pelvic Positioning Belt**

Add £106.60

- IVC-NXT-SFDI 3AP4-S  
 IVC-NXT-SFDI 3AP4-M  
 IVC-NXT-SFDI 3AP4-L


**Step 1. Choose Hardware for Planar / Flat Back**  
**Note: Price includes 1/2" wood & 1" 4# foam.**

Client \_\_\_\_\_


- |                                     |   |         |
|-------------------------------------|---|---------|
| <input type="checkbox"/> EE2        | Back - Econo-eze 2 point hardware (Add Seat Hardware)                 | £364.70 |
| <input type="checkbox"/> Fastclamp2 | Back - Fast Clamp 2 point hardware (Add seat hardware)                | £307.20 |
|                                     | <input type="checkbox"/> 1" posts <input type="checkbox"/> 7/8" posts |         |
| <input type="checkbox"/> 102-7500   | Back -Econo-eze 4 point hardware                                      | £398.10 |

**Step 2. Choose Modifications to Planar / Flat Back**


- |                                  |                  |            |
|----------------------------------|------------------|------------|
| <input type="checkbox"/> Mod-2I  | Curved Back-Wood | Add £78.40 |
| <input type="checkbox"/> Mod-2IA | Curved Back-Foam | Add £70    |




- |  |            |
|--|------------|
| <input type="checkbox"/> Mod-2H I-Modification | Add £75.20 |
|--|------------|



- |  |             |
|--|-------------|
| <input type="checkbox"/> Mod-2DD Built-In Adj.Lumbar Support | Add £177.70 |
|--|-------------|



- |  |            |
|--|------------|
| <input type="checkbox"/> Mod-2J T-Modification | Add £75.20 |
|--|------------|



- |  |            |
|--|------------|
| <input type="checkbox"/> Mod-1U Removable Seat Cover | Add £61.70 |
|--|------------|

- |                                     |                   |            |
|-------------------------------------|-------------------|------------|
| <input type="checkbox"/> IVC-NXT-6B | Strap risers pair | Add £25.10 |
|-------------------------------------|-------------------|------------|

**Medial Thigh - Abductor/Pommel**


Client \_\_\_\_\_

- |  |         |
|--|---------|
| <input type="checkbox"/> EE3A Adjustable/Removable Pelvic Stabilizer | £173.50 |
| <input type="checkbox"/> EE3S Push Button Flip Down Abductor         | £289.50 |

Standard Sizes	
<input type="checkbox"/> Small	3L x 3H x 3-2 Taper
<input type="checkbox"/> Medium	4L x 3H x 3-2 Taper
<input type="checkbox"/> Large	5L x 4H x 4-2 Taper
<input type="checkbox"/> X-Large	6L x 5H x 5-3 Taper

**Hip Guides**

- |                                    |                                   |            |
|------------------------------------|-----------------------------------|------------|
| <input type="checkbox"/> EE4A-450A | Laterals-L Bracket-Std 4#         | £164.10/pr |
| <input type="checkbox"/> EE4B-442  | Laterals-Offset- 1" O/S4#         | £164.10/pr |
| <input type="checkbox"/> EE4B-443  | Laterals-Offset- 2" O/S4#         | £164.10/pr |
| <input type="checkbox"/> EE4B-442B | Laterals-Offset- 1" REV O/S4#     | £164.10/pr |
| <input type="checkbox"/> EE4D-STD  | Hanger Mounted Adductor-Std 4"    | £437.90/pr |
| <input type="checkbox"/> EE4D-78   | Hanger Mounted Adductor—7/8"      | £437.90/pr |
| <input type="checkbox"/> EE4E-450A | Laterals Fix/Arm Hgt. Adj.-Std 4# | £232/pr    |
| <input type="checkbox"/> EE4E-442  | Laterals Fix/Arm Hgt. Adj.-1" O/S | £232/pr    |
| <input type="checkbox"/> EE4E-443  | Laterals Fix/Arm Hgt. Adj.-2" O/S | £232/pr    |
| <input type="checkbox"/> Mod-4A    | Custom Contour                    | £61.70/pr  |



A = 4" high or \_\_\_\_\_"  
 B = Set up seat depth or \_\_\_\_\_"

**Trunk Supports**

- |                                    |   |            |  |  |
|------------------------------------|---|------------|--|--|
| <input type="checkbox"/> EE5A      | Fixed Thoracic Support - STD 4#         | £164.10/pr |  |  |
| <input type="checkbox"/> EE5B-442  | Thoracic Support-1"O/S-Std 4#           | £164.10/pr |  |  |
| <input type="checkbox"/> EE5B-443  | Thoracic Support-2" O/S Std 4#          | £164.10/pr |  |  |
| <input type="checkbox"/> EE5B-442B | Thoracic Support-1" Reverse O/S         | £164.10/pr |  |  |
| <input type="checkbox"/> EE5I      | Tracks-Thoracic Support (Chose size)    | £35.50/pr  |  |  |
| <input type="checkbox"/> EE5J-450C | Fixed Thoracic-Summer Winter-Std 4#     | £275.90/pr |  |  |
| <input type="checkbox"/> EE5J-442C | Fixed Thoracic-Summer Winter-1"O/S 4#   | £275.90/pr |  |  |
| <input type="checkbox"/> EE5J-443C | Fixed Thoracic-Summer Winter-2"O/S 4#   | £275.90/pr |  |  |
| <input type="checkbox"/> EE5P-510A | S/A Lever Release Thoracic-Std-4#       | £362.60/pr |  |  |
| <input type="checkbox"/> EE5P-510B | S/A Lever Release Thoracic- 1"O/S 4#    | £362.60/pr |  |  |
| <input type="checkbox"/> EE5P-510C | S/A Lever Release Thoracic- 2"O/S 4#    | £362.60/pr |  |  |
| <input type="checkbox"/> EE5Q-497A | S/A Trunks lever w/Summer Winter        | £468.20/pr |  |  |
| <input type="checkbox"/> EE5Q-497C | S/A Trunks lever w/Summer Winter 1" O/S | £468.20/pr |  |  |
| <input type="checkbox"/> EE5Q-497D | S/A Trunk levers w/Summer Winter 2" O/S | £468.20/pr |  |  |
| <input type="checkbox"/> EE5W      | S/A Thoracic Support push button        | £332.30/pr |  |  |
| <input type="checkbox"/> Mod-5F    | Protective Boots                        | £36.60/pr  |  |  |
- 
- |                                |                                 |                                |                                     |                                    |
|--------------------------------|---------------------------------|--------------------------------|-------------------------------------|------------------------------------|
|                                |                                 |                                | <b>Pads</b>                         | <b>Tracks</b>                      |
| <input type="checkbox"/> Small | <input type="checkbox"/> Medium | <input type="checkbox"/> Large | <input type="checkbox"/> 3" Single  | <input type="checkbox"/> 5" Single |
|                                |                                 |                                | <input type="checkbox"/> 7" Single  | <input type="checkbox"/> 9" Single |
|                                |                                 |                                | <input type="checkbox"/> 11" Single | <input type="checkbox"/> 9" Double |
|                                |                                 |                                | <input type="checkbox"/> 11" Double |                                    |
|                                |                                 |                                | <b>Trunk Supports</b>               |                                    |
|                                |                                 |                                | Flat <input type="checkbox"/>       |                                    |
|                                |                                 |                                | Curved <input type="checkbox"/>     |                                    |

**Headrest Hardware**

- |                                      |   |         |
|--------------------------------------|---|---------|
| <input type="checkbox"/> 7700        | Standard headrest                       | £103.50 |
| <input type="checkbox"/> 7708        | Rock N Roll Headrest                    | £117.70 |
| <input type="checkbox"/> 7709        | Rock N Roll Flip Back Headrest          | £219.50 |
| <input type="checkbox"/> 7733        | Adjustable, Universal Headrest Hardware | £43.90  |
| <input type="checkbox"/> 77Platemate | Platemate                               | £27.20  |

**Headrest Pads**



- |  |   |  |  |   |         |
|--|---|--|--|---|---------|
| <input type="checkbox"/> 7750 Curved pad- 5#             | <input type="checkbox"/> Sm. 11W* x 3.75H | <input type="checkbox"/> Med. 11W* x 3.75H | <input type="checkbox"/> Lg. 12W* x 4H     | <input type="checkbox"/> Xlg.12W*x4H    | £113.90 |
| <input type="checkbox"/> 7753 Tripad Lock 3-Piece Pad-5# | <input type="checkbox"/> Sm. 4W x 4H      | <input type="checkbox"/> Med. 5W x 4H      | <input type="checkbox"/> Lg. 5W x 5H       |   | £314.50 |
| <input type="checkbox"/> 7756 Soft Healdrest Pad         | <input type="checkbox"/> Sm. 10W x 4H     | <input type="checkbox"/> Med. 12W x 5H     | <input type="checkbox"/> Lg. 13.5W x 4.5 H |   | £103.50 |
| <input type="checkbox"/> 7760 Slimline Headrest Pad      | <input type="checkbox"/> Sm. 8W x 4.5H    | <input type="checkbox"/> Med. 10W x 4.5H   | <input type="checkbox"/> Lg. 12W x 4.5H    | <input type="checkbox"/> Xlg. 14W x4.5H | £103.50 |
| <input type="checkbox"/> 7757 Occipital low profile      | <input type="checkbox"/> Sm. 8W x 4.5H    | <input type="checkbox"/> Med. 10W x 4.5H   | <input type="checkbox"/> Lg. 12W x 4.5H    | <input type="checkbox"/> Xlg. 14W x4.5H | £103.5  |

\* Outside Width

Client \_\_\_\_\_

**Foot Boxes**

- |  |         |  |
|--|---------|--|
| <input type="checkbox"/> 08-10000 Leg Cradles      |         | <input type="checkbox"/> 08-11000 Foot Box                 |
| <input type="checkbox"/> Footplate/Footboard Mount | £376.20 | <input type="checkbox"/> Footplate/Footboard Mount £320.80 |
| <input type="checkbox"/> Hanger Mount              |         | <input type="checkbox"/> Hanger Mount                      |



**Shoe Holders**

- SFDI-9BB Select Line High Tops      £89.90



- Small 4.5" x 7.5"  
 Medium 5" x 9"  
 Large 5.5" x 10"

**Trays**

What is the width of the frame \_\_\_\_\_ ?

- |                                   |   |         |
|-----------------------------------|---|---------|
| <input type="checkbox"/> 08-12000 | EZ-Lite 1/4" thick Lexan Tray           | £122.30 |
| <input type="checkbox"/> 08-13000 | EZ-Lite 1/4" thick Lexan Tray with Ribs | £131.70 |

Tray Sizes	
Size	Description
<input type="checkbox"/> X- Small	10" - 11" wide chair, Tray dimensions—21"wide x 18" deep, body cut out— 10" w x 7" deep
<input type="checkbox"/> Small	12" - 13" wide chair, Tray dimensions—21"wide x 18" deep, body cut out— 12" w x 7" deep
<input type="checkbox"/> Medium	14" - 15" wide chair, Tray dimensions—24"wide x 22" deep, body cut out— 14" w x 9" deep
<input type="checkbox"/> Large	16" - 17" wide chair, Tray dimensions—26"wide x 24" deep, body cut out— 16" w x10" deep
<input type="checkbox"/> X- Large	18" - 19" wide chair, Tray dimensions—28"wide x 24" deep, body cut out— 18" w x10" deep