

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
www.invacare.co.uk

Matrx Flo-tech back support
April 2020

Account Name :	Account No :
Address :	Order No :
Tel No:	Order Qty :
Delivery address :	Order Date :

Matrx Flotech Back Support

Retail order form
Maximum user weight 26 stone (165kg)


The Matrx *Flo-tech Back Support* provides a high degree of stability and provides improved posture resulting in enhanced protection from pressure ulcer development.



- Moisture resistant and vapour permeable cover
- Easily fixed to armchair or wheelchair with extendable straps
- Built in lumbar option available for increased protection
- One common size available

Back Support Dimensions:

- (A) Outside Width 43cm
- (B) Inside Width 36cm
- (C) Height 46cm
- (D) Depth 10cm



Code	X	Description	Price
Back supports			



Back support



Back support with lumbar



Back support with lumbar and handle



Lumbar pad

FBB1		Back support	Black Smooth (ingress resistant) cover	£155
FBB1/LUM		Back support with lumbar	Black Smooth (ingress resistant) cover	£170
FBB1/LUM/Handle		Back support with lumbar and handle	Black Smooth (ingress resistant) cover	£171
1538975		Back support	Spacer fabric cover	£158
LUMPAD/VIS		Lumbar pad		£43
FBB1B		Spare smooth cover		£50

Special Instructions - Please supply the following information:



Name:

Email:

Phone:

Description of special request (please include relevant option code):