



RETAIL PRESCRIPTION FORM

Yes, you can:

Tel: 01656 776222


Email: ordersuk@invacare.comOnline spares available at www.invacare.co.uk

Flo-tech Contour Lo-Back cushion

June 2025

Customer name:	Account No :
Contact no:	Order No :
Delivery address:	Order Qty :
	Order Date :

Matrx Flo-tech Contour Lo-Back
Retail order form
Maximum User Weight 152kg (24 stone)



Contoured into sections offering excellent stability and positioning

Lower rear edge which will allow the person to sit further back in the seat

4" height

Contoured surface optimizes the support surface contact area for pressure distribution

**Features**

At Risk/Increased Vulnerability

Tapered front to back 4" to 3"

Two way stretch vapour permeable cover

Accessories available

Anti Slip Base with handle

3 Year Warranty

**matrx**
SEATING SERIES**Flo-tech**

Code	X	Description	Price
Flo-tech Contour Lo Back			
HCL0002		Complete Cushion	£165.00
HCL0003		Cover Only	£62.00
HCL0004		Foam Only	£110.00
HCL0005		Pack of 3 Replacement Platilon Liners	£25.00
Foam Type			
HCL0006		Standard Foam (nc 3, 5)	◆
Base Type			
HCL0007		Flat Base (nc 5)	◆
Seat Width			
HCL0008		36cm (14 Inches) (nc 23, 24, 25, 26, 27, 28)	○
HCL0009		38cm (15 Inches) (nc 28)	○
HCL0010		41cm (16 Inches) (nc 28)	○
HCL0011		43cm (17 Inches) (nc 28)	○
HCL0012		46cm (18 Inches) (nc 28)	○
HCL0013		48cm (19 Inches) (nc 19, 20, 28)	○
HCL0014		51cm (20 Inches) (nc 19, 20)	○

HCL0015	53cm (21 Inches) (nc 19, 20, 21, 22)	○
HCL0016	55cm (22 Inches) (nc 19, 20, 21, 22)	○
HCL0017	58cm (23 Inches) (nc 19, 20, 21, 22)	○
HCL0018	60cm (24 Inches) (nc 19, 20, 21, 22)	○
Seat Depth		
HCL0019	33cm (13 Inches)	○
HCL0020	36cm (14 Inches)	○
HCL0021	38cm (15 Inches)	○
HCL0022	41cm (16 Inches)	○
HCL0023	43cm (17 Inches)	○
HCL0024	46cm (18 Inches)	○
HCL0025	48cm (19 Inches)	○
HCL0026	51cm (20 Inches)	○
HCL0027	53cm (21 Inches)	○
HCL0028	55cm (22 Inches)	○
Fabric Type		
HCL0029	Black Smooth (nc 4, 5)	◆
HCL0038	Black Spacer (nc 4, 5)	○
HCL0039	Black Wipedown (nc 4, 5)	○
Cover Branding		
HCL0030	Standard Invacare Branded Cover (nc 4, 5)	◆
Liner		
HCL0034	Platilon Liner Covering Moulded Foam (nc 3, 4, 5)	○
HCL0035	No Platilon Liner Required	○
Loop Band		
HCL0036	Loop Band Pieces Sewn On Cover Base (nc 4, 5)	£25.00
HCL0037	No Loop Band Required	○
		Total Price
		Quantity
Accessories		
1494494	Rigidiser W14" x D16"	£38.00
1487165	Rigidiser W16" x D16"	£38.00
1494495	Rigidiser W16" x D18"	£38.00
1494496	Rigidiser W16" x D20"	£76.00
1494497	Rigidiser W18" x D16"	£39.00
1494498	Rigidiser W18" x D18"	£39.00
1494499	Rigidiser W18" x D20"	£76.00
1494500	Rigidiser W20" x D16"	£39.00
1494501	Rigidiser W20" x D18"	£39.00
1494502	Rigidiser W20" x D20"	£76.00
1494503	Rigidiser W22" x D18"	£76.00
1494504	Rigidiser W22" x D20"	£76.00
Office use only: Process cushion on configurator and accessories as line items.		