



Tel: 01656 776222  
 Email: ordersuk@invacare.com  
 www.invacare.co.uk


Matrx Libra cushion  
 June 2023 V2

<b>Customer name:</b>	<b>Account No :</b>
<b>Contact no:</b>	<b>Order No :</b>
<b>Delivery address:</b>	<b>Order Qty :</b>
	<b>Order Date :</b>

**Matrx Libra**  
 Retail order Form  
 Maximum User Weight 21 stone (136kg) / HD 36 stone (227kg)

		WIDTH											
		36	38	41	43	46	48	51	53	56			
		inch	14	15	16	17	18	19	20	21	22		
DEPTH	30	12											
	36	14											
	41	16											
	43	17										HD SIZE*	HD SIZE*
	46	18										HD SIZE*	HD SIZE*
	51	20										HD SIZE*	HD SIZE*
56	22										HD SIZE*	HD SIZE*	

\* PART NUMBERS AND PRICE LISTED UNDER HEAVY DUTY SIZES



**Ordering the correct size**

Please note the complete product dimensions.  
 The example shown is of a 16"x16" cushion part number LC1616.  
 Please allow an additional 1 1/2" to the required seat depth for the rear cushion radius.

The rear cushion radius is a curved extension designed to maximise surface area and minimise the gap between seat and backrest; while stabilising fluid in tilted positions, with coccyx cut out for added relief.

**Features:**  
 Very High Risk  
 3" height  
 3 Year Warranty  
 Anatomically contoured foam base Innovative dual layer fluid sac  
 Choice of 3 covers  
 Positioning accessories available

Code	X	Description	Price
<b>Libra Cushion</b>			
LC _____ (Width x Depth)		Matrx Libra standard sizes with black smooth cover	£410.00
LCW _____ (Width x Depth)		Matrx Libra standard sizes with black wipedown cover	£410.00
LCS _____ (Width x Depth)		Matrx Libra standard sizes with black spacer cover	£410.00
<b>Heavy Duty sizes</b>			
LC2117		Heavy Duty Libra 21x17" (53x43cm)	Black Smooth cover £563.00
LC2118		Heavy Duty Libra 21x18" (53x46cm)	Black Smooth cover £563.00
LC2120		Heavy Duty Libra 21x20" (53x51cm)	Black Smooth cover £563.00
LC2122		Heavy Duty Libra 21x22" (53x56cm)	Black Smooth cover £563.00
LC2217		Heavy Duty Libra 22x17" (56x43cm)	Black Smooth cover £563.00
LC2218		Heavy Duty Libra 22x18" (56x46cm)	Black Smooth cover £563.00
LC2220		Heavy Duty Libra 22x20" (56x51cm)	Black Smooth cover £563.00
LC2222		Heavy Duty Libra 22x22" (56x56cm)	Black Smooth cover £563.00
LCW2117		Heavy Duty Libra 21x17" (53x43cm)	Black Wipedown Cover £563.00
LCW2118		Heavy Duty Libra 21x18" (53x46cm)	Black Wipedown Cover £563.00
LCW2120		Heavy Duty Libra 21x20" (53x51cm)	Black Wipedown Cover £563.00
LCW2122		Heavy Duty Libra 21x22" (53x56cm)	Black Wipedown Cover £563.00
LCW2217		Heavy Duty Libra 22x17" (56x43cm)	Black Wipedown Cover £563.00
LCW2218		Heavy Duty Libra 22x18" (56x46cm)	Black Wipedown Cover £563.00
LCW2220		Heavy Duty Libra 22x20" (56x51cm)	Black Wipedown Cover £563.00
LCW2222		Heavy Duty Libra 22x22" (56x56cm)	Black Wipedown Cover £563.00
LCS2117		Heavy Duty Libra 21x17" (53x43cm)	Spacer Fabric Cover £563.00
LCS2118		Heavy Duty Libra 21x18" (53x46cm)	Spacer Fabric Cover £563.00
LCS2120		Heavy Duty Libra 21x20" (53x51cm)	Spacer Fabric Cover £563.00
LCS2122		Heavy Duty Libra 21x22" (53x56cm)	Spacer Fabric Cover £563.00
LCS2217		Heavy Duty Libra 22x17" (56x43cm)	Spacer Fabric Cover £563.00
LCS2218		Heavy Duty Libra 22x18" (56x46cm)	Spacer Fabric Cover £563.00
LCS2220		Heavy Duty Libra 22x20" (56x51cm)	Spacer Fabric Cover £563.00
LCS2222		Heavy Duty Libra 22x22" (56x56cm)	Spacer Fabric Cover £563.00



**Accessories**



LABD	<b>Matrix Abductor</b> For comfortable mid-line positioning of thighs and hips (Provides 1 1/8" (29mm) increased height abductor contour)	<b>£39.00</b>
1487157	<b>Anterior Seat Wedge</b> Increases slope, prevents sliding and reduces shear (1" Anterior wedge increases pre-ischial contour - Underside mounted)	<b>£27.00</b>
1487158	<b>Tapered Wedges (pair)</b> For obliquity and pelvic rotation when used in rear (3/4" Tapered wedge - Underside mounted)	<b>£20.00</b>
1487155	<b>Lateral Wedge</b> Reinforces lateral stability for enhanced midline positioning (1" Lateral wedge increases hip/thigh loading, - Underside mounted)	<b>£30.00</b>
1578357	<b>Libra Fluid Supplement</b> For targeted protection/obliquity accommodation (May be used to accommodate mild pelvic obliquity/additional pressure reduction)	<b>£54.00</b>
1589366	<b>Libra Accessory Kit</b> Comprises of all the above items in an assessment presentation box with instruction guide. Parts included: LABD, 1487157, 1487158, 1487155, 1578357	<b>£120.00</b>
	<b>Spare Platilon Inner Liner</b> Provides added protection to the foam cushion from infection control issues.	<b>£13.00</b>
	<b>Additional Cover (Widths 14" to 20")</b> <input checked="" type="checkbox"/> Black Smooth <input type="checkbox"/> Black Spacer <input type="checkbox"/> Black Wipedown	<b>£96.00</b>
	<b>Additional Cover (Widths 21" and over )</b> <input checked="" type="checkbox"/> Black Smooth <input type="checkbox"/> Black Spacer <input type="checkbox"/> Black Wipedown	<b>£112.00</b>
1494494	14" x 16" Rigidiser	<b>£39.00</b>
1487165	16" x 16" Rigidiser	<b>£39.00</b>
1494495	16" x 18" Rigidiser	<b>£39.00</b>
1494496	16" x 20" Rigidiser	<b>£76.00</b>
1494497	18" x 16" Rigidiser	<b>£39.00</b>
1494498	18" x 18" Rigidiser	<b>£39.00</b>
1494499	18" x 20" Rigidiser	<b>£76.00</b>
1494500	20" x 16" Rigidiser	<b>£39.00</b>
1494501	20" x 18" Rigidiser	<b>£39.00</b>
1494502	20" x 20" Rigidiser	<b>£76.00</b>
1494503	22" x 18" Rigidiser	<b>£76.00</b>
1494504	22" x 20" Rigidiser	<b>£76.00</b>
<b>TOTAL PRICE:</b>		

**Special Instructions - Extended lead times will apply**

**Name:** \_\_\_\_\_ **Email:** \_\_\_\_\_ **Phone:** \_\_\_\_\_

**Description of special request (please include relevant option code):**