



Yes, you can.®

Tel: 01656 776222
Email: ordersuk@invacare.com
www.invacare.co.uk

RETAIL PRESCRIPTION FORM

Matrx Posture back Deep
June 2023

Account Name :	Account No :
Address :	Order No :
Tel No:	Order Qty :
Delivery address :	Order Date :

Matrx Posture Back Deep

Retail order form
Maximum User Weight: 227 kg

Features

5" Contour Depth
40° angle adjustment
4" depth adjustment
20° rotation
136kg (21 stone) user weight
Hardware not included



Deep MatRX Posture Back Backrest

Chest Width	Backrest Width	31cm	41cm	51cm	Fits W/ch Widths		With Extended Mounting Pins EMP*
		(12")	(16")	(20")			
≤26cm (10")	36cm (14")	PBD1412NH	PBD1416NH		33cm - 38cm	(13" - 15")	41cm (16")
≤31cm (12")	41cm (16")	PBD1612NH	PBD1616NH	PBD1620NH	38cm - 43cm	(15" - 17")	46cm (18")
≤33cm (13")	43cm (17")	PBD1712NH	PBD1716NH	PBD1720NH	41cm - 46cm	(16" - 18")	48cm (19")
≤36cm (14")	46cm (18")	PBD1812NH	PBD1816NH	PBD1820NH	43cm - 48cm	(17" - 19")	51cm (20")
≤41cm (16")	51cm (20")	PBD2012NH	PBD2016NH	PBD2020NH	48cm - 53cm	(19" - 21")	56cm (22")

*Not suitable for HD users (max. 136kg/300lb weight capacity)

Code	X	Description	Fits wheelchairs		Price
PB Deep			Width	Height	
1545001		Posture Back Deep PBD1412NH	330mm-380mm	310mm	£410.00
1545002		Posture Back Deep PBD1416NH	330mm-380mm	410mm	£410.00
1545003		Posture Back Deep PBD1612NH	380mm-430mm	310mm	£410.00
1545058		Posture Back Deep PBD1616NH	380mm-430mm	410mm	£410.00
1545059		Posture Back Deep PBD1620NH	380mm-430mm	510mm	£410.00
1545060		Posture Back Deep PBD1712NH	410mm-460mm	310mm	£410.00
1545061		Posture Back Deep PBD1716NH	410mm-460mm	410mm	£410.00
1545062		Posture Back Deep PBD1720NH	410mm-460mm	510mm	£410.00
1545063		Posture Back Deep PBD1812NH	430mm-480mm	310mm	£410.00
1545064		Posture Back Deep PBD1816NH	430mm-480mm	410mm	£410.00
1545065		Posture Back Deep PBD1820NH	430mm-480mm	510mm	£410.00
1545066		Posture Back Deep PBD2012NH	480mm-530mm	310mm	£410.00
1545067		Posture Back Deep PBD2016NH	480mm-530mm	410mm	£410.00
1545068		Posture Back Deep PBD2020NH	480mm-530mm	510mm	£410.00
Hardware					
1545275		Easy Set mounting hardware (pair) PBMH	Fits 19mm, 22mm and 25mm back canes		£182.00
1545276		Easy Set mounting hardware (pair) PBMHQ	Fits 29mm back canes		£182.00
1545277		Fixed Space Saver Mounting Hardware FSSH			£200.00
1545278		Split Mounting hardware SMH			£182.00
1544776		Heavy Duty Easy Set Mounting Hardware HDMH	Fits 19mm, 22mm and 25mm back canes		£260.00

Order No :

Customer Ref :

Posture back Deep
1 of 2

Subject to Invacare's standard Terms and Conditions.
Invacare reserve the right to alter product specification without prior notice.

