



Yes, you can.®

Tel: 01656 776222
 Email: ordersuk@invacare.com
 www.invacare.co.uk

Matrx MX2 Back
 June 2023

Account Name :	Account No :
Address :	Order No :
Tel No:	Order Qty :
Delivery address :	Order Date :

Matrx MX2
 Retail order form
 Max.user weight: 113kg / 17 stone



		WIDTH					
		355	380	405	430	455	
		inch	14	15	16	17	18
HEIGHT	mm						
	150	6					
	230	9					
	310	12					
410	16						



Standard Features:
 Lightweight carbon-fibre back shell
 Anatomical contouring with 76 mm depth for enhanced lateral stability
 Angle and depth adjustable
 Cutouts for breathability to help manage microclimate
 Slotted design to allow height adjustments even when in the chair
 Breathable Spacer cover
 Backrests 230mm or higher come with leatherette side panels
 Oversized foam for enhanced comfort and protection

Code	X	Description	Price
MX2 Back shell			
1629002		MX2 Back (35.5w x 15h)	14x6" £519.00
1629003		MX2 Back (35.5w x 23h)	14x9" £519.00
1629004		MX2 Back (35.5w x 31h)	14x12" £519.00
1629005		MX2 Back (35.5w x 41h)	14x16" £519.00
1629006		MX2 Back (38w x 15h)	15x6" £519.00
1629007		MX2 Back (38w x 23h)	15x9" £519.00
1629008		MX2 Back (38w x 31h)	15x12" £519.00
1629009		MX2 Back (38w x 41h)	15x16" £519.00
1629010		MX2 Back (40.5w x 15h)	16x6" £519.00
1629011		MX2 Back (40.5w x 23h)	16x9" £519.00
1629012		MX2 Back (40.5w x 31h)	16x12" £519.00
1629013		MX2 Back (40.5w x 41h)	16x16" £519.00
1629014		MX2 Back (43w x 15h)	17x6" £519.00
1629015		MX2 Back (43w x 23h)	17x9" £519.00
1629016		MX2 Back (43w x 31h)	17x12" £519.00
1629017		MX2 Back (43w x 41h)	17x16" £519.00
1629018		MX2 Back (45.5w x 15h)	18x6" £519.00
1629019		MX2 Back (45.5w x 23h)	18x9" £519.00
1629020		MX2 Back (45.5w x 31h)	18x12" £519.00
1629027		MX2 Back (45.5w x 41h)	18x16" £519.00

Hardware

Hardware is width-adjustable allowing the backrest to be fitted to wheelchairs that are up to 25mm narrower or wider than the actual backrest width for a precise fit and positioning in rigid or foldable wheelchairs with backrest tube diameter of 19 mm, 22 mm and 25 mm.



Fixed hardware



Mini set hardware

1629027	Fixed carbon fibre hardware	£153.00
1629028	Mini set (quick release) hardware	£206.00

Privacy flap

1551799	Extra Small	Fits 35.5cm width back	£49.00
1551800	Small	Fits 40.5cm width back	£49.00
1551903	Medium	Fits 43 and 45.5cm width back	£49.00

Special Instructions - Extended lead times will apply



Name:

Email:

Phone:

Description of special request (please include relevant option code):